
# Avance Anesthesia Machine Technical Reference Manual



Datex-Ohmeda products have unit serial numbers with coded logic which indicates a product group code, the year of manufacture, and a sequential unit number for identification. The serial number can be in one of two formats.

AAA <b>X</b> 11111	AAA <b>XX</b> 111111AA
The <b>X</b> represents an alpha character indicating the year the product was manufactured; <b>H</b> = 2004, <b>J</b> = 2005, etc. <b>I</b> and <b>0</b> are not used.	The <b>XX</b> represents a number indicating the year the product was manufactured; <b>04</b> = 2004, <b>05</b> = 2005, etc.

**Avance** is a registered trademark of Datex-Ohmeda, Inc.

Other brand names or product names used in this manual are trademarks or registered trademarks of their respective holders.

### **Avance Anesthesia Machine**

Machines with the original display unit (DU) with System Software 3.2 or earlier

Machines with the high performance display unit (HPDU) with System Software 4.X or greater

This document is not to be reproduced in any manner, nor are the contents to be disclosed to anyone, without the express authorization of the product service department, Datex-Ohmeda, Ohmeda Drive, PO Box 7550, Madison, Wisconsin, 53707.

1009-0357-000 09/07 i

 $<sup>^{\</sup>hbox{\scriptsize (C)}}$  2007 Datex-Ohmeda Inc.

# **Important**

The information contained in this Technical Reference manual pertains only to those models of products which are marketed by Datex-Ohmeda as of the effective date of this manual or the latest revision thereof. This Technical Reference manual was prepared for exclusive use by Datex-Ohmeda service personnel in light of their training and experience as well as the availability to them of parts, proper tools and test equipment. Consequently, Datex-Ohmeda provides this Technical Reference manual to its customers purely as a business convenience and for the customer's general information only without warranty of the results with respect to any application of such information. Furthermore, because of the wide variety of circumstances under which maintenance and repair activities may be performed and the unique nature of each individual's own experience, capacity, and qualifications, the fact that customer has received such information from Datex-Ohmeda does not imply in anyway that Datex-Ohmeda deems said individual to be qualified to perform any such maintenance or repair service. Moreover, it should not be assumed that every acceptable test and safety procedure or method, precaution, tool, equipment or device is referred to within, or that abnormal or unusual circumstances, may not warrant or suggest different or additional procedures or requirements.

This manual is subject to periodic review, update and revision. Customers are cautioned to obtain and consult the latest revision before undertaking any service of the equipment. Comments and suggestions on this manual are invited from our customers. Send your comments and suggestions to the Manager of Technical Communications, Datex-Ohmeda, Ohmeda Drive, PO Box 7550, Madison, Wisconsin 53707.

### **⚠** CAUTION

Servicing of this product in accordance with this Technical Reference manual should never be undertaken in the absence of proper tools, test equipment and the most recent revision to this service manual which is clearly and thoroughly understood.

# **Technical Competence**

The procedures described in this Technical Reference manual should be performed by trained and authorized personnel only. Maintenance should only be undertaken by competent individuals who have a general knowledge of and experience with devices of this nature. No repairs should ever be undertaken or attempted by anyone not having such qualifications.

Datex-Ohmeda strongly recommends using only genuine replacement parts, manufactured or sold by Datex-Ohmeda for all repair parts replacements.

Read completely through each step in every procedure before starting the procedure; any exceptions may result in a failure to properly and safely complete the attempted procedure.

ii 09/07 1009-0357-000

# **Table of Contents**

	Important	ii
	Technical Competence	il
Introduction		
	1.1 What this manual includes	1-2
	1.2 User's Reference manuals	1-2
	1.3 What is an Avance anesthesia machine?	1-3
	1.4 Anesthesia system components	1-4
	1.5 Breathing system components	1-6
	1.5.1 Optional ABS components	
	1.6 Display controls	1-8
	1.7 Anesthesia system display	1-9
	1.7.1 Using menus	1-11
	1.8 Symbols used in the manual or on the equipment	1-12

1

1009-0357-000 09/07 iii

### 2 Theory of Operation

2.1 Electrical system	2-2
2.2 Power subsystem (original)	2-4
2.2.1 Power Controller board	2-5
2.3 Power subsystem (with U-Frame power supply)	2-6
2.3.1 Power Controller board	2-7
2.3.2 Power distribution	
2.4 System communications	2-10
2.4.1 Software Power On Self Tests (POST)	2-12
2.5 Display Unit	
2.5.1 Software requirements	2-14
2.6 System connections	
2.6.1 Display Unit	2-15
2.6.2 High Performance Display Unit	
2.6.3 Display Connector board	
2.7 Power Controller and Anesthesia Control board connections	
2.8 Anesthesia Control board	2-18
2.9 Electronic Gas Mixer	
2.10 Ventilator Interface board	
2.11 Gas flow through the anesthesia machine	
2.11.1 Overview	
2.11.2 Physical connections	
2.11.3 Suction regulators	
2.12 Flow through the breathing system	
2.12.1 Overview of flow paths	
2.12.2 Manual ventilation	
2.12.3 Mechanical ventilation	
2.12.4 Fresh gas and O <sub>2</sub> flush flow (with SCGO)	2-41
$2.12.5$ Fresh gas and $0_2$ flush flow (with ACGO)	2-43
2.13 Ventilator mechanical subsystems	2-45
2.13.1 Drive gas filter and Gas Inlet Valve	2-45
2.13.2 Pressure regulator	2-46
2.13.3 Flow control valve	2-46
2.13.4 Drive Gas Check Valve (DGCV)	2-47
2.13.5 Bellows Pressure Relief Valve	
2.13.6 Exhalation valve	
2.13.7 Mechanical Overpressure Valve	
2.13.8 Reservoir and bleed resistor	
2.13.9 Free breathing valve	
2.13.10 Breathing circuit flow sensors	2-51

iv 09/07 1009-0357-000

### **3 Checkout Procedure**

3.1 Inspect the system	
3.2 System checkout (for System software 3.X or greater)	3-2
3.2.1 Leak < 250 ml	
3.2.2 Machine Check	
3.2.3 Machine Check - System (Ventilator Circuit Testing)	
3.2.4 Machine Check - Circuit (Bag Circuit Testing)	
3.2.5 Machine Check - Circuit O2	3-4
3.3 Individual Checks (for System software 3.X or greater)	3-4
3.3.1 System	
3.3.2 Circuit	
3.3.3 Circuit 02 Cell	
3.3.4 Low P Leak	
3.3.5 Low P Leak (machines with ACGO)	
3.4 System "All checks" (for System software 2.X)	
3.4.1 Low P leak check	
3.4.2 Quick check	
3.4.3 Vent check	
3.4.4 Circuit O <sub>2</sub> cell check	
3.5 Bellows drop test	
3.6 Backlight test	
3.7 Vaporizer back pressure test	
3.8 Pipeline and cylinder tests	
3.8.1 O <sub>2</sub> supply alarm test	
3.9 Pressure relief tests	3-10
3.10 Flush Flow Test	3-11
3.11 Alarm tests	3-12
3.12 Alternate O2 flowmeter tests	3-13
3.13 Auxiliary O2 flowmeter tests	3-13
3.14 Integrated Suction Regulator tests	3-13
3.15 Power failure test	3-14
3.16 Electrical safety tests	3-14

1009-0357-000 09/07 v

### 4a Installation and Service Menus (DU)

4a.1	Service and Installation menu structure	.4a-2
4a.2	Install/Service Menu (Super User)	.4a-3
	4a.2.1 Colors Menu	.4a-4
	4a.2.2 Units Menu	.4a-4
	4a.2.3 Factory Defaults	.4a-5
4a.3	Installation Menu	.4a-6
	4a.3.1 Configuration	.4a-7
	4a.3.2 Units Menu	.4a-8
	4a.3.3 Options Key	.4a-8
	4a.3.4 Copy Configuration	.4a-9
4a.4	Service Menu	la-10
	4a.4.1 Software/Hardware Ver Menu	la-11
	4a.4.2 Service Log Menu	la-12
4a.5	Calibration	la-13
	4a.5.1 User Calibration menu	la-13
	4a.5.2 Manifold P Span	la-14
	4a.5.3 Insp Flow Zero	la-15
	4a.5.4 Inspiratory Flow Valve	la-16
	4a.5.5 Bleed Resistor	la-18
	4a.5.6 Paw Span	la-19
	4a.5.7 Zero Gas Xducrs	la-20
	4a.5.8 Cal Config	la-21
	4a.5.9 Mixer P Zero	la-21

vi 09/07 1009-0357-000

# 4b Install/Service Menus (HPDU)

**5 Calibration** 

4b.1 Service and Installation menu structure
4b.2 Install/Service Menu (Super User)
4b.2.1 Install/Service - Page 1       4b-3         4b.2.2 Install/Service - Page 2       4b-9         4b.3 Installation Menu       4b-10
4b.3.1 Configuration
4b.3.2 Configuration Units
4b.3.3 Options Key
4b.3.4 Copy Configuration
4b.4 Service Menu
4b.4.1 Software/Hardware Ver Menu
4b.4.2 Service Log Menu
4b.5 Calibration
4b.5.1 Spiro Calibration
4b.5.2 User Calibration menu
4b.5.3 Manifold P Span
4b.5.4 Inspiratory Flow Valve       4b-21         4b.5.5 Insp Flow Zero       4b-22
4b.5.6 Bleed Resistor
4b.5.7 Paw Span
4b.5.8 Zero Gas Xducrs
4b.5.9 Cal Config
4b.5.10 Mixer P Zero
5.1 Primary Regulators5-2
5.1.1 Test setup5-3
5.1.2 Testing Primary Regulators5-3
5.1.3 Adjusting Primary Regulators
5.2 O <sub>2</sub> Flush Regulator5-9
5.3 Adjust Drive Gas Regulator
5.4 Ventilator Calibrations
5.4.1 Cal Config 5-12
5.4.2 Manifold P Span
5.4.3 Inspiratory Flow Valve Cal 5-14
5.4.4 Insp Flow Zero
5.4.5 Bleed Resistor Cal
5.4.6 Paw Span

1009-0357-000 09/07 vii

### **6 Installation and Maintenance**

7 Troubleshooting

6.1 Avance Installation Checklist6-2
6.2 Avance Planned Maintenance6-4
6.3 Free breathing valve maintenance
6.4 MOPV pressure relief valve test6-7
6.4.1 Test setup6-7
6.4.2 Test procedure6-7
6.5 Pressure Limit Circuit test6-8
6.6 Mixer test 6-10
6.6.1 Mixer outlet check valve leak test
6.6.2 Mixer flow verification
6.7 Alternate O2 flowmeter tests 6-11
6.8 Auxiliary O2 flowmeter tests
6.9 Integrated Suction Regulator tests
6.10 Battery capacity test         6-15
7.1 Troubleshooting Guidelines
7.2 Troubleshooting high pressure and low pressure leaks $\dots$
7.3 Troubleshooting Startup Screen (POST) messages — for DU
7.4 Troubleshooting Startup Screen (POST) messages — for HPDU7-5
7.5 Troubleshooting the HPDU Display7-6
7.6 Troubleshooting System Malfunction (safe-state) screen7-7
7.7 Breathing System Leak Test Guide7-8
7.7.1 Check Valves
2.7.2 Breathing System Troubleshooting Flowcharts
7.7.3 Leak Isolation Tests
7.8 System Troubleshooting Flowcharts
7.9 System Malfunction and Alt O2 Flowchart Table $\hdots$ 7-42
7.10 Technical Alarms
$7.11 \ \text{Steps and Messages displayed during the System Checkout} - \\ \text{(for System software 3.X or greater)} \ \dots \ 7-71$
$7.12~{ m Steps}$ and Messages displayed during the System Checkout $-$
(for System software 2.X)
7.12.1 Steps for the Quick Check
7.12.2 Steps for the Vent Check

viii 09/07 1009-0357-000

# 8a Service Diagnostics and Software Download (DU)

8a.1 Avance Service Application	8a-2
8a.1.1 Main Menu and System Information	8a-2
8a.1.2 Power Diagnostics	8a-3
8a.1.3 Power Controller Power Diagnostics	8a-4
8a.1.4 Anesthesia Control Board Power Diagnostics	8a-6
8a.1.5 Electronic Mixer Power Diagnostics	8a-8
8a.1.6 Ventilator Interface Board Power Diagnostics	8a-9
8a.1.7 Display Unit Power Diagnostics	8a-10
8a.2 Gas Diagnostics	8a-11
8a.2.1 Gas Supplies	8a-12
8a.2.2 Mixer Output	8a-13
8a.2.3 Mixer Tests and Pressure	8a-14
8a.2.4 Mixer Temperature	8a-15
8a.2.5 Setting Gas Flow	8a-16
8a.2.6 Breathing System Leak Test	8a-17
8a.3 Ventilation Diagnostics	8a-18
8a.3.1 Status	8a-19
8a.3.2 Vent Flow and Pressure	8a-20
8a.4 Display Diagnostics	8a-21
8a.5 Special Functions	8a-22
8a.5.1 Mixer Service Menu	8a-23
8a.5.2 View Revision Log	8a-24
8a.5.3 View PC Card Install Log	8a-24
8a.6 Software Download	8a-25

1009-0357-000 09/07 ix

# 8b Software Download and Special Functions (HPDU)

8b.1	Overview	. 8b-2
	8b.1.1 Main Menu and System Information	. 8b-2
8b.2	Software Download	. 8b-3
8b.3	Special Functions	. 8b-5
	8b.3.1 Power Diagnostics	. 8b-5
	8b.3.2 Power Controller Power Diagnostics	. 8b-6
	8b.3.3 Display Unit Power Diagnostics	. 8b-8
	8b.3.4 Anesthesia Control Board Power Diagnostics	. 8b-9
	8b.3.5 Electronic Mixer Power Diagnostics	8b-11
	8b.3.6 Ventilator Interface Board Power Diagnostics	8b-12
8b.4	Ventilation Diagnostics	8b-13
	8b.4.1 Vent Status	8b-14
	8b.4.2 Vent Flow and Pressure	8b-15
8b.5	Gas Diagnostics	8b-16
	8b.5.1 Setting Gas Flow	8b-17
	8b.5.2 Gas Supplies	8b-18
	8b.5.3 Mixer Output	8b-19
	8b.5.4 Mixer Tests and Pressure	8b-20
	8b.5.5 Mixer Temperature	8b-21
	8b.5.6 Breathing System Leak Test	8b-22
8b.6	Mixer Service Menu	8b-23
8b.7	Display Diagnostics	8b-25
	8b.7.1 Test Keys and Battery	8b-26
8b.8	Compatibility Table	8b-27
	8b.8.1 System Download Log	8b-28
	8b.8.2 CF Card Install Log	8b-28
	8b.8.3 View Install Errors	8b-29

x 09/07 1009-0357-000

### **9 Repair Procedures**

9.1 Circuit board replacement precautions	9-4
9.2 How to bleed gas pressure from the machine	9-5
9.3 How to remove the rear panels	9-6
9.3.1 To remove the rear upper panel	9-6
9.3.2 To remove the lower access panels	9-6
9.4 How to remove the tabletop	9-7
9.5 Servicing the Display Unit (DU)	9-8
9.5.1 Remove the Display Unit	9-8
9.5.2 Disassemble the Display Unit	9-9
9.5.3 To replace the CPU board	9-10
9.5.4 To replace the LCD display	9-11
9.5.5 To replace the backlights	9-13
9.5.6 To replace the Inverters	9-14
9.5.7 To replace the front enclosure or components	9-15
9.6 Servicing the High Performance Display Unit (HPDU)	9-17
9.6.1 Remove the Display Unit	9-17
9.6.2 Disassemble the Display Unit	9-18
9.6.3 CPU Fan	9-18
9.6.4 To replace the CPU board	9-19
9.6.5 To replace the LCD display	9-20
9.6.6 To replace the backlights	9-21
9.6.7 To replace the Inverters	9-22
9.6.8 To replace the front enclosure or components	9-23
9.7 Replacing the Display and MGAS cables	9-25
9.7.1 Remove the MGAS oxygen partition	9-25
9.8 Servicing the lower electrical enclosure components	9-26
9.8.1 Power Controller board (original)	9-26
9.8.2 Power Controller board and Universal Power Supply	9-27
9.8.3 Anesthesia Control board	9-28
9.8.4 Backup batteries	9-29
9.8.5 Fan	9-30
9.8.6 Display Connector board	9-30
9.9 Servicing the pan electrical enclosure components	9-31
9.9.1 Electronic Gas Mixer assembly	9-31
9.9.2 Ventilator Interface board	9-32
9.9.3 Filter board	9-33
9.9.4 Pan Connector board	9-34
9.9.5 Pan enclosure fan	9-34

1009-0357-000 09/07 xi

9.10 Servicing the Vent Engine	9-35
9.10.1 To remove the Vent Engine	9-36
9.10.2 Replacing Vent Engine components	9-37
9.10.3 Replacing GIV components	9-38
9.11 Servicing the pipeline inlet manifold components	9-39
9.11.1 Replace pipeline inlet filter	9-39
9.11.2 Replace pipeline inlet check valve	9-39
9.11.3 Replace the inlet manifold	
9.12 Service the cylinder supply modules	9-41
9.12.1 Replace primary regulator module (complete replacement)	9-41
9.12.2 Replace cylinder inlet filter	9-42
9.12.3 Replace cylinder check valve	
9.13 Replace gas-supply pressure transducers	9-43
9.14 Service vaporizer manifold parts	9-44
9.14.1 Repair manifold port valve	9-44
9.14.2 Checkout procedure for manifold port valve	9-45
9.14.3 Replace vaporizer manifold check valve	9-46
9.14.4 Replace vaporizer pressure relief valve	
9.14.5 Replace vaporizer manifold	9-49
4.15 Replace ACGO selector switch	9-50
9.16 Clean or replace ACGO port flapper valve	9-52
9.17 Replace the APL valve	9-53
9.18 Replace the bag support arm	9-54
9.18.1 Servicing the bag support arm	9-55
9.18.2 Replace friction pad in lower bag arm assembly	9-56
9.18.3 Replace bag port housing	9-57
9.19 Replace system switch assembly	9-58
9.20 Replace Alt O2 components	9-60
9.21 Replace auxiliary O <sub>2</sub> flowmeter	9-61
9.22 Replace the suction regulator	9-62
9.23 Replace task light components	9-63
9.23.1 To replace the task light switch	9-63
9.23.2 To replace the upper task light	9-63
9.23.3 To replace the lower task light	9-64
9.24 Replace ABS breathing system components	9-65
9.24.1 Replace Bag/Vent switch assembly	9-65
9.24.2 Replace bellows base latch assembly	9-66
9.24.3 EZchange Canister spring replacement	9-67
9.25 Replace casters	9-68
9.26 Reconfigure sample gas return line	9-69
9.27 Change drive gas	9-70

xii 09/07 1009-0357-000

### **10 Illustrated Parts**

10.1 Service tools
10.1.1 Software tools
10.1.2 Manifold pressure test adapter
10.1.3 Test Devices
10.1.4 Lubricants and Adhesives
10.1.5 Test Tools
10.2 External components - front view
10.3 External components - front view references
10.4 External Components - rear view
10.5 AC Power cords and AC Inlet filter
10.6 AC Inlet/Outlet Components
10.7 Display Unit (DU)
10.8 High Performance Display Unit (HPDU)
10.9 Lower electronic enclosure components
10.9.1 Anesthesia Control and Display Connector board
10.10 Pan electronic enclosure components
10.11 Electronic Gas Mixer
10.12 Pipeline inlet fittings
10.13 Cylinder Gas Supplies
10.13.1 Cylinder inlet fittings
10.14 Vaporizer manifold
10.15 Vent Engine Housing
10.16 Vent Engine
10.16.1 Vent Engine - under side
10.17 ABS to machine Interface Components (SCGO)
10.18 ABS to machine Interface Components (ACGO)
10.19 Flush Regulator and Flush Valve
10.20 Front panel, Alt 02, and system switch
10.21 Breathing system interface

1009-0357-000 09/07 xiii

10.22 Breathing System	-32
10.22.1 APL Valve	-32
10.22.2 Bag/Vent Switch	-33
10.22.3 Absorber canister	-34
10.22.4 Flow Sensor Module	-35
10.22.5 Breathing Circuit Module	-36
10.22.6 Exhalation valve	
10.22.7 Bellows	
10.22.8 Bellow base	
10.22.9 Bag Arms	
10.22.10 EZchange Canister system (CO <sub>2</sub> Bypass)	
10.22.11 Condenser	
10.22.12 Anesthetic Gas Scavenging System – AGSS	
10.22.13 Passive AGSS	
10.22.14 Adjustable AGSS       10-         10.22.15 Active AGSS       10-	
10.22.16 AGSS gauge, and sample return	
10.22.17 Airway module (MGAS) components	
10.23 Integrated Suction Regulator	
10.23.1 Major Components (Continuous and Venturi suction)	
10.23.3 Venturi assembly	
10.24 Auxiliary O <sub>2</sub> Flowmeter	
<del>-</del>	
10.25 Rear panel components	
10.26 Tabletop components	
10.27 Right-side Components	-58
10.28 External components - lower assembly	-59
10.29 Drawer	-60
10.30 Legris quick-release fittings	-61
10.31 Vent Drive and low-pressure tubing	-62
10.32 Tubing for use with Legris fittings	-64
10.33 Cables and harnesses (power supply)	-66
10.34 Cables and harnesses in lower electronic enclosure	-68
10.34.1 Machines with original Power Controller board	-68
10.34.2 Machines with Universal Power Supply (U-Frame)	
10.35 Cables and harnesses (Filter Board interface)	
10.36 Cables and harnesses in Pan enclosure	
10.37 Optional Monitor Display mounts	
10.38 Display arm mounting kits for optional equipment	
20.00 2.0p.0, ann mounting the for optional oquipmont in this interest in 10	

xiv 09/07 1009-0357-000

# **11 Schematics and Diagrams**

### **12 Service Application**

12.1 Avance Service Application (PC based)	12-2
12.1.1 PC Requirements	12-2
12.2 Startup screen — System Status	12-3
12.3 System Schematics	12-4
12.3.1 Power Schematic	12-4
12.3.2 Gas Delivery Schematic	12-5
12.3.3 Vent Schematic	12-6
12.4 Menu Items	12-7
12.5 File menu	12-8
12.5.1 File — Preferences	12-8
12.6 Tools menu	12-9
12.6.1 Tools — Communication Status	12-9
12.6.2 Tools — System Calibrations	12-10
12.6.3 Tools — Transfer Logs	12-11
12.7 Power Diagnostics menu	12-12
12.7.1 Power Diagnostics — Power Board	12-12
12.7.2 Power Diagnostics — Anesthesia Control Board Power	12-13
12.7.3 Power Diagnostics — Mixer Board Power	12-14
12.7.4 Power Diagnostics — Vent Interface Board Power	12-15
12.7.5 Power Diagnostics — Display Unit Power	12-16
12.8 Gas Delivery Subsystem menu	12-17
12.8.1 Gas Delivery Subsystem — Gas Supply Status	12-17
12.8.2 Gas Delivery Subsystem — Mixer Output	12-18
12.8.3 Gas Delivery Subsystem — Mixer Pressure and Temperature	12-19
12.8.4 Gas Delivery Subsystem — Gas Delivery Status	12-20
12.8.5 Gas Delivery Subsystem — Mixer Post/Checkout Test Results	12-21
12.8.6 Gas Delivery Subsystem — Perform Mixer Checkout Tests	12-22
12.9 Vent Subsystem menu	12-24
12.9.1 Vent Subsystem — Vent Status	12-24
12.9.2 Vent Subsystem — Vent Flow and Pressure	12-25
12.10 Window menu	12-26
12 11 Help menu	12-26

1009-0357-000 09/07 xv

Notes

xvi 09/07 1009-0357-000

# 1 Introduction

# In this sectionThis section provides a general overview of the Avance anesthesia machine.1.1 What this manual includes1-21.2 User's Reference manuals1-21.3 What is an Avance anesthesia machine?1-31.4 Anesthesia system components1-41.5 Breathing system components1-61.5.1 Optional ABS components1-71.6 Display controls1-81.7 Anesthesia system display1-91.7.1 Using menus1-11

1009-0357-000 09/07 1-1

### 1.1 What this manual includes

This manual covers the service information for the S/5 Avance line of anesthesia machines. It covers the following components:

- Display Unit
- Integral electronics
- Gas delivery components
- Breathing system components
- Frame component
- Optional suction regulator
- Optional auxiliary O<sub>2</sub> flowmeter

### Other equipment

Other equipment may be attached to the system on a display mount, the top shelf, or on the side dovetail rails. Consult separate documentation relative to these items for details.

# 1.2 User's Reference manuals

Some sections of this manual refer you to the User's Reference manual for the S/5 Avance. To expedite repairs, you must have, and be familiar with, the User's Reference manuals for this product.

Refer to the S/5 Avance User's Reference manual if you need further information about the operation of the system.

1-2 09/07 1009-0357-000

### 1.3 What is an Avance anesthesia machine?

The Avance anesthesia system combines advanced anesthesia delivery, patient monitoring, and care information management in a contemporary, compact design. The system features electronic gas mixing for up to three gases.

Optional, integrated features include auxiliary 0<sub>2</sub>, suction control, and gas monitoring (E- and M-Series respiratory gas module capable).

The Avance system uses SmartVent ventilation technology offering Volume Control Ventilation with tidal volume compensation, Pressure Control Ventilation, and electronic PEEP. The proven SmartVent also features optional Pressure Support Ventilation with an Apnea Backup (PSVPro) that is used for spontaneously breathing patients, Synchronized Intermittent Mandatory Ventilation (SIMV) modes, Pressure control ventilation-volume guarantee (PCV-VG), and VCV cardiac bypass.

The Avance system is not suitable for use in a MRI environment.

Note

Configurations available for this product depend on local market and standards requirements. Illustrations in this manual may not represent all configurations of the product.





Figure 1-1 ■ S/5 Avance system

1009-0357-000 09/07 1-3


# 1.4 Anesthesia system components



- 1. Light switch
- 2. Dovetail/GCX rails
- 3. Vaporizer
- 4. Alternate O<sub>2</sub> control
- 5. System switch
- 6. Mains indicator
- 7. Integrated suction (optional)
- 8. Brake
- 9. 0<sub>2</sub> flush button
- 10. Advanced breathing system
- 11. Auxiliary O<sub>2</sub> flow control (optional)
- 12. Anesthesia display

Figure 1-2 • Front view

1-4 09/07 1009-0357-000




- 1. Outlet circuit breaker
- 2. Isolated electrical outlet (optional)
- 3. Cable access door
- 4. Vacuum connection
- 5. Collection bottle connection
- 6. Cylinder yoke
- 7. AGSS (Anesthesia Gas Scavenging System)
- 8. Equipotential stud
- 9. Mains inlet
- 10. System circuit breaker
- 11. Pipeline connections

Figure 1-3 • Rear view

1009-0357-000 09/07 1-5

# 1.5 Breathing system components




- 1. Expiratory check valve
- 2. Inspiratory check valve
- 3. Inspiratory flow sensor
- 4. Expiratory flow sensor
- 5. Absorber canister
- 6. Absorber canister release
- 7. Leak test plug
- 8. Breathing system release


- 9. Manual bag port
- 10. Adjustable pressure-limiting (APL) valve
- 11. Bag/Vent switch
- 12. Bellows assembly
- 13. Sample gas return port
- 14. AGSS indicator (only present on some AGSS versions)
- 15. Airway module (optional)

Figure 1-4 ■ Breathing system

1-6 09/07 1009-0357-000

# 1.5.1 Optional ABS components






- 1. Bag support arm
- 2. Auxiliary Common Gas Outlet (ACGO) switch
- 3. ACGO port
- 4. Airway module
- 5. EZchange Canister system (CO<sub>2</sub> Bypass)
- 6. EZchange Canister release
- 7. Condenser drain button
- 8. Condenser

Figure 1-5 • Breathing system options

1009-0357-000 09/07 1-7

# 1.6 Display controls



1. Silence Alarms key Push to silence any active, silenceable high and medium priority alarms or to suspend/acknowledge any non-active medium or high priority alarms.

Alarm is silenced for 120 seconds or alarm is suspended for 90 seconds.

2. Menu keys Push to show corresponding menu.

3. ComWheel Push to select a menu item or confirm a setting. Turn clockwise or counterclockwise to scroll menu items or change settings.

4. Normal Screen key Push to remove all menus from the screen.

Quick keys
 Push to change corresponding gas setting or ventilator setting. Turn the
 ComWheel to make a change. Push the ComWheel to activate the change.

### **HPDU** only

6. Timer keys Push to start or stop the timer. Push to reset the timer back to zero.


7. MV/TV Alarms key Push to turn off the MV and TV alarms. Push again to turn the MV and TV

alarms back on.

Figure 1-6 • Ventilator controls

1-8 09/07 1009-0357-000

# 1.7 Anesthesia system display




- 1. Electronic gas flow tubes
- 2. Alarm silence countdown
- 3. Alarm message fields
- 4. Waveform field
- 5. Clock
- 6. Number field
- 7. Free number display
- 8. Ventilator settings
- 9. Gas settings

Figure 1-7 • Normal view

1009-0357-000 09/07 1-9

When a menu key is selected, the menu field overlays the gas flow tubes and the waveform fields start at the right edge of the menu.




- 1. Menu
- 2. Waveform fields

Figure 1-8 • Menu view

1-10 09/07 1009-0357-000

### 1.7.1 Using menus

Push a menu key to display the corresponding menu. Use the ComWheel to navigate through the menu.



- 1. Menu title
- 2. Present selection
- 3. Adjustment window
- 4. Indicates submenu
- 5. Short instructions
- 6. Menu selections

Figure 1-9 • Example menu

- 1. Push the menu key to display the corresponding menu.
- 2. Turn the ComWheel counterclockwise to highlight the next menu item. (Turn the ComWheel clockwise to highlight the previous menu item.)
- 3. Push the ComWheel to enter the adjustment window or a submenu.
- 4. Turn the ComWheel clockwise or counterclockwise to highlight the desired selection.
- 5. Push the ComWheel to confirm the selection.
- 6. Select **Normal Screen** or push the **Normal Screen** key to exit the menu and return to the normal monitoring display. (Select **Previous Menu** to return to the last displayed menu, if available.)

1009-0357-000 09/07 1-11

# 1.8 Symbols used in the manual or on the equipment

Symbols replace words on the equipment, on the display, or in Datex-Ohmeda manuals. No one device or manual uses all of the symbols.

Warnings and Cautions tell you about dangerous conditions that can occur if you do not follow all instructions in this manual:


- Warnings tell about a condition that can cause injury to the operator or the patient.
- Cautions tell about a condition that can cause damage to the equipment.

Read and follow all warnings and cautions.

1	On (power)	×	Alarm silence
0	Off (power)	$\bowtie$	Alarm silence
0	Standby	<b>†</b>	Type B equipment
Ċ	Standby or preparatory state for part of the equipment	<b>†</b>	Type BF equipment
$\odot$	"ON" only for part of the equipment	•	Type CF equipment
Ċ	"OFF" only for part of the equipment	Ţ	Caution, ISO 7000-0434
===	Direct current	$\triangle \mathbf{A}$	Attention, refer to product instructions, IEC 60601-1
$\sim$	Alternating current	4	Dangerous voltage
	Protective earth ground	<u></u>	Electrical input
Ţ	Earth ground	$\Longrightarrow$	Electrical output
,,,	Frame or chassis ground		Pneumatic input
$\Diamond$	Equipotential	$\qquad \qquad \longrightarrow$	Pneumatic output


1-12 09/07 1009-0357-000

### 1 Introduction



1009-0357-000 09/07 1-13

### Avance



1-14 09/07 1009-0357-000

# 2 Theory of Operation

In this section	2.1 Electrical system	2-2
	2.2 Power subsystem (original)	2-4
	2.2.1 Power Controller board	2-5
	2.3 Power subsystem (with U-Frame power supply)	2-6
	2.3.1 Power Controller board	2-7
	2.3.2 Power distribution	2-8
	2.4 System communications	. 2-10
	2.4.1 Software Power On Self Tests (POST)	. 2-12
	2.5 Display Unit	. 2-14
	2.5.1 Software requirements	. 2-14
	2.6 System connections	. 2-15
	2.6.1 Display Unit	. 2-15
	2.6.2 High Performance Display Unit	. 2-15
	2.6.3 Display Connector board	. 2-16
	2.7 Power Controller and Anesthesia Control board connections	. 2-17
	2.8 Anesthesia Control board	. 2-18
	2.9 Electronic Gas Mixer	. 2-20
	2.10 Ventilator Interface board	. 2-22
	2.11 Gas flow through the anesthesia machine	. 2-24
	2.11.1 Overview	. 2-24
	2.11.2 Physical connections	. 2-28
	2.11.3 Suction regulators	. 2-29
	2.12 Flow through the breathing system	. 2-30
	2.12.1 Overview of flow paths	. 2-30
	2.12.2 Manual ventilation	
	2.12.3 Mechanical ventilation	
	$2.12.4$ Fresh gas and $0_2$ flush flow (with SCGO)	
	$2.12.5$ Fresh gas and $0_2$ flush flow (with ACGO)	
	2.13 Ventilator mechanical subsystems	
	2.13.1 Drive gas filter and Gas Inlet Valve	
	2.13.2 Pressure regulator     2.13.3 Flow control valve	
	2.13.4 Drive Gas Check Valve (DGCV)	
	2.13.5 Bellows Pressure Relief Valve	
	2.13.6 Exhalation valve	
	2.13.7 Mechanical Overpressure Valve	
	2.13.8 Reservoir and bleed resistor	
	2.13.9 Free breathing valve	. 2-50
	2.13.10 Breathing circuit flow sensors	. 2-51

1009-0357-000 09/07 2-1

## 2.1 Electrical system

The electrical system consists of two main computing units: the Display Unit and the Anesthesia Control board. Additional subsystems interact with these computing hosts to perform various gas delivery, ventilation, and monitoring functions.

The Display Unit handles the main user interface functions and connections to external devices. The Display Unit software runs on the Windows CE operating system.

Therapy functions are handled by the Anesthesia Control board. The Anesthesia Control board is based on the Motorola Coldfire processor with a Nucleus operating system.


Embedded controllers are used to perform specific machine functions on subsystems like the Power Controller board and the Mixer board.

The processors communicate through serial bus channels.

The various functions of the electrical system are accomplished on the following:

- Display Unit CPU (A)
- Display Unit System Interconnect assembly (B)
- Display Connector board (C)
- Power Control assembly (D\*)
- Anesthesia Control board (E)
- Pan Connector board (F)
- Electronic Mixer board (G)
- Ventilator Interface board (H)
- ABS Filter board (I)
- Vent Engine Connector board (J)
- MGAS Power Supply board (K)
- Light Strip board (L)
- Inrush board (M)
- Toroid (N)
- \* Early production machines use a full-size Power Controller board that includes an on-board AC/DC power supply. Late production machines (or machines serviced with replacement assemblies) use a universal power supply (U-Frame) that feeds into a half-width Power Controller board.

2-2 09/07 1009-0357-000



1009-0357-000 09/07 2-3

# 2.2 Power subsystem (original)

Mains power enters the system through the AC Inlet module ( $\bf A$ ), which includes a line filter and the system circuit breaker. Mains power is routed through the Inrush ( $\bf B$ ) circuit board to the isolation transformer ( $\bf C$ ).

The isolated secondary output of the transformer is routed through fuses (**D**) and a second line filter (**E**) to the input of the Power Controller board (**F**). If the system is equipped with electrical power outlets, the transformer (larger size) also supplies isolated power to the electrical outlets through individual circuit breakers.

The Power Controller board interfaces with the system through:

- the Anesthesia Control board connector (G),
- the Display Connector board connector (H),
- the battery connector (I) and fan connector (J).




Figure 2-10 • Power subsystem (original)

2-4 09/07 1009-0357-000

# 2.2.1 Power Controller board

The system uses a distributed power bus. The Power Controller board contains:

- an AC/DC converter that converts line voltage to high voltage DC.
- a DC/DC converter that converts the high voltage DC to battery voltage.
- a DC/DC converter that converts battery voltage to the 12.5 VDC system bus voltage.

The Power Controller contains supervisory circuitry that performs:

- battery charge control (battery switch circuits provide a minimum of 30 minutes of system power in the event of AC power failure).
- current, voltage, and temperature monitoring.
- AC sensing.
- fan control.

Two 12-volt batteries, wired in series, provide the back-up power.

The Power Controller communicates with the Display Unit through a RS-422, 9.6 kB channel. It receives the On/Standby signal from the system switch through the Anesthesia Control board.




Figure 2-11 • Power subsystem (original)

### 2.3 Power subsystem (with U-Frame power supply)

Mains power enters the system through the AC Inlet module ( $\bf A$ ), which includes a line filter and the system circuit breaker. Mains power is routed through the Inrush ( $\bf B$ ) circuit board to the isolation transformer ( $\bf C$ ).

The isolated secondary output of the transformer is routed through fuses ( $\mathbf{D}$ ) and a second line filter ( $\mathbf{E}$ ) to the universal power supply ( $\mathbf{F}$ ). The DC output of the power supply feeds into the Power Controller board ( $\mathbf{G}$ ). If the system is equipped with electrical power outlets, the transformer (larger size) also supplies isolated power to the electrical outlets through individual circuit breakers.

The Power Controller board interfaces with the system through:

- the Anesthesia Control board connector (H),
- the Display Connector board connector (I),
- the battery connector (J) and fan connector (K-Fan1).




Figure 2-12 • Power subsystem (with U-Frame Power Supply)

2-6 09/07 1009-0357-000

# 2.3.1 Power Controller board

The system uses a distributed power bus. The Power Controller board contains:

• a DC/DC converter that converts the input from the universal power supply to the 12.5 VDC system bus voltage.

The Power Controller contains supervisory circuitry that performs:

- battery charge control (battery switch circuits provide a minimum of 30 minutes of system power in the event of AC power failure);
- current, voltage, and temperature monitoring;
- AC sensing;
- fan control.

Two 12-volt batteries, wired in series, provide the back-up power.

The Power Controller communicates with the Display Unit through a RS-422, 9.6 kB channel. It receives the On/Standby signal from the system switch through the Anesthesia Control board.




Figure 2-13 • Power subsystem diagram (with U-Frame Power Supply)

# 2.3.2 Power distribution

The Power Controller board provides outputs to the Anesthesia Control board and the Display Connector board. These boards provide distribution of power supplies required by the system.

The Anesthesia Control board interfaces with the Mixer board and the Ventilator Interface board through the Pan Connector board.

The Display Connector board interfaces with the Display Unit and the Module assembly.




Figure 2-14 • Power distribution (original)

2-8 09/07 1009-0357-000




Figure 2-15 • Power distribution (with U-Frame Power Supply and an HPDU)

### 2.4 System communications

RS-422 serial communication is used between the two main processors — Display Unit and Anesthesia Computer — and the subsystem processors. Various baud rates accommodate data requirements between subsystem and host. External communication uses the standard RS-232 interface.




Figure 2-16 • System communications (with DU)

2-10 09/07 1009-0357-000




Figure 2-17 • System communications (with HPDU)

### 2.4.1 Software Power On Self Tests (POST)

When the machine is in the "Off" state (No AC Power connected to the machine and the On/Standby Switch in the "Standby" position, there is no machine activity. Circuitry on the Power Controller Board and on the Anesthesia Control Board monitors the On/Standby Switch for movement.

If the machine state changes from the "Off" state to the "On" state (AC Power connected to the machine while the On/Standby Switch in the "Standby" position), the hospital AC (stepped to approximately 180v AC by the isolation transformer) enters the 150w Power Supply (or PCB). The AC voltage is converted to +12v dc. Standard CPU tests are performed (including but not limited to RAM, ROM, Watchdog, and application CRC) and the PCB application starts. Once the PCB has passed its CPU tests and the application is running, the +12 v circuit to power the DU and the ACB is activated. The DU and ACB simultaneously begin powering on. The PCB evaluates the battery capacity and charges the batteries, if necessary.

Both the ACB and the DU begin by converting the incoming +12 v dc to local needed power (+3.3 v, +4 v, +5 v, +8 v, etc.) and perform standard CPU tests. Each board loads their software (that resides locally) and begins their appropriate self tests described below:

- Once the ACB passes all CPU tests and application is loaded, independent circuitry turns on 10 VA limited (+12 v dc) power to the Gas Mixer, the Ventilator Interface Board, and the M-Gas Power Supply Board. These boards convert the incoming +12 v dc into locally need power supplies and simultaneously begin to power up. The ACB energizes the Alt O2 Selector Valve, which closes the valve. The ACB begins testing the GIV.
  - The processor on the Gas Mixer board performs standard CPU tests, checks communication link with the ACB, reports to the ACB that it has begun to perform Power-On Self Tests (POST). The Gas Mixer tests the O2 and Balance Gas channels for leaks and flow delivery from both channels. Once completed, reports to the ACB that all Self-Tests have been completed and the Gas Mixer compatibility information (serial number, hardware revision, and software revision).
  - The processor on the VIB performs standard CPU tests and begins to download application software and ventilator calibration constants from the ACB. The VIB works in conjunction with the ACB to perform the GIV test (ACB provides valve status information to the VIB, the VIB calculates the voltage necessary to open the valve to the desired value). Once completed, the VIB reports to the ACB Self-Tests have been completed and the VIB compatibility information (serial number, hardware revision, and software revision).
  - There is no processor on the M-Gas Power Supply. Local voltages are produced for the M-Gas and available if a M-Gas is installed. The M-Gas module performs its on POST.

2-12 09/07 1009-0357-000

• Once the DU passes all the CPU tests, the DU application is started. Part of this application is software that enables the DU CPU to communicate with the M-Gas module. These applications take longer to start than any other system in the machine. When the applications have completely loaded and communication has been established with the ACB, all the systems compatibility information is transferred to the DU for comparison with the compatibility table created during the last software download. If the compatibility does not match, the machine enters the System Malfunction.

If the On /Standby switch is moved to the "On" position while the machine is in the "Off" state, the circuitry on the Power Controller Board engages the batteries to power the CPU on the PCB. Standard CPU tests are performed and the PCB application starts. Once the PCB has passed its CPU tests and the application is running, it activates the +12 v circuit to power the DU and the ACB. The DU and ACB simultaneously begin powering on. The PCB evaluates the battery capacity.

### 2.5 Display Unit

The Avance anesthesia machine can use one of two display units:

- the original Display Unit (DU),
- or the High Performance Display Unit (HPDU).

The display unit handles most of the machine's user interface functions through the front panel controls and the LCD screen. It is the primary interface to external peripherals.

The main components of the display unit include:

- An active matrix thin film transistor liquid crystal display (A)
- The CPU board (B)
- The System Interconnect assembly (C)

The CPU board includes a host processor and three coprocessors to handle display, front panel, and monitoring interfaces.

Both the DU and HPDU use a lithium battery to power the real time clock when the machine is in Standby or Off states. The HPDU also uses the lithium battery to retain the BIOS setup. If the battery is removed or becomes disconnected, the BIOS setup needs to be restored by booting the system off the Software Download Card.

# 2.5.1 Software requirements

The DU uses a PCMCIA interface to handle software upgrades and to load the diagnostics Service Application. The DU is compatible with system software version 3.2 or earlier.

The HPDU uses a Compact Flash interface to handle software upgrades and to load the diagnostics Service Application. The HPDU requires system software version 4.X or greater.




Figure 2-18 • Display Unit


2-14 09/07 1009-0357-000

### 2.6 System connections

### 2.6.1 Display Unit

The DU accommodates the following connections:


- System Power Interface (1).
- System Signal Interface (2).
- Serial Port standard interface for external communication (3).
- DIS connector supports D-O Device Interface Solution (DIS) (4).
- Network connection Standard Ethernet port for network connectivity (5).
- Network ID accept D-O proprietary network identification plug (6).
- USB port standard USB 1.1 interface (7).



### 2.6.2 High Performance Display Unit

The HPDU accommodates the following connections:


- System Power Interface (1).
- System Signal Interface (2).
- Serial Port standard interface for external communication (3).
- Remote monitor On/Standby (4).
- Network connection Standard Ethernet port for network connectivity (5).
- USB port standard USB 2.0 interface (6).
- USB port standard USB 2.0 interface (7).



# 2.6.3 Display Connector board

The top side of the Display Connector board accepts the following cables:

- System Power Interface to Display Unit (1).
- System Signal Interface to Display Unit (2).
- Airway Module (MGAS) Power Supply board (8).
- Not used (9).
- The under side of the Display Connector board accepts the following cables:
- Power Controller board (10).
- Anesthesia Control board (MGAS power) connector (11).
- Anesthesia Control board (signal) connector (12).
- Not used (13).



2-16 09/07 1009-0357-000


### 2.7 Power Controller and Anesthesia Control board connections

The Power Controller:

- Distributes 12.5 VDC power and communicates with the Display Unit (by way of the Display Connector board) through connector (**10**).
- Distributes 12.5 VDC power to the Anesthesia Control board through connector (14).


The Anesthesia Control board:

- Receives power from the Power Controller board through connector (14).
- Distributes 10VA power supplies to the Pan Connector board through connector (15).
- Communicates with Pan assemblies through connector (16).
- Communicates with Display Unit through connector (12).
- Distributes 10VA power supplies to the Display Unit through connector (11).



Display Connector board (topside)

**Power Controller board** 



Display Connector board (underside)

**Anesthesia Control board** 

### 2.8 Anesthesia Control board

The Anesthesia Control board (**A**) uses a Motorola MCF5307 Coldfire microcontroller with 4M Flash and 16M error correcting DRAM. The Anesthesia Control board includes 6 UARTs with a 64 byte FIFO and RS-422 communications to interface with the Display Unit, an accessory port, and anesthesia delivery subsystems located in the pan electronic enclosure. These include the Gas Mixer and the Ventilator Interface board.



Figure 2-19 • Anesthesia Control board

2-18 09/07 1009-0357-000




Figure 2-20 • Anesthesia Control board block diagram

### 2.9 Electronic Gas Mixer

The Gas Mixer receives its pneumatic inputs from the pipeline and cylinder supplies and sends mixed gas to the vaporizer manifold. The Gas Mixer interfaces to the Anesthesia Control board for power and communications.

The Gas Mixer consists of the following subassemblies and main components:

- Gas Mixer board (A)
- Control Manifold (**B**) manifold, selector valves, proportional valves
- Flow sensor assembly (C)
- Mixed gas manifold and exit check valve (D)




Figure 2-21 • Electronic Gas Mixer

2-20 09/07 1009-0357-000

Desired gas flows are sent from the Anesthesia Control board to the Gas Mixer.

Gas Mixer operation is controlled through a microcontroller which:

- Sends requests for the Anesthesia Control board to open and close selector valves for O<sub>2</sub>, N<sub>2</sub>O and Air.
- Regulates flow control valves for O<sub>2</sub> and balance gas (N<sub>2</sub>O or Air).

Closed-loop flow control is accomplished through a hot-wire anemometer in concert with the flow control valves. Gas flow, based on a calibration table, is on target when the reference measurement equals the flow measurement.

Pressure measurements across each of the flow sensor channels are used as checks on the flow measurement for hazard mitigation, ambient pressure compensation, and compensation for back pressure downstream of the Mixer.

In case of certain failures or errors, Alternate  $\rm O_2$  control activates automatically to delivers  $\rm O_2$  (and agent, if turned on) through an alternate pneumatic path to the patient circuit. Alternate  $\rm O_2$  can be activated manually through a front panel control.




Figure 2-22 • Electronic Gas Mixer block diagram

### 2.10 Ventilator Interface board

The Ventilator Interface board (A) provides the electrical and/or pneumatic interface to the following:

- Inspiratory (**B**) and expiratory (**C**) flow sensors (transducers)
- Patient airway (**D**) and manifold (**E**) pressures (transducers)
- Oxygen sensor (in breathing system)
- ABS On switch
- Canister Release switch
- ACGO position switch (if ACGO installed)
- SCGO solenoid, SCGO/CGO position switches (if SCGO installed)
- Bag/Vent switch
- 0<sub>2</sub> Flush switch
- Gas Inlet Valve
- Inspiratory Flow Valve
- Accessory Power (for task lights)

The Ventilator Interface board functions are managed locally by a microcontroller. The microcontroller communicates data values to the controlling CPU via an RS-422 serial interface.




Figure 2-23 • Ventilator Interface board

2-22 09/07 1009-0357-000




Figure 2-24 • Ventilator Interface board block diagram

### 2.11 Gas flow through the anesthesia machine

### 2.11.1 Overview

Refer to Figure 2-25.

#### Gas supplies

Gas comes into the system through a pipeline (1) or cylinder (6) connection. All connections have indexed fittings, filters, and check valves (one-way valves). Pressure transducers monitor the pipeline (2) and cylinder (7) pressures.

The  $\rm O_2$  supply failure alarm is derived from the  $\rm O_2$  pipeline and the  $\rm O_2$  cylinder pressure transducer inputs.

A primary regulator (8) decreases the cylinder pressures to approximately pipeline levels. A pressure relief valve (3) helps protect the system from high pressures.

To help prevent problems with the gas supplies:

- Install yoke plugs on all empty cylinder connections.
- When a pipeline supply is adequate, keep the cylinder valve closed.

#### **Gas flow**

Pipeline or regulated cylinder pressure supplies  $O_2$  or Air directly to the ventilator engine (**4a** or **4b**) and as pilot pressure (**4**) for the SCGO assembly (**E**). Connection points are also available for venturi suction (**5a** or **5b**) drive gas supply. An additional  $O_2$  regulator (**18**) decreases the pressure for the  $O_2$  Flush valve (**19**) and the auxiliary  $O_2$  flowmeter (**24**).

The  $O_2$  Flush valve supplies high flows of  $O_2$  to the fresh gas outlet (**22** or **23**) through the SCGO/ACGO assembly (**E/F**). The flush pressure switch (**20**) monitors activation of the flush valve.

#### **Gas mixing**

Under normal conditions, with the system switch ( $\mathbf{10}$ ) in the On position, the Alternate  $O_2$  Disable valve ( $\mathbf{13}$ ) is energized to block alternate O2 flow. Normal gas flows are enabled through their respective selector valves ( $\mathbf{11}$ ). The system controls gas flow through the flow control valves ( $\mathbf{12}$ ) and derives the individual flow rates through the hot-wire anemometers ( $\mathbf{14}$ ).

Under system failure conditions (or if Alt  $O_2$  is selected), the normally-open Alternate  $O_2$  Disable valve (13) allows delivery of  $O_2$  through the Alternate  $O_2$  Flowmeter when the system switch is in the On position.

#### Mixed gas

The mixed gas flows through the vaporizer manifold ( $\mathbf{D}$ ), and vaporizer ( $\mathbf{16}$ ) that is On, to the SCGO/ACGO assembly ( $\mathbf{E}/\mathbf{F}$ ). A pressure relief valve ( $\mathbf{17}$ ) on the vaporizer manifold limits the maximum outlet pressure.

The SCGO assembly (**E**) directs the mixed gas to the selected circuit: **22** (ABS-circle) or **23** (to Inspiratory port of ABS). On SCGO assemblies, a relief valve (**21**) limits pressure in the breathing system to approximately 150 cmH<sub>2</sub>O.

The ACGO assembly (**F**) directs the mixed gas to the selected circuit: **22** (ABS-circle) or **23** (external ACGO port).

2-24 09/07 1009-0357-000




Figure 2-25 • Pneumatic circuit




Figure 2-26 • Pneumatic circuit

2-26 09/07 1009-0357-000

#### Refer to Figure 2-26.

## Key to Numbered Components

- 1. Pipeline inlet
- 2. Pipeline pressure transducer
- 3. High-pressure relief valve (758 kPa / 110 psi)\*
- 4. Supply connections for the ventilator and pilot pressure for SCGO
  - a. O<sub>2</sub> drive gas
  - b. Air drive gas
- 5. Venturi suction supply connection
  - a.  $0_2$  drive gas
  - b. Air drive gas
- 6. Cylinder inlet
- 7. Cylinder pressure transducer
- 8. Primary regulator (cylinder pressure)
- 9. Test port (primary regulator)
- 10. System switch
- 11. Selector valve

$$a = O_2$$
;  $b = Air$ ;  $c = N_2O$ 

12. Flow controller

$$a = 0_2$$
;  $b = balance gas$ 

- 13. Alternate 02 disable valve
- 14. Hot-wire anemometer

a = O<sub>2</sub> flow sensor channel; b = balance gas flow sensor channel

- 15. Vaporizer port valve
- 16. Vaporizer
- 17. Low-pressure relief valve (38 kPa / 5.5 psi)\*
- 18. O<sub>2</sub> flush and auxiliary flowmeter regulator (241 kPa / 35 psi)\*
- 19. 02 Flush valve
- 20. Pressure switch (used with the ventilator)
- 21. Breathing system pressure relief valve (SCGO only 150 cmH<sub>2</sub>O)\*
- 22. To Port 3 of ABS interface (circle)
- 23. For SCGO, to Port 2 of ABS interface (non-circle Inspiratory port) For ACGO, to external 22-mm ACGO connector
- 24. Auxiliary O<sub>2</sub> flowmeter (optional)
- \* Approximate values

### **Key to Symbols**

→ → Pneumatic Connection


← Filter

Direction of Flow

Check Valve

# 2.11.2 Physical connections

Figure 2-27 shows the physical path that the gas takes. The item numbers are described on the previous page.



2-28 09/07 1009-0357-000

Figure 2-27 • Typical tubing connections - pictorial

# 2.11.3 Suction regulators

### **Pipeline vacuum**

The suction regulator (shown in Figure 2-27) uses an external vacuum source.

#### **Venturi Drive vacuum**

The suction regulator (shown in Figure 2-28) uses an internal, venturi derived vacuum source.

Drive gas (internally plumbed  $\operatorname{Air}$  or  $\operatorname{O}_2$ ) enters the Venturi Module (VM) at the drive port (A). As the drive gas passes through the venturi module, a vacuum is created at port B. The drive gas exits the venturi module at port C and is exhausted outside the machine through the muffler (D).

The control port (**E**) on the venturi module responds to pneumatic signals from the front panel switch on the Suction Control Module (**SCM**) to turn the venturi vacuum drive gas on or off. The check valve (**CV**) helps prevent pressurization of the suction circuitry if the exhaust is occluded or the venturi unit fails.




Figure 2-28 • Venturi suction

### 2.12 Flow through the breathing system

# 2.12.1 Overview of flow paths

This section looks at four types of flow paths.

- **Ventilation paths:** How gas flows from the drive source (bag or bellows) to and from the patient.
- Fresh gas paths: Fresh gas can flow from the machine interface directly to the patient through the inspiratory check valve, or through the absorber into the expiratory flow, or directly to an external circuit through the optional auxiliary common gas outlet.
- Scavenged gas paths: APL or Pop-off.
- Flow through the optional EZchange Canister and Condenser: EZchange ON and EZchange OFF (CO<sub>2</sub> bypass).

2-30 09/07 1009-0357-000


### 2.12.2 Manual ventilation

### **Manual inspiration**

The Bag/Vent switch closes the ventilator path (B).

Gas flows from the bag (1), through the absorber (2), into the breathing circuit module, and through a unidirectional valve (inspiratory check valve) to the patient (3).

During inspiration, fresh gas (**FG**) flows from the machine into the inspiratory limb, upstream of the inspiratory check valve.



- **AP** Airway Pressure
- **B** Bag/Vent switch to Bag
- **FG** Fresh Gas
- 1 Flow to absorber
- **2** Flow from absorber
- **3** Inspiratory flow

Figure 2-29 • Gas flow during manual inspiration

#### **Manual expiration**

The Bag/Vent switch keeps the ventilator path closed (**B**).

Gas flows from the patient (4), through a unidirectional valve (expiratory check valve), and into the bag (5).

During exhalation, fresh gas flows backwards through the absorber (**FG**) into the expiratory limb, downstream of the expiratory check valve.

For machines that are plumbed to return sample gas to the breathing system, the returned gas (**SGR**) enters the breathing system after the expiratory check valve.





Figure 2-30 • Flow during manual expiration

2-32 09/07 1009-0357-000

**APL Valve** The APL valve sets a pressure limit for manual ventilation.

As you turn the APL knob, it puts more or less force on the APL disc and seat (**D/S**). If the circuit pressure is too high (**6**), the disc and seat inside the diaphragm opens and vents gas to the scavenging system (**7**).



**D/S** APL disc and seat

**6** APL flow

**7** To scavenging

Figure 2-31 • Flow through the APL Valve


### 2.12.3 Mechanical ventilation

### **Mechanical inspiration**

The Bag/Vent switch closes the manual path  $(\mathbf{V})$ . Pilot pressure  $(\mathbf{P})$  closes the exhalation valve.

Drive gas (**D**) pushes down on the bellows. Gas flows from the bellows (**1**), through the absorber (**2**), and through a unidirectional valve (inspiratory check valve) to the patient (**3**).

During inspiration, fresh gas flows into the inspiratory limb, upstream of the inspiratory check valve.



- **AP** Airway Pressure
- **D** Drive gas
- **FG** Fresh Gas
- P Pilot pressure
- **V** Bag/Vent switch to Vent
- 1 Flow to absorber
- **2** Flow from absorber
- 3 Inspiratory flow

Figure 2-32 • Mechanical inspiration

2-34 09/07 1009-0357-000


#### **Mechanical expiration**

Drive-gas flow stops and the exhalation valve opens. Exhaled gas flows from the patient (4), through a unidirectional valve (expiratory check valve) and into the bellows (5). Residual drive gas (**D**) flows out of the bellows to the scavenging system (6).

If PEEP is selected, static pressure on the pilot port of the exhalation valve sets the PEEP level.

During exhalation, fresh gas flows backwards through the absorber (**FG**) into the expiratory limb, downstream of the expiratory check valve.

For machines that are plumbed to return sample gas to the breathing system, the returned gas (**SGR**) enters the breathing system after the expiratory check valve.



**AP** Airway Pressure

**D** Drive gas

**FG** Fresh Gas

**SGR** Sample Gas Return

- **4** Expiratory flow
- **5** Flow to bellows
- **6** To scavenging

Figure 2-33 • Mechanical expiration

### Mechanical inspiration (EZchange and condenser ON)

The Bag/Vent switch closes the manual path  $(\mathbf{V})$ . Pilot pressure  $(\mathbf{P})$  closes the exhalation valve.

Drive gas (**D**) pushes down on the bellows. Gas flows from the bellows (**1**), through the absorber (**2a**), Condenser (**2b**), and through a unidirectional valve (inspiratory check valve) to the patient (**3**).

During inspiration, fresh gas flows into the inspiratory limb, upstream of the inspiratory check valve.




Figure 2-34 • Mechanical inspiration through Condenser with EZchange Canister and Condenser ON

2-36 09/07 1009-0357-000

### Mechanical expiration (EZchange and condenser ON)

Drive-gas flow stops and the exhalation valve opens. Exhaled gas flows from the patient (4), through a unidirectional valve (expiratory check valve) and into the bellows (5). Residual drive gas (D) flows out of the bellows to the scavenging system (6).

If PEEP is selected, static pressure on the pilot port of the exhalation valve sets the PEEP level.

During exhalation, fresh gas flows backwards through the Condenser and absorber (**FG**) into the expiratory limb, downstream of the expiratory check valve.

For machines that are plumbed to return sample gas to the breathing system, the returned gas (**SGR**) enters the breathing system after the expiratory check valve.




Figure 2-35 • Mechanical expiration through Condenser with EZchange Canister and Condenser ON

### Mechanical inspiration (EZchange and condenser OFF)

The Bag/Vent switch closes the manual path  $(\mathbf{V})$ . Pilot pressure  $(\mathbf{P})$  closes the exhalation valve.

Drive gas ( $\mathbf{D}$ ) pushes down on the bellows. Gas flows from the bellows ( $\mathbf{1}$ ), through the EZchange module bypassing the absorber ( $\mathbf{2}$ ), and through a unidirectional valve (inspiratory check valve) to the patient ( $\mathbf{3}$ ).

During inspiration, fresh gas flows into the inspiratory limb, upstream of the inspiratory check valve.




Figure 2-36 • Mechanical inspiration with EZchange Canister and Condenser OFF

2-38 09/07 1009-0357-000

# Mechanical expiration (EZchange and condenser OFF)

Drive-gas flow stops and the exhalation valve opens. Exhaled gas flows from the patient (4), through a unidirectional valve (expiratory check valve) and into the bellows (5). Residual drive gas (**D**) flows out of the bellows to the scavenging system (6).

If PEEP is selected, static pressure on the pilot port of the exhalation valve sets the PEEP level.

During exhalation, fresh gas flows backwards through the EZchange module (**FG**) into the expiratory limb, downstream of the expiratory check valve.

For machines that are plumbed to return sample gas to the breathing system, the returned gas (**SGR**) enters the breathing system after the expiratory check valve.





Figure 2-37 • Mechanical expiration with EZchange Canister and Condenser OFF

### Pop-off valve

The pop-off valve limits the pressure inside the bellows to  $2.5~{\rm cm}~{\rm H}_2{\rm O}$  above the drive gas pressure. This normally occurs when the bellows reaches the top of the housing at the end of exhalation (5).

Excess gas (6) vents to the scavenging system (7) through the pop-off valve and the exhalation valve.



- 5 Flow to bellows
- **6** Pop-off flow
- **7** To scavenging

Figure 2-38 • Flow through the pop-off valve

2-40 09/07 1009-0357-000

### 2.12.4 Fresh gas and 0<sub>2</sub> flush flow (with SCGO)

## To ABS (Circle) breathing system

Fresh gas (1) flows from the vaporizer manifold outlet to the SCGO assembly.

With the Circle system selected, fresh gas flow is channeled to Port 3 of the breathing system (before the inspiratory check valve).

The output of the  $O_2$  Flush regulator (2) is channeled to the  $O_2$  Flush valve. When activated,  $O_2$  flush flow joins the fresh gas flow in the SCGO assembly.




Figure 2-39 • Fresh gas and O<sub>2</sub> flush flow (to ABS)

1009-0357-000 09/07 2-41

#### Switched (Non-circle) Common Gas Outlet

Fresh gas (1) flows from the vaporizer manifold outlet to the SCGO assembly.

With the Non-Circle system selected, fresh gas flow is channeled to Port 2 of the breathing system (after the inspiratory check valve - to an external patient circuit through the Inspiratory port).

The output of the  $\rm O_2$  Flush regulator (2) is channeled to the  $\rm O_2$  Flush valve. When activated,  $\rm O_2$  flush flow joins the fresh gas flow in the SCGO assembly.




Figure 2-40 • Fresh gas and O<sub>2</sub> flush flow (to Insp port)

2-42 09/07 1009-0357-000

## $\bf 2.12.5$ Fresh gas and $\bf 0_2$ flush flow (with ACGO)

## To ABS (Circle) breathing system

Fresh gas (1) flows from the vaporizer manifold outlet to the ACGO Selector Switch.

With the ACGO Selector Switch in the ABS position, fresh gas flow is channeled to the breathing system.

The output of the  $\rm O_2$  Flush regulator (2) is channeled to the  $\rm O_2$  Flush valve. When activated,  $\rm O_2$  flush flow joins the fresh gas flow in the ACGO Selector Switch.




Figure 2-41 • Fresh gas and O<sub>2</sub> flush flow (to ABS)

1009-0357-000 09/07 2-43

## Auxiliary (Non-circle) Common Gas Outlet

Fresh gas (1) flows from the vaporizer manifold outlet to the ACGO Selector Switch.

With the ACGO Selector Switch in the ACGO position, fresh gas flow is channeled to the ACGO outlet.

At the ACGO outlet, a small sample is diverted to the  $\rm O_2$  Cell in the ABS for  $\rm O_2$  monitoring.

The output of the  $\rm O_2$  Flush regulator (2) is channeled to the  $\rm O_2$  Flush valve. When activated,  $\rm O_2$  flush flow joins the fresh gas flow in the ACGO Selector Switch.




Figure 2-42  $\bullet$  Fresh gas and  $O_2$  flush flow (to ACGO)

2-44 09/07 1009-0357-000

### 2.13 Ventilator mechanical subsystems

Refer to Figure 11-1, "System circuit diagram" in Section 11, for the complete pneumatic/mechanical subsystem diagram.

The mechanical subsystems for the ventilator include:

**Pneumatic Vent Engine** 

- Drive gas inlet filter
- Gas inlet valve
- Supply gas pressure regulator
- Flow control valve
- Drive gas check valve
- Mechanical Overpressure Valve (MOPV)
- Bleed resistor
- Free breathing valve

Exhalation valve


Bellows assembly

Breathing circuit flow sensors

# 2.13.1 Drive gas filter and Gas Inlet Valve

Drive gas (can be selected from  $O_2$  or Air) enters the Vent Engine (1) at a pressure of 241 to 690 kPa (35 to 100 psi) through a 2-micron filter (2) that is located under the Gas Inlet Valve (3).

During normal operation the Gas Inlet Valve (GIV) is open to let supply gas flow. The GIV shuts off supply gas to the ventilator under failure conditions detected by the CPU or over-pressure switch. The output from the GIV stays at the filtered supply gas pressure.







Figure 2-43 • Inlet filter and Gas Inlet Valve (GIV)

1009-0357-000 09/07 2-45

## 2.13.2 Pressure regulator

The pressure regulator (4) is a non-relieving pressure regulator that regulates high pressure filtered supply gas down to 172 kPa (25 psi).







Figure 2-44 • Pressure regulator

## 2.13.3 Flow control valve

The flow control valve ( $\mathbf{5}$ ) is controlled by the CPU. Signals are sent to the flow control valve of the necessary flow determined by ventilator settings and sensor signals. The flow control valve modulates the incoming 172 kPa (25 psi) drive gases to an output from 0 to 120 liters per minute at pressures ranging from 0 to 100 cm  $H_2O$ .







Figure 2-45 • Flow control valve

2-46 09/07 1009-0357-000

#### 2.13.4 Drive Gas Check Valve (DGCV)

The Drive Gas Check Valve ( $\mathbf{6}$ ) is used downstream of the flow control valve to create the pilot pressure for closing the exhalation valve during inspiratory phases. The DGCV is biased shut by an integral weight that supplies approximately 3.5 cm H $_2$ O of bias pressure before permitting flow downstream to the bellows assembly. When the ventilator is exhausting flow from the breathing circuit, the DGCV permits the exhalation valve pilot pressure to be de-coupled from the circuit pressure. This permits the exhalation valve to open and lets gas flow to the exhaust and the gas scavenging system.






Figure 2-46 • Drive Gas Check Valve

#### 2.13.5 Bellows Pressure Relief Valve

The Bellows assembly is the interface between drive gas and patient gas in the breathing system. The pressure relief valve (or pop-off valve) in the bellows assembly (7) controls the pressure in the breathing circuit and exhausts excess patient gas through the exhalation valve.

The pressure relief valve is normally closed, maintaining approximately 1.5 cm  $H_2O$  in the breathing circuit in a no-flow condition, enough to keep the bellows inflated. It is piloted closed during inspiration and remains closed until the bellows is refilled during exhalation. It will exhaust  $\leq 4$  L/min excess fresh gas flow at  $\leq 4$  cm  $H_2O$ .






Figure 2-47 • Bellows pressure relief valve

1009-0357-000 09/07 2-47

## 2.13.6 Exhalation valve

The exhalation valve contains an elastomeric diaphragm that is used along with the flow valve to control the pressures in the breathing circuit. The exhalation valve includes two male ports on the bottom for:


- Bellows drive gas (8)
- Exhalation valve pilot (9) (manifold pressure)

The exhalation valve includes three ports on top that connect to the bellows base manifold:

- Drive gas pass through (10)
- Drive gas return and pop-off valve flow (11)
- APL exhaust flow to scavenging (12)

A port at the back of the exhalation valve (13) connects to the down tube that directs all the exhaust flows to the scavenging receiver.

The exhalation valve is normally open. Approximately  $2 \text{ cm H}_2\text{O}$  of pilot pressure is necessary to close the valve. When the exhalation port is open, gas flows from the bellows housing to the scavenging port.






Figure 2-48 • Exhalation valve

2-48 09/07 1009-0357-000

## 2.13.7 Mechanical Overpressure Valve

The Mechanical Overpressure Valve (MOPV) is a mechanical valve ( $\mathbf{14}$ ) that operates regardless of electrical power. It functions as a third level of redundancy to the ventilator's pressure limit control functions, supplying pressure relief at approximately  $110 \text{ cm H}_2O$ .




Figure 2-49 • Mechanical overpressure valve

## 2.13.8 Reservoir and bleed resistor

The reservoir (**15**) is a 200 ml chamber that dampens the manifold (pilot) pressure pulses to the exhalation valve.

The bleed resistor (16) is a "controlled leak" from 0 to 12 l/min in response to circuit pressures from 0 to 100 cm  $\rm H_2O$ . The small quantity of pneumatic flow exhausting through the bleed resistor permits control of the exhalation valve's pilot pressure by modulation of the valve output. The bleed resistor exhausts only clean drive gas and must not be connected to a waste gas scavenging circuit. The output is routed away from the electrical components to make sure that systems using oxygen drive gas meet the 10VA limitation requirement for oxygen enrichment.





Figure 2-50 • Reservoir and bleed resistor

1009-0357-000 09/07 2-49

## 2.13.9 Free breathing valve

The free breathing valve (17) helps assure the patient can spontaneously breathe. The ventilator is programmed to supply a specified number of breaths per minute to the patient. If, in between one of these programmed cycles, the patient needs a breath (spontaneous), the free breathing valve permits the patient to inhale. The free breathing valve is closed on mechanical inspiration.






Figure 2-51 • Free breathing valve

2-50 09/07 1009-0357-000

# 2.13.10 Breathing circuit flow sensors

Two flow sensors are used to monitor inspiratory and expiratory gas flow:

- The inspiratory flow sensor is downstream of the breathing system inspiratory check valve.
- The expiratory flow sensor is located at the input to the breathin system expiratory check valve.

#### For System software prior to 5.0

Feedback from the inspiratory transducer is used to supply tidal volumes that make allowances for the effects of fresh gas flow and circuit compressibility.

Feedback from the expiratory flow sensor is used to supply signals for expiratory tidal volume monitoring and the breath rate.

#### For System software 5.0 or greater

Feedback from both the inspiratory and expiratory transducers is used to:

- supply tidal volumes that make allowances for the effects of fresh gas flow and circuit compressibility.
- supply signals for expiratory tidal volume monitoring and the breath rate.




Figure 2-52 • Flow sensors

1009-0357-000 09/07 2-51

Notes

2-52 09/07 1009-0357-000

## **3 Checkout Procedure**

In this section	3.1 Inspect the system	3-2
	3.2 System checkout (for System software 3.X or greater)	3-2
	3.2.1 Leak < 250 ml	3-2
	3.2.2 Machine Check	
	3.2.3 Machine Check - System (Ventilator Circuit Testing)	
	3.2.4 Machine Check - Circuit (Bag Circuit Testing)	
	3.2.5 Machine Check - Circuit O2	
	3.3 Individual Checks (for System software 3.X or greater)	
	3.3.1 System	
	3.3.2 Circuit	
	3.3.3 Circuit 02 Cell	
	3.3.4 Low P Leak	
	3.3.5 Low P Leak (machines with ACGO)	
	3.4 System "All checks" (for System software 2.X)	
	3.4.1 Low P leak check	
	3.4.2 Quick check	
	3.4.3 Vent check	
	3.4.4 Circuit O <sub>2</sub> cell check	
	3.5 Bellows drop test	
	3.6 Backlight test	
	3.7 Vaporizer back pressure test	3-8
	3.8 Pipeline and cylinder tests	
	3.8.1 O <sub>2</sub> supply alarm test	
	3.9 Pressure relief tests	3-10
	3.10 Flush Flow Test	3-11
	3.11 Alarm tests	3-12
	3.12 Alternate O2 flowmeter tests	3-13
	3.13 Auxiliary O2 flowmeter tests	3-13
	3.14 Integrated Suction Regulator tests	3-13
	3.15 Power failure test	3-14
	3 16 Flectrical safety tests	3-1/

#### **⚠ WARNINGS**

After any repair or service of the Avance system, complete all tests in this section.

Before you do the tests in this section:

- Complete all necessary calibrations and subassembly tests. Refer to the individual procedures for a list of necessary calibrations.
- Completely reassemble the system.

If a test failure occurs, make appropriate repairs and test for correct operation.

1009-0357-000 09/07 3-1

### 3.1 Inspect the system

#### **⚠** CAUTION

The upper shelf weight limit is 34 kg (75 lb).

#### **⚠** WARNING

Do not leave gas cylinder valves open if the pipeline supply is in use. Cylinder supplies could be depleted, leaving an insufficient reserve supply in case of pipeline failure.

Before testing the system, ensure that:

- The equipment is not damaged.
- Components are correctly attached.
- The breathing circuit is correctly connected, not damaged.
- Pipeline gas supplies are connected.
- Cylinder valves are closed.
- Models with cylinder supplies have a cylinder wrench attached to the system.
- Models with cylinder supplies have a reserve supply of O<sub>2</sub> connected to the machine during system checkout.
- The casters are not loose and the brakes are set and prevent movement.
- The power cord is connected to a wall outlet. The mains indicator comes on when AC Power is connected.

## 3.2 System checkout (for System software 3.X or greater)

For System software 2.X refer to Section 3.4.

#### 3.2.1 Leak < 250 ml

The **Leak < 250 ml** setting is used during the circuit leak check portion of the checkout procedures. This check tests for leaks in the machine, breathing circuit, patient circuit, and manual bag. The default setting is **No**.

**Note** 

Extraction of gas by external gas monitors may cause failure of the leak checks during tests.

When No is selected, the leak test will pass for leaks below 250 ml at 3 kPa (30 cmH $_2$ 0) pressure with no user interaction required. For leaks between 250 ml and 750 ml, the user can fix the leak and rerun the test or accept the leak and continue. For leaks above 750 ml, the test will fail and the user must fix the leak and rerun the test.

Set to **Yes** to quantify small leaks above 100 ml during the checkout procedures. Selecting **Yes** will display the measured leak at 3 kPa (30 cmH<sub>2</sub>0) pressure and result in the test taking somewhat longer.

3-2 09/07 1009-0357-000

#### 3.2.2 Machine Check

The machine check runs automatically and beeps to indicate when it is finished or if interaction is required.

The Machine Check does a:

- Machine Check System check (Ventilator Circuit),
- Machine Check- Circuit check (Bag Circuit),
- and a *Machine Check Circuit 02* cell check (if circuit 0<sub>2</sub> cell is present).

When one of the checks is completed, the system will transition to the next check.

- 1. Turn the System switch to On.
- 2. Select Machine Check and follow the instructions.
- 3. If a check fails, follow the instructions to perform a recheck or accept the results.

### 3.2.3 Machine Check -System (Ventilator Circuit Testing)

The *Machine Check-System* checks the Bag/Vent switch, proper gas supply pressures, ventilator operation and leak, battery and electrical power, circuit compliance, flow control operation, and vaporizer operation. This is a two-step check.

- 1. Set the Bag/Vent switch to Vent.
- 2. Open the patient Y.
- 3. (ACGO option only.) Set the ACGO switch to Circle.
- 4. Select **Start**. The display shows the checks being run.
  - The system beeps when this portion of the check is done.
  - The results are shown on the display.
- 5. Make sure the bellows is fully collapsed.
- 6. Occlude the patient Y.
- 7. Select **Continue**. The display shows the checks being run.
- 8. When the check passes, the system will transition to the next step.

### 3.2.4 Machine Check -Circuit (Bag Circuit Testing)

The **Machine Check-Circuit** checks the Bag/Vent switch, proper gas supply pressures, airway pressure measurement transducer, APL valve, and manual circuit leak.

- 1. Occlude the patient Y.
- Set Bag/Vent switch to Bag.
- 3. Set the APL valve halfway between 30 and 70.
- 4. (ACGO option only.) Set the ACGO switch to Circle.
- 5. Select **Start**. The display shows the checks being run.
  - The system beeps when the check is done.
  - The results are shown on the display.
- 6. When the check passes, the system will transition to the next step.

1009-0357-000 09/07 3-3

## 3.2.5 Machine Check - Circuit 02

The *Machine Check-Circuit 02* check measures the  $0_2$ %.

- 1. Open the patient Y.
- 2. Set the Bag/Vent switch to Vent.
- 3. (ACGO option only.) Set the ACGO switch to Circle.
- 4. The display will show the  $O_2$ %. Do not select **Done** when 21 is first displayed. Allow the reading to stabilize, then select **Done**. Calibrate the  $O_2$  cell if necessary (measured reading outside 21%  $\pm$ 3%).

## 3.3 Individual Checks (for System software 3.X or greater)

Individual checks allow you to perform any combination of single checks. These checks are helpful if there is a specific problem/alarm and you want to test only that portion of the system.

The checks do not automatically move on to the next check.

#### **3.3.1 System**

The **System** check checks the Bag/Vent switch, proper gas supply pressures, ventilator operation and leak, battery and electrical power, circuit compliance, flow control operation, and vaporizer operation. This is a two-step check.

- 1. Set the Bag/Vent switch to Vent.
- 2. Open the patient Y.
- 3. (ACGO option only.) Set the ACGO switch to Circle.
- 4. Select **Start**. The display shows the checks being run.
  - The results are shown on the display.
- 5. Make sure the bellows is fully collapsed.
- 6. Occlude the patient Y.
- 7. Select **Continue**. The display shows the checks being run.
- 8. When the check passes, select **Back**.
- 9. Select another check or select **Start Case** to go to the **Start Case** menu.

3-4 09/07 1009-0357-000

#### 3.3.2 Circuit

The *Circuit* check checks the Bag/Vent switch, proper gas supply pressures, airway pressure measurement transducer, APL valve, and manual circuit leak.

- 1. Occlude the patient Y.
- Set Bag/Vent switch to Bag.
- 3. Set the APL valve halfway between 30 and 70.
- 4. (ACGO option only.) Set the ACGO switch to Circle.
- 5. Select **Start**. The display shows the checks being run.
- 6. When the check passes, select Back.
- 7. Select another check or select **Start Case** to go to the **Start Case** menu.

#### 3.3.3 Circuit 02 Cell

The **Circuit O2 cell** check measures the  $0_2$ %.

- 1. Open the patient Y.
- 2. Set the Bag/Vent switch to Vent.
- 3. (ACGO option only.) Set the ACGO switch to Circle.
- 4. The display will show the  $O_2$ %. Do not select **Done** when 21 is first displayed. Allow the reading to stabilize, then select **Done**. Calibrate the  $O_2$  cell if necessary (measured reading outside 21%  $\pm$ 3%).
- 5. Select another check or select **Start Case** to go to the **Start Case** menu.

#### 3.3.4 Low P Leak

The positive pressure **Low P Leak** check measures machine leaks before the breathing system and between the gas mixer and the common gas outlet. It measures low pressure pneumatic leaks with a pass/fail limit of 50 ml.

- 1. Occlude the inspiratory (right-hand) port.
- 2. Select Start.
- 3. The display shows the checks being run.
- 4. Open the inspiratory port and reconnect the breathing circuit.
- 5. Select another check or select **Start Case** to go to the **Start Case** menu.

## 3.3.5 Low P Leak (machines with ACGO)

The negative low P leak check measures machine leaks before the breathing system and between the gas mixer and the common gas outlet.

- 1. Make sure the ACGO switch is set to ACGO.
- 2. Attach the squeeze bulb to the ACGO outlet.
- 3. Squeeze (collapse) the bulb.
- 4. If the bulb inflates in less than 30 seconds, select **Fail**.
- 5. If the bulb remains collapsed, select Pass.
- 6. Remove the squeeze bulb from the ACGO outlet.

1009-0357-000 09/07 3-5

## 3.4 System "All checks" (for System software 2.X)

On the system "Checkout" menu, select **All Checks** and follow the instructions for "Low P leak check", "Quick check", "Vent check", and "Circuit 02 cell check".

If a check fails, follow the instructions on the display to perform a recheck or accept the results.

#### 3.4.1 Low P leak check

The low P leak check looks for leaks between the mixer, vaporizer, and the inspiratory side of the breathing circuit.

#### For machines with SCGO:

- 1. Ensure the vaporizers are turned off.
- 2. Plug the inspiratory (right-hand) port.
- 3. Select **Start**. The display shows the checks being run.
- 4. Repeat the check for each vaporizer with the vaporizer turned on.
- 5. When the checks pass, turn the vaporizer off.
- 6. Remove the plug from the inspiratory port.
- 7. Select **Next** to go to the next check.

#### For machines with ACGO:

- 1. Ensure the vaporizers are turned off.
- Set the ACGO switch to ACGO.
- 3. Attach the "negative low-pressure leak test" device to the ACGO outlet.
- 4. Collapse the bulb (squeeze).
- 5. If the bulb inflates in < 30 seconds, there is a leak.
- 6. If the bulb remains collapsed, repeat the check for each vaporizer with the vaporizer turned on.
- 7. When the checks pass, turn the vaporizer off. Select **Pass** to go the next check.

3-6 09/07 1009-0357-000

#### 3.4.2 Quick check

The quick check checks that the:

- Bag/Vent switch works in Bag position.
- Gas supply pressures are OK.
- Power cord is connected and the mains power is OK.
- Battery is fully charged.
- Manual circuit leak is OK.
- Flow controls operate correctly.
- 1. Occlude the patient Y piece.
- 2. Set the Bag/Vent switch to Bag.
- 3. Set the ACGO switch to Circle (ACGO option only).
- 4. Set the APL valve to approximately 50 (between the 30 and 70 marks).
- 5. Select **Start**. The display shows the checks being run.
- 6. When the checks pass, select **Next** to go to the next check.

#### 3.4.3 Vent check

The vent check measures circuit compliance and checks that the:

- Bag/Vent switch works in Vent position.
- Ventilator drive gas and O<sub>2</sub> pressure are OK.
- Ventilator circuit leak is OK.
- Ventilator delivers correctly.
- Alarms for ventilator failure or problems do not occur.
- 1. Set the Bag/Vent switch to Vent.
- 2. Open the patient Y piece.
- 3. Set the ACGO switch to Circle (ACGO option only).
- 4. Select **Start**. The display shows time remaining for check.
- 5. When check passes, follow the instructions on the display.
- 6. Make sure the bellows is collapsed.
- 7. Occlude the patient Y piece.
- 8. Select **Continue**. The display shows the checks being run.
- When the checks pass, the circuit compliance is displayed.
- 10. Select *Next* to go to the next check.

#### 3.4.4 Circuit O<sub>2</sub> cell check

The circuit  $O_2$  cell check measures the  $O_2$ %.

- 1. Unplug the patient Y piece.
- 2. Set the Bag/Vent switch to Vent.
- 3. Set the ACGO switch to Circle (ACGO option only).
- 4. The display will show the  $0_2$ %.
- Select *Exit* when measured O<sub>2</sub> is stable.
   O<sub>2</sub> cell calibration is recommended if displayed value is less than 21%.

1009-0357-000 09/07 3-7

### 3.5 Bellows drop test

- 1. End a case.
- 2. Set the Bag/Vent switch to Vent.
- 3. Occlude the patient Y piece.
- 4. Push the **02 Flush** button until the bellows is full.
- 5. After the initial drop, if the bellows falls more than 100 ml/min, it has a leak.

## 3.6 Backlight test

- 1. Push the **Main Menu** key.
- 2. Select Calibration.
- 3. Select Backlight Test.
- 4. Select Start Test.
- 5. The display will show the test running on light 1 and then on light 2. If the display goes completely blank or flickers during the test, one of the lights has failed.

### 3.7 Vaporizer back pressure test

#### WARNING

Anesthetic agent comes out of the circuit during this test. Use a safe, approved procedure to collect and remove the agent.

- 1. Set the System switch to On.
- 2. Start a case.
- 3. Set the  $O_2$  flow to 6 l/min.
- 4. Slowly adjust the vaporizer concentration from 0 to 1%.
  - Make sure that the O<sub>2</sub> flow stays constant.
  - Verify that the system continues to operate without issuing any related alarms.
- 5. Repeat the test for both vaporizer positions.

3-8 09/07 1009-0357-000

### 3.8 Pipeline and cylinder tests

- 1. Connect the pipeline supplies one at a time and ensure that the corresponding display indicates pipeline pressure.
- 2. Disconnect all pipeline supplies.
  - a. Open each cylinder valve.
  - b. Make sure that each cylinder has sufficient pressure. If not, close the applicable cylinder valve and install a full cylinder.
- 3. Test the cylinder supplies for a high pressure leak. Make sure that each cylinder has sufficient pressure:
  - a. If equipped, turn the auxiliary O<sub>2</sub> flow control fully clockwise (no flow).
  - b. If equipped, turn off venturi derived suction.
  - c. Open each cylinder.
  - d. Record the cylinder pressure.
  - e. Close each cylinder valve.
  - f. Record the cylinder pressure after one minute. If the pressure decreases more than indicated below, there is a leak.

5000 kPa (725 psig) for ventilator drive gas.

690 kPa (100 psig) for non ventilator drive gas.

If a cylinder supply fails this test, install a new cylinder gasket and do this step again.

4. Close all cylinder valves.

#### **⚠** WARNING

Do not leave gas cylinder valves open if the pipeline supply is in use. Cylinder supplies could be depleted, leaving an insufficient reserve supply in case of pipeline failure.

### $3.8.10_2$ supply alarm test

- 1. Establish O<sub>2</sub>, Air, and (if equipped) N<sub>2</sub>O gas supplies.
- 2. Set  $O_2$  to 25% and (if equipped)  $N_2O$  as balance gas. For machines without  $N_2O$ , set Air as balance gas.
- 3. Set total flow to 3 L/min.
- Stop the O<sub>2</sub> supply. (Disconnect the pipeline supply or close the cylinder valve.)
- 5. Make sure that:
  - a. The low "O<sub>2</sub> supply pressure low" alarm occurs.
  - b. The  $N_2O$  (if equipped) and  $O_2$  flows stop.
  - c. Air (if selected) flow continues or an Air selection prompt appears.
- 6. Reconnect the  $O_2$  supply.

1009-0357-000 09/07 3-9

#### 3.9 Pressure relief tests

To check the pressure relief valve in the vaporizer manifold outlet.

#### For machines with SCGO:

- 1. Remove the back cover to access the vaporizer manifold.
- 2. Remove the outlet tubing and connect a test device (pressure gauge or a digital manometer) to the vaporizer manifold outlet.
- 3. Adjust the  $O_2$  flow to 0.5 L/min.
- 4. Verify that the test device reading stabilizes within the following range: **31–60 kPa (230–450 mm Hg) (4.5–8.5 psi)**.
- 5. Remove the test device and reconnect the outlet tubing.
- 6. Replace the back cover.

#### For machines with ACGO:

- 1. Set the ACGO selector switch to ACGO.
- 2. Connect a test device (pressure gauge or a digital manometer) to the ACGO outlet using the positive pressure leak test adapter.
- 3. Adjust the  $O_2$  flow to 0.5 L/min.
- 4. Verify that the test device reading stabilizes within the following range:
  - 31-60 kPa (230-450 mm Hg) (4.5-8.5 psi).
- 5. Remove the test device and the adapter.




3-10 09/07 1009-0357-000

#### 3.10 Flush Flow Test

- 1. With Bag/Vent switch in Bag, verify case has ended.
- 2. Set the Bag/Vent switch to Vent.
- 3. Attach a patient circuit and plug the patient port.
- 4. For ACGO equipped machines, set the ACGO selector switch to Circle.
- 5. Ensure that the bellows is completely collapsed.
- 6. Measure the amount of time it takes to fill the bellows when the O<sub>2</sub> Flush button is fully and continuously depressed.
- 7. Repeat the above measurement two more times (deflate bellows by removing the plug from the patient port).
  - The bellows should fill in 1.8 to 2.3 seconds.

#### **Possible Causes of Failure**

- Large leak (if long filling time).
- Flush regulator setting (Section 5.2).
- Flush regulator cross-connection (if long filling time).
- SCGO/ACGO selector valve inlet cross-connection (if short filling time).



1009-0357-000 09/07 3-11

### 3.11 Alarm tests

4. S

**NOTE**: If an Airway Gas Module is installed, *FiO2* readings are taken from the module instead of the O<sub>2</sub> sensor in the breathing circuit. When using an Airway Gas Module, a sample line must be connected to the patient circuit for testing the O<sub>2</sub> alarms.

- 1. Connect a test lung to the patient connection.
- Start a case.
- 3. Set the Bag/Vent switch to Vent.
- 4. Set the  $O_2$  concentration to 30%, and allow the  $O_2$  reading to stabilize.
- 5. Test the  $0_2$  alarms:
  - Set the Fi02 low alarm limit to 50%. Make sure an Fi02 low alarm occurs.
  - Set the Fi02 low alarm limit back to 21% and make sure that the Fi02 low alarm cancels.
  - Set the Fi02 high alarm limit to 50%.
  - Push the O<sub>2</sub> flush button.
  - Make sure the Fi02 high alarm occurs.
  - Set the Fi02 high alarm limit back to 100%. Make sure that the Fi02 high alarm cancels.
- 6. Test the **MVexp low** alarm:
  - Go to the *Alarm Setup* menu.
  - Set the **MV low** alarm limit to greater than the measured minute volume.
  - Make sure that a **MVexp low** alarm occurs.
  - Set the MV low alarm limit to off.
- 7. Test the **Ppeak high** alarm:
  - Set the **Pmax** to less than the peak airway pressure.
  - Make sure that the **Ppeak high** alarm occurs.
  - Set the **Pmax** to the desired level.
- 8. Test the **PEEP high. Blockage?** alarm:
  - Close the APL valve.
  - Set the Bag/Vent switch to Bag. Mechanical ventilation stops.
  - Block the patient connection and push the O<sub>2</sub> flush button.
  - Make sure that the **PEEP high. Blockage?** alarm occurs after approximately 15 seconds.
- 9. Test the **Ppeak low. Leak?** alarms:
  - Unblock the patient connection.
  - Set the Bag/Vent switch to Vent.
  - Set the tidal volume and total flow to minimum.
  - Other alarms such as MVexp low can occur.
  - Make sure that the **Ppeak low. Leak?** alarms occur.

10. Set all alarm limits to approved clinical values.

3-12 09/07 1009-0357-000

## 3.12 Alternate 0<sub>2</sub> flowmeter tests

- 1. Open the  $O_2$  cylinder valve or connect an  $O_2$  pipeline.
- 2. Rotate the Alt O<sub>2</sub> flow control fully clockwise to minimum flow.
- Press the Alternate O<sub>2</sub> switch to turn on Alternate O<sub>2</sub> flow.
   The flowmeter should indicate 0.5 to 0.7 L/min.
- 4. Rotate the flow control counterclockwise (increase). The ball should rise immediately after rotation is begun. It should rise smoothly and steadily with continued counterclockwise rotation. When a desired flow is set, the ball should maintain in a steady position.
- 5. Rotate the flow control clockwise to minimum flow.
- 6. Press the Alternate  $O_2$  switch to turn off Alternate  $O_2$  flow; confirm yes.

## 3.13 Auxiliary 0<sub>2</sub> flowmeter tests

- 1. Open the  $O_2$  cylinder valve or connect an  $O_2$  pipeline.
- 2. Rotate the flow control clockwise (decrease) to shut off the flow. The ball should rest at the bottom of the flow tube and not move.
- Rotate the flow control counterclockwise (increase). The ball should rise immediately after rotation is begun. It should rise smoothly and steadily with continued counterclockwise rotation. When a desired flow is set, the ball should maintain in a steady position.
- 4. Occlude the auxiliary O<sub>2</sub> outlet. The ball should rest at the bottom of the flow tube and not move. A ball that does not rest at the bottom of the flow tube indicates a leak and requires service.
- 5. Rotate the flow control clockwise to shut off the flow.

### 3.14 Integrated Suction Regulator tests

The gauge needle should come to rest within the zero range bracket when no suction is being supplied. Gauges which do not comply may be out of calibration.

- 1. Adjust the regulator setting to minimum.
- 2. Turn the mode selector to I (On).
- 3. Ensure the gauge remains less than 200 mmHg (26 kPa, 0.26 Bar).
- 4. Occlude the inlet.
- 5. Ensure the gauge remains less than 200 mmHg (26 kPa, 0.26 Bar).
- 6. Adjust the regulator in an increasing vacuum level.
- 7. The gauge should rise after rotation has begun. The gauge should rise with continued rotation of the regulator adjustment.
- 8. Adjust the regulator setting to minimum.
- 9. Turn the Mode selector to O (Off).

1009-0357-000 09/07 3-13

#### 3.15 Power failure test

1. Connect the power cord to a wall outlet. The mains indicator on the front panel comes on when AC Power is connected.



- 2. Set the system switch to On and Start a case.
- 3. Unplug the power cord with the system turned on.
- 4. Make sure that the power failure alarm comes on.
- 5. Make sure the following message is displayed:
  - Plug in power cable. On battery
- 6. Connect the power cable again.
- 7. Make sure the alarm cancels.

## 3.16 Electrical safety tests

Make sure the system is completely assembled and all accessory devices are connected to electrical outlets.

1. Connect an approved test device (e.g. UL, CSA, or AAMI) and verify that the leakage current is less than:

Voltage	Max. Leakage Current
120/100 Vac	300 μAmps
220/240 Vac	500 μAmps

2. Make sure that the resistance to ground is less than  $0.2\Omega$  between an exposed metal surface and the ground pin on the power cord.

3-14 09/07 1009-0357-000

## **4a Installation and Service Menus (DU)**

In this section	4a.1 Service and Installation menu structure	4a-2
	4a.2 Install/Service Menu (Super User)	4a-3
	4a.2.1 Colors Menu	4a-4
	4a.2.2 Units Menu	4a-4
	4a.2.3 Factory Defaults	4a-5
	4a.3 Installation Menu	4a-6
	4a.3.1 Configuration	4a-7
	4a.3.2 Units Menu	4a-8
	4a.3.3 Options Key	4a-8
	4a.3.4 Copy Configuration	4a-9
	4a.4 Service Menu	4a-10
	4a.4.1 Software/Hardware Ver Menu	4a-11
	4a.4.2 Service Log Menu	4a-12
	4a.5 Calibration	4a-13
	4a.5.1 User Calibration menu	4a-13
	4a.5.2 Manifold P Span	4a-14
	4a.5.3 Insp Flow Zero	4a-15
	4a.5.4 Inspiratory Flow Valve	4a-16
	4a.5.5 Bleed Resistor	4a-18
	4a.5.6 Paw Span	4a-19
	4a.5.7 Zero Gas Xducrs	4a-20
	4a.5.8 Cal Config	4a-21
	4a 5 9 Mixer P 7ero	4a-21

1009-0357-000 09/07 (DU) 4a-1

#### 4a.1 Service and Installation menu structure

This section describes the Service level functions that are part of the main software installed in the anesthesia machine.

Section 8a, "Service Diagnostics and Software Download (DU)," covers a separate service application that loads from a PCMCIA card and is used to download system software and run service diagnostics and other service tests.

#### Menu structure

The Service menu structure has three levels which are password protected:

- Install/Service (super-user)
- Installation
- Service


The **Install/Service** level (super-user password) supports standard hospital preferences: choosing units; setting ventilator, alarm, and gas delivery defaults.

The **Installation** level requires the service password and supports language, gas color codes, flow tube position, country, hardware flags for system components (acgo or scgo etc.), enabling software options, and cloning a system.

The **Service** level requires the service password and supports diagnostic tools and automated component tests.

Follow the menu structure to access the various service screens:

- on the Main Menu, select Screen Setup;
- on the **Screen Setup** menu, select **Install/Service** to access the Install/Service (with super-user password) menu;
- on the Install/Service menu, select Installation (with service password) to access the Installation menu.
- to access the Service menu, select Service (with service password) on the Install/Service menu; or, from the Installation menu, select Service to access the same Service menu without having to enter the service password.




(DU) 4a-2 09/07 1009-0357-000

## 4a.2 Install/Service Menu (Super User)

Use the super-user password to access the Install/Service menu: "16-4-34"

Menu Item	Message text	Comments
Colors	Set colors of parameters.	
Units	Set units of weight, CO2, gas supply pressure, and Paw.	
Show Alarm Limits	Select yes to show alarm limits in digit fields.	Default is Yes
Save Default Case	Save normal screen, air/ N20, circuit type, ventilator settings, and alarm limits from the last case as defaults	Last used alarm settings (including hide/show alarm limits), screen layout (middle waveform selection, sweep speed, scaling), ventilator mode and setting, balance gas, and start case gas outlet selection are saved as facility defaults.  Note: The Pmax alarm limit shall not be saved higher than 40 cmH20. The low FiO2 alarm limit shall not be saved lower than 21%.
Factory Defaults	Return to default factory settings.  After selecting Factory Defaults:  "Reset machine for defaults to take effect."	Action: All facility defaults get replaced with factory defaults. Super User settings also get set to Factory Defaults. No Service level configuration settings are changed.
Installation	Set language, gas colors, hardware, and enable options.	Navigate with password to Installation menu. Password is "26-23-8"
Service	Show technical data for troubleshooting and calibration.	Navigate with password to Service menu. Password is "26-23-8"
Exit	Turn power off to exit the service and super user menus.	



1009-0357-000 09/07 (DU) 4a-3

#### 4a.2.1 Colors Menu

Menu Item	Message text	Values
Paw	Change color of Paw waveform, digits and trend.	Yellow, White, Green, Red, or Blue
Flow	Change color of Flow waveform, Flow and Volume digits and trends.	Yellow, White, Green, Red, or Blue
Resp	Change color of respiration, digits and trend.	Yellow, White, Green, Red, or Blue
C02	Change color of CO2 waveform, digits and trend.	Yellow, White, Green, Red, or Blue
Previous Menu	Return to previous menu.	

### 4a.2.2 Units Menu

The Units menu can be accessed here in the super-user level to change individual preferences, or if required during installation, in the service level Installation menu.

Menu Item	Message text	Values
Weight	Change weight unit: kg or lb.	kg or lb
C02	Change CO2 unit: %, kPa, or mmHg.	%, kPa, or mmHg;
Gas Supply Pressure	Change gas supply pressure unit: kPa, psi, or bar.	psi, kPa, or bar
Paw	Change Paw unit: kPa, hPa, cmH2O, mmHg, mbar.	kPa, hPa, cmH2O, mmHg, or mbar
Previous Menu	Return to previous menu.	

(DU) 4a-4 09/07 1009-0357-000

## **4a.2.3 Factory Defaults** The following table lists the factory defaults for parameters and alarm limits:

Parameter	Value
Vent Mode	VCV
TV (tidal volume)	500 ml
Pinsp	5 cmH20
	(5 hPa, 0.5 kPa, 5 mbar, 3.6
	mmHg)
RR	12 /min
Mech RR	12 /min
Tinsp	1.70
I:E	1:2.0
Trig. Window	25%
Flow Trig.	2 I/min
End Breath	25%
Psupport	Off
PEEP	Off
Tpause	Off
Backup Time	30 s
02%	100%
Balance Gas	Air
Circuit	Circle
Gas Outlet installed	SCG0
Paw Color	Yellow
Flow Color	Green
Resp Color	White
CO2 Color	White
Paw Units	cmH2O
Weight Units	kg
Temperature Units	С
CO2 Units	%
Altitude	300 m
Gas Supply Pressure	kPa
Units	- Alak
Decimal marker	. <dot></dot>
Language	English
Gas Supply Colors	180 (02 - white N20 - blue
	(02 = white, N20 = blue, Air = blk/wht)
02 flow tube	Right side
Vent drive gas	02
PSV Pro	Enabled
SIMV/PSV	Enabled
PCV	Enabled
SIMV-PC	Enabled
N20 enabled	No

Alarm Limit	Value
Pmax High	40 cmH20
	(40 hPa, 4 kPa, 40 mbar,
	30 mmHg)
MV High	10 l/min
MV Low	2 I/min
TVIII d	4000
TV High	1000 ml
TV Low	Off
RR High	Off
RR Low	Off
Et CO2 High	6.5%
Ft 000 I	(50 mmHg or 6.5kPa)
Et CO2 Low	Off
Fi CO2 High	Off
Fi O2High	Off
Fi O2 Low	21%
EtO2 High	Off
EtO2 Low	Off
Fi Iso High	5%
Fi Iso Low	Off
Et Iso High	Off
Et Iso Low	Off
Fi Sev High	8%
Fi Sev Low	Off
Et Sev High	Off
Et Sev Low	Off
Fi Des High	15%
Fi Des Low	Off
Et Des High	Off
Et Des Low	Off
Fi Enf High	5%
Fi Enf Low	Off
Et Enf High	Off
Et Enf Low	Off
Fi Hal High	5%
FI Hal Low	Off
Et Hal High	Off
Et Hal Low	Off
_	

1009-0357-000 09/07 (DU) 4a-5

## 4a.3 Installation Menu

Use the service-level password to access the Installation menu:  $\mbox{``26-23-8.''}$ 

Whenever the installation menu is entered, "Enter Service dd-mmm-yyyy hh:mm:ss" is recorded in the Event log.

Menu Item	Message text
Configuration	Set language, gas color code, and O2 flowmeter position.
Units	Set units.
Options Key	Enable software options.
Options List	Display software options.
Copy Config	Normal Message "Save or install configuration and default settings using memory card."  Blocked Message "Please insert memory card."
Comitos	,
Service	Show error, event, and alarm logs. (Accessing the Service menu from the Installation menu does not require second use of service password.)
Previous Menu	Return to previous menu.

Configuration
Decimal Marker
Language
Gas supply Colors
02 Flowtube
Ventilator Drive Gas
Altitude
Gas Outlet
N20 Enabled

Units Menu
Weight
CO2
Gas Supply Pressure
Paw

Options Key
Current Key
Entry 1
Entry 2
Entry 3
Entry 4
Entry 5
Entry 6
Entry 7
Save New Key
Control Board ID

Options List
Available Options SIMV/PSV PCV PSV Pro

Copy Configuration
Save to Card Copy from Card

(DU) 4a-6 09/07 1009-0357-000

### 4a.3.1 Configuration

Menu Item	Message text	Values	Comments
Decimal Marker	Select decimal delineator.	0.01, 0 01 or 0,01	
Language	Change language translation of screen texts.	English, French, German, Spanish, Italian, Japanese, Portuguese, Dutch, Chinese (simplified), Finnish, Norwegian, Hungarian, Polish, Greek, Czech, Turkish, and Russian.	Default: English
Gas supply Colors	Change color of O2, N2O, and Air.	ANSI, ISO, Neutral	ANSI: 02 green, Air yellow, N20 blue; ISO: 02 white, Air black/white, N20 blue; Neutral: All gases white.
02 Flowtube	O2 on left or right- hand side.	Left, Right	
Ventilator Drive Gas	Change drive gas to match machine configuration.	Air, 02	
Altitude	Change altitude used for gas calculations.	-400 to 3000 m in 100-m increments	
Gas Outlet*	Change type of fresh gas outlet.	SCGO, ACGO	SCGO: Use insp port. ACGO: Use auxiliary port.
N20 Enabled	Change to match machine configuration.	Yes; No	

<sup>\*</sup> For machines without a separate auxiliary common gas outlet and selector switch, set **Gas Outlet** to SCGO: Selectable Common Gas Outlet.

1009-0357-000 09/07 (DU) 4a-7

<sup>\*</sup> For machines with an external auxiliary common gas outlet and selector switch, set **Gas Outlet** to ACGO: Auxiliary Common Gas Outlet

#### 4a.3.2 Units Menu

Menu Item	Message text	Values
Weight	Change weight unit: kg or lb.	kg or lb
CO2	Change CO2 unit: %, kPa, or mmHg.	%, kPa, or mmHg
Gas Supply Pressure	Change gas supply pressure unit: kPa, psi, or bar.	psi, kPa, bar
Paw	Change Paw unit: kPa, hPa, cmH2O, mmHg, mbar.	kPa, hPa, cmH2O, mmHg, or mbar

#### 4a.3.3 Options Key

The Options Key menu is used to configure the software to include the features that the customer has purchased. The included features are shown in the Options List menu.

#### **Options Key menu**

Menu Item	Message text	Values
Current Key	Enter key code to enable options.	XXXAXBC
Entry 1	Enter first entry of key-code.	0 to 9, A to Z, ~,
Entry 2	Enter second entry of key-code.	!, @, #, \$, %, ^, *, (,), ?
Entry 3	Enter third entry of key-code.	
Entry 4	Enter fourth entry of key-code.	
Entry 5	Enter fifth entry of key-code.	
Entry 6	Enter sixth entry of key-code.	
Entry 7	Enter seventh entry of key-code.	
Save New Key	Confirm entries for key-code.	
Control Board ID	Control number used by key-code.	XXX

When options are added, "Add <option> dd-MMM-yyy hh:mm:ss" is written to the event log.

If more than one option is added, each option is be listed separately.

(DU) 4a-8 09/07 1009-0357-000

#### **Options List menu**

The options list shows which options are enabled.

Menu Item	Message text	Values
Available Options		
SIMV/PSV	SIMV vent w/pressure support.	On, Off
PCV	Pressure controlled ventilation.	On, Off
PSV Pro	Pressure support ventilation w/backup.	On, Off

## 4a.3.4 Copy Configuration

#### **Copy Configuration menu**

Menu Item	Message text	Values	Comments
Save to Card	Save Configuration and defaults to card.	   The field is blank until the data has either been written to the card (OK) or the system determines it cannot write to the card (Fail).	Saves all settings that are not hardware dependent, including facility defaults, colors, units, 02 flow tube position, decimal marker, and altitude.
Copy from Card	Copy Configuration and defaults from card.  When completed: Copy from card complete. Please reboot system.	    The field is blank until the data has either been read from the card (OK) or the system determines it cannot read the card or the card does not have the required data (Fail).	

Systems cannot accept configuration files from a different product model.

The software version is stored with the saved configuration. A system will reject any configurations from other than the current version of software.

Selecting Save to Card overwrites any configuration on the card.

1009-0357-000 09/07 (DU) 4a-9

### 4a.4 Service Menu

Use the service-level password to access the Service menu:  $^{\circ}26-23-8.^{\circ}$ 

Whenever service menu is entered, "Enter Service dd-mmm-yyyy hh:mm:ss" is recorded in the Event log.

Menu Item	Message text
SW HW versions	Scroll through system information.
Service Log	Show error, event, and alarm histories.
Calibration	Push ComWheel to perform service calibrations.
Previous Menu	Return to previous menu.

Service Log Men
Scroll Recent
Error History Event History Alarm History Copy Logs Reset Logs Previous Menu

Calibration
Instructions
User Calibration
Manifold P Span
Insp Flow Zero
Insp Flow Valve
Bleed Resistor
Paw Span
Zero Gas Xducer
Cal Config
Previous Menu

(DU) 4a-10 09/07 1009-0357-000

# 4a.4.1 Software/ Hardware Ver Menu

Turn the ComWheel to scroll through the list box.

Push the ComWheel to return to the Service menu.

# **System Information menu**

List box text with X=Number, A, B, C = letter
Total Time: XXXXX (Minutes)
Software Release: XX.XX
Model Code: XXX
Machine Serial Number: ABCDXXXXX
Option Package: XXX
Options Code: XXXXX
Anes Software Version: XX.XX
Anes Hardware Version: XXXX-XXXX REV A
Anes Board Serial Number: ABCXXXXX
Disp Software Version: XX.XX
Disp Hardware Version: XXXX-XXXX REV A
Disp Hardware Serial Number: ABCXXXXX
Mixer Software Version: XX.XX
Mixer Hardware Version: XXXX-XXXX-XXX REV A
Mixer Board Serial Number: ABCXXXXX
Mixer O2 Flow Sensor Serial Number: XXXXXXXXX
Mixer Balance Gas Flow Sensor Serial Number: XXXXXXXXX
PCA Serial Number: ABCXXXXX
Vent Software Version: XX.XX
Vent Hardware Version: XXXX-XXXX REV A
Vent Intf Board Serial Number: ABCXXXXX
Power Software Version: XX.XX
Power Hardware Version: XXXX-XXXX REV A
Power Board Serial Number: ABCXXXXX
MGas Software Version: X.X
MGas Hardware Version: GAS SW Pr. XXXXXXX-X
MGas Hardware Serial Number: ABCXXXXX

The MGas information is only displayed when an Airway module is present.

#### 4a.4.2 Service Log Menu

The Service log menu is an organized listing of stored events.

Menu Item	Message text
Scroll Recent	Scroll through newest entries.
Error History	Show error history.
Event History	Show event history.
Alarm History	Show alarm history.
Copy Logs	Copy logs to PCMCIA card. Takes about 1 minute.
Reset Logs	Erase Error and Alarm log entries
Previous Menu	Return to previous menu.

Each history log shows at the top of the screen the total "Running Hours" and the date when the logs were last reset. The running hours number is the same number as the operating hours shown on the system status page during normal operation.

Whenever logs are reset, "Reset Logs dd-MMM-yyy hh:mm:ss" is recorded in the Event log.

If the logs are saved to a memory card, the machine's serial number is saved along with the current contents of the logs and the date and time.

#### **Error History** ◀

The Error History log lists the last 200 errors logged since the last log reset, starting with the most recent. The system stores the last 1,000 errors logged since the last log reset.

#### **Event History ◀**

The Event History log records the service history of the device. This includes: service calibrations, entry into the service mode, options enabled, and software installation. In the event of a board replacement, it is understood that this log like all others could be lost.

The Event History menu lists the last 200 events logged starting with the most recent. The Event History log stores the last 1000 events.

The Event History log cannot be reset.

#### Alarm History ◀

The Alarm History log lists the last 200 medium and high priority parameter alarms since the last log reset starting with the most recent. The Alarm History log store the last 1000 entries.

#### **Copy Logs**

The Copy Logs function copies Error, Event, and Alarm logs along with the software/hardware configuration to a text file on a PCMCIA card. The copying takes about one minute.

**Note**: Do not remove the Flash Card until the screen shows copy is complete.

(DU) 4a-12 09/07 1009-0357-000

# 4a.5 Calibration

For step-by-step instruction, refer to Section 5.4, "Ventilator Calibrations."

Menu Item	Message text
Instructions	
User Calibration	Show the normal user calibration menu.
Manifold P Span	Calibrate manifold pressure transducer.
Insp Flow Zero	Zero inspiratory flow valve.
Insp Flow Valve	Calibrate inspiratory flow valve.
Bleed Resistor	Calibrate bleed resister flow.
Paw Span	Calibrate the airway pressure transducer.
Zero Gas Xducrs	Calibrate the gas supply transducers.
Cal Config	Set vent drive gas and altitude.
Previous Menu	Return to previous menu.

# 4a.5.1 User Calibration menu

Menu Item	Message text
Flow and Pressure	Calibrate the flow and pressure sensors.
Circuit 02 Cell	Calibrate Circuit O2 Cell.
Airway Gas	Start Gas Calibration. Calibrate CO2, O2, N2O, and agent measurements.
Backlight Test	Push ComWheel to test back lights. Test every month.
Previous Menu	Return to the previous menu.

# 4a.5.2 Manifold P Span

The Manifold P Span instructions appear when the focus is on Manifold P Span menu item.

Refer to Section 5.4.2, "Manifold P Span."

#### Instructions

Read all steps before you start:

- 1. Remove the breathing system, the exhalation valve, and the metal plate.
- 2. Put #2 plugs in the manifold and the drive gas ports of the vent engine.
- 3. Connect a pressure gauge in line with the manifold pressure transducer.
- 4. Push the ComWheel to continue.
- 5. Select Start Manifold P Span.
- 6. Increase the Flow valve setting until the gauge shows 100 cmH20.
- 7. When the gauge shows 100 cmH2O, select Save Calibration.

#### **Manifold P Span menu**

Menu Item	Message text	Values/ Comments
Insp Flow Valve (DAC)	Increase setting until test gauge shows 100 cmH20 (approx 1020 counts). Then save calibration.	O to 4095  Disabled until user selects Start Manifold P Span.
Start Manifold P Span	Start Calibration. Increase flow valve setting until test gauge = 100 cmH20 (approximately 1020 counts). Then save calibration.  Blocking message: "Connect a supply of the drive gas to continue."	Blocked when the ventilator drive gas supply pressure would cause a gas supply failure alarm during normal operation.
Save calibration	Save Manifold P Span calibration.	
Previous Menu	Return to the previous menu.  During calibration:  "Calibration in progress. Push ComWheel to cancel."	

(DU) 4a-14 09/07 1009-0357-000

# 4a.5.3 Insp Flow Zero

The Insp Flow Zero instructions appear when the focus is on the Insp Flow Zero menu item.

Refer to Section 5.4.4, "Insp Flow Zero."

#### Instructions

Read all steps before you start:

- 1. Push the ComWheel to start the zero check
- 2. No disassembly is required.
- If the outcome of the calibration is Pass, the new calibration data is saved.
- If the outcome is Fail, the old calibration data is retained.
- The result of the calibration is saved to the Event Log.

Selecting Previous Menu before the calibration is done aborts the calibration in progress and keeps the old calibration constants.

#### **Insp Flow Zero menu**

Menu Item	Message text	Values
Start	If the result is failed, do the insp flow valve calibration.	Pass or Fail
Previous Menu	Return to previous menu.  During Calibration:  "Calibration in progress. Push ComWheel to cancel."	

# 4a.5.4 Inspiratory Flow Valve

The Inspiratory Flow Valve instructions appear when the focus is on the Insp Flow Valve menu item.

Refer to Section 5.4.3, "Inspiratory Flow Valve Cal."

#### Instructions

Read all steps before you start:

- 1. Complete the Manifold P Span calibration.
- 2. Put #2 plugs in the manifold and the drive gas ports of the vent engine.
- 3. Push the ComWheel to show the next menu.
- 4. Select Stage 1 calibration.
- 5. After Pass, replace the manifold port plug with the calibration orifice.
- 6. Select Stage 2 calibration.
- 7. You MUST do both stages for the calibration to be saved.

During calibration, a separate menu shows the counts and corresponding flow at each step.

- If the outcome of both stages of the calibration is Pass, the new calibration data is saved.
- If the outcome of either stage is Fail, the old calibration data is retained.
- The results of each stage of the calibration are saved to the Event Log.

Selecting Previous Menu before the calibration is done, aborts the calibration in progress and keep the old calibration constants.

(DU) 4a-16 09/07 1009-0357-000

## **Insp Flow Valve Menu**

Menu Item	Message text	Comments
Stage 1	Calibrate the insp flow valve at low flows.  Blocked text:  "Connect a supply of the drive gas to continue."	Blocked if the ventilator drive gas supply pressure would cause a gas supply alarm during normal operation.
Stage 2	Calibrate the insp flow valve at high flows.  Blocked text: "Stage 1 calibration is required first."	Blocked if Stage 1 has not been completed.
Insp Flow Valve Data	Show insp flow valve calibration table.	
Previous Menu	Return to previous menu.  During Calibration:  "Calibration in progress. Push ComWheel to cancel."	

## **Insp Flow Valve Data menu**

The Insp Flow Valve Data menu contains a table of 24 entries from the previous calibration. The table is erased at the start of Stage 1. The table is updated in real time during the calibration.

#### 4a.5.5 Bleed Resistor

The Bleed Resistor instructions appear when the focus is on the Bleed Resistor menu item.

Refer to Section 5.4.5, "Bleed Resistor Cal."

#### Instructions

Read all steps before you start:

- 1. Complete the Insp Flow Valve calibration.
- 2. Put #2 plugs in the manifold and the drive gas ports of the vent engine.
- 3. Push the ComWheel to show the next menu.
- 4. Select Start.

The calibration fails if the flow required to reach 91 cmH20 is > 16 l/min.

- If the outcome of the calibration is Pass, the new calibration data is saved.
- If the outcome is Fail, the old calibration data is retained.
- The result of the calibration is saved to the Event Log.

Selecting Previous Menu before the calibration is done aborts the calibration in progress and keep the old calibration constants.

#### **Bleed Resistor Menu**

Menu Item	Message text	Comments
Start	Calibrate manifold pressure to bleed resistor flow.	Blocked if the ventilator drive gas supply pressure would
	Blocked text: "Connect a supply of the drive gas to continue."	cause a gas supply failure alarm during normal operation.
Bleed Resistor Data	Show bleed resistor calibration table.	
Previous Menu	Return to the previous menu.	
	During Calibration: "Calibration in progress. Push ComWheel to cancel."	

#### **Bleed Resistor Data menu**

The Bleed Resistor Data menu contains a table of 17 entries from a previous calibration. The table is erased at the start of the calibration. The table is updated in real time during the calibration.

(DU) 4a-18 09/07 1009-0357-000

# 4a.5.6 Paw Span

The Airway P Span instructions appear when the focus is on Paw Span menu item.

Refer to Section 5.4.6, "Paw Span."

#### Instructions

Read all steps before you start:

- 1. Complete the Bleed Resistor calibration.
- 2. Install the flow sensor and circuit module
- 3. Put a #2 plug in the drive gas port of the vent engine.
- 4. Put the calibrated orifice in the manifold port of the vent engine.
- 5. Connect the pressure tee to the insp port. Connect the tee to the calibrated orifice with a 22 mm tube.
- 6. Connect a pressure gauge to the pressure tee.
- 7. Select Start Paw Span.
- 8. Increase the flow valve setting until the gauge shows 100 cmH20.
- 9. Select Save Calibration

#### **Paw Span menu**

Menu Item	Message text Values/Commo	
Insp Flow Valve (DAC)	Increase setting until test gauge shows 100 cmH2O (approximately 1020 counts). Then save calibration.	0 to 4095 (initially set to 800)
Start Paw Span	Start Calibration. Increase flow valve setting until test gauge = 100 cmH20 (approximately 1020 counts). Then save calibration.  Blocked text: Connect a supply of the drive gas to continue.	Blocked if the ventilator drive gas supply pressure would cause a gas supply failure alarm during normal operation.
Save calibration	Save Paw Span calibration.	Saves new calibration data. Writes calibration result, date and time to the event log.
Previous Menu	Return to the previous menu.  During Calibration:  "Calibration in progress. Push ComWheel to cancel."	

## 4a.5.7 Zero Gas Xducrs

The Zero Gas Xducrs instructions appear when the focus is on the Zero Gas Xducer menu item.

#### Instructions

Read all steps before you start:

- 1. Remove all cylinders.
- 2. Disconnect all pipeline supplies.
- 3. Select Zero Gas Xducers.
- 4. Select Start Zero on the next menu.

This page also shows:

- Gas supply pressures
- Gas supply ID

A failed test is usually the result of a pipeline or cylinder still connected to the system.

- If the outcome of the calibration is Pass, the new calibration data is saved.
- If the outcome is Fail, the old calibration data is retained.
- The result of the calibration is saved to the Event Log.

Selecting Previous Menu before the calibration is done aborts the calibration in progress and keep the old calibration constants.

#### Zero Gas Xducrs menu

The Zero Gas Xducrs menu shows only transducers that are installed. If not installed, the menu row is blank.

Menu Item	Message text	Values
02 Pipeline		0-4095 Counts
02 Cylinder 1		0-4095 Counts
02 Cylinder 2		0-4095 Counts
N20 Pipeline		0-4095 Counts
N20 Cylinder		0-4095 Counts
Air Pipeline		0-4095 Counts
Air Cylinder		0-4095 Counts
Start Zero	Disconnect all pipelines and remove cylinders. Then select Start Zero.	
Previous Menu	Return to the previous menu. During Calibration: "Calibration in progress. Push ComWheel to cancel."	

(DU) 4a-20 09/07 1009-0357-000

#### 4a.5.8 Cal Config

Before calibration, you must verify that the Ventilator Drive Gas and the Altitude settings are set appropriately to match the current drive gas configuration and machine location.

If you change any of the settings in the Cal Config menu, you must restart the system.

#### **Cal Config menu**

Menu Item	Message text	Values
Ventilator Drive Gas	Change drive gas to match machine configuration.	Air, 02
Altitude	Change altitude used for gas calculations.	-400 to 3000 m (in 100-m increments)

#### 4a.5.9 Mixer P Zero

(for System software 3.X; previous versions do not include this functiin)

The Mixer P Zero instructions appear on the Mixer P Zero menu.

#### Instructions

Back	Default	Start	
To go back to fa	actory defaults, select C	Defaults (above).	
To Zero Pres Se	ensors:		
1. Disconnect p	oipeline gas supplies.		
2. Close the ga	s cylinders.		
3. Remove the	flow sensors.		
4. Insert a non-	-Des cassette.		
5. Push the 02	Flush button for 3 seco	inds.	
6. Let the syste	em sit WITHOUT gas flow	for > 5 min.	
7. Select Start.			
Do not disturb	the system while waitin	g for results.	

#### Note

If repeated zero attempts fail, follow the procedure below:

- 1. Gain access to the components in the pan electrical enclosure (Section 9.9).
- 2. Disconnect the Alt O2 inlet tubing elbow fitting from the Mixer manifold.
- 3. Disconnect the tubing from the outlet elbow fitting.
- 4. Repeat the Mixer P Zero following the instructions as they appear on the screen.
- 5. Reassemble in reverse order.

Notes

(DU) 4a-22 09/07 1009-0357-000

# 4b Install/Service Menus (HPDU)

In this section	4b.1 Service and Installation menu structure	4b-2
	4b.2 Install/Service Menu (Super User)	4b-3
	4b.2.1 Install/Service - Page 1	4b-3
	4b.2.2 Install/Service - Page 2	4b-9
	4b.3 Installation Menu	4b-10
	4b.3.1 Configuration	4b-11
	4b.3.2 Configuration Units	4b-12
	4b.3.3 Options Key	4b-12
	4b.3.4 Copy Configuration	4b-13
	4b.4 Service Menu	4b-14
	4b.4.1 Software/Hardware Ver Menu	4b-15
	4b.4.2 Service Log Menu	4b-16
	4b.5 Calibration	4b-17
	4b.5.1 Spiro Calibration	4b-18
	4b.5.2 User Calibration menu	4b-19
	4b.5.3 Manifold P Span	4b-20
	4b.5.4 Inspiratory Flow Valve	4b-21
	4b.5.5 Insp Flow Zero	4b-22
	4b.5.6 Bleed Resistor	
	4b.5.7 Paw Span	
	4b.5.8 Zero Gas Xducrs	
	4b.5.9 Cal Config	
	Ah 5 10 Miyor D Zoro	1h 26

# 4b.1 Service and Installation menu structure

This section describes the Service level functions that are part of the main software installed in the anesthesia machine.

Section 8a, "Service Diagnostics and Software Download (DU)," covers the functions of the Compact Flash card used to download system software.

Section 12, "Service Application," covers a separate, Windows based service application used to run service diagnostics and other service tests.

#### Menu structure

The Service menu structure has three levels which are password protected:

- Install/Service (super-user)
- Installation
- Service


The **Install/Service** level (super-user password) supports standard hospital preferences: choosing units; setting ventilator, alarm, and gas delivery defaults.


The **Installation** level requires the service password and supports language, gas color codes, flow tube position, country, hardware flags for system components (acgo or scgo etc.), enabling software options, and cloning a system.

The **Service** level requires the service password and supports diagnostic tools and automated component tests.

Follow the menu structure to access the various service screens:

- on the Main Menu, select Screen Setup;
- on the **Screen Setup** menu, select **Install/Service** to access the Install/Service (with super-user password) menu;
- on the Install/Service menu, select Installation (with service password) to access the Installation menu.
- to access the Service menu, select Service (with service password) on the Install/Service menu; or, from the Installation menu, select Service to access the same Service menu without having to enter the service password.





(HPDU) 4b-2 09/07 1009-0357-000

# 4b.2 Install/Service Menu (Super User)

Use the super-user password to access the Install/Service menu: "16-4-34."

# 4b.2.1 Install/Service - Page 1

Menu Item	Message text	Comments
Trends Setup	Configure graphical trend pages.	
Colors and Units	Set colors and units of parameters.	Refer to section
Show Alarm Limits	Select yes to show alarm limits in digit fields.	Default is Yes.
External Gas Monitor	Yes disables O2 limit alarms, the "No O2 sensor alarm", and the "No CO2 or AA monitor" alarm.	Default is No. Select yes only if system is using external monitor for O2, AA, and CO2.
Save Default Case	Save normal screen, gas settings, circuit type, ventilator, and alarm settings from the last case or presettings as defaults.	Last used alarm settings (including hide/show alarm limits, Auto MV Limit, alarm volume), screen layout (middle waveform selection, sweep speed), ventilator mode and setting, balance gas, and start case gas outlet selection are saved as facility defaults.  Note: The Pmax alarm limit shall not be saved higher than 40 cmH2O. The low FiO2 alarm limit shall not be saved lower than 21%.
Factory Defaults	Return to default factory settings.  After selecting Factory Defaults:  "Reset machine for defaults to take effect."	Action: All facility defaults get replaced with factory defaults. Super User settings also get set to Factory Defaults. No Service level configuration settings are changed.
Parameter Settings	Set volume conditions and CO2 humidity compensation.	
Install/Service - Page 2	Show page 2 of the Install/Service Menu.	
Exit	Turn power off to exit the Install/Service menu.	

# **Trends Setup**

Menu Item	Message text	Values
Default Trend	Change default trend type: graphical, numerical, or settings.	Num (default), Graph, or Set
Graphical Trends	Configure graphical trend pages.	
Previous Menu	Return to previous menu.	

# **Graphical Trends**

Menu Item	Message text	Values
Page 1 (Page 2 to Page 5)	Configure first graphical trend page (second, third, fourth, fifth)	
Previous Menu	Return to previous menu.	

# Page Menus

Menu Item	Options	Page 1 Default	Page 2 Default	Page 3 Default	Page 4 Default	Page 5 Default
Field 1	Off—Select Off to clear trend field	Pres	02	AA2	Bal	rr+CO2
Field 2	rr+co2—respiration rate and CO2 Pres—Ppeak, Pplat, and PEEP	TVexp	N20	N20	MAC	Compl
Field 3	MVexp—expired minute volume and respiration rate CO2—CO2	C02	AA1	MAC	MVexp	Off
	02–02  Bal—balance gas  AA1—current anesthetic  AA2—previous anesthetic agent if used  N20—N20  MAC—minimum alveolar concentration  TVexp—tidal volume and respiratory rate  Pmean—Pmean  Spont—spontaneous MVexp and respiration rate  Compl—compliance and Raw					
Previous Menu	Return to previous menu.					

(HPDU) 4b-4 09/07 1009-0357-000

#### **Colors and Units Menu**

The Units menu can be accessed here in the super-user level to change individual preferences, or if required during installation, in the service level Installation menu.

Menu Item	Message text	Values
Colors	Set colors of parameters.	
Weight	Change weight unit: kg or lb.	kg or lb
CO2	Change CO2 unit: %, kPa, or mmHg.	%, kPa, or mmHg;
Gas Supply Pressure	Change gas supply pressure unit: kPa, psi, or bar.	psi, kPa, or bar
Paw	Change Paw unit: kPa, hPa, cmH2O, mmHg, mbar.	kPa, hPa, cmH2O, mmHg, or mbar
Previous Menu	Return to previous menu.	

#### **Colors Menu**

Menu Item	Message text	Values
Paw	Change color of Paw waveform, digits and trend.	Yellow, White, Green, Red, or Blue
Flow	Change color of Flow waveform, Flow and Volume digits and trends.	Yellow, White, Green, Red, or Blue
Resp	Change color of respiration, digits and trend.	Yellow, White, Green, Red, or Blue
C02	Change color of CO2 waveform, digits and trend.	Yellow, White, Green, Red, or Blue
Previous Menu	Return to previous menu.	

# **Factory Defaults**

The following table lists the factory defaults for parameters. The table on the next page lists the factory defaults for alarm limits.

# Factory defaults — Parameters

Parameter	Value
Vent Mode	VCV
TV (tidal volume)	500 ml
Pinsp	5 cmH20 (5 hPa, 0.5 kPa, 5 mbar, 4.0 mmHg)
RR	12 /min
Mech RR	12 /min
Tinsp	1.70
I:E	1:2.0
Trig. Window	25%
Flow Trig.	2 I/min
End Breath	30%
Psupport	2 cmH20
Tpause	Off
Backup Time	30 s
Total Gas Flow	6 I/min with Circle and 10 I/min with Non-Circle
02 Flow I/min	6 I/min with Circle and 10 I/min with Non-Circle
Control Style used for default setting	02% when Install/Service Fresh Gas Controls menu item is set to "02% or "User."
	Flow when Install/Service Fresh Gas Controls menu item is set to "Flow."
Balance Gas	Air
Circuit	Circle
Gas Outlet installed	SCGO
Paw Color	Yellow
Flow Color	Green
Resp Color	White
CO2 Color	White
Paw Units	cmH2O
Weight Units	kg
Temperature Units	С

%
300 m
kPa
. <dot></dot>
English
SO (O2 = white, N2O = blue, Air = blk/wht)
Right side
02
Auto
5
No
Vent
Adult
Off
3
On
On
Fast
Gas
Paw
Flow
CO2

(HPDU) 4b-6 09/07 1009-0357-000

# Factory defaults — Alarm limits

Alarm Limit	Value
Pmax High	40 cmH20 (40 hPa, 4 kPa, 40 mbar, 30 mmHg)
MV High	10 l/min
MV Low	2 I/min
TV High	1000 ml
TV Low	Off
RR High	Off
RR Low	Off
Et CO2 High	8.0% (60 mmHg or 8.0 kPa)
Et CO2 Low	3.0% (23 mmHg or 2.0 kPa)
Fi CO2 High	Off
Fi O2High	Off
Fi 02 Low	21%
EtO2 High	Off
EtO2 Low	Off
Fi Iso High	5%
Fi Iso Low	Off
Et Iso High	Off
Et Iso Low	Off
Fi Sev High	8%
Fi Sev Low	Off
Et Sev High	Off
Et Sev Low	Off
Fi Des High	15%
Fi Des Low	Off
Et Des High	Off
Et Des Low	Off
Fi Enf High	5%
Fi Enf Low	Off
Et Enf High	Off
Et Enf Low	Off
Fi Hal High	5%
FI Hal Low	Off
Et Hal High	Off
Et Hal Low	Off

# **Parameter Settings**

Menu Item	Message text	Values
TV Based on	Change volume calculation conditions: ATPD or BTPS.	ATPD - default (Ambient temperature and pressure, dry humidity condition)  BTPS (Body temperature, ambient pressure, saturated humidity condition)
CO2 Numbers	Change humidity compensation type in CO2 partial pressure values.	Dry - default Wet
Previous Menu	Return to previous menu.	

(HPDU) 4b-8 09/07 1009-0357-000

# 4b.2.2 Install/Service - Page 2

Menu Item	Message text	Comments
Fresh Gas Controls	Select style for Fresh Gas Controls: 02% and Total Flow or Individual Gas Flows.	Default is 02%.
VCV Cardiac Bypass	Allow machine breaths during cardiac bypass.	Default is No.
Installation	Set language, gas colors, hardware, and install options.	Navigate with password to Installation menu. Password is "26-23-8"
Service	Show technical data for troubleshooting and calibration.	Navigate with password to Service menu. Password is "26-23-8"
Install/Service - Page 1	Show page 1 of the Install/Service Menu.	
Exit	Turn power off to exit the Install/Service menu.	

#### **Fresh Gas Controls**

Menu Item	Message text	
User	Allow flow control styles to change between cases.	
Flow	Select Individual Gas Flow controls.	
02%	Select 02% and Total Flow controls.	

# **VCV Cardiac Bypass**

Menu Item	Message text	
No	Disable alveolar support in cardiac bypass.	
Yes	Enable alveolar support in cardiac bypass.	

# 4b.3 Installation Menu

Use the service-level password to access the Installation menu: "26-23-8."

Whenever the installation menu is entered, "Enter Service dd-mmm-yyyy hh:mm:ss" is recorded in the Event log.

Menu Item	Message text	
Configuration	Set language, gas color code, O2 flowmeter position.	
Units	Set units.	
Options Key	Enable software options.	
Options List	Display software options.	
Copy Config	Normal Message "Save or install configuration and default settings using memory card."	
	Blocked Message "Please insert memory card."	
Service	Show error, event, and alarm logs. (Accessing the Service menu from the Installation menu does not require second use of service password.)	
Previous Menu	Return to previous menu.	

C	Configuration
D	ecimal Marker
L	anguage
G	as supply Colors
C	2 Flowtube
۷	entilator Drive Gas
Α	ltitude
G	as Outlet
Ν	I20 Enabled

<b>Configuration Units</b>
Weight CO2 Gas Supply Pressure Paw

Options Key
Current Key
Entry 1
Entry 2
Entry 3
Entry 4
Entry 5
Entry 6
Entry 7
Save New Key
Control Board ID

Options List		
Available Options		
SIMV/PSV		
PCV		
PSV Pro		
PCV-VG		
VCV Cardiac Bypass		

Copy Configuration
Save to Card Copy from Card

(HPDU) 4b-10 09/07 1009-0357-000

# 4b.3.1 Configuration

Menu Item	Message text	Values	Comments
Decimal Marker	Select decimal delineator.	0.01, 0 01 or 0,01	
Language	Change language translation of screen texts.	Chinese (simplified) Czech Danish Dutch English Finnish French German Greek Hungarian Italian Japanese Norwegian Polish Portuguese Russian Spanish Swedish Turkish	Default: English
Gas supply Colors	Change color of 02, N2O, and Air.	ANSI, ISO, Neutral	ANSI: 02 green, Air yellow, N20 blue; ISO: 02 white, Air black/white, N20 blue; Neutral: All gases white.
02 Flowtube	O2 on left or right- hand side.	Left, Right	
Ventilator Drive Gas	Change drive gas to match machine configuration.	Air, 02	
Altitude	Change altitude used for gas calculations.	-400 to 3000 m in 100-m increments	
Gas Outlet*	Change type of fresh gas outlet.	SCGO, ACGO	SCGO: Use insp port. ACGO: Use auxiliary port.
N20 Enabled	Change to match machine configuration.	Yes; No	

<sup>\*</sup> For machines without a separate auxiliary common gas outlet and selector switch, set *Gas Outlet* to SCGO: Selectable Common Gas Outlet.

<sup>\*</sup> For machines with an external auxiliary common gas outlet and selector switch, set *Gas Outlet* to ACGO: Auxiliary Common Gas Outlet

# 4b.3.2 Configuration Units

Menu Item	Message text	Values
Weight	Change weight unit: kg or lb.	kg or lb
CO2	Change CO2 unit: %, kPa, or mmHg.	%, kPa, or mmHg
Gas Supply Pressure	Change gas supply pressure unit: kPa, psi, or bar.	psi, kPa, bar
Paw	Change Paw unit: kPa, hPa, cmH2O, mmHg, mbar.	kPa, hPa, cmH2O, mmHg, or mbar

## 4b.3.3 Options Key

The Options Key menu is used to configure the software to include the features that the customer has purchased. The included features are shown in the Options List menu.

#### **Options Key menu**

Menu Item	Message text	Values
Current Key	Enter key code to enable options.	XXXAXBC
Entry 1	Enter first entry of key-code.	0 to 9, A to Z, ~,
Entry 2	Enter second entry of key-code.	!, @, #, \$, %, ^, *, (,), ?
Entry 3	Enter third entry of key-code.	
Entry 4	Enter fourth entry of key-code.	
Entry 5	Enter fifth entry of key-code.	
Entry 6	Enter sixth entry of key-code.	
Entry 7	Enter seventh entry of key-code.	
Save New Key	Confirm entries for key-code.	
Control Board ID	Control number used by key-code.	XXX

When options are added, "Add <option> dd-MMM-yyy hh:mm:ss" is written to the event log.

If more than one option is added, each option is be listed separately.

(HPDU) 4b-12 09/07 1009-0357-000

#### **Options List menu**

The options list shows which options are enabled.

Menu Item	Message text	Values *
Available Options		
SIMV/PSV	SIMV vent w/pressure support.	On, Off
PCV	Pressure controlled ventilation.	On, Off
PSV Pro	Pressure support ventilation w/backup.	On, Off
PCV-VG	Pressure controlled volume guaranteed ventilation.	On, Off
VCV Cardiac Bypass	Allow VCV during cardiac bypass.	On, Off

<sup>\*</sup> On if option enabled. Off if disabled.

# 4b.3.4 Copy Configuration

#### **Copy Configuration menu**

Menu Item	Message text	Values	Comments
Save to Card	Save Configuration and defaults to card.	   The field is blank until the data has either been written to the card (OK) or the system determines it cannot write to the card (Fail).	Saves all settings that are not hardware dependent, including facility defaults, colors, units, 02 flow tube position, decimal marker, and altitude.
Copy from Card	Copy Configuration and defaults from card.  When completed: Copy from card complete. Please reboot system.	    The field is blank until the data has either been read from the card (OK) or the system determines it cannot read the card or the card does not have the required data (Fail).	

Systems cannot accept configuration files from a different product model.

The software version is stored with the saved configuration. A system will reject any configurations from other than the current version of software.

Selecting Save to Card overwrites any configuration on the card.

# 4b.4 Service Menu

Use the service-level password to access the Service menu: "26-23-8."

Whenever service menu is entered, "Enter Service dd-mmm-yyyy hh:mm:ss" is recorded in the Event log.

Menu Item	Message text
SW HW versions	Scroll through system information.
Service Log	Show error, event, and alarm histories.
Calibration	Push ComWheel to perform service Calibrations.
Previous Menu	Return to previous menu.

SW HW Versions
Total Time: Software Release: Model Code: Machine Serial Number: Option Package: Option Code: Anes Disp Mixer Vent
Power MGas EVap

# Scrvice Log Menu Scroll Recent Error History Event History Alarm History Copy Logs Reset Logs Previous Menu

Calibration	
Instructions	
Spiro Calibration	
User Calibration	
Manifold P Span	
Insp Flow Valve	
Insp Flow Zero	
Bleed Resistor	
Paw Span	
Zero Gas Xducrs	
Cal Config	
Mixer P Zero	
Previous Menu	

(HPDU) 4b-14 09/07 1009-0357-000

# 4b.4.1 Software/ Hardware Ver Menu

Turn the ComWheel to scroll through the list box.

Push the ComWheel to return to the Service menu.

# **System Information menu**

List box text with X=Number, A, B, C = letter
Total Time: XXXXX (Minutes)
Software Release: XX.XX
Model Code: XXX
Machine Serial Number: ABCDXXXXX
Option Package: XXX
Options Code: XXXXX
Anes Software Version: XX.XX
Anes Hardware Version: XXXX-XXXX REV A
Anes Board Serial Number: ABCXXXXX
Disp Software Version: XX.XX
Disp BIOS Ver: XX.XX
Disp Hardware Version: XXXX-XXXX REV A
Disp Hardware Serial Number: ABCXXXXX
<u>'</u>
Mixer Software Version: XX.XX
Mixer Hardware Version: XXXX-XXXX REV A
Mixer Board Serial Number: ABCXXXXX
Mixer O2 Flow Sensor Serial Number: XXXXXXXXX
Mixer Balance Gas Flow Sensor Serial Number: XXXXXXXXX
Vent Software Version: XX.XX
Vent Hardware Version: XXXX-XXXX REV A
Vent Intf Board Serial Number: ABCXXXXX
Power Software Version: XX.XX
Power Hardware Version: XXXX-XXXX REV A
Power Board Serial Number: ABCXXXXX
MGas Software Version: X.X
MGas Hardware Version: <module type=""></module>
MGas Hardware Serial Number: XXXXXXXX

The MGas information is only displayed when an Airway module is present.

#### 4b.4.2 Service Log Menu

The Service log menu is an organized listing of stored events.

Menu Item	Message text
Scroll Recent	Scroll through newest entries.
Error History	Show error history.
Event History	Show event history.
Alarm History	Show alarm history.
Copy Logs	Save HW/SW info and all logs to memory card.
Reset Logs	Erase Error and Alarm log entries
Previous Menu	Return to previous menu.

Each history log shows at the top of the screen the total "Running Hours" and the date when the logs were last reset.

Whenever logs are reset, "Reset Logs dd-MMM-yyy hh:mm:ss" is recorded in the Event log.

If the logs are saved to a memory card, the machine's serial number is saved along with the current contents of the logs and the date and time.

#### **Error History** ◀

The Error History log lists the last 200 errors logged since the last log reset, starting with the most recent. The system stores the last 1,000 errors logged since the last log reset.

#### **Event History ◀**

The Event History log records the service history of the device. This includes: service calibrations, entry into the service mode, options enabled, and software installation. In the event of a board replacement, it is understood that this log like all others could be lost.

The Event History menu lists the last 200 events logged starting with the most recent. The Event History log stores the last 1000 events.

The Event History log cannot be reset.

#### **Alarm History ◀**

The Alarm History log lists the last 200 medium and high priority parameter alarms since the last log reset starting with the most recent. The Alarm History log store the last 1000 entries.

#### Copy Logs

The Copy Logs function copies Error, Event, and Alarm logs along with the software/hardware configuration to a text file on a PCMCIA card. The copying takes about one minute.

**Note**: Do not remove the Flash Card until the screen shows copy is complete.

(HPDU) 4b-16 09/07 1009-0357-000

# 4b.5 Calibration

For step-by-step instruction, refer to Section 5.4, "Ventilator Calibrations."

Menu Item	Message text
Instructions	These values are used for calibration: Ventilator drive gas - Air or O2 Altitude - XXXX m Change these values on the Cal Config menu.
Spiro Calibration	Check gas module spirometry gains.  Blocking message: Insert gas module with spirometry
User Calibration	Show the normal user calibration menu.
Manifold P Span	Calibrate manifold pressure transducer.
Insp Flow Zero	Zero inspiratory flow valve.
Insp Flow Valve	Calibrate inspiratory flow valve.
Bleed Resistor	Calibrate bleed resister flow.
Paw Span	Calibrate the airway pressure transducer.
Zero Gas Xducrs	Calibrate the gas supply transducers.
Cal Config	Set vent drive gas and altitude.
Mixer P Zero	Zero mixer pres transducer.
Previous Menu	Return to previous menu.

# 4b.5.1 Spiro Calibration

The Spiro Calibration instructions appear when the focus is on Spiro Calibration menu item.

#### **Instructions**

To display TV data:

- 1. Connect a spirometry sensor to MGAS.
- 2. Push the ComWheel to continue.
- 3. Select the correct spirometry sensor type (Adult or Pedi).
- 4. Use the test device to deliver a known TV through the sensor.
- 5. If necessary, repeat steps 1-4 for the other type of sensor.

### **Spiro Calibration menu**

Menu Item	Values	Comments
Sensor Type Sensor Type	Adult or Pedi	Select Sensor type: Adult (Dlite) or Pedi (Pedilite).
Insp Gain	1000	
TVinsp		ml
Exp Gain	1000	
TVexp		ml

(HPDU) 4b-18 09/07 1009-0357-000

# 4b.5.2 User Calibration menu

Menu Item	Message text	
Flow and Pressure	Calibrate the flow and pressure sensors.	
	Remove flow sensor module to start. Replace when Pass or Fail message appears.	
Circuit O2 Cell	Calibrate Circuit O2 Cell.	
	21% 02 Remove Flow Sensor Module. Expose O2 Cell to Room Air. Start 21% calibration. May take 3 min.  100% 02 1. Reconnect flow sensor module. 2. Set Bag/Vent switch to Vent. 3. Select 100% 02.	
Airway Gas	Start Gas Calibration. Calibrate CO2, O2, N2O, and agent measurements. (MGAS module must be installed).  Gas calibration is not available during gas sampling warm-up and certain alarms.	
Backlight Test	Push ComWheel to test back lights. Test every month.  This test turns off one backlight to test the other light. Screen brightness may change during test.	
Previous Menu	Return to the previous menu.	

# 4b.5.3 Manifold P Span

The Manifold P Span instructions appear when the focus is on Manifold P Span menu item.

Refer to Section 5.4.2, "Manifold P Span."

#### Instructions

Read all steps before you start:

- 1. Remove the breathing system, the exhalation valve, and the metal plate.
- 2. Put #2 plugs in the manifold and the drive gas ports of the vent engine.
- 3. Connect a pressure gauge in line with the manifold pressure transducer.
- 4. Push the ComWheel to continue.
- 5. Select Start Manifold P Span.
- 6. Increase the Flow valve setting until the gauge shows 100 cmH20.
- 7. When the gauge shows 100 cmH2O, select Save Calibration.

#### **Manifold P Span menu**

Menu Item	Message text	Values/ Comments
Insp Flow Valve (DAC)	Increase setting until test gauge shows 100 cmH2O (approx 1020 counts). Then save calibration.	O to 4095  Disabled until user selects "Start Manifold P Span".
Start Manifold P Span	Start Calibration. Increase flow valve setting until test gauge = 100 cmH20 (approximately 1020 counts). Then save calibration.  Blocking message: "Connect a supply of the drive gas to continue."	Blocked when the ventilator drive gas supply pressure would cause a gas supply failure alarm during normal operation.
Save calibration	Save Manifold P Span calibration.	
Previous Menu	Return to the previous menu.  During calibration:  "Calibration in progress. Push ComWheel to cancel."	

(HPDU) 4b-20 09/07 1009-0357-000

# 4b.5.4 Inspiratory Flow Valve

The Inspiratory Flow Valve instructions appear when the focus is on the Insp Flow Valve menu item.

Refer to Section 5.4.3, "Inspiratory Flow Valve Cal."

#### Instructions

Read all steps before you start:

- 1. Complete the Manifold P Span calibration.
- 2. Put #2 plugs in the manifold and the drive gas ports of the vent engine.
- 3. Push the ComWheel to show the next menu.
- 4. Select Stage 1 calibration.
- 5. After Pass, replace the manifold port plug with the calibration orifice.
- 6. Select Stage 2 calibration.
- 7. You MUST do both stages for the calibration to be saved.

During calibration, a separate menu shows the counts and corresponding flow at each step.

- If the outcome of both stages of the calibration is Pass, the new calibration data is saved.
- If the outcome of either stage is Fail, the old calibration data is retained.
- The results of each stage of the calibration are saved to the Event Log.

Selecting Previous Menu before the calibration is done, aborts the calibration in progress and keep the old calibration constants.

#### **Insp Flow Valve Menu**

Menu Item	Message text	Comments
Stage 1	Calibrate the insp flow valve at low flows.	Blocked if the ventilator
		drive gas supply pressure
	Blocked text:	would cause a gas supply
	"Connect a supply of the drive gas to continue."	alarm during normal
		operation.
Stage 2	Calibrate the insp flow valve at high flows.	Blocked if Stage 1 has not
		been completed.
	Blocked text:	
	"Stage 1 calibration is required first."	
Insp Flow Valve Data	Show insp flow valve calibration table.	
Previous Menu	Return to previous menu.	
	During Calibration:	
	"Calibration in progress. Push ComWheel to	
	cancel."	

#### **Insp Flow Valve Data menu**

The Insp Flow Valve Data menu contains a table of 24 entries from the previous calibration. The table is erased at the start of Stage 1. The table is updated in real time during the calibration.

# 4b.5.5 Insp Flow Zero

The Insp Flow Zero instructions appear when the focus is on the Insp Flow Zero menu item.

Refer to Section 5.4.4, "Insp Flow Zero."

#### Instructions

Read all steps before you start:

- 1. Push the ComWheel to start the zero check
- 2. No disassembly is required.
- If the outcome of the calibration is Pass, the new calibration data is saved.
- If the outcome is Fail, the old calibration data is retained.
- The result of the calibration is saved to the Event Log.

Selecting Previous Menu before the calibration is done aborts the calibration in progress and keeps the old calibration constants.

#### **Insp Flow Zero menu**

Menu Item	Message text	Values
Start	If the result is failed, do the insp flow valve calibration.	Pass or Fail
Previous Menu	Return to previous menu.  During Calibration:  "Calibration in progress. Push ComWheel to cancel."	

(HPDU) 4b-22 09/07 1009-0357-000

#### 4b.5.6 Bleed Resistor

The Bleed Resistor instructions appear when the focus is on the Bleed Resistor menu item.

Refer to Section 5.4.5, "Bleed Resistor Cal."

#### Instructions

Read all steps before you start:

- 1. Complete the Insp Flow Valve calibration.
- 2. Put #2 plugs in the manifold and the drive gas ports of the vent engine.
- 3. Push the ComWheel to show the next menu.
- 4. Select Start.

The calibration fails if the flow required to reach 91 cmH20 is > 16 l/min.

- If the outcome of the calibration is Pass, the new calibration data is saved.
- If the outcome is Fail, the old calibration data is retained.
- The result of the calibration is saved to the Event Log.

Selecting Previous Menu before the calibration is done aborts the calibration in progress and keep the old calibration constants.

#### **Bleed Resistor Menu**

Menu Item	Message text	Comments
Start	Calibrate manifold pressure to bleed resistor flow.  Blocked text: "Connect a supply of the drive gas to continue."	Blocked if the ventilator drive gas supply pressure would cause a gas supply failure alarm during normal operation.
Bleed Resistor Data	Show bleed resistor calibration table.	
Previous Menu	Return to the previous menu.  During Calibration:  "Calibration in progress. Push ComWheel to cancel."	

#### **Bleed Resistor Data menu**

The Bleed Resistor Data menu contains a table of 17 entries from a previous calibration. The table is erased at the start of the calibration. The table is updated in real time during the calibration.

## 4b.5.7 Paw Span

The Airway P Span instructions appear when the focus is on Paw Span menu item.

Refer to Section 5.4.6, "Paw Span."

#### Instructions

Read all steps before you start:

- 1. Complete the Bleed Resistor calibration.
- 2. Install the flow sensor and circuit module
- 3. Put a #2 plug in the drive gas port of the vent engine.
- 4. Put the calibrated orifice in the manifold port of the vent engine.
- 5. Connect the pressure tee to the insp port. Connect the tee to the calibrated orifice with a 22 mm tube.
- 6. Connect a pressure gauge to the pressure tee.
- 7. Select Start Paw Span.
- 8. Increase the flow valve setting until the gauge shows 100 cmH20.
- 9. Select Save Calibration

#### **Paw Span menu**

Menu Item	Message text	Values/Comments
Insp Flow Valve (DAC)	Increase setting until test gauge shows 100 cmH2O (approximately 1020 counts). Then save calibration.	0 to 4095 (initially set to 800) Disabled until user selects "Start Paw Span".
Start Paw Span	Start Calibration. Increase flow valve setting until test gauge = 100 cmH20 (approximately 1020 counts). Then save calibration.  Blocked text: Connect a supply of the drive gas to continue.	Blocked if the ventilator drive gas supply pressure would cause a gas supply failure alarm during normal operation.
Save calibration	Save Paw Span calibration.	Saves new calibration data. Writes calibration result, date and time to the event log.
Previous Menu	Return to the previous menu.  During Calibration:  "Calibration in progress. Push ComWheel to cancel."	

(HPDU) 4b-24 09/07 1009-0357-000

#### 4b.5.8 Zero Gas Xducrs

The Zero Gas Xducers instructions appear when the focus is on the Zero Gas Xducer menu item.

#### Instructions

Read all steps before you start:

- 1. Remove all cylinders.
- 2. Disconnect all pipeline supplies.
- 3. Select Zero Gas Xducrs.
- 4. Select Start Zero on the next menu.

This page also shows:

Gas supply counts

Gas supply ID

A failed test is usually the result of a pipeline or cylinder still connected to the system.

- If the outcome of the calibration is Pass, the new calibration data is saved.
- If the outcome is Fail, the old calibration data is retained.
- The result of the calibration is saved to the Event Log.

Selecting Previous Menu before the calibration is done aborts the calibration in progress and keep the old calibration constants.

#### Zero Gas Xducrs menu

The Zero Gas Xducrs menu shows only transducers that are installed. If not installed, the menu row is blank.

Menu Item	Message text	Values
O2 Pipeline		0-4095 Counts
02 Cylinder 1		0-4095 Counts
02 Cylinder 2		0-4095 Counts
N20 Pipeline		0-4095 Counts
N20 Cylinder		0-4095 Counts
Air Pipeline		0-4095 Counts
Air Cylinder		0-4095 Counts
Start Zero	Disconnect all pipelines and remove cylinders. Then select Start Zero.	
Previous Menu	Return to the previous menu. During Calibration: "Calibration in progress. Push ComWheel to cancel."	

1009-0357-000 09/07 (HPDU) 4b-25

#### 4b.5.9 Cal Config

Before calibration, you must verify that the Ventilator Drive Gas and the Altitude settings are set appropriately to match the current drive gas configuration and machine location.

If you change any of the settings in the Cal Config menu, you must restart the system.

#### **Cal Config menu**

Menu Item	Message text	Values
Ventilator Drive Gas	Change drive gas to match machine configuration.	Air, 02
Altitude	Change altitude used for gas calculations.	-400 to 3000 m (in 100-m increments)

#### 4b.5.10 Mixer P Zero

The Mixer P Zero instructions appear on the Mixer P Zero menu.

#### Instructions

)efault	Start	
lefaults, select	Defaults (above).	
gas supplies.		
lers.		
nsors.		
ssette.		
outton for 3 sec	onds.	
/ITHOUT gas flo	w for > 5 min.	
	e gas supplies. ders. nsors. assette. outton for 3 sec	defaults, select Defaults (above). e gas supplies. ders. nsors.

#### **Note**

If repeated zero attempts fail, follow the procedure below:

Do not disturb the system while waiting for results.

- 1. Gain access to the components in the pan electrical enclosure (Section 9.9).
- 2. Disconnect the Alt O2 inlet tubing elbow fitting from the Mixer manifold.
- 3. Disconnect the tubing from the outlet elbow fitting.
- 4. Repeat the Mixer P Zero following the instructions as they appear on the screen.
- 5. Reassemble in reverse order.

(HPDU) 4b-26 09/07 1009-0357-000

## **5 Calibration**

<b>⚠ WARNING</b>	After adjustments and calibration are completed, always per procedure. Refer to Section 3 of this manual.	form the checkout
In this section	5.1 Primary Regulators	5-2
	5.1.1 Test setup	5-3
	5.1.2 Testing Primary Regulators	5-3
	5.1.3 Adjusting Primary Regulators	5-8
	5.2 O <sub>2</sub> Flush Regulator	5-9
	5.3 Adjust Drive Gas Regulator	5-10
	5.4 Ventilator Calibrations	5-11
	5.4.1 Cal Config	5-12
	5.4.2 Manifold P Span	5-13
	5.4.3 Inspiratory Flow Valve Cal	5-14
	5.4.4 Insp Flow Zero	5-16
	5.4.5 Bleed Resistor Cal	5-17
	5.4.6 Paw Span	5-18

## **5.1 Primary Regulators**

Follow the procedure in Section 5.1.1 to gain access to the regulators. Then, in Section 5.1.2, select the test that is appropriate for the regulator you are testing.

#### **⚠** WARNING

When testing/adjusting  $N_2O$  regulators, nitrous oxide flows through the system. Use a safe and approved procedure to collect and remove it.

To test or calibrate the primary regulators, you must boot the system with the PCMCIA Service Application and access the Gas Diagnostics function as detailed in Section 8.

5-2 09/07 1009-0357-000

#### 5.1.1 Test setup


#### **⚠** WARNING

Wear safety glasses while test device is connected to the test port.

#### **⚠** CAUTION

Be careful not to plug the output of the primary regulator without having a pressure relief valve in the output circuit.

- 1. Set the system switch to Standby.
- 2. Disconnect all pipeline supplies.
- 3. Remove the upper rear panel (Section 9.3).
- 4. If equipped, turn the auxiliary O<sub>2</sub> flowmeter control fully clockwise (no flow).
- 5. Install a full cylinder in the cylinder supply to be tested. It is essential that the cylinder be within 10% of its full pressure.
- 6. Remove the plug from the test port and connect a test device capable of measuring 689 kPa (100 psi).



# 5.1.2 Testing Primary Regulators

There are two variations of the test procedure for the primary regulators:

- Test A For primary regulators that supply drive gas to the ventilator.
- Test B For all gases not used to supply drive gas to the ventilator.

#### **Test A** For primary regulators that supply drive gas to the ventilator $(0_2$ or Air):

Under low flow conditions, the output pressure of a properly adjusted/functioning regulator should fall within specifications listed in step 4e.

Under high flow conditions, the output pressure should not drop below the specifications listed in step 5g.

- 1. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 2. On the Gas Diagnostics menu, select Gas Supplies.
- 3. Slowly open the cylinder valve for the regulator being tested and observe the pressure reading for the cylinder.

#### 4. Low Flow Test:

- a. On the Gas Diagnostics menu, set the flow of the tested gas to 0.5 l/min.
- b. On the Gas Diagnostics menu, select Gas Supplies.
- c. Close the cylinder valve and allow the pressure to decay to 2068 kPa (300 psi) as indicated on the cylinder pressure display.
- d. At the time that the cylinder pressure reaches 2068 kPa (300 psi), select Main Menu on the Gas Diagnostics page to turn off gas flow.
- e. Within one minute, the test device reading must stabilize between:
  - (60) DIN 372-400 kPa (54-58 psi)
  - (50) Pin Indexed 310-341 kPa (45.0-49.5 psi).
  - If the test device pressure does not stabilize within one minute, replace the cylinder supply.
  - If the test device stabilizes within one minute, but the readings are not within specifications, readjust the regulator (Section 5.1.3).

#### 5. High Flow Test:

- a. Slowly open the cylinder valve.
- b. Remove the ABS breathing system from the machine to allow continuous Insp Valve flow through the exhalation valve.
- c. On the Main Menu of the Service Application, select Vent Diagnostics.
- d. On the Vent Diagnostics page, select **Gas Inlet VIv ON/OFF** to Open the gas inlet valve.
  - Press the ComWheel to return focus to the selection menu.
- e. Select Flow Valve Control.

5-4 09/07 1009-0357-000

- f. Adjust the Flow Valve counts until the **Flow Valve Setting** at the bottom of the screen reads 65 I/min.
- g. While watching the test device, toggle the Gas Inlet Valve several times (Closed, Open, Closed):
  - Press the ComWheel to return focus to the selection menu.
  - Set Gas Inlet VIv to Closed.
  - Press the ComWheel to return focus to the selection menu.
  - Set Gas Inlet VIv to Open.
  - Press the ComWheel to return focus to the selection menu.
  - Set Gas Inlet VIv to Closed.
  - The minimum test device reading observed must be greater than:
    - (60) DIN 221 kPa (32 psi)
    - (50) Pin Indexed 207 kPa (30 psi)
  - Repeat this step (4g) three times.

If the test device reading under "high flow" conditions is less than specified, readjust the regulator per the procedure in Section 5.1.3; however, set the regulated pressure higher by the difference you noted in this step plus 7 kPa (1 psi). This adjusts the "low flow" regulated output to the high side of the specification so that the "high flow" regulated pressure can fall within the specification.

If the regulator subsequently fails the "low flow" specification (step 3d) because the reading is too high, replace the cylinder supply.

- 6. Set the system switch to Standby.
- 7. Close the cylinder valve.
- 8. Bleed the system of all pressure.
- 9. Disconnect the test device and plug the test port (pull on the plug to ensure it is locked in the fitting).
- 10. Replace the ABS breathing system.
- 11. Replace the rear panel.
- 12. Perform the checkout procedure (Section 3).

#### Test B For all gases not used to supply drive gas to the ventilator:

Under low flow conditions, the output pressure of a properly adjusted/functioning regulator should fall within specifications listed in step 4e.

Under high flow conditions, the output pressure should not drop below the specifications in step 5b.

- 1. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 2. On the Gas Diagnostics menu, select Gas Supplies.
- 3. Slowly open the cylinder valve for the regulator being tested and observe the pressure reading for the cylinder.

#### 4. Low Flow Test:

- a. On the Gas Diagnostics menu, set the flow of the tested gas to 0.5 l/min.
- b. On the Gas Diagnostics menu, select Gas Supplies.
- c. Close the cylinder valve and allow the pressure to decay to 2068 kPa (300 psi) as indicated on the cylinder pressure display.
- d. At the time that the cylinder pressure reaches 2068 kPa (300 psi), select Main Menu on the Gas Diagnostics page to turn off gas flow.
- e. Within one minute, the test device reading must stabilize between:
  - (60) DIN 372-400 kPa (54-58 psi)
  - (50) Pin Indexed 310-341 kPa (45.0-49.5 psi).
  - If the test device pressure does not stabilize within one minute, replace the cylinder supply.
  - If the test device stabilizes within one minute, but the readings are not within specifications, readjust the regulator (Section 5.1.3).

#### 5. High Flow Test:

- a. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- b. On the Gas Diagnostics menu, set the flow of the tested gas to 10.0 l/min.

The test device reading must be greater than:

(60) DIN 221 kPa (32 psi)

(50) Pin Indexed 221 kPa (32 psi)

- If the test device reading under "high flow" conditions is less than specified, readjust the regulator per the procedure in Section 5.1.3; however, set the regulated pressure higher by the difference you noted in this step plus 7 kPa (1 psi). This adjusts the "low flow" regulated output to the high side of the specification so that the "high flow" regulated pressure can fall within the specification.
- If the regulator subsequently fails the "low flow" specification (step 3b) because the reading is too high, replace the cylinder supply.

5-6 09/07 1009-0357-000

- 6. Set the system switch to Standby.
- 7. Close the cylinder valve.
- 8. Bleed the system of all pressure.
- 9. Disconnect the test device and plug the test port (pull on the plug to ensure it is locked in the fitting).
- 10. Replace the rear panel.
- 11. Perform the checkout procedure (Section 3).

# 5.1.3 Adjusting Primary Regulators

**Important:** Cylinder supplies in an S/5 Avance machine must have all primary regulators set to the same pressure range: (50) Pin Indexed or (60) DIN.

If a regulator is replaced, the replacement regulator must be set (as required) to the same specification as the one removed.

**Important:** Install a full cylinder in the cylinder supply to be adjusted. It is essential that the cylinder be within 10% of its full pressure.

To adjust the primary regulators, follow the procedure in Section 5.1.1 to gain access to the regulators.

Do not attempt to adjust without flow.

- 1. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 2. On the Gas Diagnostics menu, select Gas Supplies.
- 3. Slowly open the cylinder valve for the regulator being tested and observe the pressure reading for the cylinder.
- 4. On the Gas Diagnostics menu, set the flow of the tested gas to 0.5 l/min.
- 5. On the Gas Diagnostics menu, select Gas Supplies.
- 6. Close the cylinder valve and allow the pressure to decay to 2068 kPa (300 psi) as indicated on the cylinder pressure display.
- 7. When the cylinder gauge reaches 2068 kPa (300 psi), adjust the regulator output pressure to:
  - (60) DIN 386-400 kPa (56-58 psi)
  - (50) Pin Indexed 327-341 kPa (47.5-49.5 psi).

**Note:** It may be necessary to open the cylinder valve and repeat steps 6 and 7 a number of times to achieve the above setting.

- 8. Test the regulator settings per the appropriate test in Section 5.1.2:
  - **Test A** For primary regulators that supply drive gas to the ventilator.
  - **Test B** For all gases not used to supply drive gas to the ventilator.





Adjust clockwise to increase setting

5-8 09/07 1009-0357-000

## 5.2 0<sub>2</sub> Flush Regulator

- 1. Bleed all gas pressure for the machine (Section 9.2).
- 2. Remove the upper rear panel (Section 9.3).
- 3. Remove the  $O_2$  Flush Regulator output tubing. Attach a 6-mm tee and a test device to the open port.





- 4. Connect an  $O_2$  pipeline supply or slowly open the  $O_2$  cylinder valve.
- 5. Push the flush button just enough to achieve a slight flow or open the auxiliary flowmeter if equipped with this option. Read the pressure shown on the test device.

The pressure should be  $241 \pm 7 \text{ kPa}$  ( $35 \pm 1.0 \text{ psi}$ ).

- 6. If adjustment is required:
  - a. Loosen the adjustment screw's jam nut.
  - b. Adjust the regulator (in small steps) to the above specification.
  - c. Tighten the jam nut.
  - d. Verify the reading.
- 7. Disconnect the pipeline supply or close the cylinder valve.
- 8. Bleed gas pressure by pushing the flush button; then, disconnect the tee and test device.
- 9. Reattach the output tubing to the regulator.
- 10. Install the rear panel.


## **5.3 Adjust Drive Gas Regulator**

The drive gas regulator must be adjusted while maintaining a flow of 15 l/min. To adjust the flow, you must boot the system with the PCMCIA Service Application and access the Vent Diagnostics function as detailed in Section 8.

The drive gas regulator should provide a constant gas input pressure of 172 kPa (25 psi).

#### **Calibration setup**

- 1. Attach a pressure test device to the regulator pressure port (shown below)
  - Remove the 6.35-mm (1/4 inch) plug.
  - Attach test device to the open port.
- 2. Remove the ABS breathing system from the machine to allow continuous Insp Valve flow through the exhalation valve.



#### **Calibration procedure:**

- 1. On the Main Menu of the Service Application, select **Vent Diagnostics**.
- 2. On the Vent Diagnostics page, select **Gas Inlet VIv ON/OFF** to Open the gas inlet valve.

Press the ComWheel to return focus to the selection menu.

- 3. Select Flow Valve Control.
- 4. Adjust the Flow Valve counts until the *Flow Valve Setting* at the bottom of the screen reads 15 l/min.
- 5. If required, adjust the regulator to  $172 \pm 1.72$  kPa ( $25 \pm 0.25$  psi) through the access hole in the Vent Engine cover (Section 9.10.1).

5-10 09/07 1009-0357-000

#### 5.4 Ventilator Calibrations

Before performing the ventilator calibrations, verify that the drive gas regulator is adjusted to specifications (Section 5.3).

The Service menu structure is detailed in Section 4. To access the Ventilator Calibrations menu:

- 1. Turn on the system.
- 2. Navigate the menu selections to the Calibration menu.
  - On the Checkout menu, select Bypass Checks.
  - On the **Start Case** menu, press the **Main Menu** button.
  - On the Main Menu, select Screen Setup.
  - On the **Screen Setup** menu, select **Install/Service** (dial in 16 4 34).
  - On the Install/Service menu, select Service (dial in 26 23 8).
  - On the **Service** menu, select **Calibration**.

Unless otherwise specified, perform the ventilator calibrations in the order that they appear on the Calibration menu.

- User Calibration
- Manifold P Span
- Insp Flow Zero
- Insp Flow Valve
- Bleed Resistor
- Paw Span

The following calibrations should be performed as required:

- Zero Gas Xducer:
  - The pipeline and cylinder pressure transducer should be "zeroed" at least once a year.
  - Whenever a pipeline or cylinder pressure transducer is replaced.
- Cal Config:
  - Reset the **Ventilator Drive Gas** to match the machine configuration.
  - Reset the *Altitude* whenever the machine is moved to a new location that differs by more than 100 meters.

#### 5.4.1 Cal Config


Before calibration, you must verify that the Ventilator Drive Gas and the Altitude settings are set appropriately to match the current drive gas configuration and machine location.

If you change any of the settings in the Cal Config menu, you must restart the system.

- 3. On the **Installation** menu, select **Configuration**.
- 4. On the **Configuration** menu, verify the **Ventilator Drive Gas** and the **Altitude** setting; adjust as necessary.
- 5. When done, reboot the system (System switch to Standby; then On).

5-12 09/07 1009-0357-000

#### 5.4.2 Manifold P Span





#### **Calibration setup:**

- 1. Remove the ABS breathing system from the machine.
- 2. Remove the Exhalation Valve.
- 3. Remove the Vent Engine cover.
- 4. Plug the Drive Port (**A**) and the Manifold Port (**B**) on the Vent Engine interface valve.
- 5. Connect the manifold pressure tee adapter (**C**) refer to Section 10.1.2 to the Manifold Pressure Transducer tubing (white inline connectors).
- 6. Connect a manometer to the open port of the tee adapter.

#### **Calibration procedure:**

- 1. On the Calibration menu, select *Manifold P Span*.
- 2. Select Start Manifold P Span.
- 3. Adjust the Insp Flow Valve (DAC) setting until the manometer reading equals 100 cmH<sub>2</sub>O:
  - start at approximately 950 counts (press the ComWheel to activate).
  - continue to increment the count until the manometer reading equals 100 cmH<sub>2</sub>0.
- 4. Select Save Calibration.
- 5. Select Previous Menu.
- 6. Disconnect the manometer from the tee adapter.
- 7. Remove the tee adapter and reconnect the Manifold Pressure Transducer tubing.

#### **Troubleshooting**

#### **Manifold P Span Calibration Failure**


The Calibration will fail if the:


• ADC value calculated for span is outside the range of 21000-27000 counts.

Possible causes for calibration failure:

- Occlusion or moisture in bulkhead or tubing to VIB transducers.
- Pressure transducer outside of range limits Check Service Application for A/D value.

#### **5.4.3 Inspiratory Flow Valve Cal**





#### **Calibration setup**

Leave the Drive Port (**A**) and the Manifold Port (**B**) on the interface valve plugged.

#### **Calibration procedure:**

- 1. On the Calibration menu, select *Insp Flow Valve*.
- 2. Push the ComWheel to enable the **Stage 1** calibration.
- 3. When Stage 1 is completed, remove the plug from the Manifold port and insert the calibrated orifice (**C**)
- Push the ComWheel to enable the **Stage 2** calibration.
   (May take two minutes before you see any effects of the test on the screen.)
- 5. When Stage 2 is completed, select **Previous Menu**.



#### **Troubleshooting**

#### **Stage 1 Calibration Failures**

The Calibration will fail if the:

- Flow valve DAC counts are ≥ 1000 counts while finding points 2 through 6.
- Previously found DAC value is ≥ to the current DAC value while finding points 2 through 6.
- Previously recorded flow for a previous DAC is > the previous flow for a previous DAC.
- Points 2, 3, and 4 have the same value stored for flow (this would cause a divide by zero when extrapolating).

Possible causes for calibration failure:

- Check Altitude and Drive Gas selection
- Leaks around the test plugs
- Leaks in Vent Engine Interface Manifold Inspect for leaks
- Insp Flow Valve not closing completely (leaky) Replace Insp Flow Valve
- Drive gas regulator not adjusted / stable Check regulator calibration
- Insp Flow Valve not linear Replace Insp Flow Valve

5-14 09/07 1009-0357-000

#### **Troubleshooting** Stage 2 Calibration Failures


The Calibration will fail if the:

- Flow valve DAC reaches 4095 before determining the Lift-Off Point.
- Previously found DAC value is ≥ to the current DAC value while finding points 7 through 24.
- Previously recorded flow for a previous DAC is > the previous flow for a previous DAC.
- End point DAC of 4095 does not give a flow > 100 l/min.

Possible causes for calibration failure:

- Check Altitude and Drive Gas selection.
- Leaks around the test plug or Calibration Flow Orifice.
- Leaks in Vent Engine Interface Manifold Inspect for leaks.
- Inadequate drive gas supply (cannot deliver > 100 l/min).
- Drive gas regulator not adjusted / stable Check regulator calibration.
- Insp Flow Valve not linear Replace Insp Flow Valve.

#### 5.4.4 Insp Flow Zero Calibration setup



- 1. Remove the Calibration Orifice from the Manifold port.
- 2. Plug the Manifold (B) port.
- 3. Leave the Drive Gas (A) port plugged.

#### **Calibration procedure:**

- 1. On the Calibration menu, select *Insp Flow Zero*.
- 2. Select Start.
- 3. Select Previous Menu.

#### **Troubleshooting**

#### **Insp Flow Zero Calibration Failures**

The Calibration will fail if the:


- Inspiratory Valve Calibration (stage 1 and stage 2) has not been performed.
- Corresponding flow is not between the 2nd and 5th point in the Inspiratory valve calibration table.

Possible causes for calibration failure:

- Inspiratory Valve Calibration not complete Perform Calibration.
- Insp Flow Valve not closing completely (leaky) Replace Insp Flow Valve.
- Insp Flow Valve not linear Replace Insp Flow Valve.

5-16 09/07 1009-0357-000

#### 5.4.5 Bleed Resistor Cal



#### **Calibration setup**

Leave the Drive Port (**A**) and the Manifold Port (**B**) on the interface valve plugged.

#### **Calibration procedure**

- 1. On the Calibration menu, select **Bleed Resistor**.
- 2. Select Start.
- 3. When the test is completed, select **Previous Menu**.

#### **Troubleshooting**

#### **Bleed Resistor Calibration Failures**

The Calibration will fail if the:

- Flow is greater than 50 I/min before 105 cmH20 of pressure is reached.
- Table created has a pressure or flow that is greater than or equal to the next flow or pressure point in the table.


Possible causes for calibration failure:


- Check Altitude and Drive Gas selection.
- Leaks around the test plugs.
- Drive gas regulator not adjusted / stable Check regulator calibration.
- Inspiratory Valve Calibration not complete Perform Calibration.
- Insp Flow Valve not closing completely (leaky) Replace Insp Flow Valve.
- Insp Flow Valve not linear Replace Insp Flow Valve.

#### 5.4.6 Paw Span

#### **Calibration setup**

- 1. Leave the Drive port (A) port plugged.
- 2. Remove the plug from the Manifold port.
- 3. Attach a patient circuit tube to the Calibrated Flow Orifice.
- 4. Insert the Calibrated Flow Orifice into the Manifold port (B).
- 5. Separate the Circuit module from the ABS Bellows module.
- 6. Install only the Circuit module (**C**) on to the machine.
- 7. Connect a pressure sensing tee (**D**) to the inspiratory flow patient connection.
- 8. Connect the open end of the patient circuit tube to the flow port of the pressure sensing tee.
- 9. Connect a manometer to the pressure sensing port of the tee connector.





#### **Calibration procedure**

- 1. On the Calibration menu, select *Paw Span*.
- 2. Select Start Paw Span.
- 3. Adjust the Insp Flow Valve (DAC) setting until the manometer reading equals 100 cm ${\rm H}_2{\rm O}$ :
  - start at approximately 950 counts (press the ComWheel to activate).
  - continue to increment the count until the manometer reading equals 100 cmH<sub>2</sub>0.
- 4. Select Save Calibration.
- 5. Select Previous Menu.

5-18 09/07 1009-0357-000

#### **Troubleshooting** Paw Span Calibration Failure

The Calibration will fail if the:

ADC value calculated for span is outside the range of 21000-27000 counts.
 Possible causes for calibration failure:

- Occlusion or moisture in bulkhead or tubing to VIB transducers.
- Pressure transducer outside of range limits Check Service Application for A/D value.

Notes

5-20 09/07 1009-0357-000

## **6 Installation and Maintenance**

#### In this section

This section covers the regular maintenance procedures (minimum requirements) needed to make sure that the S/5 Avance anesthesia machine operates to specifications.

6.1 Avance Installation Checklist6-2
6.2 Avance Planned Maintenance6-4
6.3 Free breathing valve maintenance
6.4 MOPV pressure relief valve test6-7
6.4.1 Test setup6-7
6.4.2 Test procedure6-7
6.5 Pressure Limit Circuit test6-8
6.6 Mixer test
6.6.1 Mixer outlet check valve leak test 6-10
6.6.2 Mixer flow verification
6.7 Alternate O2 flowmeter tests
6.8 Auxiliary O2 flowmeter tests
6.9 Integrated Suction Regulator tests
6.10 Battery capacity test

#### **⚠** WARNINGS

Do not perform testing or maintenance on the S/5 Avance anesthesia machine while it is being used on a patient. Possible injury can result.

Items can be contaminated due to infectious patients. Wear sterile rubber gloves. Contamination can spread to you and others.

Obey infection control and safety procedures. Used equipment may contain blood and body fluids.

## **6.1 Avance Installation Checklist**

Serial Number:	Date: (YY/MM/DD) / /
Hospital:	Performed by:
	Unpack and assemble the Avance System.
	2. Reconfigure the sample gas return line as required (Section 9.26).
	3. Access the Install/Service menu and configure the parameter colors as required (Section 4a.2.1).
	4. Access the Installation menu from the Install/Service menu and change the following as required:
	<ul> <li>a. Configuration (Section 4a.3.1)</li> <li>Decimal Marker</li> <li>Language</li> <li>Gas Supply Colors</li> <li>O<sub>2</sub> Flowtube</li> <li>Ventilator Drive Gas</li> <li>Altitude</li> <li>Gas Outlet</li> <li>N<sub>2</sub>O Enabled</li> <li>b. Units Menu (Section 4a.3.2)</li> <li>Weight</li> <li>CO<sub>2</sub></li> </ul>
	<ul><li>Gas Supply Pressure</li><li>Paw</li><li>Options List (Section 4a.3.3)</li></ul>
	<ul> <li>Check that the factory installed ventilation options match the configuration purchased with the machine.</li> </ul>
	d. Copy Configuration Menu (Section 4a.3.4)
	<ul> <li>Can be used to save a configuration to a Compact Flash card and then copy the configuration to additional machines.</li> </ul>
	e. From the Service Menu select the Service Log Menu (Section 4a.4.2)
	<ul> <li>Review and reset the error and alarm log entries.</li> </ul>

6-2 09/07 1009-0357-000

- f. From the Service Menu select Calibration and perform the following calibrations (Section 4a.5):
  User Calibration
  Manifold P Span
  Insp Flow Zero
  Insp Flow Valve
  Bleed Resistor
  Paw Span
  Zero Gas Xducers
  5. Verify the "Schedule Service Calibration" message is not present in the
- normal display.
- 6. Complete the System Checkout by performing the following steps:
  - a. Inspect the system (Section 3.1)
  - b. System checkout (Section 3.4)
  - c. Pipeline and Cylinder tests (Section 3.8)
  - d. Flush Flow test (Section 3.10)
  - e. Alarm tests (Section 3.11)
  - f. Alternate O<sub>2</sub> flowmeter tests (Section 3.12)
  - g. Auxiliary  $O_2$  flowmeter tests, if equipped with option (Section 3.13)
  - h. Integrated suction regulator tests, if equipped with option (Section 3.14)
  - i. Power failure test (Section 3.15)
  - j. Electrical safety tests (Section 3.16)


## **6.2 Avance Planned Maintenance**

Serial Number:			Date: (YY/MM/DD) / /	
Hospital:			Performed by:	
☐ 12 months ☐ 24	4 month	☐ 48 month		
Every twelve (12) months			ving steps every 12 months. o the sections listed.	
		Avance anesthe	ed URM are found in the User's Reference manuals for the esia system.  Ed TRM are found in this Technical Reference manual.	
Parts Replac	cement			
			rizer port o-rings (Section 9.14.1) r 1102-3016-000)	
Checks an	d Tests			
		1. AGSS Mainter	nance (URM - Part 2, Section 2):	
		<ul><li>Empty any of</li></ul>	condensate from the reservoir (disposable item).	
•		· · · · · · · · · · · · · · · · · · ·	orake for occlusion on active AGSS.	
		•	an or replace filter on active AGSS.	
			tem Maintenance (URM - Part 2, Section 2)	
		3. Bellows Asser	nbly Maintenance ( <i>URM - Part 2, Section 2</i> )	
		4. Bellows Asser	nbly Tests (URM - Part 2, Section 2)	
		5. O <sub>2</sub> Cell Calibra	ation (URM - Part 2, Section 3)	
		6. Flow Sensors	Calibration (URM - Part 1, Section 5)	
		7. Perform the ch	neckout procedures in Section 3.	
			system (TRM - Section 3.1)	
		_	st (TRM - Section 3.6)	
			ack pressure test (TRM - Section 3.7)	
		•	d cylinder tests (TRM - Section 3.8)  larm test (TRM - Section 3.8.1)	
		-	ief tests (TRM - Section 3.9)	
			Fest (TRM - Section 3.10)	
		<ul><li>Alarm tests</li></ul>	(TRM - Section 3.11)	
		<ul><li>Power failur</li></ul>	e test (TRM - Section 3.15)	
		8. Alternate 0 <sub>2</sub> fl	owmeter tests (TRM - Section 6.7)	
		9. Auxiliary O <sub>2</sub> flo	owmeter tests (TRM - Section 6.8)	
		10.Integrated Suc	ction Regulator tests (TRM - Section 6.9)	
		11. Mixer outlet o	heck valve leak test (TRM - Section 6.6.1)	

6-4 09/07 1009-0357-000

	12. Zero the Mixer (TRM - Section 8a.5.1)
	13. Mixer flow verification test (TRM - Section 6.6.2)
	14. Perform the following diagnostics using the PCMCIA Service Application.
	<ul> <li>Display Diagnostics (TRM - Section 8a.4).</li> <li>MOPV pressure relief valve test (TRM - Section 6.4).</li> <li>Pressure Limit Circuit test (TRM - Section 6.5).</li> <li>Adjust Drive Gas Regulator (TRM - Section 5.3).</li> </ul>
	15. From the Service Calibration menu ( <i>TRM - Section 4a.5</i> ), perform the following (refer to <i>TRM - Section 5.4</i> for details):
	<ul> <li>User Cals</li> <li>Manifold P Span</li> <li>Inspiratory flow valve</li> <li>Inspiratory flow zero</li> <li>Bleed resistor</li> <li>Paw Span</li> <li>Zero Gas Transducers</li> </ul>
	16. From the Service Log menu (TRM - Section 4a.4.2), perform the following:
	<ul> <li>Access the Error History log. If any error codes have been logged, follow the appropriate troubleshooting procedures. Clear the error log.</li> </ul>
	17. Perform the system "All Checks" (TRM - Section 3.4).
	<ul> <li>Low P leak check (TRM - Section 3.4.1)</li> <li>Quick check (TRM - Section 3.4.2)</li> <li>Vent check (TRM - Section 3.4.3)</li> <li>Circuit O<sub>2</sub> cell check (TRM - Section 3.4.4)</li> </ul>
	18. Electrical safety tests (TRM - Section 3.16).
Every twenty-four (24) months	In addition to the 12-month requirements, replace the following parts every 24 months. All parts should be replaced before performing the checks, tests, and calibrations.
<b>Parts Replacement</b>	Refer to TRM - Section 6.3.
	Perform the following step:
	$1. \ \ \text{Replace the free breathing flapper valve (Stock  \text{Number}  0211\text{-}1454\text{-}100)}.$
	$2. \ \ \mbox{Replace the free breathing valve o-ring (Stock Number 1503-3208-000)}.$
Every forty-eight (48) months	In addition to the 24-month requirements, replace the following parts every 48 months. All parts should be replaced before performing the checks, tests, and calibrations.
<b>Parts Replacement</b>	Refer to TRM - Section 9.8.4.
	Perform the following step:
	1. Replace the system batteries* (Stock Number 1009-5682-000).
	*Note: Refer to the "Battery capacity test" in TRM - Section 6.10.

## 6.3 Free breathing valve maintenance



Refer to Section 9.10 to access the Pneumatic Vent Engine.

- 1. Unscrew the valve seat (A) from the side of the interface manifold.
- 2. Inspect the flapper (B) and valve seat for nicks, debris and cleanliness.

#### To replace the flapper valve

- 3. If necessary, clean the new flapper valve with alcohol.
- 4. Pull the tail (**C**) of the new free breathing valve flapper through the center of the valve seat until it locks in place.
- 5. Trim the tail flush with outside surface of the valve seat (refer to the removed flapper).
- 6. Replace the O-ring (**D**). Lubricate with a thin film of Krytox.
- 7. Hand screw the assembly into the interface manifold.
- 8. Reassemble the system.
- 9. Perform the Preoperative Checkout Procedure (refer to the User's Reference manual).

6-6 09/07 1009-0357-000

### **6.4 MOPV pressure relief valve test**


#### **⚠** WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Make sure that there are no test plugs or other objects caught in the breathing system.

#### 6.4.1 Test setup

- 1. Remove the ABS breathing system.
- 2. Plug the inspiratory flow (drive pressure) port of the exhalation valve with a stopper.



#### **6.4.2 Test procedure**

To test the pressure relief valve, you must establish a flow (blocked by setup above) of 30 I/min through the Inspiratory Flow Control valve. To adjust the flow, you must boot the system with the PCMCIA Service Application and access the Vent Diagnostics function as detailed in Section 8.

- 1. On the Main Menu of the Service Application, select **Vent Diagnostics**.
- 2. On the Vent Diagnostics page, select **Gas Inlet VIv ON/OFF** to Open the gas inlet valve.

Press the ComWheel to return focus to the selection menu.

- 3. Select Flow Valve Control.
- 4. Adjust the Flow Valve counts until the *Flow Valve Setting* at the bottom of the screen reads approximately 30 I/min.
- Carefully listen for the MOPV relief weight to be relieving and "popping off" from its seat (a purring sound). This indicates the valve is functioning correctly.
- 4. Set the system switch to Standby.
- 5. Remove the stopper from the inspiratory flow port.
- 6. Reassemble the system.
- 7. Perform the Preoperative Checkout Procedure (refer to the User's Reference manual).


#### **6.5 Pressure Limit Circuit test**

To perform the test:


- establish a closed patient airway circuit.
- increment the pressure in the airway circuit.
- observe the output of the airway pressure transducer.
- note that the "pressure limit circuit" trips at approximately 109 cmH<sub>2</sub>0.

#### **Test setup**

- 1. Remove the ABS breathing system from the machine.
- 2. Remove the Exhalation Valve.
- 3. Remove the Vent Engine cover.
- 4. Separate the Circuit Module from the ABS Bellows Module.
- 5. Install the Circuit Module only.
- 6. Plug the Drive Port (A) on the Vent Engine interface valve.
- 7. Attach a patient circuit tube to the Calibrated Flow Orifice test tool.
- 8. Insert the Calibrated Flow Orifice into the Manifold (pilot) Port (B).
- 9. Connect the open end of the patient circuit tube to the inspiratory flow patient connection (**C**).



Calibrated Flow Orifice 1504-3016-000






6-8 09/07 1009-0357-000

#### **Test Procedure**

- 1. On the Main Menu of the Service Application, select **Vent Diagnostics**.
- On the Vent Diagnostics menu, select **Status** and verify that "Over Pressure Circuit" reads OK.
- 3. Press the ComWheel to return focus to the selection menu.
- Select Gas Inlet VIv ON/OFF to Open the gas inlet valve.
   Press the ComWheel to return focus to the selection menu.
- 5. Select Flow Valve Control.
- 6. Adjust the Flow Valve counts to approximately 1000 counts and observe the "Airway Pressure" reading.
- 7. Increase the flow count slowly until the "Airway Pressure" reading reaches approximately 109 cmH<sub>2</sub>O.
- 8. Continue to increase the flow by one encoder click and observe the airway pressure until gas flow stops.
- 9. Select the Status page and verify that:
  - "Over Pressure Circuit" reads High Pressure.
  - "Gas Inlet Valve Feedback" reads Closed.
- 8. Reassemble the system.
- 9. Perform the Preoperative Checkout Procedure (refer to the User's Reference manual).

#### 6.6 Mixer test



To perform the mixer tests, you must gain access to the mixer outlet tubing which is connected to the inlet of the vaporizer manifold.

- 1. Remove the rear upper panel.
- 2. Disconnect the mixer outlet tube at the inlet to the vaporizer manifold.

# 6.6.1 Mixer outlet check valve leak test

To test the mixer outlet check valve you must apply back pressure to the check valve through the mixer outlet tubing and time the leak down rate of the pressure.

- 1. Tee in a pressure gauge and a syringe to the mixer outlet tube.
- 2. Slowly pressurize the mixer outlet check valve to 200 mmHg.
- 3. The pressure shown on the test gauge should not decrease by more than 10 mmHg in 30 seconds.

# 6.6.2 Mixer flow verification

To perform the flow verification test, you must attach a flowmeter to the mixer outlet tubing and access the Gas Diagnostics function of the Service Diagnostics application on the Software Download card.

- 1. Connect a flowmeter to mixer outlet tubing.
- 2. If the system includes an N<sub>2</sub>O supply, connect the output of the flowmeter to the input of the vaporizer manifold to scavenge the flow.
- 3. On the Gas Diagnostics menu
  - (Section 8a.2 for machines with a DU)
  - $\bullet$  (Section 8b.5.1 for machines with an HPDU),

select the following and verify the readings on the test flowmeter.

Note: On machines with an
<b>HPDU and System Software</b>
5.X, you can also set up gas
flows on the Gas Delivery
Schematic (Section 12.3.2)
of the PC Service
Application.

	Verify Flowmeter Reading		
Select	Lower Limit (I/min)	Upper Limit (I/min)	
100% 02 at 10 I/min	9.0	11.0	
100% 02 at 0.5 I/min	0.45	0.55	
Air at 10 l/min	9.0	11.0	
Air at 0.5 I/min	0.45	0.55	
100% N20 at 10 l/min	9.0	11.0	
100% N20 at 0.5 I/min	0.45	0.55	

#### Note

If you will be testing the Auxiliary  $O_2$  flowmeter (Section 6.7), you can proceed to the Alternate  $O_2$  "Flow Accuracy Test" at this point without reassembling the machine.

- 4. Remove the test device.
- 5. Connect the mixer outlet tubing to the vaporizer manifold.
- 6. Replace the back cover.
- 7. Perform the Preoperative Checkout Procedure (refer to the User's Reference manual).

6-10 09/07 1009-0357-000

## **6.7 Alternate 0<sub>2</sub> flowmeter tests**

- 1. Open the  $O_2$  cylinder valve or connect an  $O_2$  pipeline.
- 2. Rotate the Alt  $O_2$  flow control fully clockwise to minimum flow.
- 3. Press the Alternate  $O_2$  switch to turn on Alternate  $O_2$  flow. The flowmeter should indicate 0.5 to 0.7 L/min.
- 4. Rotate the flow control counterclockwise (increase). The ball should rise immediately after rotation is begun. It should rise smoothly and steadily with continued counterclockwise rotation. When a desired flow is set, the ball should maintain in a steady position.
- 5. Rotate the flow control clockwise to minimum flow.
- 6. Press the Alternate  $O_2$  switch to turn off Alternate  $O_2$  flow; confirm yes.

#### **Flow Accuracy Test**

**Note:** To check flow accuracy, be sure that the flow test device is capable of measuring 0-15 l/min with an accuracy of  $\pm 2\%$  of reading.

To perform the test, you must gain access to the mixer outlet tubing which is connected to the inlet of the vaporizer manifold.

- 1. Remove the rear cover.
- 2. Disconnect the mixer outlet tube at the inlet to the vaporizer manifold.
- 3. Connect a flowmeter to mixer outlet tubing.
- 4. Press the Alternate  $O_2$  switch to turn on Alternate  $O_2$  flow.
- 5. Adjust the flowmeter so the **center** of the ball aligns with the selected test point (observe that the ball maintains a steady position for 10 seconds).
- 6. The test device reading should be between the limits shown for each of the selected settings in the table below.

#### **Flow Tester Reading**

Flowmeter Setting L/min	Lower Limit I/min	Upper Limit I/min
minimum (valve fully closed)	0.5	0.7
1	0.5	1.5
3	2.5	3.5
5	4.5	5.5
10	9.0	11.0
maximum (valve fully open)	10.0	13.0

- 7. Rotate the flow control clockwise to minimum flow.
- 8. Close the  $O_2$  cylinder valve or disconnect the  $O_2$  pipeline.
- 9. Remove the test device.
- 10. Connect the mixer outlet tubing to the vaporizer manifold.
- 11. Replace the back cover.
- 10. Perform the Preoperative Checkout Procedure (refer to the User's Reference manual).

## 6.8 Auxiliary 02 flowmeter tests

- 1. Open the O<sub>2</sub> cylinder valve or connect an O<sub>2</sub> pipeline.
- 2. Rotate the flow control clockwise (decrease) to shut off the flow. The ball should rest at the bottom of the flow tube and not move.
- Rotate the flow control counterclockwise (increase). The ball should rise immediately after rotation is begun. It should rise smoothly and steadily with continued counterclockwise rotation. When a desired flow is set, the ball should maintain in a steady position.
- 4. Rotate the flow control clockwise to shut off the flow.

#### **Flow Accuracy Test**

**Note:** To check flow accuracy, be sure that the flow test device is capable of measuring 0 to 15 L/min with an accuracy of  $\pm 2\%$  of reading.

- 1. Connect the flowmeter outlet to the flow test device.
- 2. Adjust the flowmeter so the **center** of the ball aligns with the selected test point (observe that the ball maintains a steady position for 10 seconds).
- 3. The test device reading should be between the limits shown for each of the selected settings in the table below.

#### **Flow Tester Reading**

Flowmeter Setting L/min	Lower Limit L/min	Upper Limit L/min
1	0.5	1.5
3	2.5	3.5
5	4.5	5.5
10	9.0	11.0
maximum (valve fully open)	12.0	

- 4. Rotate the flow control clockwise to shut off the flow.
- 5. Close the  $O_2$  cylinder valve or disconnect the  $O_2$  pipeline.

6-12 09/07 1009-0357-000

## **6.9 Integrated Suction Regulator tests**

#### Note

There are two types of integrated suction systems for the Avance anesthesia machine:

- Continuous Vacuum Regulator, Three-Mode, Pipeline Vacuum
- Continuous Vacuum Regulator, Three-Mode, Venturi Derived Vacuum

#### For Pipeline Vacuum systems,

a vacuum source of at least 500 mm Hg (67 kPa or 20 in Hg) is required for testing. The supply open flow must be a minimum of 50 L/min.

#### For Venturi Derived Vacuum systems,

an O2 or Air source of at least 282 kPa (41 psi) is required for testing.

#### **Gauge Accuracy**

The gauge needle should come to rest within the zero range bracket when no suction is being supplied. Gauges which do not comply may be out of tolerance.

#### Note

To check gauge accuracy, be sure that the test gauge is capable of measuring 0 to 550 mm Hg with an accuracy of  $\pm 1\%$  of reading.

- 1. Connect the suction patient port to the test gauge.
- 2. Turn the mode selector switch to I (ON).
- 3. Ensure that the vacuum test gauge is in agreement with the suction vacuum gauge  $\pm$  38 mm Hg/5 kPa at the following test points.

Test points	
Suction vacuum gauge	Test gauge tolerance
100 mm Hg (13.3 kPa)	62-138 mm Hg (8.3-18.4 kPa)
300 mm Hg (40 kPa)	262-338 mm Hg (35-45 kPa)
500 mm Hg (66.7 kPa)	462-538 mm Hg (61.6-71.7) kPa)

#### **Flow Test**

**Note:** To check flow accuracy, be sure that the flow test device is capable of measuring 0-30 L/min.

- 1. Connect the patient port of the suction regulator to the flow test device.
- 2. Rotate the suction control knob fully clockwise (increase).
- 3. Turn the mode selector switch to I (ON) and verify that the flow rate is:
  - at least 20 L/min.
- 4. Disconnect the test flowmeter.

#### (Tests continue on next page.)

#### **Regulation Test**

- 1. Turn the mode selector switch to I (ON).
- 2. Occlude the patient port of the suction regulator.
- 3. Set the vacuum regulator gauge to 100 mm Hg/13 kPa.
- 4. Open and close the patient port several times.
- 5. With the patient port occluded, the gauge should return to 100 mm Hg/13 kPa within a tolerance of  $\pm 10 \text{ mm Hg}/1.3 \text{ kPa}$ .

#### **Vacuum Bleed Test**

- 1. Occlude the patient port of the suction regulator.
- 2. Set the vacuum regulator gauge to 100 mm Hg/13 kPa.
- 3. Turn the mode selector switch to 0 (OFF) and observe the gauge needle. It must return to the zero range bracket or stop pin within 10 seconds.

#### **Vacuum Leak Test**

- 1. Turn the mode selector switch to O (OFF).
- Rotate the suction control knob a minimum of two full turns in the clockwise direction (increase suction) to ensure its setting is not at the off position.
- 3. Occlude the patient port of the suction regulator.
- 4. Observe the suction gauge, the needle should not move.
- 5. Rotate the suction control knob fully counterclockwise to ensure its setting is at the fully off position.
- 6. Turn the mode selector switch to I (ON).
- 7. Observe the suction gauge, the needle should not move.

6-14 09/07 1009-0357-000

### **6.10 Battery capacity test**

Although replacement of the backup batteries is recommended at the end of 4 years, batteries that pass the capacity test can be considered viable for battery backup of the system for up to 6 years at the discretion of the hospital.

Before testing the batteries, ensure that they are fully charged.

### **Test procedure**

- 1. Turn the system on and start a case (simulated).
- 2. Turn off the mains system breaker on the AC Inlet.
- 3. Allow the system to run on battery until it does an orderly shutdown and powers off (can be in excess of 90 minutes).
- 4. Set the system switch to Standby and turn on the mains system breaker.
- 5. Start up the system with the appropriate Software Download card and access Page 1 of the Power Controller Diagnostics function.

Note: On machines with an HPDU and System Software 5.X, you can also access the Power Board Diagnostics data (Section 12.7.1) with the PC Service Application.

DU (Section 8a.1.3)	HPDU (Section 8b.3.2)
Power Diagnostics	Special Functions
Power Control	Power Diagnostics
	Power Controller

- 6. Page 1 of the Power Control Power Diagnostics screen shows the "Date battery Tested" (the date of the last full discharge) and the "Last Full Discharge Time".
  - If the "Last Full Discharge Time" is greater than 45 minutes, the batteries can be left in service for one more year.
  - If the "Last Full Discharge Time" is less than 45 minutes, both batteries should be replaced.

Notes

6-16 09/07 1009-0357-000

# 7 Troubleshooting

In this section	7.1 Troubleshooting Guidelines	7-2
	7.2 Troubleshooting high pressure and low pressure leaks	7-2
	7.3 Troubleshooting Startup Screen (POST) messages – for DU	7-4
	7.4 Troubleshooting Startup Screen (POST) messages – for HPDU	7-5
	7.5 Troubleshooting the HPDU Display	7-6
	7.6 Troubleshooting System Malfunction (safe-state) screen	7-7
	7.7 Breathing System Leak Test Guide	7-8
	7.7.1 Check Valves	7-9
	2.7.2 Breathing System Troubleshooting Flowcharts	7-10
	7.7.3 Leak Isolation Tests	7-15
	7.8 System Troubleshooting Flowcharts	7-30
	7.9 System Malfunction and Alt O2 Flowchart Table	7-42
	7.10 Technical Alarms	7-44
	7.11 Steps and Messages displayed during the System Checkout —  (for System software 3.X or greater)	7-71
	7.12 Steps and Messages displayed during the System Checkout — (for System software 2.X)	7-73
	7.12.1 Steps for the Quick Check	7-73
	7.12.2 Steps for the Vent Check	7-75

### 7.1 Troubleshooting Guidelines

Review system error logs using the Special Functions feature of the Software Download card (Section 8a.5). Review the logs to identify issues and follow the appropriate subsystem troubleshooting procedures.

Troubleshooting high pressure and low pressure leaks	Section 7.2 on page 7-2
Troubleshooting Startup Screen (POST) messages — for DU	Section 7.3 on page 7-4
${\it Troubleshooting Startup Screen (POST) messages-for HPDU}$	Section 7.4 on page 7-5
Troubleshooting the HPDU Display	Section 7.5 on page 7-6
Troubleshooting System Malfunction (safe-state) screen	Section 7.6 on page 7-7
Breathing System Leak Test Guide	Section 7.7 on page 7-8
System Troubleshooting Flowcharts	Section 7.8 on page 7-30
System Malfunction and Alt 02 Flowchart Table	Section 7.9 on page 7-42
Technical Alarms	Section 7.10 on page 7-44

### **⚠ WARNING**

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Make sure that there are no test plugs or other objects caught in the breathing system.

### 7.2 Troubleshooting high pressure and low pressure leaks

Problem	Possible Cause	Action
High Pressure Leak	Pipeline leak	Use a leak detector or Snoop to check for source of leak. Repair or replace defective parts.
	0 <sub>2</sub> flush valve	Use a leak detector or Snoop to check for source of leak. Make sure tubing connections are tight. Replace valve if defective.
	System switch	Use a leak detector or Snoop to check for source of leak. Make sure tubing connections are tight. Replace switch if defective.
	Cylinder not installed properly	Make sure cylinder is correctly aligned. Verify that tee handles are tight.
	Cylinder transducer	Use a leak detector or Snoop to check for source of leak. Tighten/replace transducer if defective.
	Cylinder gaskets	Use a leak detector or Snoop to check for source of leak. Replace gasket if defective.
	Relief valves	Use a leak detector or Snoop to check for source of leak. Replace valve if defective.

7-2 09/07 1009-0357-000

Problem	Possible Cause	Action
Low Pressure Leak (with vaporizer mounted)	Vaporizer not installed properly	Reseat vaporizer if not installed properly.  Have vaporizer serviced at vaporizer center if vaporizer leaks.
	Missing or damaged o-ring on vaporizer manifold	Check condition of o-ring. Replace if missing or damaged.
	Loose fill port	Check fill port. Tighten if loose.
Low Pressure Leak (with or without vaporizer)	Leaking port valve on vaporizer manifold	Use the Vaporizer Manifold Valve Tester to check for leak. See Section 9.14.2 for instructions. If test fails, tighten, repair, or replace as needed.
	Leak in mixer	If vaporizer manifold passed previous tests: Remove tubing from inlet port of vaporizer manifold (mixer outlet tube) and perform leak test of mixer.
	Leaking relief valve on vaporizer manifold	Remove relief valve. Occlude opening. Perform leak test. If test passes, replace valve.
	Leaking flush valve	Attach pressure measuring device on CGO. Replace valve if device shows increased pressure.
	Leaking system switch	Attach pressure measuring device on CGO. Replace switch if device shows increased pressure.
Bellows leak	Pop-off valve diaphragm not sealing properly	Disassemble pop-off valve; inspect and clean seats; reseat; reassemble.
	Bellows mounting rim loose	Remove rim and pop-off valve diaphragm; reseat diaphragm; snap rim (2) into place.
	Bellows improperly mounted or has a hole or tear	Check that only the last bellows convolute is mounted to the rim and that the ring roll is in the groove under the rim. Inspect the bellows for damage; replace.
Breathing System Leak	Absorber canister open or missing	Install canister properly.
	Damaged/missing canister o-ring	Check/replace o-rings.
Breathing System Leak (Intermittent)	ACGO O <sub>2</sub> sense check valve	Replace.
Unable to begin mechanical	ABS not fully engaged	Remount ABS.
ventilation	No O <sub>2</sub> supply	Check O <sub>2</sub> supply.
	Defective Bag/Vent switch	Check Bag/Vent switch.

## 7.3 Troubleshooting Startup Screen (POST) messages — for DU

If the Avance system encounters a problem at startup to where it cannot initiate system software, a BIOS error message indicating the failure will be displayed.

Message	What it indicates	Troubleshooting Action Required
***NOTE: RTC failed. xx	This indicates an illegal value in the CPU's Real Time Clock date or time registers.	Set the system date and time manually to the correct values.
***NOTE: CMOS defaults reloaded.	This is an informational message.	No action needed.
***NOTE: CMOS battery failed. Replace soon.	Service is required to replace the CPU battery.	Replace the battery on the DU CPU board. Reload software and check out the system.
***NOTE: Alarm tone detection failed.	Service is required to correct a faulty connection to the speaker.	Reconnect the speaker if possible. Replace DU CPU board if speaker connection can not be corrected.
***NOTE: ECC failed. xx	The RAM memory ECC circuitry is faulty.	Replace the DU CPU board.
***NOTE: System reset. EC xx	This usually indicates a software (BIOS) failure.	Report this error, along with the machine logs, to Technical Support. Reload software. If problem persists, replace the DU CPU board.
***NOTE: Critical hardware failure. xx	Some of the DU hardware could not be initialized.	Replace the DU CPU board.
***NOTE: OS/App ROM CRC failed.	The software must be reloaded.	Reload software and check out the system.
***NOTE: System reset.	The software must be reloaded.	Reload Software and check out the system.

7-4 09/07 1009-0357-000

## 7.4 Troubleshooting Startup Screen (POST) messages — for HPDU

If the Avance system encounters a problem at startup to where it cannot initiate system software, a BIOS error message indicating the failure will be displayed.

Message	What it indicates	Troubleshooting Action Required
***ERROR: CPU data cache fault.	This indicates a hardware failure.	Replace the HPDU CPU board.
***ERROR: No bootable device available.	This indicates a problem with the internal CF card.	Check or replace the internal compact flash card.
***ERROR: Program load failed - CRC.	This usually indicates a software file corruption.	Reload the software and check out the system.
*** ERROR: RAM memory error.	This indicates a hardware failure.	Replace the HPDU CPU board.
***ERROR: System reset: ECxx xx xx	This usually indicates a software failure.	Report this error, along with the machine logs, to Technical Support. Reload software. If problem persists, replace the internal flash card and reload software. If problem persists, replace the HPDU CPU board.
*** ERROR: System reset: FFFF FF FF	Indicates the HPDU lithium battery has lost contact with the holder or the battery is below voltage.	Replace the lithium battery if below rated voltage. Install the software downloader card to restore system setups.
***ERROR: watchdog circuit failed.	This indicates a hardware failure.	Replace the HPDU CPU board.
***NOTE: Alarm speaker not detected. Check connection.	Service is required to correct a faulty connection to the speaker.	Reconnect the speaker if possible. Replace the HPDU CPU board if speaker connection can not be corrected.
***NOTE: CMOS battery is weak. Please replace.	Service is required to replace the CPU battery.	Replace the battery on the HPDU CPU board. Reload software and check out the system.
***NOTE: RTC date/time error. Battery may be weak.	Service is required to replace the CPU battery.	Replace the battery on the HPDU CPU board. Reload software and check out the system.

# **7.5 Troubleshooting the HPDU Display**

Symptom	Resolution
System will not boot from external Compact Flash card during software installation process	Verify that the Compact Flash card is properly inserted.     Insert a backup Compact Flash card.     Open the HPDU and verify that the external Compact Flash card carrier socket (1009-5961-000) is properly seated.     Replace external Compact Flash card carrier socket (1009-5961-000).     Replace HPDU main PCB.
Display appears mostly white and the green LED is on	<ol> <li>Verify that the cable connecting the HPDU to the system's rear panel is properly seated.</li> <li>Open the HPDU and verify that the cable connecting the main PCB to the display at J28 is properly seated within the mating housing and check for damage and/or wear.</li> <li>Replace the HPDU.</li> </ol>
Rotary encoder fails to work	Open the HPDU and verify that the cable connecting the main PCB to the rotary encoder at J38 is properly seated within the mating connector.      Verify that the revision of the flex cables connecting the main PCB to the keypads at J23 and J26 are at revision 101 or greater by inspecting the labels directly adhered to the flex cables.      Replace the rotary encoder,      Replace the HPDU main PCB.
Unit fails to boot and the green LED is on	<ol> <li>Open the HPDU and ensure that the lithium coin cell at J17 on the main PCB is properly installed.</li> <li>Replace the lithium coin cell at J17 on the main PCB.</li> <li>Verify that the internal Flash card at J10 on the main PCB is properly seated.</li> <li>Attempt to boot the system using a spare compact flash card and replace the card if necessary.</li> <li>Attempt to boot the system from the external Compact Flash card.</li> <li>Replace the HPDU main PCB.</li> </ol>
Excessive fan noise	Clean the HPDU fan inlet filter.     Check for obstructions within the external fan and verify source of fan noise.     Open the HPDU and verify internal CPU fan noise.     Replace internal or external fan assembly if causing excessive noise.

7-6 09/07 1009-0357-000

### 7.6 Troubleshooting System Malfunction (safe-state) screen

### **Machine logs**

- Compatibility incomplete: No versions received from Vent SIB
- System Self-tests failed

The above error log entries are due to the Ventilator failing its Power-On Self Test (Post).

When the ventilator fails its self-test, the system enters the safe-state.

The likely cause is a Gas Inlet Valve Solenoid that is not transitioning to the de-energized state when valve drive power is removed.

### **Troubleshooting above stated condition:**

- Replace the Gas Inlet Valve Solenoid (Section 9.10.2).
- If problem continues, replace the Gas Inlet Valve components (Section 9.10.3).

### 7.7 Breathing System Leak Test Guide

#### Note

Always do the **System "All Checks"** (Section 3.4) on the machine before proceeding with these breathing system leak tests.

- The **Low P leak check** looks for leaks between the mixer, vaporizer, and the inspiratory side of the breathing circuit.
- A failure in **Quick check** indicates a leak in Bag Mode.
- A failure in **Vent check** indicates a leak in Vent Mode.

Follow the troubleshooting flowcharts in Section 2.7.2 to determine the best sequence of tests for locating a breathing system leak.

The procedures in Section 7.7.3 test specific components of the breathing system for leaks.

### **⚠** WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Make sure that there are no test plugs or other objects caught in the breathing system.

7-8 09/07 1009-0357-000

# 7.7.1 Check Valves

Make sure that the check valves on the breathing circuit module work correctly: The Inspiratory check valve rises during inspiration and falls at the start of expiration. The Expiratory check valve rises during expiration and falls at the start of inspiration. A leak across one of the check valves may be great enough to cause a "reverse flow" alarm.


### Inspiratory check valve

- 1. Set the system switch to On.
- 2. Set fresh gas flow to minimum.
- 3. If equipped with an ACGO, connect a tube between the ACGO outlet and the Inspiratory port.
  - Set the ACGO switch to the ACGO position.
  - Verify that the Airway Pressure reading increases to 10 cm H<sub>2</sub>0 in 30 seconds.
- 4. If not equipped with an ACGO, select End Case and connect a tube to the Inspiratory port.
  - Stretch the tube approximately 5 cm.
  - Occlude the open end of the tube.
  - Release the tension on the tube.
  - Ensure that the Airway Pressure reading increases to between 20 and 40 cm H<sub>2</sub>0.
     If not, repeat the above steps, but stretch the tube a little further.
  - Verify that the Airway Pressure reading does not drop by more than 10 cm H<sub>2</sub>O in 30 seconds.

### **Expiratory check valve**

- 1. Set all gas flows to minimum.
- 2. Set the Bag/Vent switch to Bag.
- 3. Fully close the APL valve (70 cm H<sub>2</sub>0).
- 4. Connect a tube between the Inspiratory port and the Bag port.
- 5. Slowly increase the  $O_2$  flow to achieve 30 cm  $H_2O$ .
  - The leak rate is equal to the flow needed to maintain 30 cm H<sub>2</sub>O.
  - The leak rate should be less than 500 mL/min.


### 2.7.2 Breathing System Troubleshooting Flowcharts




	Vent Mode	Bag Mode
Software ≥ 3.00	Machine Check - System	Machine Check - Circuit
Software < 3.00	Vent Check	Quick Check

# Flowchart 1

7-10 09/07 1009-0357-000




### Flowchart 2




## Flowchart 3

7-12 09/07 1009-0357-000







# Flowchart **5**

7-14 09/07 1009-0357-000

# 7.7.3 Leak Isolation Tests

The previous flowcharts refer you to the following tests.

These tests require the use of the Low Pressure Leak Test Device and the Leak Test Tool Kit (refer to Section 10.1, "Service tools").

The Leak Test Tool Kit includes:

- the Machine Test Tool
- the Circuit Test Tool
- and various Test Plugs

When performing these tests on machines with an ACGO outlet, ensure that the ACGO selector switch is set to the ABS (Circle circuit) position.

#### Note

To perform most of these tests, you must boot the system with the PCMCIA Service Application and access the Diagnostics function as described in the test and detailed in Section 8.

Test 1:	Verifying the integrity of the test tools
Test 2:	Low-pressure leak testing the machine
Test 3:	Testing the airway pressure transducer, and Port 1 and Port 3 u-cup seals
Test 4:	Testing the bag port cover, the APL valve, the Bag/Vent switch, and the negative pressure relief valve 7-19
Test 5:	Testing the APL diaphragm
Test 6:	Testing the bellows module and the Bag/Vent switch
Test 7:	Testing the bellows, the bellows pop-off valve, the bellows base manifold, and the Bag/Vent switch 7-22
Test 8:	Testing the bellows assembly
Test 9:	Testing the flow sensor module, the circuit module, and the soda lime canister
<b>Test 10:</b>	Testing the circuit module and the canister
Test 11:	Testing the circuit module
Test 12:	Testing the inspiratory side of the circuit module
Test 13:	Testing the negative pressure relief valve
Test 14:	Testing the flow sensors only
Test 15:	Testing a flow sensor including the Ventilator Monitoring Assembly and interfacing components

### **⚠** WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:


- Do not use a test plug that is small enough to fall into the breathing system.
- Make sure that there are no test plugs or other objects caught in the breathing system.

### **⚠** CAUTION


Do not use  $\rm O_2$  Flush for leak isolation tests. Do not leave pressurized systems unattended. High pressure and equipment damage may result.

### Test 1 Verifying the integrity of the test tools

#### Machine Test Tool Front View



**Back View** 





- 1. Verify integrity of low-pressure leak test device.
  - Put your hand on the inlet of the leak test device. Push hard for a good seal.
  - Squeeze the bulb to remove all air from the bulb.
  - If the bulb completely inflates in less than 60 seconds, replace the leak test device.



- 2. Attach the low-pressure leak test device to the Machine Test Tool.
- 3. Plug the two pressure orifices.
- 4. Repeatedly squeeze and release the hand bulb until it remains collapsed.
- 5. If the bulb inflates in less than 30 seconds, locate and correct the leak.

7-16 09/07 1009-0357-000

### Test 2 Low-pressure leak testing the machine




- 1. Remove the breathing system from the machine.
- 2. Attach the Machine Test Tool (using only the Thru Port) and the low-pressure leak test device to **Port 3** of the breathing system interface as shown above.

**Note**: To prevent damage to the airway pressure transducer, ensure that the gauge port (**Port 1**) is not connected to the Test Tool.

- 3. On the Main Menu of the Service Application, select **Vent Diagnostics**.
- 4. On the Vent Diagnostics page, select **Status**:
  - Ensure that the *Circuit Setting* shows "Circle". If not, select *Toggle Circuit* in the Vent Diagnostics menu selections.
  - For machines with an ACGO outlet, ensure that the ACGO selector switch is set to the ABS (circle breathing circuit).
- 5. Turn off all vaporizers.
- 6. Compress and release the bulb until it is empty.
- 7. If the bulb completely inflates in 30 seconds or less, there is a leak in the low-pressure circuit.


# Test 3 Testing the airway pressure transducer, and Port 1 and Port 3 u-cup seals



- 1. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 2. On the Gas Diagnostics page, select **02 Flow**.
- 3. Set  $O_2$  Flow to **0.2 l/min**.
- 4. Attach the Machine Test Tool to the breathing system interface ports (using the alignment post) as shown above.
- 5. Occlude the tapered plug.
- 6. On the Gas Diagnostics page, reselect the  ${\it 02\,Flow}$  screen.
  - the Airway Pressure reading should increase.
  - If not, there is a leak in the tested circuit.
- 7. Set **02 Flow** to OFF.

7-18 09/07 1009-0357-000

Test 4 Testing the bag port cover, the APL valve, the Bag/Vent switch, and the negative pressure relief valve




- 1. Separate the Bellows Module from the Circuit Module and re-install the Bellows Module.
- 2. Occlude the Bag Port connector.
- 3. Connect the Machine Test Tool to the interface ports as shown above.
- 4. Set the Bag / Vent switch to Bag and close the APL Valve (70 cm H<sub>2</sub>0).
- 5. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 6. On the Gas Diagnostics page, select **02 Flow**.
- 7. Set  $O_2$  Flow to **0.2 l/min**.
- 8. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>O.

**Note**: If the bellows rises, it indicates a leak in the Bag / Vent Switch.

9. Set O2 Flow to OFF.

Test 5 Testing the APL diaphragm




**Note** If required, set up the Machine Test Tool and breathing system as shown in Test 4.

- 1. Slide the Bellows Module away from the machine.
- 2. Remove the APL ramp and diaphragm.
- 3. Insert a Test Plug into the APL scavenging port, as shown above.
- 4. Slide the Bellows Module partially back onto the machine casting.
- 5. Ensure that the Bag Port is plugged and that the Bag/Vent switch is set to Bag.
- 6. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 7. On the Gas Diagnostics page, select **02 Flow**.
- 8. Set  $O_2$  Flow to **0.2 l/min**.
- 9. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.

10. Set **02 Flow** to OFF.

7-20 09/07 1009-0357-000






- 1. Separate the Bellows Module from the Circuit Module and re-install the Bellows Module.
- 2. Connect the Machine Test Tool to the interface ports as shown above.
- 3. Set the Bag/Vent switch to the Vent position.
- 4. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 5. On the Gas Diagnostics page, select **Breathing Sys Leak** to access the Breathing Sys Leak screen (Section 8a.2.6).
- 6. Set  $O_2$  Flow to **0.2 l/min**.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.
- 7. Press the ComWheel to Exit the Breathing Sys Leak test.
- 8. Set **02 Flow** to OFF.

Test 7 Testing the bellows, the bellows pop-off valve, the bellows base manifold, and the Bag/Vent switch





- 1. Separate the Bellows Module from the Circuit Module.
- 2. Insert appropriate test plugs into the bellows base manifold as shown to the left.

**Note**: Position the bellows assembly so that the bellows remain collapsed as you plug the ports.

- 3. Set Bag/Vent switch to Vent.
- 4. Position the bellows upright with the bellows collapsed.
- 5. Connect the Machine Test Tool to the interface ports as shown above.
- 6. On the Main Menu of the Service Application, select Gas Diagnostics.
- 7. On the Gas Diagnostics page, select **02 Flow**.
- 8. Set  $O_2$  Flow to **0.2 l/min**.
- 9. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.

10. Set **02 Flow** to OFF.

7-22 09/07 1009-0357-000

Test 8 Testing the bellows assembly



**Note** 


If required, set up the Machine Test Tool and breathing system as shown in Test 7.



2. Insert appropriate test plugs into the bellows base manifold as shown to the left.

**Note**: Position the bellows assembly so that the bellows remain collapsed as you plug the ports.

- 3. Connect the tapered plug of the Machine Test Tool to the bellows base inlet as shown to the left.
- 4. Position the bellows upright with the bellows collapsed.
- 5. On the Main Menu of the Service Application, select **Gas Diagnostics**.
- 6. On the Gas Diagnostics page, select **02 Flow**.
- 7. Set  $O_2$  Flow to **0.2 l/min**.
- 8. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.
- 9. Set **02 Flow** to OFF.



Test 9 Testing the flow sensor module, the circuit module, and the soda lime canister



- 1. Separate the Bellows Module from the Circuit Module and re-install the Circuit/Flow Sensor Module.
- 2. Connect short tubing between the inhalation and exhalation ports of the breathing system.
- 3. Insert an appropriate test plug in the outlet port of the Circuit Module.
- 4. On the Main Menu of the Service Application, select Gas Diagnostics.
- 5. On the Gas Diagnostics page, select **02 Flow**.
- 6. Set  $O_2$  Flow to **0.2 I/min**.
- 7. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.
- 8. Set O2 Flow to OFF.
- 9. Remove the plug to release pressure.


7-24 09/07 1009-0357-000

### Test 10 Testing the circuit module and the canister



- 1. Remove the Flow Sensor module.
- 2. Connect the Circuit Test Tool to the Circuit Module as shown above.
- 3. Set  $O_2$  Flow to **0.2 I/min**.
- 4. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.
- 5. Set O2 Flow to OFF.


### Test 11 Testing the circuit module



**Note**: If required, set up the machine as in Test 10.

- 1. Remove the Soda Lime Canister.
- 2. Using appropriate Test Plugs, plug the three canister ports in the Circuit Module as shown above.
- 3. Set  $O_2$  Flow to **0.2 l/min**.
- 4. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.
- 5. Set O2 Flow to OFF.

Test 12 Testing the inspiratory side of the circuit module



Note: If required, set up the machine as in Test 10 and 11.


- 1. Connect the Circuit Test Tool to the Circuit Module as shown above.
- 2. Insert an appropriate test plug in the inspiratory outlet to the canister as shown above.
- 3. Set  $O_2$  Flow to **0.2 I/min**.
- 4. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.

5. Set O2 Flow to OFF.


7-26 09/07 1009-0357-000

### Test 13 Testing the negative pressure relief valve

- 1. Separate the Bellows Module from the Circuit Module.
- 2. Remove the Bellows Interface Manifold.
- 3. Insert test plug (recessed end) into the rear Bag/Vent switch port as shown.




- 4. Install the Bellows Module.
- 5. Connect the Machine Test Tool to the interface ports and the Bellows Module as shown above.



- 6. Set the Bag/Vent Switch to Vent.
- 7. Set  $O_2$  Flow to **0.2 I/min**.
- 8. On the Gas Diagnostics page, reselect the **02 Flow** screen.
  - Ensure that the Airway Pressure rises to  $\geq$  30 cm H<sub>2</sub>0.
- 9. Set *02 Flow* to OFF.

Test 14 Testing the flow sensors only

Note: To ensure a air-tight seal, use the corresponding plug as illustrated for the original flow sensor (A) or the new, moisture resistant (offset) flow sensor (B).




- 1. Remove the Flow Sensor Module.
- 2. Plug each Flow Sensor as shown above.
- 3. Connect the low-pressure leak test device to the open end of the Flow Sensor.
- 4. Block the connector end of the Flow Sensor with your hand.
- 5. Compress and release the bulb until it is empty.
- 6. If the bulb inflates in 30 seconds or less, there is a leak in the flow sensor.
- 7. If there are no leaks in the flow sensors, go to Test 15.

7-28 09/07 1009-0357-000


Test 15 Testing a flow sensor including the Ventilator Monitoring Assembly and interfacing components

Note: To ensure a air-tight seal, use the corresponding plug as illustrated for the original flow sensor (A) or the new, moisture resistant (offset) flow sensor (B).



- 1. Remove Flow Sensors from the Flow Sensor Module.
- 2. Attach the Flow Sensor to the bulkhead connector.
- 3. Plug each Flow Sensor as shown.
- 4. Connect the low-pressure leak test device to the open end of the Flow Sensor.
- 5. Compress and release the bulb until it is empty.
- 6. If the bulb inflates in 30 seconds or less, there is a leak. The leak may be through the connector o-rings, in the internal tubing, or in the Transducer on the VIB.


### 7.8 System Troubleshooting Flowcharts






7-30 09/07 1009-0357-000

### **Display Troubleshooting**



## Flowchart 7


### **Inaccurate Volume Ventilation Troubleshooting**



Flowchart 8


7-32 09/07 1009-0357-000

### **No Ventilation Troubleshooting**



Flowchart **9** 


### **High Intrinsic PEEP Troubleshooting**



Flowchart 10


7-34 09/07 1009-0357-000

### Alternate 02 / System Malfunction Screen Troubleshooting



# Flowchart 11


#### **Anesthesia Control Board Troubleshooting**



## Flowchart 12


7-36 09/07 1009-0357-000

### **ACB - Mixer Communication Troubleshooting**



# Flowchart 13


### **DU - ACB Communication Troubleshooting**



Flowchart 14


7-38 09/07 1009-0357-000

### **Mixer Specific Troubleshooting**



## Flowchart 15


#### **Power - Valves - Mixer Troubleshooting**



## Flowchart 16

7-40 09/07 1009-0357-000

#### **Valves - Mixer Troubleshooting**



# Flowchart 17

# 7.9 System Malfunction and Alt $\mathbf{0_2}$ Flowchart Table

Error	Display Type	Flow Chart
+ 12.0V H AMPS Gas SEL Valves +12.0V H AMPS ALT 02 +12.0V H AMPS MIXER +12.5V TO ACB	Alternate 02 Screen	Power - Valves - Mixer Troubleshooting (Flowchart 16)
ACB CLOCK SPEED	System Malfunction	Anesthesia Control Board Troubleshooting (Flowchart 12)
ACB COM FAIL	System Malfunction	DU - ACB Communication Troubleshooting (Flowchart 14)
ACB CPU TEST FAIL	System Malfunction	Anesthesia Control Board Troubleshooting (Flowchart 12)
ACB DCB COM FAIL	System Malfunction	DU - ACB Communication Troubleshooting (Flowchart 14)
ACB EEPROM FAIL ACB FLASH FAIL ACB HW WATCHDOG ACB MICROPROC ERROR ACB RAM ERROR ACB RAM MEMORY TEST FAILURE ACB REDUNDANT MEMORY FAIL ACB SW ERROR ACB SW WATCHDOG ACB UNEXPECTED RESET	System Malfunction	Anesthesia Control Board Troubleshooting (Flowchart 12)
ALT 02 SWITCH FAIL	Alternate 02 Screen	Check Alt O2 Switch Harness and Connections
DCB RAM ERROR	System Malfunction	DU - ACB Communication Troubleshooting (Flowchart 14)
FLOW ATTAIN FAIL BAL GAS CHANNEL FLOW ATTAINMENT FAILURE CH1	Alternate 02 Screen	Power - Valves - Mixer Troubleshooting (Flowchart 16)
FLOW SENSOR FAIL BALANCE GAS CH FLOW SENSOR FAILURE CH1 FLOW VERIF FAIL (dP) XX CH	Alternate 02 Screen	Mixer Specific Troubleshooting (Flowchart 15)
LOSS OF GAS DELIVERY USER SETTINGS LOSS OF VENT DELIVERY USER SETTINGS	System Malfunction	DU - ACB Communication Troubleshooting (Flowchart 14)
MIXER ACB COM FAIL MIXER ACB COM FAILLOST CMD	Alternate 02 Screen	ACB - Mixer Communication Troubleshooting (Flowchart 13)
MIXER BAL GAS CHANGE FAIL	Alternate 02 Screen	Power - Valves - Mixer Troubleshooting (Flowchart 16)

7-42 09/07 1009-0357-000

Error	Display Type	Flow Chart
MIXER CRC EEPROM FAIL MIXER CRC FLASH FAIL MIXER CRC RAM FAIL MIXER P SENSOR P2 FAIL	Alternate O2 Screen	Mixer Specific Troubleshooting (Flowchart 15)
MIXER POST FAIL	Alternate 02 Screen	See Related Errors in Error Logs
MIXER PRES SENSOR P1P3 FAIL MIXER SW WDOG FAIL	Alternate 02 Screen	Mixer Specific Troubleshooting (Flowchart 15)
MIXER VOLTS FAIL	Alternate 02 Screen	Power - Valves - Mixer Troubleshooting (Flowchart 16)
MIXER XX F SENSOR FAIL	Alternate 02 Screen	Mixer Specific Troubleshooting (Flowchart 15)
MIXER XX FLOW CHECK FAIL	Alternate 02 Screen	Power - Valves - Mixer Troubleshooting (Flowchart 16)
MIXER XX FLOW FAIL	Alternate 02 Screen	Mixer Specific Troubleshooting (Flowchart 15)
MIXER XX SELECTION VLV FAIL	Alternate 02 Screen	Power - Valves - Mixer Troubleshooting (Flowchart 16)
MIXER XX TEMP LIMIT MIXER XX TSENSOR FAIL MIXER XX ZERO FAILLEAK TEMP MEASUREMENT CH1 HIGH TEMP MEASUREMENT CH2 HIGH TEMP SENSOR FAIL BALANCE GAS CH TEMP SENSOR FAILURE CH1	Alternate O2 Screen	Mixer Specific Troubleshooting (Flowchart 15)
XX CHK VALVE LEAK	None	Reference TB ADV MSN 04 011
XX FLOW CTRL FAIL XX PROP VALVE FAIL	Alternate 02 Screen	Power - Valves - Mixer Troubleshooting (Flowchart 13)

#### 7.10 Technical Alarms

The Error Log includes technical alarms and other error conditions reported by the system.

A technical alarm, as apposed to a parameter alarm, is an alarm condition that exists whether or not a patient is connected to the machine. Technical alarms include:

- Failed state alarms internal problem prevents normal operation
- Ventilator failure alarms
- Vent Fail. Monitoring Only alarms
- Alternate O<sub>2</sub> state alarms can be caused by multile issues. The error logs will indicate the failure.

Alarms that do not fit into any particular category but are technical in nature are referred to as a Status alarms in this table.

Source table: AC = Anesthesia Computer

DC = Display Controller
Mixer = Electronic Gas Mixer
PC = Power Controller
Vent = Ventilator Interface

Log entries may be preceded with the subsystem that reported the error condition (example: "ACB: +12.0V H AMPS ALT 02"). The Error Log entry listed below does not include the subsystem prefix that reported the alarm (ACB:).

Error Log Entry		Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria	
	Ac	tion/Troubleshooting					
+12.0V H AMPS GAS SEL VALVES		Alternate O <sub>2</sub> Screen.	AC detected high current to the Gas Select Valves.	Medium	AC	Fresh gas select valves +10VA is turned On.	
		sconnect the flex cable from the from the from the resistance of $\epsilon$	om each three-way and NC gas s	select valv	es.		
		should be approximately					
+12.0V H AMPS ALT 02		Alternate O <sub>2</sub> Screen.	AC detected high current.	Medium	AC	Alternate O <sub>2</sub> valve +10VA turned On.	
	Dis	sconnect the flex cable fr	om the NO Alternate $O_2$ valve.				
	Мє	easure the resistance of t	he NO Alternate O <sub>2</sub> Bypass Valve	e:			
	•	should be approximately	75Ω				
+12.0V H AMPS MIXER		Alternate O <sub>2</sub> Screen.	Status bit shows current high.	Medium	AC	Mixer +10VA turned On.	
	In Service Software, under "Anesthesia Control Board Power Diagnostics" (page 2 of Section 8a.1.4), observe that <b>Mixer 10VA Amps</b> is reported as <b>Fail</b> .						
	Turn off power to the machine and disconnect the system interface harness from the Mixer.						
	If the <i>Mixer 10VA Amps</i> is now reported as <i>OK</i> ,						
	• replace the Mixer.						
	If the <i>Mixer 10VA Amps</i> is still reported as <i>Fail</i> ,  • inspect the harnesses from the ACB to PCB and PCB to Mixer for cross connections or damaged pins.  • replace the ACB.						

7-44 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria			
	Action/Troubleshooting							
+12.0V H AMPS VENTSIB	Ventilator failure!	Status bit shows current high.	High	AC	Ventilator Interface board 10VA is turned on.			
	Reboot system.  If problem continues, re	place VIB.						
+12.0V H AMPS MGAS	Gas monitoring not available	Status bit shows current high.	Medium	AC	MGAS 10 VA is turned on after 3 consecutive "ACB: +12.0V H AMPS MGAS" error log messages			
	Note: Single occurrence Remove Gas Module fro If the problem continued If the message disappear	MGAS" is not the same as "ACB: +1 es of "ACB: +12.0V H AMPS MGAS" m the Module Bay. s, replace the M-Gas Monitoring boa ars when module is removed, repair al Reference Manual for repair instru	require no ard. the M-Ga	action.				
+12.0V H AMPS VENT&OUTLET VALVES	Vent Fail. Monitoring Only.	Status bit shows current high.	Medium	AC Vent	Ventilator valves +10VA is On.			
	Disconnect the GIV and Insp Flow Valve. Measure the resistance of each valve:  • should be approximately $25\Omega$ for the GIV and $75\Omega$ for the Insp Flow Valve.							
+12.0V L VOLT ALT 02		Indicates Low Volts to the Mixer Alt O2 Valve.		ACB				
	This error can be falsely triggered. If no machine issues exist, ignore this entry.  If machine issues exist, disconnect the flex cable from the NO Alternate O2 valve. Use a Multimeter to measure the resistance of the NO Alternate O2 Bypass Valve. It should be approximately 75Ω.							
+12.0V L VOLT VENT&OUTLET VALVES		Indicates the Gas Inlet Valve Solenoid or Inspiratory Flow Control Valve is drawing too much power.		ACB				
	Disconnect the GIV and Insp Flow Valve.  Measure the resistance of each valve:  • should be approximately 25Ω for each valve.							
+12.0V L VOLT VENTSIB		Indicates low voltage to the VIB.		ACB				
	Use the Service Application to isolate VIB from cable.							
+12.5V TO ACB	Alternate O <sub>2</sub>	<11.9 or > 12.9 Vdc	High	AC - DC checks the service state.				
	Reboot system. If problem continues, replace the Power Controller board.							

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>		
	Action/Troubleshooting						
+5V H AMP GAS SUPPLY XDUCERS	Cannot read gas supply pressures	Status bit shows current high.	Medium	AC	Pressure transducer +10VA turned On.		
	<ul><li>2. Reboot system.</li><li>If problem continues, re</li><li>If error is no longer pres</li></ul>	, disconnect all gas supply presson place the ACB. ent, set system to Standby and re ck for error with each transducer	econnect o	ne pressu	re transducer at a time.		
12 HR TEST	Turn power Off and On for self tests	System has been operating for longer than 12 hours without a power up self test.	Low	AC - Vent DC checks enable criteria	System state is in Checkout.		
	At next available time, mov On position.	ve the system switch from the On	position to	the Off p	osition, then back to the		
ACB 4.096V ADC REF	Cannot monitor gas supplies	<4.018 or > 4.176 Vdc	Low	AC			
	Reboot System. If problem continues, replace the ACB.						
ACB ADC FAIL	Cannot monitor gas supplies	ADC timeout on any MUX channel.	Low	AC			
	Reboot System. If problem	continues, replace the ACB.					
ACB CLOCK SPEED	System Malfunction	AC clock frequency is > 1.1*(expected value) or <0.9*(expected value). AC clock frequency incorrect.	High	AC			
	Reboot System. If problem	continues, replace the ACB.	1				
ACB COM FAIL	System Malfunction	After establishing initial communication, the DC does not receive any messages from AC in 10 sec.	High	DC			
	Reboot System. If problem	continues, replace the ACB.	•				
ACB CPU TEST FAIL	System Malfunction	CPU instruction Test Failure	High	AC			
	Reboot System. If problem continues, replace the ACB						
ACB DCB COM FAIL	System Malfunction	The Anesthesia Computer receives no system state messages from the Display Computer for 10 seconds.	High	AC	Initial communications established.		
	Reboot System. If problem continues, check the ACB to Mixer communication LED's (or VIB communication LED's). If the RCV and XMT (or TXD and RXD) LED's indicate activity, check DU cable connections, replace Display Controller PCB if problem continues. If the RCV and XMT (or TXD and RXD) LED's indicate no activity, check the Anesthesia Control board connection, replace the Anesthesia Control board if problem continues.						

7-46 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria
	Action/Troubleshooting				
ACB EEPROM FAIL	Memory (EEPROM) failure	Read/Write failure or CRC failure of the EEPROM located on the Anesthesia Control Board.	Low	AC	
	Reboot System. If problem	n continues, replace the ACB.			
ACB FLASH FAIL	System Malfunction	CRC Failure in code space.	High	AC	
	Reboot System. If problem	n continues, replace the ACB.			
ACB HW WATCHDOG	System Malfunction	Hardware watchdog fails boot up test, times out, or detects an incorrect code sequence	High	AC	
	Reboot System. If problem	n continues, replace the ACB.			1
ACB MICROPROC ERROR	System Malfunction	Unexpected microcontroller exception (bus error, address error, etc.)	High	AC	
	Reboot System. If problem	n continues, replace the ACB.			
ACB RAM ERROR	System Malfunction	Memory Test Failure, Multiple bit errors detected.	High	AC	
	Reboot System. If problem	n continues, replace the ACB.			
ACB REDUNDANT MEMORY FAIL	System Malfunction	A redundantly stored parameter could not be stored properly or was corrupted.	High	AC	
	Reboot System. If problem	n continues, replace the ACB.			
ACB SW ERROR	System Malfunction	Unexpected software error	High	AC	
	Reboot System. If problem If problem continues, repl	n continues, reload ACB Software ace the ACB.			
ACB SW WATCHDOG	System Malfunction	Software watchdog failed power-up test, timed out, or a software function was delinquent for too long.	High	AC	
	Reboot System. If problem	n continues, replace the ACB.			
ACB UNEXPECTED RESET	System Malfunction	Unexpected reset of AC	High	AC	
	Reboot system. If problem	continues, replace the ACB.			
ACG0	Vol and Apnea monitoring off	Non Circle (ACGO) selected	Low	DC	System has ACGO
	No Service Action Require	d.			
ACMains POWER FAIL	Plug in power cable. On battery	ACMains_GOOD goes and stays low for at least 300 msec (3 software loops)	Medium	PC	30 minutes of battery power available.
	No Service Action Require	d.	!		1
AIR PIPE INVALID	Cannot monitor Air pipeline	Air Pipeline pressure is invalid.	Medium	DC	

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria		
	Action/Troubleshooting						
	Check Air Pipeline Supply. Check/Replace Air Pipelin	e Pressure Transducer.					
AIR PRESS LOW	Air supply pressure low	Air pipeline pressure is less than 252 kPa and the air cylinder has a pressure less than 2633 kPa for one second.	Medium	AC, DC	Air is selected as the balance gas with a non zero flow of air or the ventilator uses air as the drive gas and mechanical ventilation is ON		
	Check Air Supply. Check/Replace Air Pipelin	e/Cylinder Pressure Transducer					
AIR PRESS LOW DURING 21% 02	Air pressure low. Increase O2%.	Air pipeline pressure is less than 252 kPa and the air cylinder has a pressure less than 2633 kPa for one second.	High	AC DC	21% O2 (Air) is selected for fresh gas flow		
	Check Air Supply. Check/Replace Air Pipelin	e/Cylinder Pressure Transducer					
AIRWAY SENSOR CAL ERROR	Calibrate flow sensors	Airway Pressure Sensor zero offset out of range	Low	AC, Vent	Flow sensor detected		
	In Service Software, under "Vent Flow and Pressure Diagnostics" (Section 8a.3.2), verify the Airway Pressure counts is 800 ± 250.  Disconnect the Black in-line connector in the Patient Airway. If the counts return within specified range, check for occlusions in the Bulkhead harness.  If the counts do not return within the specified range, replace the VIB.						
ALT 02 SWITCH FAIL		Alternate O <sub>2</sub> switch status indicates Alt O <sub>2</sub> switch fault. The fault detection condition must persist for 1 second.	Medium	AC			
	Replace the Alt $O_2$ Switch.	<u> </u>					
AUX OUTLET FAIL	No fresh gas flow?	The measured SCGO position does not match commanded position.	High	AC, Vent			
	In Service Software, under "Vent Status" (Section 8a.3.1), view the Circuit Feedback status.  If the feedback indicates "Fault", toggle the Circuit. If the Status changes to match the Circuit setting, check/replace the SCGO/ACGO microswitches.						
BACKUP MODE ENTERED	Backup Mode active	No spontaneous breaths in set period of time (Backup Time (sec)) and 30 seconds has elapsed since starting PSVPro mode.	Low	DC			
	No spontaneous breaths in PSVPro mode. No Service Action Require	n set period of time (Apnea time) d.	and 30 se	conds has	s elapsed since starting		

7-48 09/07 1009-0357-000

Error Log Entry		Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria			
	Ac	Action/Troubleshooting							
BAL CHANNEL PROP VALVE LEAK FAIL		Alternate 02	Likely caused by a leaky Balance Proportional Valve.		Mixer				
		place the proportional v place Mixer if issue cont							
BAL FLOW CTRL FAIL			Mixer status bit STS_FLOW_CTRL_CH2_FAIL indicates flow attainment failure.	Medium	AC, Mixer	Balance gas supply pressure OK			
	Re	boot System. If problem	continues, replace the Mixer.		•				
BAL PROP VALVE FAIL			Mixer status bit STS_CH2_PROPN VALVE FAIL indicates proportional valve failure (over current, etc.)	Medium	AC, Mixer				
	Re	boot System. If problem	continues, replace the Mixer.						
BATT V VERY LOW		Plug in power cable. On battery	Available battery power decreases to between 10 and 5 min	Medium	PC	AC Mains Power Failure in progress.			
	Leave the system plugged in to charge the battery.  If problem continues, check the battery charge circuit in Service Software.  Replace Battery.								
BATTERY < 1MIN		System shutdown in <1min	Available battery power is <1min	High	PC	AC Mains Power Failure in progress.			
	Leave the system plugged in to charge the battery.  If problem continues, check the battery charge circuit in Service Software.  Replace Battery.								
BATTERY CHARGE FAIL		No battery backup	The system is in standby and the battery charge current is >4.0 amps.  or The system is powered on with a battery current >1.7	Medium	PC				
			amps.						
		eck the battery charge of place Battery.	ircuit in Service Software.						
BATTERY EMPTY		System shutdown in <5 min	Available battery power is between 1 and 5 minutes	High	PC	AC Mains Power Failure in progress.			
	lf p	ave the system plugged problem continues, chec place Battery.	in to charge the battery. k the battery charge circuit in Se	rvice Softv	vare.				

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria		
	Action/Troubleshooting						
BATTERY FAIL	No Battery Backup.	Battery voltage <10.5 V or While in bulk, over, or float charging battery is <10.5VDC or Battery has been bulk charging for >12 h in Standby or 24 h while powered on. or Voltage > 16.5V during bulk or over charging and normal current >0.25 Amps	Medium	PC			
	Leave the system plugged in If problem continues, check Replace Battery.	the battery charge circuit in Se	rvice Softw	vare.			
BATTERY LOW	Plug in power cable. On battery	Available battery power decreases to between 20 and 30 min	Medium	PC	Mains AC Mains Power Failure in progress.		
	Leave the system plugged in to charge the battery.  If problem continues, check the battery charge circuit in Service Software.  Replace Battery.						
BATTERY MISSING	No battery backup	Any battery voltage is between ±1.0 VDC.	Medium	PC	POST state		
	Connect Battery Leave the system plugged in to charge the battery. If problem continues, check the battery charge circuit in Service Software. Replace Battery.						
BATTERY REVERSED CONNECTIONS	No battery backup	Any battery voltage is less than –1.0 VDC	Medium	PC			
	Check Battery Connections						
BATTERY V LOW	Plug in power cable. On battery	Available battery power decreases to between 10 and 20 minutes	Medium	PC	AC Mains Power Failure in progress.		
	Leave the system plugged in to charge the battery.  If problem continues, check the battery charge circuit in Service Software.  Replace Battery.						
BELLOWS COLLAPSED	Unable to drive bellows	Manifold pressure is > Paw + 10 + [0.25*(Insp valve flow)]	Low	AC, Vent	In range Paw and manifold pressure data available and mechanical ventilation On.		
	Check the breathing circuit for leaks or hose occlusions. Perform flow sensor calibration. Check drive gas check valve. Check VIB cabling. Replace VIB.						

7-50 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria		
	Action/Troubleshooting						
BREATHING SYSTEM NOT LATCHED	Breathing system loose	Breathing system detection switch indicates breathing system not latched.	Low	AC, Vent			
	Check/replace ABS On sw Check/replace harness (A	ritch. BS switches to Filter board).					
CAL DATA FAILURE IN EEPROM	Service calibration advised	Default cal data is being used due to corrupt data in cal region.	Low	AC			
	Perform complete service	level calibrations (ventilator).					
CHECK FLOW SENSOR	Check flow sensors	During Mechanical breaths, the measured flow for 6 consecutive breaths, to and from the patient, does not meet certain criteria.  No or negative flow on Insp flow sensor during inspiration or negative flow on Exp flow sensor.	Medium	AC, Vent	In-range flow data available during mechanical ventilation		
	Check flow sensor connections. Check the breathing circuit. Check VIB sensor tubing for leaks. Perform flow sensor calibration. Check Insp/Exp check valves. Check/Replace flow sensors.						
Circuit check failed.		A message or failure displayed during the system checkout.		DU			
	Perform suggested action	and repeat the system Checkout					
Circuit O2 check skipped.		A message or failure displayed during the system checkout.		DU			
	Perform suggested action and repeat the system Checkout.						
COM ERROR VENT TO ACB	System Malfunction	After regular communications has been established between the Ventilator boundary object and the Vent SIB CPU, a total loss of communications shall be declared if the Ventilator boundary object receives no messages from the Vent SIB CPU for 35 milliseconds.	High	AC Vent			
	Reboot System. If problem 1. Check cabling. 2. Replace VIB. 3. Replace ACB.	n continues:					

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria
	Action/Troubleshoot	ting			
Compatibility failure: No version info in file for subsystem 0.		Indicates a subsystem did not report compatibility information to the Display Unit.		DU	
	Look for other entries Perform Software Dov	in the Error Logs. i.e. "Self-tests Faile vnload.	ď".		
Compatibility failure: Software Error		Indicates a subsystem did not report compatibility information to the Display Unit.		DU	
	Look for other entries Perform Software Dov	in the Error Logs. i.e. "Self-tests Faile vnload.	ď".		
Compatibility incomplete: No versions from Vent SIB		Indicates the Compatibility information for the Ventilator Interface Board does not match the Compatibility Table created during the last software download or the GIV Test did not pass.		DU	
	Check for other errors Perform Download Ne If persists, replace VII	ew to rebuild Compatibility Table.	1		
Compatibility incomplete: No versions received from Vent SIB		Indicates the Compatibility information for the Ventilator Interface Board does not match the Compatibility Table created during the last software download or the GIV Test did not pass.		DU	
	Replace the GIV Sole	of the Gas Inlet Valve operation.	d".		
COOLING FAN CURRENT LOW FAILURE		Indicates the Pan Fan is drawing too little current.		Mixer	
	Replace Fan. Replace Mixer.		1	I	1

7-52 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>				
	Action/Troubleshooting								
CPU FAN SPEED FAIL	Cooling fans failed. May overheat.	CPU fan speed less than 50% of nominal speed	Medium	DC					
	_	e CPU heatsink fan in the HPDU. over removed and verify CPU fan lugged in.	is not worl	king.					
CPU OVERHEAT	Cooling fans failed. May overheat.	Temperature reading of either DU thermistor > 60 degrees C	Medium	DC					
		e case fan in the HPDU. ve the fan filter from back of unit aat fan connector is plugged in.	and feel if	fan is work	king.				
DCB RAM ERROR	System Malfunction	Self test failure or multi bit error detected.	High	DC					
	Reboot System. If problem	Reboot System. If problem continues, replace the Display Controller Board.							
DC COMMANDED AC TO FAILURE		Indicates the Display Controller detected issues and commanded the Anesthesia Controller to safe state.		DC					
	Check for other errors in the error logs (Compatibility failure, Compatibility incomplete, etc.).								
DC: FRONT PANEL KEY STUCK		Indicates a stuck Keypad or encoder on the Display.		DU					
	Replace the Keypad / Encoder.								
DRIVE GAS LOST	Ventilator has no drive gas	O <sub>2</sub> supply low if O <sub>2</sub> is selected drive gas OR AIR supply low if Air is selected drive gas.	High	AC DC checks enable criteria	Mechanical Ventilation is ON.				
	Connect O <sub>2</sub> or AIR supply.								
	See Action/Troubleshooting	ng for O <sub>2</sub> PRESS LOW or AIR PRES	S LOW.						
EXP FLOW SENSOR CAL ERROR	Calibrate flow sensors	Exp Flow Sensor zero offset out of range	Low	AC, Vent	Flow sensor detected				
	In Service Software, under "Vent Flow and Pressure Diagnostics" (Section 8a.3.2), verify the Expiratory Flow counts is 2050 ± 250.  Disconnect the Blue and Yellow in-line connectors. If the counts return within specified range, check for occlusions in the Bulkhead harness.  If the counts do not return within the specified range, replace the VIB.								
EXP FLOW SENSOR EEPROM FAILURE	Replace exp flow sensor	EEPROM cal data read failure	Low	AC, Vent					
	Replace Exp Flow Sensor.	1	1	1	1				

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria			
	Action/Troubleshooting							
FAN FAIL	Cooling fan needs service. System OK	Fan Power Status Bit is Low (FAN1_GOOD).	Medium	PC	Communication between Power Controller and Display Computer.			
	Connect cooling fan. Replace cooling fan.							
FANS FAIL	Cooling fans failed. May overheat.	Both of the Fan Power Status Bits are Low (FAN1_GOOD, FAN2_GOOD)	Medium	PC	Communication between Power Controller and Display Computer.			
	Connect cooling fans. Replace cooling fans. Replace PCB.							
FLOW ATTAINMENT FAILURE CH1		Indicates the commanded flow through the O2 Gas channel does not match the measured flow via the flow sensor and the differential pressure transducers.		Mixer				
	Replace the O2 Proportional Valve. Replace the Mixer.							
FLOW SENSOR CAL ERROR	Calibrate flow sensors	Insp or Exp flow sensor or the airway or manifold pressure sensor zero offset out of range (flow calibration failure)	Low	AC, Vent	Flow sensor detected			
	See associated Errors. i.e.	EXP FLOW SENSOR CAL ERROR	or AIRWAY	SENSOR (	CAL ERROR			
FLOW VALVE CURRENT FAILURE		Indicates the current feedback from the Insp Flow Valve was incorrect for seven consecutive readings.		ACB				
	In the Service Software / Vo Flow Valve Current mA and	ent Flow & Pressure Diagnosis, i Counts.	ncrease th	e Flow Val	ve counts and view the			
FLOW_VALVE_DAC_F AILURE		Indicates the current feedback from the Insp Flow Valve was incorrect for seven consecutive readings.		ACB				
	In the Service Software / Vent Flow & Pressure Diagnosis, increase the Flow Valve counts and view the Flow Valve Current mA and Counts.							

7-54 09/07 1009-0357-000

Error Log Entry	Alarm	Text	Condition (Basic info)	Priority	Source	Enabling Criteria			
	Action/Troubleshooting								
FLOW VERIFICATION FAILURE (dP) CH1			Indicates the commanded flow through the O2 Gas Channel and the flow measured by the Hot-wire anemometer agrees but the flow as measured by the pressure transducers does not agree.		Mixer				
	Perform t Replace N	he Mixer zero. Mixer							
FRONT PANEL COM FAIL	Displa failure	y panel controls	Key pad controller fails to send "life tick" for greater than 10 Sec.	Medium (Yellow)	DC				
	Reboot system. If problem continues, replace Display Controller Board.								
GAS INLET VALVE BOOTUP TEST FAIL	Vent F Only	ail. Monitoring	Boot-up test failed.	High	AC, Vent				
	2. Replac	d the GIV compon							
INSP FLOW SENSOR CAL ERROR	Calibr	ate flow sensors	Insp Flow Sensor zero offset out of range.	Low	AC, Vent	Flow sensor detected			
	verify the Disconne occlusion	Inspiratory Pressuct the Black and V s in the Bulkhead	"Vent Flow and Pressure Diagnoure counts is 2050 ± 250. White in-line connectors. If the conharness. White the specified range, replace	ounts retur	n within s <sub>i</sub>	,			
INSP FLOW SENSOR EEPROM FAILURE	Repla senso	ce insp flow r	EEPROM cal data read failure	Low	AC, Vent				
	Replace t	he Inspiratory Flo	w Sensor.	•	•				

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>
	Action/Troubleshootin	g			
LOSS OF GAS DELIVERY USER SETTINGS		After regular communications has been established between the AC and the		AC	
LOSS OF VENT PARAMETER SETTINGS	Vent Fail. Monitoring Only	Display Computer, this alarm is declared if the system is in the Therapy State and the AC determines the Gas Delivery User Setting (Ventilator Parameter Settings) from the Display Computer arrived more than 10 seconds ago.			
	communication LED's).  1. If the RCV and XMT (o    check DU cable conn    replace Display Contr  2. If the RCV and XMT (o    check the Anesthesia	em continues, check the ACB to Mix r TXD and RXD) LED's indicate activi- nections. roller board if problem continues. r TXD and RXD) LED's indicate no ac Control board connection. ia Control board if problem continue	ity, ctivity,	nication L	ED's (or VIB
Low Pressure Leak check failed.		Indicates the LowP Leak section of the System Checkout Failed.		DU	
	Troubleshoot the Low Pr	essure Leak.			1
Low Pressure Leak check fails.		A message or failure displayed during the system checkout. Low Pressure Leak Check with SCGO failed automated check. Leak measured is greater than 50 ml/min.		DU	
	Check Vaporizer for leak Check integrity of low-pr	ressure circuit (Mixer outlet to SCGC	) / Insp Flo	w Sensor	).
MANIFOLD PAW SENSOR FAIL	Vent Fail. Monitoring Only	Calibration failure at bootup.	Medium	AC, Vent	
	verify the Manifold Flow Disconnect the White in range, check for occlusi	der "Vent Flow and Pressure Diagno counts is 800 ± 250. -line connector in the Manifold Pres ons in the Bulkhead harness. rn within the specified range, replac	ssure. If the	e counts re	

7-56 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria
	Action/Troubleshooting	g		-	
MANIFOLD PRESSURE SENSR FAILURE		Indicates a calibration failure at bootup.		ACB	
	800 ± 250. Disconnect the White in- range, check for occlusion	Vent Flow & Pressure Diagnostics line connector in the Manifold Presons in the Bulkhead harness. In within the specified range, replace	ssure. If the	e counts re	
MANIFOLD SENSOR CAL ERROR		Indicates the Manifold Pressure zero failed.		ACB	
	Could be caused by bad	span calibration or leaky vent Insp	Flow Valve	9.	
MGAS CHECK SAMPLE GAS OUT >20 SEC	Check sample gas ou	MGAS SPEC. Continuous Occlusion Bit set.	Medium	MGAS	MGAS present and MGAS communicates continuous occlusion for 20 seconds
	Replace sample line. Se	e AM TRM for further Troubleshooti	ng.		
MGAS INLET FILTER RESIDUE >40 SEC	Replace D-Fend	MGAS SPEC (Residue build- up on the water trap membrane. This decreases air flow).	Medium	MGAS	MGAS present and MGAS communicates this the Replace Trap alarm bit for 40 seconds
	Replace D-Fend. See AM	TRM for further Troubleshooting.			
MGAS LINE BLOCKED >20 SEC	Sample line blocked	MGAS SPEC states The sample tubing inside or outside the monitor blocked, or the water trap is occluded.	Medium	MGAS	MGAS present and MGAS communicates this the continuous occlusion alarm for 20 seconds
	Replace sample line. Se	e AM TRM for further Troubleshooti	ng.		
MGAS SAMPLE LINE NOT CONNECTED >40 S	Check D-Fend	MGAS SPEC states The sample tubing or the D-Fend module is not installed.	Medium	MGAS	MGAS present and MGAS communicates this the OpenGasCircuit alarm for 40 seconds
	Replace D-Fend. See AM	TRM for further Troubleshooting.		ļ.	
MGAS SENSOR INOP > XX	Module fail. No CO <sub>2</sub> , AA, O <sub>2</sub> data	MGAS SPEC Mgas communicates hardware failure (RAM failure; ROM checksum error; Error in CPU eeprom; Error O <sub>2</sub> preamp eeprom; Error in SSS board eeprom; Voltage error; Lamp control failure.) or UPI does not initialize.	Medium	MGAS	
	See AM TRM for further T	roubleshooting.	•	•	

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>
	Action/Troubleshooting	ng			
MIXER BAL GAS CHANGE FAIL	Alternate O <sub>2</sub> Screen	Mixer Status Bit: STSBALGAS_CHANGE_OV ER_FAIL After the mixer commanded a change to the balance gas, the status of the selector valve shows the old balance gas is still connected.	Medium	AC, Mixer	
	Reboot System. If prob	em continues, replace the Mixer.	•		
MIXER BALGAS Flow FAIL		Mixer error bit STS_CH2_DELTAP_FLOW_FA IL Pressure difference between P3 and P2 differs from the drop expected at the measured flow for channel 2 (Balance Gas).	Medium	AC, Mixer	
	Reboot System. If prob	em continues, replace the Mixer.			1
MIXER 02 Flow FAIL	Alternate O <sub>2</sub> Screen	Mixer error bit STS_CH1_DELTAP_FLOW_FA IL Pressure difference differs from the drop expected at the measured flow for Channel 1 $(0_2)$ .	Medium	AC, Mixer	
	Reboot System. If prob	em continues, replace the Mixer.			
MIXER ACB COM FAIL	Alternate O <sub>2</sub> Screen  Reboot System. If prob	measured flow data from the mixer. em continues,	Medium	AC	Communication has been established between mixer and AC.
	<ul> <li>replace the Mixer.</li> </ul>	onnector to Mixer cable.			
MIXER ACB COM FAILLOST CMD	Alternate O <sub>2</sub> Screen	. Mixer status Bit STS_LOSS_OF_SETFLOW_C MD.  Mixer has lost AC flow commands for 5 sec or received "illegal" commands.(hypoxic mix, settings not allowed)	Medium	AC, Mixer	
	Reboot System. If problem check/replace Pan Control replace the Mixer.	em continues, connector to Mixer cable.			1

7-58 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria			
	Action/Troubleshooting							
MIXER AIR SELECTION VLV FAIL	Alternate O <sub>2</sub> Screen.	Mixer Status Bit: STS_SELV_VAIR_NOTIFY_FAI L The status of the air selector valve does not match the commanded state.	Medium	AC, Mixer				
	Reboot System. If probler	n continues, replace the Mixer.						
MIXER BAL GAS F SENSOR FAIL	Alternate O <sub>2</sub> Screen.	Mixer error bit STS_F2_SENSOR_FAIL Balance gas flow sensor failure.	Medium	AC, Mixer				
	Reboot System. If probler	n continues, replace the Mixer.	'					
MIXER BAL GAS TSENSOR FAIL		Mixer error bit STS_T2_SENSOR_FAIL (Balance Gas). Balance gas temperature sensor failure.	Medium	AC, Mixer				
	Reboot System. If probler	n continues, replace the Mixer.	•					
MIXER BAL GAS FLOW CHECK FAIL		Mixer status bit  1LPM_FLOW_TEST_FAIL. Bal gas proportional valve fails flow check STS_FLOW_TEST_BAL_CHAN _FAILshows balance gas proportional valve failed self test.	Medium	AC, Mixer				
	Reboot System. If probler	m continues, replace the Mixer.	1					
MIXER BAL GAS TEMP LIMIT		Mixer error bit STS_CH2_TEMP_LIMIT (Balance Gas). Balance gas temperature exceeds 50 °C.	Medium	AC, Mixer				
	Reboot System. If probler	n continues, replace the Mixer.						
MIXER BAL GAS ZERO FAILLEAK		Mixer status bit STS_CH2_ZERO_FLOWPROP N_V_CH2_LEAK_FAIL_TEST_ FAIL. Bal gas proportional valve fails zero flow check shows flow while closed.	Medium	AC, Mixer				
	Reboot System. If probler	n continues, replace the Mixer.	-					
MIXER CRC EEPROM FAIL	Alternate 0 <sub>2</sub> Screen	Runtime CRC check on EEPROM failed. Mixer Status Bit STS_EEPROM_CRC_FAIL.	Medium	AC, Mixer				
	Reboot System. If probler	n continues, replace the Mixer.						

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria
	Action/Troubleshooting				·
MIXER CRC FLASH FAIL	Alternate 0 <sub>2</sub> Screen	Runtime CRC check on Flash failed. Mixer Status Bit STS_FLASH_CRC_FAIL	Medium	AC, Mixer	
	Reboot System. If probler	n continues, replace the Mixer.			
MIXER CRC RAM FAIL	Alternate O <sub>2</sub> Screen	Runtime CRC walking pattern check on RAM failed. Mixer Status Bit STS_RAMCRC_FAIL.	Medium	AC, Mixer	
	Reboot System. If probler	n continues, replace the Mixer.			
MIXER N20 SELECTION VLV FAIL		Mixer Status Bit: STS_SELV_VN2O_NOTIFY_FA IL The status of the N <sub>2</sub> O selector valve does not match the commanded state.	Medium	AC, Mixer	
	Reboot System. If probler	n continues, replace the Mixer.	'		
MIXER O2 TSENSOR FAIL	Alternate 0 <sub>2</sub> Screen.	Mixer error bit STS_T1_SENSOR_FAIL O <sub>2</sub> temperature sensor failure	Medium	AC, Mixer	
	Reboot System. If probler	n continues, Replace the Mixer.			
MIXER 02 F SENSOR FAIL	Alternate 0 <sub>2</sub> Screen	Mixer error bit STS_F2F1_SENSOR_FAIL $(O_2)$ . $O_2$ flow sensor fail.	Medium	AC, Mixer	
	Reboot System. If probler	n continues, replace the Mixer.	1		
MIXER 02 FLOW CHECK FAIL	Alternate O <sub>2</sub> Screen	Mixer status bit STS_FLOW_TEST_CH1_FAIL1 LPM_FLOW_TEST_FAIL. O <sub>2</sub> proportional valve fails flow check.	Medium	AC, Mixer	
	Reboot System. If probler	n continues, replace the Mixer.			
MIXER O2 SELECTION VLV FAIL	Alternate 0 <sub>2</sub> Screen	Mixer Status Bit: STS_SELV_VOXY_NOTIFY_FAI L The status of the O <sub>2</sub> selector valve does not match the commanded state.	Medium	AC, Mixer	
	Rehoot System If probler	n continues, replace the Mixer.			
MIXER O2 TEMP LIMIT	Alternate O <sub>2</sub> Screen	Mixer error bit STS_CH1_TEMP_LIMIT (O <sub>2</sub> ). O <sub>2</sub> temperature exceeds 50 °C.	Medium	AC, Mixer	
	Reboot System. If probler	n continues, replace the Mixer.			I

7-60 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>			
	Action/Troubleshooting							
MIXER O2 ZERO LEAK FAIL		Mixer status bit STS_CH1_ZERO_FLOW_TESP ROPN_V_LEAK_FAILT_FAIL. O <sub>2</sub> proportional valve fails zero flow checks for leaks when it should be closed.	Medium	AC, Mixer				
	Reboot System. If problem	n continues, replace the Mixer.						
MIXER P SENSOR P2 FAIL		Mixer error bit STS_PRESS_SENSOR_FAIL_ P2 Pressure sensor 2 in the mixer has failed.	Medium	AC, Mixer				
	Reboot System. If problem	n continues, replace the Mixer.						
MIXER POST FAIL	Alternate O <sub>2</sub> Screen	Mixer tells AC that Power Up Self Test Fail	Medium	AC, Mixer				
	See associated Error in Er FAILURE (dP) CH1"	ror Log. i.e. "MIXER O <sub>2</sub> FLOW CHE	CK FAIL" o	r "Mix: FLC	OW VERIFICATION			
MIXER PRES SENSOR P1P3 FAIL		Mixer error bit STS_PRESS_SENSOR_FAIL_ P1P3 One of the pressure sensors in the Mixer has failed (P1 or P3).	Medium	AC, Mixer				
	Reboot System. If problem	n continues, replace the Mixer.	•	•				
MIXER P SENSOR FAIL		Indicates one of the Mixer Pressure Transducers have failed		Mixer				
	See associated Error in Error Log. i.e. "MIXER P Sensor P2 Fail" or "Mixer Pres Sensor P1P3 Fail".							
MIXER SW WDOG		Indicates the Mixer Watchdog has been activated.		Mixer				
	See associated Error in Er	ror Log. i.e. "MIXER O2 FLOW CHE	CK FAIL" o	r"Mix: FL0	OW VERIFICATION			
MIXER SW WDOG FAIL	Alternate O <sub>2</sub> Screen	Mixer status Bit STS_SW_WDOG_FAIL.	Medium	AC, Mixer				
	Reboot System. If problem	n continues, replace the Mixer.	l	l	1			
MIXER VOLTS FAIL	Alternate O <sub>2</sub> Screen	Mixer power supply (on board) is out of tolerance. Status bit STS_VOLT_REF_FAIL.	Medium	AC, Mixer	+12.5 V (10 VA) to mixer OK.			
	view the "Mixer 10VA Volta	s "OK, and +12.5 Vdc reads "Fail" s "Fail",			,			

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria		
	Action/Troubleshooting				'		
MIX O2 BYPASS VLV 10VA OVER CURR		Indicates the current feedback from the O2 Bypass Selector was incorrect for seven consecutive readings.		ACB			
	Measure the O2 Bypass Se Replace ACB.	elector Valve. Should be approxin	nately 75 (	2.			
MODULE NOT COMPATIBLE	Module not compatible	The Monitoring Module detected is not compatible with system software. System is designed to work with the following Compact Airway Module versions: M-CaiO (HW rev 00 and above, SW rev 3.2 and above) and M-CaiOV (HW rev 00 and above, SW Rev 3.2 and above).	Low	DC			
	Replace M-Gas module wi	th compatible module.					
N20 PRESS LOW	N <sub>2</sub> O supply pressure low	$\rm N_2O$ pipeline pressure is less than 252 kPa and the $\rm N_2O$ cylinder pressure is less than 2633 kPa.	Low	AC	N <sub>2</sub> O is selected as the balance gas with a non zero flow of N <sub>2</sub> O		
	Check N <sub>2</sub> O Supply. Check / Replace N <sub>2</sub> O Pipe	line/Cylinder Pressure Transduce	er.				
NO EXPIRATORY FLOW SENSOR	No exp flow sensor	No Expiratory sensor connected and not calibrating	Medium	A, Vent			
	Connect Expiratory flow sensor. Check/Replace Bulkhead harness. Replace VIB Board.						
NO INSPIRATORY FLOW SENSOR	No insp flow sensor	No inspiratory sensor connected and not calibrating.	Medium	AC, Vent	AC -Vent		
	Connect Inspiratory flow so Check/Replace Bulkhead Replace VIB Board.						
02 PROP VALVE FAIL		Mixer status bit STS_CH1_PROPN VALVE FAIL indicates proportional valve failure.	Medium	AC, Mixer			
	Reboot System. If problem	continues, Replace the Mixer.		ı			
02 FLOW CTRL FAIL		Mixer status bit STS_FLOW_CTRL_CH1_FAIL indicates flow control failure.	Medium	AC, Mixer	O <sub>2</sub> gas supply pressure OK		
	Reboot System. If problem	continues, Replace the Mixer.					

7-62 09/07 1009-0357-000

Error Log Entry		Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>				
	Ac	tion/Troubleshooting								
02 FLUSH FAILURE		O <sub>2</sub> flush stuck on?	Switch is detected "on" continuously > 30 sec.	Low	AC, Vent					
	Ala	Alarm condition becomes false for 2 consecutive switch readings.								
02 PIPE INVALID		Cannot monitor O <sub>2</sub> pipeline	O <sub>2</sub> Pipeline pressure is invalid.	Medium	DC					
		eck O <sub>2</sub> Pipeline Supply. eck / Replace O <sub>2</sub> Pipelin	e Pressure Transducer							
02 PRESS LOW		O <sub>2</sub> supply pressure low	$\rm O_2$ pipeline pressure is less than 252 kPa and the $\rm O_2$ cylinder has a pressure less than 2633 kPa for one second.	High	AC, DC	N <sub>2</sub> O flow stops on threshold detection and Air continues to flow if selected.				
		eck O <sub>2</sub> Supply.	e/Cylinder Pressure Transducer							
02 SENSOR FAILURE	-	Replace O <sub>2</sub> sensor	0 <sub>2</sub> < 5%	Low	AC, Vent	Galvanic O <sub>2</sub> sensor connected				
	lf (	calibration fails, replace ( calibration continues to fa calibration fails after 90 r	ail, wait 90 minutes and repeat o	calibration						
O2CAL ERROR	If (	calibration fails after 90 r Calibrate O <sub>2</sub> sensor	Offset, slope, or cell voltage	Low	AC,	Galvanic O <sub>2</sub> sensor				
	_		not in range or 0 <sub>2</sub> > 110%		Vent	connected				
	Calibrate O <sub>2</sub> Sensor.									
	If calibration fails, replace $\rm O_2$ Sensor.  If calibration continues to fail, wait 90 minutes and repeat calibration.  If calibration fails after 90 minute, replace VIB.									
ON/STANDBY SWITCH TO STANDBY		Turn switch on to continue use	On/Standby switch transitions from On to Standby.	High	PC	System state is Therapy and Power Controller is communicating with DC				
	No	Service Action Required		!						
PATIENT VOLUME MISMATCH OCCURRED		Calibrate, dry, or replace flow sensors	PATIENT VOLUME MISMATCH alarm occurred.	Low	AC, Vent, DC checks enable criteria	System state is in Checkout.				
		Check flow sensor conne Replace flow sensors. Check the VIB tubing for Replace VIB.		1	1					

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria			
	Action/Troubleshoot	ing						
PCB Alarm Off: AC		Indicates AC Power		PCB				
lains failure		Removed.						
	No Service Action Req	uired.						
CB Alarm Off: DC-		Indicates the AC supply is		PCB				
OC power module		OK (AC GOOD HIGH) but the		1 05				
ailure		system reports using the						
		battery (BATT STAT 1 and 2						
		LOW).						
	Check U-frame wiring.				I			
	=	Standby; remove mains; wait 20 seco	nds; powe	r up systei	m.			
	If problem continues,	replace PCB.						
PCB Alarm Off:		Indicates the Battery		PCB				
Monitor Current		backup current to monitor is						
Active w/AC Mains		active while AC supply is OK.						
	Check AC power connection to anesthesia monitor and circuit breaker.							
	If problem continues,	replace anesthesia monitor.						
CB Alarm On: 1 min		A message or failure		PCB				
ime left		displayed during the system						
		checkout.						
	Leave the system plugged into AC Mains for 24 hours. If issue persists, replace batteries.							
PCB Alarm On: AC		Indicates AC Power		PCB				
nains failure		Removed.						
	No Service Action Required.							
PCB Alarm On: Bulk		Indicates the batteries were		PCB				
charge time		being Bulk Charged for a						
exceeded in stdby		period greater than 12						
		hours while the system was						
		in the Standby state.						
	Replace Batteries.	ı	1	1	1			
	Replace the Power Co	ntroller Board.						
		Indicates the Battery		PCB				
PCB Alarm On:		The state of the s	1		1			
Monitor Current		backup current to monitor is						
Monitor Current		backup current to monitor is active while AC supply is OK.						
Monitor Current	Check AC power conne		uit breake	r.				
Monitor Current	•	active while AC supply is OK.	uit breake	r.				
PCB Alarm On: Monitor Current Active w/AC Mains PCB Saved Alarm Off:	•	active while AC supply is OK. ection to anesthesia monitor and circ	uit breake	r.				
Monitor Current Active w/AC Mains	•	active while AC supply is OK. ection to anesthesia monitor and circ replace anesthesia monitor.	uit breake					
Monitor Current Active w/AC Mains PCB Saved Alarm Off:	•	active while AC supply is OK. ection to anesthesia monitor and circ replace anesthesia monitor.  A message or failure	uit breake					

7-64 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>			
	Action/Troubleshooting							
PCB Saved Alarm Off: AC Mains failure		Indicates AC Power Removed.		PCB				
	No Service Action Requir	ed.						
PCB Saved Alarm Off: Blk chrg time exceed in stdby		Indicates the batteries were being Bulk Charged for a period greater than 12 hours while the system was in the Standby state.		PCB				
	Replace Batteries. Replace the Power Contr	oller Board.						
PCB Saved Alarm Off: DC-DC power module failure		Indicates greater than 0.25 amps current draw out of the batteries for 2 minutes while AC Mains is connected to the machine.		PCB				
	Replace Batteries. Replace the Power Controller Board.							
PCB Saved Alarm On: 1 min time left		A message or failure displayed during the system checkout.		PCB				
	Leave the system plugged into AC Mains for 24 hours. If issue persists, replace batteries.							
PCB Saved Alarm On: AC Mains failure		Indicates AC Power Removed.		PCB				
	No Service Action Required.							
PCB Saved Alarm On: DC-DC power module failure		Indicates greater than 0.25 amps current draw out of the batteries for 2 minutes while AC Mains is connected to the machine.		PCB				
	Replace Batteries. Replace the Power Contr	oller Board.						
PCB Saved Alarm On: DU to PSC Comm Error		Indicates a the Power Controller to Display Unit communication was lost. The Power Controller "saved" the error and communicated the error to the Display Unit when communication was next established.		PCB				
	No Service Action Requir	red.						

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria
	Action/Troubleshootin	ng			
PCB Saved Error: POWER CNTRL COM FAIL	Check the Anesthesia C	Indicates communication between the Display Unit computer and the Anesthesia Controller Board (ACB) was lost and a Power Controller error occurred. The error was stored until communication could be reestablished and written to the Display Unit computer.	it commun	PCB	ble.
PCSELF TEST	Internal failure. Syst	em PC failed self tests (memory, voltages, or CPU).	High	PC	
	Reboot system. If proble Replace Power Controlle	em continues, check power supplies er Board if continues.	s in the Sei	rvice Soft\	ware.
PEEP PCV NOT AVAILABLE	Vol vent only. No PEI or PSV	Paw data is in range but the Pmanifold <= -15 cmH <sub>2</sub> 0	Medium Or Low	DC	None
DOWED CONTROLED	the transducer precision the Paw transducer.	bration. If calibration fails, use the S	pare linea	rity of the	_
POWER CONTROLER COM FAIL	Internal failure. Syst	em Communications with PC and DC cannot be established for ten seconds.	Medium	DC	
	Reboot system. If problem continues:  1. Check DU cable connections.  2. Check the Display Connector board cable connections.  3. Replace the Power Controller board.				
POWER SUPPLY 75C	Circuitry > 75C shutdown possible	Power supply temperature exceeds 75C.	Medium	PC	
	Check / Clean cooling f	an.			
PRESS SNSOR1 FAILURE		Indicates P1 pressure transducer on the Mixer is out of specification.		Mixer	
	Zero the Mixer pressure Replace Mixer.	transducers.			
PWR CNTRL DC-DC FAIL	Using battery. PC fai	AC supply is OK (AC GOOD HIGH) but the system reports using the battery (BATT STAT 1 and 2 LOW).	Medium	PC	
	If applicable, check U-fo Set system switch to Sta If problem continues, re	andby; remove mains; wait 20 seco	nds; powe	r up syste	m.

7-66 09/07 1009-0357-000

Error Log Entry	A	larm Text	Condition (Basic info)	Priority	Source	Enabling Criteria		
	Action/Troubleshooting							
Quick check fails.			A message or failure displayed during the system checkout.		DU			
	Chec	Check for other errors in the error logs.						
REVERSE EXPIRATORY FLOW		everse exp flow. heck valves OK?	Flow towards the patient (volume >= 20 mL) on expiratory sensor and flow towards the patient (volume >= 5 mL) on the inspiratory sensor during inspiration for 6 consecutive mechanical breaths.	Medium	AC Vent	In-range flow data available, mechanical ventilation on		
	Chec Perfo Chec Repla Chec	k flow sensor connection the breathing circuit.  It is the breathing circuit.  It is flow sensor calibrated the low sensor calibrated the low sensors.  It is the flow sensors the flow sensors the flow sensors.  It is the VIB cabling.	es.					
SCG0	m	ol and Apnea onitoring off	Non Circle SCGO selected.	Low	DC	System has SCGO		
	No se	ervice action required.						
SEVERE SUSTAINED PAW			Indicates the measures airway pressure was greater than 100 cm H20 for 10 seconds.		ACB			
	No Service Action. Reboot system. If problem continues, check Airway Pressure signal in Service Mode.							
STANDBY PATIENT DETECTION	N	o fresh gas flow!	3 volume breaths are detected within 30 seconds or 3 CO <sub>2</sub> breaths are detected within 30 seconds	High	DC	System in Checkout: General or Checkout: Start Case		
	No service action required.							
System Self-tests failed			Indicates the Power-on tests failed. Look for other entries for clarification.		DU			
	Chec	k for other errors in the	error logs (Compatibility failure	, Compatil	bility incor	mplete, etc.).		
VALVE CH1 LEAK TESTS NOT DONE			Indicates the O2 Leak Test skipped. Can be caused by no O2 supply connected at power- up.		Mixer			
	No S	ervice Action Required				<u> </u>		

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>		
	Action/Troubleshooting						
VAP-MIXER FAN FAIL	Cooling fan needs service. System OK	Fan current < 0.09 Amps > 0.25 Amps	Medium	Mixer	Communication between Display Computer - ACB and ACB - Mixer.		
	Connect rear panel fan. Replace rear panel fan. Replace Mixer if no voltage.						
Vent check stage 1 failed.		A message or failure displayed during the system checkout.		DU			
	Checkout failed. Check for other errors in the error logs.						
Vent check stage 1 fails.		A message or failure displayed during the system checkout.		DU			
	Checkout failed.Check for	other errors in the error logs.					
Vent check stage 2 failed.		A message or failure displayed during the system checkout.		DU			
	Checkout failed.Check for other errors in the error logs.						
Vent check stage 2 fails.		A message or failure displayed during the system checkout.		DU			
	Checkout failed. Check for other errors in the error logs.						
VENT FLOW VALVE FAIL DAC	Vent Fail. Monitoring Only	Incorrect DAC feedback for 3 consecutive readings	Medium	AC, Vent			
	Reboot System. If problem continues, in the Service Software / Vent Flow & Pressure Diagnosis, increase the Flow Valve counts and view the Flow Valve Feedback mV and Counts. Verify the settings match.						
VENT +12.5V FAIL	Vent Fail. Monitoring Only	Nominal 12.5V <11.3 Vdc or >13.13 Vdc	Medium	AC, Vent			
	In Service Software, under "Vent Interface Bd Power Diagnostics" (Section 8a.1.6), view the "Vent Int Bd 10VA Voltage" from Board Supplies:  If "Vent Int Bd 10VA Voltage" reads "OK, and +12.5 Vdc reads "Fail", replace the VIB.  If "Vent Int Bd 10VA Voltage" reads "Fail" and the "Vent Int Bd 10VA Voltage" from the Anes Cntrl Bd reads "OK", Check cabling between ACB and VIB.						
VENT +6V FAIL	Vent Fail. Monitoring Only	VSIB +6V out of range<5.51 Vdc or > 6.5 Vdc	High	AC, Vent	Vent +12.5 V (10 VA) is OK		
	In Service Software, under "Vent Interface Bd Power Diagnostics" (Section 8a.1.6), view the "Vent Int Bd 10VA Voltage" from Board Supplies:  If "Vent Int Bd 10VA Voltage" reads "OK, and +6.0Vdc reads "Fail", replace the VIB.  If "Vent Int Bd 10VA Voltage" reads "Fail" and the "Vent Int Bd 10VA Voltage" from the Anes Cntrl Bd reads "OK", Check cabling between ACB and VIB.						

7-68 09/07 1009-0357-000

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	<b>Enabling Criteria</b>		
	Action/Troubleshooting						
VENT 1.22V FAIL	Vent Fail. Monitoring Only	Voltage < 1.074Vdc or Voltage > 1.367 Vdc	Medium	AC, Vent	Vent +12.5 V (10 VA) is OK		
	In Service Software, under "Vent Interface Bd Power Diagnostics" (Section 8a.1.6), view the "Vent Int Bd 10VA Voltage" from Board Supplies:  If "Vent Int Bd 10VA Voltage" reads "OK, and 1.22 Vdc reads "Fail", replace the VIB.  If "Vent Int Bd 10VA Voltage" reads "Fail" and the "Vent Int Bd 10VA Voltage" from the Anes Cntrl Bd reads "OK", Check cabling between ACB and VIB						
VENT -6V FAIL	Vent Fail. Monitoring Only	VSIB -6V out of range<-6.72 Vdc or > -5.28 Vdc	High	AC, Vent	Vent +12.5 V (10 VA) is OK		
	In Service Software, under "Vent Interface Bd Power Diagnostics" (Section 8a.1.6), view the "Vent Int Bd 10VA Voltage" from Board Supplies:  If "Vent Int Bd 10VA Voltage" reads "OK, and -6.0Vdc reads "Fail", replace the VIB.  If "Vent Int Bd 10VA Voltage" reads "Fail" and the "Vent Int Bd 10VA Voltage" from the Anes Cntrl Bd reads "OK", Check cabling between ACB and VIB.						
VENT ADC VREF FAIL	Vent Fail. Monitoring Only	VSIB ADC3.200V ref voltage out of range <3.179 or >3.221 Vdc	High	AC, Vent	Vent +12.5 V (10 VA) is OK		
	In Service Software, under "Vent Interface Bd Power Diagnostics" (Section 8a.1.6), view the "Vent Int Bd 10VA Voltage" from Board Supplies:  If "Vent Int Bd 10VA Voltage" reads "OK, and 3.2 Vdc reads "Fail", replace the VIB.  If "Vent Int Bd 10VA Voltage" reads "Fail" and the "Vent Int Bd 10VA Voltage" from the Anes Cntrl Bd reads "OK", Check cabling between ACB and VIB.						
VENT AIRWAY OVERPRESS SIGNAL	Inspiration stopped	High airway overpressure signal set.	Medium	AC, Vent	Mechanical Ventilation On		
	No Service Action. Reboot system. If problem continues, check Airway Pressure signal in Service Mode.						
VENT AIRWAY OVERPRESS SIGNAL FAIL	Vent Fail. Monitoring Only	Ventilator SIB indicates the High Airway overpressure signal was set and Paw < 90 cmH <sub>2</sub> O and Pmanifold <80 cm H <sub>2</sub> O.	Medium	AC, Vent	Mechanical Ventilation On		
	No Service Action. Reboot system. If problem continues, check Airway Pressure signal in Service Mode.						
VENT FLOW VALVE FAIL CURRENT	Vent Fail. Monitoring Only	Incorrect current feedback for 7 consecutive readings.	Medium	AC, Vent			
	Reboot System.  If problem continues, in Service Software, under "Vent Interface Bd Power Diagnostics" (Section 8a.1.6), increase the Flow Valve counts and view the Flow Valve Current mA and Counts.						
VENT SIB 10VA OVER CURRENT		Indicates the current feedback from the VIB was incorrect.		ACB			
	Disconnect VIB power harness and restart the system.  If the error does not reappear in the log,  replace the VIB and retest.  If the error persists,  replace the ACB.						

Error Log Entry	Alarm Text	Condition (Basic info)	Priority	Source	Enabling Criteria	
	Action/Troubleshooting	ξ				
VENT SIB COMMUNICATION FAILURE		Indicates a loss of communication between the Anesthesia Controller Board (ACB) and the Ventilator Interface Board (VIB) after communication has been established.		ACB		
	Check cabling. Replace Pan Connector Board. Replace VIB. Replace ACB.					
VENT SUSTAINED PAW SDOWN	Vent Fail. Monitoring Only	Paw > 100 cmH <sub>2</sub> 0 for 10 seconds.	Medium	AC, Vent	In-range Paw data available	
	No Service Action. Reboot system. If problem continues, check Airway Pressure signal in Service Mode.					
VENT VALVE 10VA OVER CURRENT		Indicates the current feedback from the Insp Flow Valve was incorrect for seven consecutive readings.		ACB		
	In the Service Software / Vent Flow & Pressure Diagnosis, increase the Flow Valve counts and view the Flow Valve Current mA and Counts.					
VENT VALVE POWER FAIL	Vent Fail. Monitoring Only	Nominal 12.5V <11.3 V or >13.13Vdc	Medium	AC, Vent		
	In Service Software, under "Vent Interface Bd Power Diagnostics" (Section 8a.1.6), view the "Vent Int Bd 10VA Voltage" from Board Supplies:  If "Vent Int Bd 10VA Voltage" reads "OK, and the Vent Valve 10VA Volts reads "Fail", disconnect the VIB to Pan connector harness. If the Vent Valve 10VA Volts continues to read "Fail", replace the VIB.  If "Vent Int Bd 10VA Voltage" reads "Fail" and the "Vent Int Bd 10VA Voltage" from the Anes Cntrl Bd reads "OK", Check cabling between ACB and VIB.					
VLV BAL GAS CH LEAK TESTS NOT DONE		Indicates the Bal Gas Leak Test skipped. Can be caused by no Bal Gas connected at power-up.		Mixer		
	No Service Action Requi	red.				

7-70 09/07 1009-0357-000

# 7.11 Steps and Messages displayed during the System Checkout — (for System software 3.X or greater)

**Stage 1:** Step 1: "Bag/Vent Switch" - Verify the Bag/Vent Switch is set correctly.

- If the switch is set to Ventilator Mode continue with next step
- If the switch is set to Manual Mode, fail with "Wrong circuit selected".

Step 2: "02 Pressure" - Is O2 available and working.

- If O<sub>2</sub> Supply is adequate and mixer passes a mixer flow test of 250 ml/min of O<sub>2</sub> continue with the next step.
- If not fail with "Low 02 supply pressure" or with mixer failure.

Step 3: "Ventilator Drive Pressure" - Make sure ventilator has drive gas pressure:

- If drive gas pressure (as measured by the Manifold Pressure Transducer) continue with the next step.
- If no drive gas pressure, fail with "Ventilator has no drive gas".
- If all of the steps above pass, the ventilator will be commanded to flow 12 L/min (Software less than 3.20) or 18 L/min (Software equal to or greater than 3.20). Before selecting "Continue" on the next menu, "Make sure the bellows are fully collapsed" before you "Occlude the Patient Y".

**Stage 2:** Step 1: "Verify the Bellows Empty" - Check to make sure bellows is collapsed.

- If Airway pressure increases to or above 30 cm H<sub>2</sub>0 in 5 seconds fail the test with, "Can not empty bellows".
- If not continue with next step.

Step 2: "Circuit Leak Test" - Attempt to find the leak of the ventilator mode system.

### If the Leak < 250 ml is set to No:

- Flow 1 I/min 02 until pressure increases to 20 cm H<sub>2</sub>0.
- If pressure does not increase to 20 cm H<sub>2</sub>0 within in 15 seconds, fail with "Cannot pressurize circuit". If it does reach 20 cm H<sub>2</sub>0, change flow to 250 ml/min.
- If a flow of 250 ml/min reaches 30 cm H<sub>2</sub>0 or greater, display "Circuit leak is less than 250 ml/min" and continue to next step. If it does not, increase flow up to 750 ml/min.
- If a flow of 750 ml/min that reaches 30 cm H<sub>2</sub>0, continue to next step but indicate the "Ventilator circuit leak is between 250 ml/min and 750 ml/min".
- If the circuit pressure does not increase to 30 cm H20 with 750 ml/min flow, fail with "Ventilator circuit leak is greater than 750 ml/min".

#### If the **Leak < 250 ml** is set to **Yes**:

- Flow 1 I/min 02 until pressure increases to 20 cm H<sub>2</sub>0.
- If pressure does not increase to 20 cm H<sub>2</sub>0 within in 15 seconds, fail with "Cannot pressurize circuit". If it does reach 20 cm H<sub>2</sub>0, change flow to 100 ml/min.

1009-0357-000 09/07 7-71

- If a flow of 100 ml/min reaches 30 cm H<sub>2</sub>0 or greater, display "Circuit leak is less than 100 ml/min" and continue to next step. If it does not, gradually increase flow until pressure reaches 30 cm H<sub>2</sub>0.
- If a flow is found that reaches 30 cm H20 and that flow is less than 750 ml/ min continue to next step but indicate that, "Ventilator circuit leak is ## ml/ min".
- If the circuit pressure does not increase to 30 cm H20 with 750 ml/min flow, fail with "Ventilator circuit leak is greater than 750 ml/min".
- Step 3: "Mechanical Ventilation" Tests the Mechanical Ventilation by delivering small Pressure Controlled breaths and look for alarms:
  - If alarm condition is detected it will be stated in final menu.
  - Continue to next step.
- Step 4: "Circuit Compliance" Tests the circuit compliance by delivering small (15 cm H<sub>2</sub>0) Pressure Controlled breaths, measures the circuit volume via the expiratory flow sensor, and looks for alarms:
  - If the measured volume is less than 15 ml (either due to small patient circuit
    or flow sensor issues) a "Check Flow Sensor" alarm may be generated that
    prohibit the calculation of circuit compliance, state that "Can not measure
    circuit compliance".
  - If alarms did not occur then calculate compliance and state, "Circuit Compliance YYY ml/ cm H20"
- Step 5: " $O_2$  Flow" Run the mixer tests on  $O_2$  channel (Check gas supply and run a 3L and 10L flow delivery test and leak test:
  - If O2 pressure is low, fail with "Low O2 Supply pressure".
  - If mixer does not fail the 3L, 10L and the leak test, continue with the next step.
  - If mixer fails the 3L, 10L or the leak test, fail with the mixer failure.
- Step 6: "AIR Flow" Run the mixer tests on AIR channel (Check gas supply and run a 3L and 10L flow delivery test and leak test:
  - If 02 is drive gas and air supply is low, continue with next step and indicate, "Could not test air".
  - If Air is the drive gas and air supply is low, continue with next step and indicate, "Ventilator has not drive gas".
  - If mixer does not fail the 3L,10L and the leak test, continue with the next step.
  - If mixer fails the 3L,10L or the leak test, fail with the mixer failure.
- Step 7: " $N_2O$  Flow" Run the mixer tests on  $N_2O$  channel (Check gas supply and run a 3L and 10L flow delivery test and leak test:
  - If N<sub>2</sub>O is disabled, continue with next step.
  - If N<sub>2</sub>O supply is low, continue with next step and indicate, "Could not test N<sub>2</sub>O".
  - If mixer does not fail the 3L,10L and the leak test, continue with the next step.
  - If mixer fails the 3L,10L or the leak test, fail with the mixer failure.

7-72 09/07 1009-0357-000

Step 8: "Battery and Electrical" – Are the AC/Mains connected and the battery charged?

- If AC/Mains failed indicate, "Power cord disconnected. Using battery".
- If Battery failed indicate, "Battery failure".
- If Battery charging with 20 30 minutes available indicate, "Battery still charging."
- If Battery charging with 10 20 minutes available indicate, "Battery still charging".
- If battery charging with 0-10 minutes available indicate, "Battery still charging".
- If battery fully charged, pass step.

# 7.12 Steps and Messages displayed during the System Checkout — (for System software 2.X)

# 7.12.1 Steps for the Quick Check

Step 1: Check to make sure bag/vent switch is set correctly.

- If switch is set to manual mode continue with the next step.
- If set to ventilator mode fail with, "Wrong circuit selected".

Step2: Is 02 available and working?

- If O<sub>2</sub> Supply is Not low and mixer passes a mixer flow test of 250 ml/min of O<sub>2</sub> continue with the next step.
- If not fail with "Low 0<sub>2</sub> supply pressure" or with mixer failure.

Step3: Is circuit (airway) pressure too high?

- Check the circuit (airway) pressure is less than 30 cm H<sub>2</sub>0; continue with the next step.
- If airway pressure greater than 30 cm H<sub>2</sub>0 fail with "Circuit pressure too high".
- Step4: Find manual circuit leak:
- Increase flow to find the manual circuit leak.
- If pressure does not increase to 30 cm H<sub>2</sub>0 fail with "Cannot pressurize circuit".
- If flow required to maintain pressure at 30 cm H<sub>2</sub>0 is greater than 250 ml/min state "Manual circuit leak is XXX ml/min at 30 cm H<sub>2</sub>0".
- If flow required to maintain pressure at 30 cm H<sub>2</sub>0 is less than or equal to 250 cm H<sub>2</sub>0 continue with the next step.

Step 5: Run mixer tests on the O<sub>2</sub> channel (check supply and run a 3L, 10L and leak test):

- If O<sub>2</sub> pressure is low, fail with "Low O<sub>2</sub> Supply pressure".
- If mixer does not fail the 3L,10L and the leak test, continue with the next step.
- If mixer fails the 3L,10L or the leak test, fail with the mixer failure.

1009-0357-000 09/07 7-73

- Step 6: Run mixer tests on the Air channel (check supply and run a 3L, 10L and leak test):
  - If O<sub>2</sub> is drive gas and air supply is low, continue with next step and indicate, "Could not test air".
  - If Air is the drive gas and air supply is low, continue with next step and indicate, "Ventilator has not drive gas".
  - If mixer does not fail the 3L,10L and the leak test, continue with the next step. If mixer fails the 3L,10L or the leak test, fail with the mixer failure.
- Step 7: Run mixer tests on the N<sub>2</sub>O channel (check supply and run a 3L, 10L and leak test):
  - If N<sub>2</sub>O is disabled, continue with next step.
  - If N<sub>2</sub>O supply is low, continue with next step and indicate, "Could not test N<sub>2</sub>O".
  - If mixer does not fail the 3L,10L and the leak test, continue with the next step.
  - If mixer fails the 3L,10L or the leak test, fail with the mixer failure.

Step 8: Are the AC/Mains connected and the battery charged?

- If AC/Mains failed indicate, "Power cord disconnected. Using battery".
- If Battery failed indicate, "Battery failure".
- If Battery charging with 20 30 minutes available indicate, "Battery still charging."
- If Battery charging with 10 20 minutes available indicate, "Battery still charging".
- If battery charging with 0-10 minutes available indicate, "Battery still charging".
- If battery fully charged, pass step.

7-74 09/07 1009-0357-000

# 7.12.2 Steps for the Vent Check

## Stage 1:

- Step 1: Check to make sure bag/vent switch is set correctly.
  - If switch is set to ventilator continue with the next step.
  - If set to manual mode fail with "Wrong Circuit Selected".

### Step 2: Is O<sub>2</sub> Pressure Not Low?

- If O<sub>2</sub> Supply is Not low and mixer passes it's 3 L test continue to next step.
- If not fail with "Low O<sub>2</sub> supply pressure" or with mixer failure.

### Step 3: Make sure ventilator has drive gas pressure:

- If drive gas pressure continue with the next step.
- If no drive gas pressure, fail with "Ventilator has no drive gas".
- If all of the steps above pass, the ventilator will be commanded to flow 12 L/min. Depending on the software you have a 10 second delay may be in place before the next menu is displayed. Either way before selecting "Continue" on the next menu, "Make sure the bellows are fully collapsed" before you "Occlude the Patient Y".

## Stage 2:

- Step 1: Check to make sure bellows are collapsed.
  - If Airway pressure increases to or above 30 cm H<sub>2</sub>0 in 5 seconds fail the test with, "Can not empty bellows".
  - If not continue with next step.

Step 2: Attempt to find the leak of the ventilator mode system.

- Flow 1 I/min 0<sub>2</sub> until pressure increases to 20 cm H<sub>2</sub>0.
- If pressure does not increase to 20 cm H<sub>2</sub>0 within in 15 seconds, fail with "Cannot pressurize circuit". If it does reach 20 cm H<sub>2</sub>0, change flow to 250 ml/min.
- If a flow of 250 ml/min reaches 30 cm H<sub>2</sub>0 or greater continue to next step.
   If it does not, increase flow up to 750 ml/min.
- If a flow is found that reaches 30 cm H<sub>2</sub>0 and that flow is less than 750 ml/min continue to next step but indicate that, "Ventilator circuit leak is ## ml/min".
- If a flow greater than 750 ml/min is required, fail with, "Ventilator circuit leak is greater than 750 ml/min".

Step 3: System delivers small breaths and looks for alarms using default alarm limits.

- If alarm condition is detected it will be stated in final menu.
- Continue to next step.

#### Step 4: Calculate circuit compliance.

- If alarms occurred that prohibit the calculation of circuit compliance, state that "Can not measure circuit compliance".
- If alarms did not occur then calculate compliance and state, "Circuit Compliance YYY ml/ cm H<sub>2</sub>O"

1009-0357-000 09/07 7-75

Notes

7-76 09/07 1009-0357-000

# 8a Service Diagnostics and Software Download (DU)

In this section	8a.1 Avance Service Application	8a-2
	8a.1.1 Main Menu and System Information	8a-2
	8a.1.2 Power Diagnostics	8a-3
	8a.1.3 Power Controller Power Diagnostics	8a-4
	8a.1.4 Anesthesia Control Board Power Diagnostics	8a-6
	8a.1.5 Electronic Mixer Power Diagnostics	8a-8
	8a.1.6 Ventilator Interface Board Power Diagnostics	8a-9
	8a.1.7 Display Unit Power Diagnostics	8a-10
	8a.2 Gas Diagnostics	8a-11
	8a.2.1 Gas Supplies	8a-12
	8a.2.2 Mixer Output	8a-13
	8a.2.3 Mixer Tests and Pressure	8a-14
	8a.2.4 Mixer Temperature	8a-15
	8a.2.5 Setting Gas Flow	8a-16
	8a.2.6 Breathing System Leak Test	8a-17
	8a.3 Ventilation Diagnostics	8a-18
	8a.3.1 Status	8a-19
	8a.3.2 Vent Flow and Pressure	8a-20
	8a.4 Display Diagnostics	8a-21
	8a.5 Special Functions	8a-22
	8a.5.1 Mixer Service Menu	8a-23
	8a.5.2 View Revision Log	8a-24
	8a.5.3 View PC Card Install Log	8a-24
	8a.6 Software Download	8a-25

# 8a.1 Avance Service Application

This section documents the Avance Service Application that loads from a PCMCIA card and is used to download software or to run various diagnostic functions.

To run the application, first set the system switch to Standby and set the AC Inlet power switch to Off. Insert the card carrier (with card facing to the rear) into the rear PCMCIA interface slot of the display unit (behind left side door), then set the AC Inlet power switch and the system switch to On. The service application will load and display the Main Menu along with the System Information page.

# 8a.1.1 Main Menu and System Information

The Main Menu appears on the left-hand side of the screen and includes the following selections as shown in the table below:

Main Menu	Remarks		
Power Diagnostics	Access to the Power Supply Diagnostics functions		
Gas Diagnostics	Access to the Gas Delivery Diagnostics functions		
Vent Diagnostics	Access to the Ventilator Diagnostics functions		
Display Diagnostics	Access to the Front Panel Controls		
Special Functions	Access to logs from the Display Unit		
Software Download	Access to the Software Download function.		

#### Note

You can not return to the Diagnostic section of the service application after entering the software download section. You must reboot the system to exit Software Download.

The System Information page appears on the right-hand side of the screen and displays the following system information as shown in the table below:

### System Information

Subsystem	HW Rev	Serial #	SW Ver #	Boot	
Front Panel Cntl			XX.XX		
Power Controller	XXX/A/XXX	ABCXXXXX	XX.XX	XX.XX	
Electronic Mixer	XXX/A/XXX	ABCXXXXX	XX.XX	XX.XX	
Vent Intface Bd	XXX/A/XXX	ABCXXXXX	XX.XX	XX.XX	
Anes Control Bd	XXX/A/XXX	ABCXXXXX	XX.XX	XX.XX	
Dsply Unit BIOS	XXX/A/XXX	ABCXXXXX	XX.XX		
Dsply Unit App	XXX/A/XXX	ABCXXXXX	XX.XX		
Machine Serial Number: ABCDXXXXX					
PC Card ID: XXXX	_XX				

(DU) 8a-2 09/07 1009-0357-000

# 8a.1.2 Power Diagnostics

The service application provides power supply diagnostics for the various circuit boards in the Avance anesthesia machine.

Selecting **Power Diagnostics** on the Main Menu brings up the following menu selections in the left-hand frame and the instructions in the right-hand frame:

Main Menu
Power Diagnostics
Gas Diagnostics
Vent Diagnostics
Display Diagnostics
Special Functions
Software Download

Power Diagnostics	Power Controller Power Diagnostics			
Power Control				
Anes Control Board	Select a menu item			
Mixer	to see the power status and measured voltages.			
Vent Interface Bd	and modeline voltages.			
Display Unit	To troubleshoot a power problem, start with the power controller and work forward.			
-> Main Menu	Problem voltages are in red.			

# 8a.1.3 Power Controller Power Diagnostics

There are two pages of diagnostics for the Power Controller.

Selecting **Power Control** brings up the first page of the Power Controller Diagnostics.

Power Diagnostics

Power Control

Anes Control Board

Mixer

Vent Interface Bd

Display Unit

->Main Menu

(Page 1 of 2)
Power Control Power Diagnostics

Label	Value Format	Units	Normal range
AC Status	OK, Fail		
12Vdc Supply	XX.XX	Vdc	11.82 to 12.18
3.3Vdc Supply	X.XXX	Vdc	3.201 to 3.399
1.5 Vdc Supply	X.XXX	Vdc	1.45 to 1.55
Battery Connected	Yes, No		
Battery Status	Fail, Bulk Chg, Over Chg, Float Chg, Trickle Chg, Discharge		
Battery Current	X.XXX	Α	
Calc Battery Time	XX	Min	0 to 30
Battery 1 Volts	XX.X < 6.0 FAIL (red) <10 T Chg (yellow) 10-15.5 (Green)	Vdc	10.0 to 15.5
Battery 2 Volts	XX.X < 6.0 FAIL (red) <10 T Chg (yellow) 10-15.5 (Green)	Vdc	10.0 to 15.5
Date battery Tested (*)	XX-ABC-XXXX		
Last Full Discharge Time (*)	XX	Min	

Turn the ComWheel to select the second page.

Push the ComWheel to return focus to the Power Diagnostics selection menu.

(\*) These entries are only listed after a full discharge test of the batteries is first performed (Section 6.10, "Battery capacity test").

(DU) 8a-4 09/07 1009-0357-000

Power Diagnostics
Power Control
Anes Control Board
Mixer
Vent Interface Bd
Display Unit
->Main Menu

# (Page 2 of 2) Power Control Power Diagnostics

Label	Value Format	Units	Normal range
Board Temperature	<65C OK (green) >65C <75C Warn (Yellow) >75C Fail (red)	Deg C	Upper 64
Fan Speed	Slow, Fast		
Fan 1 Voltage	XX.XX	Vdc	11.52 to 12.48
Fan 1 Status	Low, High, Fail, OK		

Turn the ComWheel to return to the first page.

Push the ComWheel to return focus to the Power Diagnostics selection menu.

# 8a.1.4 Anesthesia Control Board Power Diagnostics

There are two pages of diagnostics for the Anesthesia Control board.

Selecting **Anes Control Board** brings up the first page of the Anesthesia Control Board Power Diagnostics.

Power Diagnostics

Power Control

Anes Control Board

Mixer

Vent Interface Bd

Display Unit

->Main Menu

(Page 1 of 2)
Anes Control Board Power Diagnostics

Label	Value Format	Units	Normal range
12.5Vdc from Pwr Cntrl	XX.XX	Vdc	11.90 to 12.90
ADC Reference	X.XXX	Vdc	4.018 to 4.176
Gas Select 10VA Volts	OK, Fail		
Gas Select 10VA Amps	OK, Fail		
P Xducer 10VA Amps	OK, Fail		
Vent Int Bd 10VA Volts	OK, Fail		
Vent Int Bd 10VA Amps	OK, Fail		
Vent Valves 10VA Volts	OK, Fail		
Vent Valves 10VA Amps	OK, Fail		
Acces 1 10VA Volts	OK, Fail		
Acces 1 10VA Amps	OK, Fail		

Turn the ComWheel to select the second page.

Push the ComWheel to return focus to the Power Diagnostics selection menu.

(DU) 8a-6 09/07 1009-0357-000

Power Diagnostics
Power Control
Anes Control Board
Mixer
Vent Interface Bd
Display Unit
->Main Menu

(Page 2 of 2)
Anes Control Board Power Diagnostics

Label	Value Format
Gas Unit 10VA Volts	OK, Fail
Gas Unit 10VA Amps	OK, Fail
Mixer 10VA Volts	OK, Fail
Mixer 10VA Amps	OK, Fail
Alt O2 10VA Volts	OK, Fail
Alt 02 10VA Amps	OK, Fail
Periph1 10VA Volts	OK, Fail
Periph1 10VA Amps	OK, Fail
Periph2 10VA Volts	OK, Fail
Periph2 10VA Amps	OK, Fail

Turn the ComWheel to return to the first page.

Push the ComWheel to return focus to the Power Diagnostics selection menu.

# 8a.1.5 Electronic Mixer Power Diagnostics

Selecting *Mixer* brings up the Electronic Mixer Power Diagnostics page.

# Power Diagnostics Power Control Anes Control Board Mixer Vent Interface Bd Display Unit ->Main Menu

# **Mixer Power Diagnostics**

Label	Value Format	Units	Normal range
From Anes Cntrl Bd Mixer 10VA Volts	OK, Fail		
<b>Board Supplies</b>			
12.5 V	XX.XX	Vdc	11.80 to 13.00
5.5V	X.XX	Vdc	5.39 to 5.61
3.3V CPU	X.XX	Vdc	3.22 to 3.38
2.5V ADC Ref	X.XX	Vdc	2.47 to 2.53

(DU) 8a-8 09/07 1009-0357-000

# 8a.1.6 Ventilator Interface Board Power Diagnostics

Selecting **Vent Interface Board** brings up the Ventilator Interface Board Power Diagnostics page.

# Power Diagnostics Power Control Anes Control Board Mixer Vent Interface Bd Display Unit ->Main Menu

# **Vent Interface Bd Power Diagnostics**

Label	Value Format	Units	Normal range
From Anes Cntrl Bd Vent Int Bd 10VA Volts Vent Valves 10VA Volts	OK, Fail OK, Fail		
Board Supplies			
Vent Int Bd 10VA Volts	XX.XX	Vdc	11.30 to 13.13
Vent Valves 10VA Volts	XX.XX	Vdc	11.30 to 13.13
3.2Vdc (12bit Vref)	X.XXX XXXX	Vdc Counts	3.179 to 3.221
1.22Vdc (10bit Vref)	X.XXX XXXX	Vdc Counts	1.074 to 1.367
+6.0Vdc	X.XX	Vdc	5.51 to 6.50
-6.0Vdc	-X.XX	Vdc	-6.72 to -5.28

Display Unit

->Main Menu

# 8a.1.7 Display Unit Power Diagnostics

Selecting **Display Unit** brings up the Display Unit Power Diagnostics page.

# Power Diagnostics Power Control Anes Control Board Mixer Vent Interface Bd

# **Display Unit Power Diagnostics**

Label	Value Format	Units	Normal range
5.0Vdc (PCMCIA)	X.XX	Vdc	4.50 to 5.50
3.3Vdc (PCMCIA)	X.XX	Vdc	2.97 to 3.63
5.0Vdc (Fan)	X.XX	Vdc	4.50 to 5.50
5.0Vdc (USB)	X.XX	Vdc	4.50 to 5.50
8.0Vdc (DIS)	X.XX	Vdc	7.20 to 8.80
11Vdc (LCD)	XX.XX	Vdc	10.35 to 13.62

(DU) 8a-10 09/07 1009-0357-000

# 8a.2 Gas Diagnostics

Selecting **Gas Diagnostics** on the Main Menu brings up the following menu selections in the left-hand frame and the instructions in the right-hand frame.

Main Menu
Power Diagnostics
Gas Diagnostics
Vent Diagnostics
Display Diagnostics
Software Download

Gas Diagnostics	Gas Diagnostics
Gas Supplies	
Mixer Output	To troubleshoot a problem with gas supplies or gas flows, start with the gas supplies menu and
Mixer Tests and Pres	work forward.
Mixer Temperatures	To turn on gas flows, select O2, N2O, or Air in the
02 Flow: OFF	left menu.
N20 Flow: OFF	Available flows are:
Air Flow OFF	100% 02: 0.2, 0.5, 5, 10 l/min
Breathing Sys Leak	100% N20: 0.2, 0.5, 5, 10 l/min Air: 0.2, 0.5, 5, 10 l/min
-> Main Menu	
	Turning ON a different gas automatically turns OFF previous flows. All gas flows stop if you go back to the main menu.

Selecting any of the first four items in the left-hand frame, brings up the corresponding diagnostic page in the right-hand frame.

Press the ComWheel to return focus to the right-hand frame to make another selection.

## Start gas flow

Selecting one of the three gas flows brings up a page on which you can set one of the four flow values for the selected gas:

• 0.2 l/min; 0.5 l/min, 5.0 l/min, 10.0 l/min.

Press the ComWheel to start the selected flow.

Press the ComWheel again to return to the flow page to observe the Airway Pressure reading. Or select one of the four diagnostics pages to view the displayed conditions.

Selecting Breathing Sys Leak, opens the GIV valve and sets a flow through the inspiratory flow valve. While on the Breathing Sys Leak page, you can set O2 flow to the breathing system and observe the Airway Pressure reading.

### Stop gas flow

Gas will continue to flow at the set rate until you:

- set a different flow for the same gas,
- set a flow for a different gas (automatically turns the previous gas flow OFF),
- set the current gas flow to OFF (all flows OFF),
- exit Gas Diagnostics to the Main Menu.

# 8a.2.1 Gas Supplies

Each gas supply shows the derived pressure in kPa and psi along with the raw voltage from the pressure transducer.

If a supply module for a gas is not installed, or if the transducer is disconnected, the supply shows 0.00 Vdc.

If a supply module is installed but no supply is connected, the supply shows approximately 0.50 Vdc.

The remaining items show the state of the Electronic Mixer selector valves.

# Gas Supplies Gas Diagnostics

Label	psi	kPa	Vdc
O2 Cylinder 1	XXXX	XXXXX	XX.XXX
02 Cylinder 2	XXXX	XXXXX	XX.XXX
Air Cylinder	XXXX	XXXXX	XX.XXX
N20 Cylinder	XXXX	XXXXX	XX.XXX
O2 Pipeline	XXXX	XXXXX	XX.XXX
Air Pipeline	XXXX	XXXXX	XX.XXX
N20 Pipeline	XXXX	XXXXX	XX.XXX
O2 Select Valve	Open	Flow comma	nd ON
N20 Select Valve	Closed	Flow comma	nd OFF
Air Select Valve	Closed	Flow comma	nd OFF
Alt 02 Valve	Closed		

Gas Diagnostics

Gas Supplies

Mixer Output

Mixer Tests and Pres

Mixer Temperatures

02 Flow: 0.5

N20 Flow: 0FF

Air Flow: OFF

Breathing Sys Leak

-> Main Menu

Note: 0.0 Vdc = not installed or transducer disconnected.

(DU) 8a-12 09/07 1009-0357-000

# 8a.2.2 Mixer Output

Selecting *Mixer Output* brings up the Mixer Output Gas Diagnostics page.

This data comes from the Anesthesia Control board. The Flow Verify signals are rough calculations of the mixer flow based on pressure drop and temperature. The ADC reference voltage is used to convert flow signals.

# **Mixer Output Gas Diagnostics**

Gas Diagnostics
Gas Supplies
Mixer Output
Mixer Tests and Pres
Mixer Temperatures
02 Flow: 0.5
N20 Flow: OFF
Air Flow: OFF
Breathing Sys Leak
-> Main Menu
-> main wenu

Label	Value	Units
02 Flow	XX.XX	I/min
Balance Flow	XX.XX	l/min
O2 Flow Verify	XX.XX	l/min
02 Flow Signal	X.XXX	Vdc
O2 Prop Valve Drive	XXXX	mA
Balance Gas ID	None, Air, N20	
Balance Flow Verify	XX.XX	I/min
Balance Flow Signal	X.XXX	Vdc
Balance Prop Valve Drive	XXXX	mA
02 Select Valve	Open	
Air Select Valve	Closed	
N20 Select Valve	Closed	
ADC Ref Voltage	X.XXX	Vdc

# 8a.2.3 Mixer Tests and Pressure

Selecting *Mixer Tests and Pres* brings up the Mixer Tests and Pres Gas Delivery Diagnostics page.

# **Mixer Tests and Pres Gas Diagnostics**

# Gas Diagnostics Gas Supplies Mixer Output Mixer Tests and Pres Mixer Temperatures O2 Flow: 0.5 N20 Flow: OFF Air Flow: OFF Breathing Sys Leak -> Main Menu

Last Power-Up Tests			
O2 Proportional Valve Leak	Not done. Sel Pass Fail. Selector	supply pressure ector valve inco valve leaks onal valve leaks	orrect state
Alt O2 Valve Leak	Pass; Fail		
Balance Gas Prop Valve Leak	Not done. Sel Pass Fail. Selector	supply pressur ector valve inco valve leaks onal valve leaks	orrect state
Balance Gas CheckValve Leak	Pass; Fail		
O2 Flow Test			
Balance Flow Test			
Balance Gas ID	None; Air; N2	0	
Pressure Data	Value	Units	kPa
02 Pressure (P1)	XX.XX	psi	XXX.XX
O2 Pres Cal (P1)	X.XXX	Vdc	
Balance Pressure (P2)	XX.XX	psi	XXX.XX
Balance Pres Cal (P2)	X.XXX	Vdc	
Mixer Output Pres (P3)	XX.XX	psi	XXX.XX
Mixer Output Pres Cal (P3)	X.XXX	Vdc	
ADC Ref Voltage	X.XXX	Vdc	

(DU) 8a-14 09/07 1009-0357-000

# 8a.2.4 Mixer Temperature

Selecting *Mixer Temperatures* brings up the Mixer Temperatures Gas Diagnostics page.

# **Mixer Temperatures Gas Diagnostics**

Gas Diagnostics
Gas Supplies
Mixer Output
Mixer Tests and Pres
Mixer Temperatures
02 Flow: 0.5
N20 Flow: OFF
Air Flow: OFF
Breathing Sys Leak
-> Main Menu

Sensor Data	Value	Units
O2 Temp (T1)	XX.X	Deg C
O2 Temp Volts (T1)	X.XXX	Vdc
Balance Temp (T2)	XX.X	Deg C
Balance Temp Volts (T2)	X.XXX	Vdc
ADC Ref Voltage	X.XXX	Vdc

# 8a.2.5 Setting Gas Flow

Selecting **02** Flow: **0FF** brings up the O2 Flow Setting page on which you can select one of four flows (or OFF if flow previously set):

- 0.2 l/min
- 0.5 l/min
- 5.0 l/min
- 10.0 l/min

Selecting "N2O Flow: OFF" or "Air Flow: OFF" brings up a comparable page for setting the above flows for the selected gas. Selecting a flow for a new gas turns the previous gas flow to OFF.

# 02 Flow Setting

	Value	Units
O2 Flow	X.X	I/min
To choose a gas flow, tu	ırn ComWheel.	
Gas flows do not chang	e until you push Con	nWheel.
Airway pressure:	XXX	cmH2O
	Push ComWheel to	Exit and start O2 flow ->
·	·	·

After setting a gas flow, push the ComWheel again to return to the Flow Setting page to observe the Airway Pressure reading, or select one of the four diagnostics pages to view the displayed conditions.

All gas flow stops when you exit Gas Diagnostics to the Main Menu.

# Gas Diagnostics Gas Supplies Mixer Output Mixer Tests and Pres

WIINEL LESIS ALIA FLES

Mixer Temperatures

# 02 Flow: OFF

N20 Flow: OFF Air Flow: OFF

Breathing Sys Leak

-> Main Menu

(DU) 8a-16 09/07 1009-0357-000

# 8a.2.6 Breathing System Leak Test

Selecting Breathing Sys Leak brings up the Breathing Sys Leak page.

Entering the Breathing Sys Leak page, opens the Gas Inlet Valve and sets up a constant flow through the Inspiratory Flow Valve.

On the Breathing Sys Leak page you can set an O2 flow through the Gas Mixer and view the Airway Pressure.

# **Breathing Sys Leak**

Breathing Sys Leak Test: Flow valve and GIV valve are open.

Gas Diagnostics

Gas Supplies

Mixer Output

Mixer Tests and Pres

Mixer Temperatures

O2 Flow: OFF

N2O Flow: OFF

Air Flow: OFF

Breathing Sys Leak

-> Main Menu

		с с с с.
	Value	Units
O2 Flow	X.X	l/min
To choose a gas flo Gas flows do not cl		
Refer to Technical F	Reference Manual	for this procedure.
Airway pressure:	XXX	cmH20

Second ComWheel Push Ends Test ->

All gas flow stops when you exit Gas Diagnostics to the Main Menu.

Note

This procedure is detailed in Test 6, "Testing the bellows module and the Bag/Vent switch" on page 7-21.

# **8a.3 Ventilation Diagnostics**

The service application provides several pages for ventilation diagnostics.

Selecting **Vent Diagnostics** on the Main Menu brings up the following menu selections in the left-hand frame and the instructions in the right-hand frame:

Main Menu
Power Diagnostics
Gas Diagnostics
Vent Diagnostics
Display Diagnostics
Special Functions
Software Download

Vent Diagnostics	Vent Diagnostics
Status	
Vent Flow and Pres	Status data shows the position of breathing circuit switches.
Gas Inlet VIv ON/OFF	_
Flow Valve Control	Vent Flow and Pres shows vent sensor readings.
Toggle Circuit	Other commands in the menu column let you control ventilator valves and send gas to either the circle or the non-circle circuit.
-> Main Menu	

Vent Diagnostics	Action when selected
Status	Selecting <b>Status</b> displays the Ventilation Status page.
Vent Flow and Pres	Selecting <b>Vent Flow and Pres</b> displays the Ventilation Flow and Pressure page.
Gas Inlet VIv ON/OFF	Selecting <b>Gas Inlet VIv ON/OFF</b> toggles the Gas Inlet Valve to the Open or Closed position and sets the right-hand screen to the Vent Flow and Pressure Diagnostics page.
Flow Valve Control	Selecting <b>Flow Valve Control</b> sets the right-hand screen to the Vent Flow and Pressure Diagnostics page and allows the user to increase or decrease the DAC count to the flow valve using the ComWheel.
Toggle Circuit	Selecting <b>Toggle Circuit</b> toggles between commanding Circle or Non-Circle for the SCGO valve and sets the right-hand screen to the Vent Status page (does not apply to machines with ACGO).
->Main Menu	Selecting <b>Main Menu</b> closes the Gas Inlet Valve and the Flow Valve and return to the Main Menu.

The actions of the Gas Inlet, the Flow Valve Control, and the Circuit selection are reflected on both the Ventilation Status page and the Ventilation Flow and Pressure page.

The Gas Inlet Valve must be in the ON position in order for gas to flow.

(DU) 8a-18 09/07 1009-0357-000

# **8a.3.1 Status** Selecting *Status* brings up the Ventilation Status page.

# Vent Diagnostics Status Vent Flow and Pres Gas Inlet VIv ON/OFF Flow Valve Control Toggle Circuit -> Main Menu

# Menu Item Value Units Gas Inlet Valve Open or Closed XXXXX Flow Valve Counts **Circuit Command** Circle or Non-Circle Vent Drive Gas Air or 02 **ABS Installed** Yes or No Flush Valve Not Pressed or Pressed 02 Cell Status Connected or None Bag/Vent Switch Bag or Vent Circuit Feedback Circle, Non-Circle, or Fault ACGO/SCGO Configuration ACGO or SCGO

OK or High Pressure

Open or Closed

**Vent Status** 

1009-0357-000 09/07 (DU) 8a-19

**Over Pressure Circuit** 

Gas Inlet Valve Feedback

# 8a.3.2 Vent Flow and Pressure

Selecting **Vent Flow and Pres** brings up the Ventilation Flow and Pressure page.

# Vent Diagnostics Status Vent Flow and Pres Gas Inlet VIv ON/OFF Flow Valve Control Toggle Circuit -> Main Menu

# **Vent Flow and Pressure Diagnostics**

Menu Item	Value	Units	Value	Units
Gas Inlet Valve	Open or Close	ed		
Flow Valve	XXXXX	Counts		
Turn ComWheel to adju	ıst flow valve	(appears whe	n Flow Valve Co	ontrol set ON)
Circuit Command	Circle or Non-	Circle		
 				_
Inspiratory Flow	XXX.X	I/min	XXXX	Counts
Expiratory Flow	XXX.X	I/min	XXXX	Counts
Airway Pressure	XX X	cmH20	XXXX	Counts
Manifold Pressure	XXX	cmH20	XXXX	Counts
O2 Cell	XXX	%	XXXX	Counts
ADC Ref Voltage	X.XXX	Vdc		
Flow Valve Setting	XXX.X	I/min	XXXX	Counts
Flow Valve Feedback	XXXX	mV	XXXX	Counts
Flow Valve Current	XXXX	mA	XXXX	Counts

(DU) 8a-20 09/07 1009-0357-000

# 8a.4 Display Diagnostics

The service application provides several pages for display diagnostics.

Selecting **Display Diagnostics** on the Main Menu brings up the following menu selections in the left-hand frame:

Main Menu
Power Diagnostics
Gas Diagnostics
Vent Diagnostics
Display Diagnostics
Special Functions
Software Download

Display Diagnostics	Display Diagnostics Instructions
Test LEDs	
Test Speaker	Select a menu item
Test Backlight 1	To troubleshoot a display problem,
Test Backlight 2	start with Test LEDs and work forward
Test Keys	
-> Main Menu	

Display Diagnostics	Action when selected
Test LEDs	Selecting <b>Test LEDs</b> causes the red and yellow LEDS next to the <b>Silence Alarms</b> key to flash for 10 seconds.
Test Speaker	Selecting <b>Test Speaker</b> causes the speaker to sound for 5 seconds.
Test Backlight 1	Selecting <b>Test Backlight 1</b> turns backlight 2 off for 10 seconds. "If screen goes black during test, a backlight is out."
Test Backlight 2	Selecting <b>Test Backlight 2</b> turns backlight 1 off for 10 seconds. "If screen goes black during test, a backlight is out."
Test Keys	Selecting <b>Test Keys</b> brings up a representative display of the front panel controls. Pressing a softkey will cause the corresponding key text to be highlighted.
Main Menu	Selecting <b>Main Menu</b> returns to the Main Menu.

# 8a.5 Special Functions

Selecting **Special Functions** on the Main Menu brings up the following menu selections in the left-hand frame:

**Note** 

The View Error Log, View Alarm Log, and View Event Log may not be present on all versions of the Software Downloader Card.

Main Menu
Power Diagnostics
Gas Diagnostics
Vent Diagnostics
Display Diagnostics
Special Functions
Software Download

Special Functions	Special Functions
Mixer Service Menu	
View Error Log	
View Alarm Log	
View Event Log	
View Revision Log	
Compatibility Table	
PC Card Install Log	
View Install Errors	
-> Main Menu	

#### **Mixer Service Menu**

Refer to section 8a.5.1.

### **Error, Alarm, Event Logs**

In a functioning system, the Error, Alarm, and Event Logs are accessible on the system's Service Log menu (refer to *Section 4a.4.2*). If a system comes up in a "Failed State", you can download the logs to the Service Application PC Card and view them on the Special Functions screen.

- 1. With the system still in the "Failed State", insert the Service Application PC Card into the Display Unit and press the "Help" softkey.
- 2. Wait approximately 60 seconds while the logs are downloaded to the card (no apparent activity).
- 3. Restart the system with the Service Application to view the logs.

# Revision and PC Card Install Logs

Whenever a Software Download is completed, the specific software download is recorded in the Revision Log that resides on the system (Display Unit) and in the PC Card Install Log that resides on the PC Card.

# **Compatibility Table**

The Compatibility Table lists the current software components that last downloaded on to the system. In essence, it is the latest listing that appears in the Revision Log, which allows you to view the current log directly without having to scroll to it.

(DU) 8a-22 09/07 1009-0357-000

# 8a.5.1 Mixer Service Menu

Selecting *Mixer Service Menu* brings up the Mixer Service Instructions.

**Note**: This function is not valid on the PCMCIA card and has been moved to the main System Software (Section 4a.5.9) for Software versions 3.X and greater.

Mixer Service Menu	Mixer Service Instructions
Zero Pres Sensors	Select Zero Pres Sensors
Reset Defaults -> Special Functions	to recalibrate the pressure sensor zero-offsets.
Special Fullcuons	Select Reset Defaults
	to restore factory default offsets.
	If you want to select Zero Pres Sensors
	1) Disconnect gas supplies from system.
	2) Close cylinders.
	3) Confirm that all vaporizers are turned "OFF".
	4) Remove the flow sensor cover.
	5) Move the Bag-to-Vent switch to the Bag position.
	6) Press the O2 Flush button for 3 seconds.
	7) WAIT AT LEAST 5 MINUTES WITH NO GAS FLOWING.
	8) Select Zero Pres Sensors on the menu.
	If you want to select Zero Pres Sensors
	DO NOT DISTURB THE SYSTEM WHILE WAITING FOR RESULTS.

# **Zero mixer pressure sensors**

Mixer Service Menu	Zero Pres Sensors
Zero Pres Sensors Reset Defaults -> Special Functions	Pressure sensor zero-offset adjust requested.  Note: this procedure will require at least 2 minutes.

# Reset factory default zero offsets

Mixer Service Menu	Zero Pres Sensors
Zero Pres Sensors	
Reset Defaults	Reset to Factory Defaults requested.
-> Special Functions	

# 8a.5.2 View Revision Log

Selecting **View Revisions Log** brings up the Revision Log for the system. The log includes chronological entries for every Software Download that was completed to the system. Each entry includes two header lines and eight data lines in the following format:

```
# Software configuration after download on (day) (date) (time)
# SvcApp Version (XX.XX), Machine Serial Number (ABCDXXXXX)

Avance ACB, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) AnesControl B

Avance MXR, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Electronic Mix

Avance VNT, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Vent Intface B

Avance FPC, *, (Software Level), (File Name) (#-------) Front Panel CN

Avance PSC, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Power Controll

Avance DUA, *, (Software Level), (File Name) (#-------) Dsply Unit App

Avance DUB, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Dsply Unit BIO

Avance DUF, *, (Software Level), (File Name) (#-------) Dsply Unit Fla
```

#### Note

The Stock Number listed is for the board assembly and may not represent an orderable service item. Refer to the parts lists in Section 10 for service level stock numbers.

The Front Panel Control (FPC), Display Unit Application (DUA), and the Display Unit Flash (DUF) reside, along with the Display Unit BIOS (DUB), on the Display Unit CPU board.

# 8a.5.3 View PC Card Install Log

Selecting **View PC Card Install Log** brings up the PC Card Install Log for the software download card. The log includes chronological entries for every Software Download that was completed with the card. Each entry includes two header lines and eight data lines in the following format:

```
# Software configuration after download on (day) (date) (time)
# SvcApp Version (XX.XX), Machine Serial Number (ABCDXXXXX), Card # XXXXXXX/

Avance ACB, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) AnesControl B

Avance MXR, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Electronic Mix

Avance VNT, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Vent Intface B

Avance FPC, *, (Software Level), (File Name) (#------) Front Panel CN

Avance PSC, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Power Controll

Avance DUA, *, (Software Level), (File Name) (#------) Dsply Unit App

Avance DUB, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Dsply Unit BIO

Avance DUF, *, (Software Level), (File Name) (#-------) Dsply Unit Fla
```

(DU) 8a-24 09/07 1009-0357-000

# 8a.6 Software Download

Main Menu
Power Diagnostics
Gas Diagnostics
Vent Diagnostics
Display Diagnostics
Special Functions
Software Download

Selecting **Software Download** bring up the following information page:

#### ENTERING SOFTWARE DOWNLOAD MODE!

To return to Diagnostics: turn On/Standby switch to Standby, and turn off AC mains switch in rear. Wait 20 seconds, then turn on power with the AC mains switch and the On/Standby switch.

(Press ComWheel to continue with Download.)

### Note

You can not return to the Diagnostic section of the service application after entering the software download section. You must reboot the system to exit Software Download.

Entering software download brings up the Software Download menu.

Software Download	Remarks
Download All	Downloads all software subsystems.
Download New	Downloads only new software versions not found on the system and compatible with installed subsystem hardware.

Since downloading all the subsystem software can take an hour or more, you should normally choose "Download New" to install only the updated subsystem software or software required for newly installed subsystems.

#### **Software Download Status**

Software Download	
Download All	
Download New	

Subsystem	HW Rev	Current SW Ver #	New SW Ver #	Status
Front Panel Cntl	XXX/A/XX	XX.XX	XX.XX	Xxxxxxx
Power Controller	XXX/A/XX	XX.XX	XX.XX	Xxxxxxx
Electronic Mixer	XXX/A/XX	XX.XX	XX.XX	Xxxxxxx
Vent Intface Bd	XXX/A/XX	XX.XX	XX.XX	Xxxxxxx
Anes Control Brd	XXX/A/XX	XX.XX	XX.XX	Xxxxxxx
Dsply Unit BIOS	XXX/A/XX	XX.XX	XX.XX	Xxxxxxx
Dsply Unit Flash				Xxxxxxx
Dsply Unit App	XXX/A/XX	XX.XX	XX.XX	Xxxxxxx
Loading Xxxx Xxxxxxxx Xxxxx:				

### Notes about downloading software

If there is no Front Panel Control software installed in the system (as would be the case when the display units control board is replaced), the Service Application automatically downloads the Front Panel Controls software at startup. During the download the two display unit LEDs will flash and the display speaker will sound an alarm tone to indicate that Software Download is proceeding. The display will be black until the automatic download is complete.

To ensure that all software versions on the system are compatible, the end result of "Download All" or "Download New" will be the same. The software loaded on the machine will exactly match what is on the card. Be sure to have the latest/correct version of software before attempting a download to avoid inadvertent overwrites of newer software with an older version.

If, during the "Download New" process, the compatibility checker detects a newer version of software component on the system, a "Notice" appears on the screen that asks you to confirm the downgrade.

"Download All" will download all compatible software from the card to the system without issuing a notice that newer version of software component may be on the system.

#### **Download process**

The PCMCIA card includes only the latest software for each subsystem along with the diagnostic application.

As each subsystem software segment is being downloaded, the following status messages note the state of each subsystem and the result of the download:

- Busy System is running its application code; not ready for download.
- **Ready** System is in its boot code; ready for download.
- **CRCtest** System is analyzing the download CRC.
- Loading System is accepting download data.
- Done Software download has completed successfully.
- **Fail** Software download did not complete successfully. A "Fail" message will require reloading of the software; or repair of the system may be necessary.
- **Skipped** Software download was bypassed.
- **Linked** System is communicating, but status is not yet known.
- Not Compatible The software version on the PCMCIA card is not compatible with the subsystem.

If the subsystem is communicating but the HW Rev or current SW Rev are not known, the message **Unknown** will appear under the columns for those values. If the HW Rev or current SW Rev are not known, the download function will still be available.

As the software loads, an activity bar at the bottom of the screen shows the download progress for each subsystem.

#### **Download complete**

When all the required subsystem software is download, the following message appears on the screen. You must shut down the system to exit the download function.

DOWNLOAD IS COMPLETE.

Remove PCMCIA card.
Turn ON/STANBY switch to STANDBY.
Turn OFF AC mains switch in rear.

Note

After powering down the system, be sure to wait at least 20 seconds before restarting the system.

(DU) 8a-26 09/07 1009-0357-000

# **8b Software Download and Special Functions (HPDU)**

In this section	8b.1 Overview	8b-2
	8b.1.1 Main Menu and System Information	8h-2
	8b.2 Software Download	
	8b.3 Special Functions	
	8b.3.1 Power Diagnostics	8b-5
	8b.3.2 Power Controller Power Diagnostics	
	8b.3.3 Display Unit Power Diagnostics	
	8b.3.4 Anesthesia Control Board Power Diagnostics	
	8b.3.5 Electronic Mixer Power Diagnostics	8b-11
	8b.3.6 Ventilator Interface Board Power Diagnostics	8b-12
	8b.4 Ventilation Diagnostics	8b-13
	8b.4.1 Vent Status	8b-14
	8b.4.2 Vent Flow and Pressure	8b-15
	8b.5 Gas Diagnostics	8b-16
	8b.5.1 Setting Gas Flow	8b-17
	8b.5.2 Gas Supplies	
	8b.5.3 Mixer Output	
	8b.5.4 Mixer Tests and Pressure	
	8b.5.5 Mixer Temperature	8b-21
	8b.5.6 Breathing System Leak Test	8b-22
	8b.6 Mixer Service Menu	8b-23
	8b.7 Display Diagnostics	8b-25
	8b.7.1 Test Keys and Battery	8b-26
	8b.8 Compatibility Table	8b-27
	8b.8.1 System Download Log	8b-28
	8b.8.2 CF Card Install Log	8b-28
	8b.8.3 View Install Errors	8b-29

# 8b.1 Overview

This section covers the functions of the Compact Flash card used to download system software or to access Special Functions (diagnostics and logs) from the High Performance Display Unit.

To run the application, first set the system switch to Standby and set the AC Inlet power switch to Off. Insert the Compact Flash card into the interface slot of the display unit, then set the AC Inlet power switch and the system switch to On. The application will load and display the Main Menu along with the System Information page.

# 8b.1.1 Main Menu and System Information

The Main Menu appears on the left-hand side of the screen and includes the following selections as shown in the table below:

Main Menu	Remarks
Software Download	Access to the Software Download function.
Special Functions	Access to diagnostics and logs from the Display Unit.

#### **Note**

You cannot return to the Special Functions section after entering the Software Download section. You must reboot the system to exit Software Download.

The System Information page appears on the right-hand side of the screen and displays the following system information as shown in the table below:

# System Serial Number ABCDXXXXXXXXX

Main Menu	
Software Download	
Special Functions	

	tion			
	Currently Installed			
Subsystem	HW Rev	Serial #	SW Ver	BootVer
Front Panel Cntl	n/a	n/a	XX.XX	XXbXX
Power Controller	XXXX/A/XX	ABCXXXXX	XX.XX	XXbXX
Electronic Mixer	XXXX/A /XX	ABCXXXXX	XX.XX	XXbXX
Vent Intface Bd	XXXX/A/XX	ABCXXXXX	XX.XX	XXbXX
Anes Control Bd	XXXX/A /XX	ABCXXXXX	XX.XX	XX.XX
Dsply Unit BIOS	XXXX/XX/XX	ABCXXXXX	XX.XX	n/a
Dsply Unit App	XXX/XX/XX	ABCXXXXX	XX.XX	n/a

(HPDU) 8b-2 09/07 1009-0357-000

# 8b.2 Software Download

Selecting **Software Download** bring up the Software Download menu.

Software Download	Remarks
Download New	Downloads only new software versions not found on the system and compatible with installed subsystem hardware.
Download All	Downloads all software subsystems.

#### Note

You cannot return to the Special Functions section after entering the Software Download section. You must reboot the system to exit Software Download.

Since downloading all the subsystem software can take up to 30 minutes, you should normally choose "Download New" to install only the updated subsystem software or software required for newly installed subsystems.

# System Serial Number ABCDXXXXXXXXX

Software Download
Download New
Download All

Loading Avance Product Software Version XX.XX				
	Currently Installed		New	
Subsystem	HW Rev	SW Ver	SW Ver	Status
Front Panel Cntl	n/a	XX.XX	XX.XX	Xxxxxxx
Power Controller	XXXX/A /XX	XX.XX	XX.XX	Xxxxxx
Electronic Mixer	XXXX/A /XX	XX.XX	XX.XX	Xxxxxx
Vent Intface Bd	XXXX/A /XX	XX.XX	XX.XX	Xxxxxx
Anes Control Bd	XXXX/A /XX	XX.XX	XX.XX	Xxxxxxx
Dsply Unit BIOS	XXXX/XX/XX	XX.XX	XX.XX	Xxxxxxx
Dsply Unit App	XXX/XX/XX	XX.XX	XX.XX	Xxxxxxx
ModBus Controllr	n/a	XX.XX	XX.XX	Xxxxxxx
Dsply Unit FontC	n/a	XX.XX	XX.XX	Xxxxxx
Dsply Unit FontJ	n/a	XX.XX	XX.XX	Xxxxxx

#### Note

Whenever a Software Download is completed, the specific software download is recorded in the following logs:

- in the System Download Log (Section 8b.8.1) that resides on the system (Display Unit)
- and in the CF Card Install Log (Section 8b.8.2) that resides on the Compact Flash card.

#### Notes about downloading software

If there is no Front Panel Control software installed in the system (as would be the case when the display units control board is replaced), the Service Application automatically downloads the Front Panel Controls software at startup. During the download the two display unit LEDs will flash and the display speaker will sound an alarm tone to indicate that Software Download is proceeding.

To ensure that all software versions on the system are compatible, the end result of "Download All" or "Download New" will be the same. The software loaded on the machine will exactly match what is on the card. Be sure to have the latest/correct version of software before attempting a download to avoid inadvertent overwrites of newer software with an older version.

If, during the "Download New" process, the compatibility checker detects a newer version of software component on the system, a "Notice" appears on the screen that asks you to confirm the downgrade.

"Download All" will download all compatible software from the card to the system without issuing a notice that newer version of software component may be on the system.

#### **Download process**

The PCMCIA card includes only the latest software for each subsystem along with the diagnostic application.

As each subsystem software segment is being downloaded, the following status messages note the state of each subsystem and the result of the download:

- In App System is running its application code; not ready for download.
- Ready System is in its boot code; ready for download.
- Loading System is accepting download data.
- Done Software download has completed successfully.
- **Fail** Software download did not complete successfully. A "Fail" message will require reloading of the software; or repair of the system may be necessary.
- Skipped Software download was bypassed.
- No Comm The subsystem is not communicating with the HPDU.

If the subsystem is communicating but the HW Rev or current SW Rev are not known, dashes will appear for those values.

As the software loads, an activity bar at the bottom of the screen shows the download progress for each subsystem.

#### **Download complete**

When all the required subsystem software is download, the following message appears on the screen. You must shut down the system to exit the download function.

#### DOWNLOAD IS COMPLETE.

Remove AC mains power. Turn on/standby switch to Standby. Then remove external CF card. Wait 20 seconds before restoring power to the system.

#### **Note**

After powering down the system, be sure to wait at least 20 seconds before restarting the system.

(HPDU) 8b-4 09/07 1009-0357-000

### **8b.3 Special Functions**

Selecting **Special Functions** on the Main Menu brings up the following menu selections in the left-hand frame:

**Note** Some functions may not be present on all versions of the Software Downloader Card.

### **Special Functions**

**Power Diagnostics** 

**Vent Diagnostics** 

Gas Diagnostics

Mixer Service

**Display Diagnostics** 

**Compatibility Table** 

System Download Log

CF Card Install Log

View Install Errors

-> Main Menu

# 8b.3.1 Power Diagnostics

The service application provides power supply diagnostics for various circuit boards in the Avance anesthesia machine.

Selecting **Power Diagnostics** on the Special Functions menu brings up the following menu selections.

### **Power Diagnostics**

**Power Controller** 

**Display Unit Power** 

**Anes Control Board** 

Mixer

Vent Interface Bd

-> Previous Menu

### 8b.3.2 Power Controller Power Diagnostics

There are two pages of diagnostics for the Power Controller.

Selecting **Power Controller** brings up the first page of the Power Controller Readings.

### Power Controller Readings, page 1

Power Diagnostics
Power Controller
Display Unit Power
Anes Control Board
Mixer
Vent Interface Bd
-> Previous Menu

Label	Value Format	Units	Normal range
AC Status	OK, Fail		
12Vdc Supply	XX.XX	Vdc	11.70 to 12.30
3.3Vdc Supply	X.XXX	Vdc	3.201 to 3.399
1.5 Vdc Supply	X.XXX	Vdc	1.450 to 1.550
Battery Connected	Yes, No		
Battery Status	Fail Bulk Chg Over Chg Float Chg Trickle Chg Discharge		
Battery Current	X.XX	Α	
Battery 1 Volts	XX.XX < 6.0 FAIL (red) <10 T Chg (yellow) 10-16.5 (Green)	Vdc	10.0 to 16.50
Battery 2 Volts	XX.XX < 6.0 FAIL (red) <10 T Chg (yellow) 10-16.5 (Green)	Vdc	10.0 to 16.50
Calc Battery Time	XX	Min	
Date battery Tested (*)	XX-ABC-XXXX		
Last Full Discharge Time (*)	XX	Min	

Turn the ComWheel to select the second page.

Push the ComWheel to return focus to the Power Diagnostics selection menu.

(\*) These entries are only listed after a full discharge test of the batteries is first performed (Section 6.10, "Battery capacity test"). These entries are not affected by the battery test in the display diagnostics section (Section 8b.7.1, "Test Keys and Battery").

(HPDU) 8b-6 09/07 1009-0357-000

### Power Controller Readings, page 2

Power Diagnostics		
Power Controller		
Display Unit Power		
Anes Control Board		
Mixer		
Vent Interface Bd		
-> Previous Menu		

Label	Value Format	Units	Normal range
Board Temperature	<65C OK (green) >65C <75C Warn (Yellow) >75C Fail (red)	Deg C	Upper 65
Fan Speed	Slow, Fast		
Fan 1 Voltage	XX.XX	Vdc	10.08 to 10.92
Fan 1 Status	Low, High, Fail, OK		

Turn the ComWheel to return to the first page.

Push the ComWheel to return focus to the Power Diagnostics selection menu.

### 8b.3.3 Display Unit Power Diagnostics

Selecting **Display Unit Power** brings up the Display Unit Power Diagnostics page.

### **Display Unit Power Readings**

Power Diagnostics
Power Controller
Display Unit Power
Anes Control Board
Mixer
Vent Interface Bd
-> Previous Menu

Label	Value Format	Units	Normal range
12.5V Backlight	XX.XX	Vdc	10.54 to 14.87
5V Audio Power	X.XX	Vdc	4.39 to 5.72
3.3V Digital	X.XX	Vdc	2.86 to 3.81
3.3V LCD	X.XX	Vdc	2.86 to 3.81
Inverter A	X.XX	Vdc	0.00 to 1.09
Inverter B	X.XX	Vdc	0.00 to 1.09
Ground 1	X.XX	Vdc	0.00 to 0.11
Ground 2	X.XX	Vdc	0.00 to 0.11

(HPDU) 8b-8 09/07 1009-0357-000

### 8b.3.4 Anesthesia Control Board Power Diagnostics

**Power Diagnostics** 

**Power Controller** 

Mixer

Display Unit Power

Anes Control Board

Vent Interface Bd
-> Previous Menu

There are two pages of diagnostics for the Anesthesia Control board.

Selecting **Anes Control Board** brings up the first page of the Anesthesia Control Board Power Readings.

### Anesthesia Control Board Power Readings, page 1

Label	Value Format	Units	Normal range
12.5Vdc from Pwr Cntrl	XX.XX	Vdc	11.90 to 12.90
ADC Reference	X.XXX	Vdc	4.018 to 4.176
Gas Select 10VA Volts	OK, Fail		
Gas Select 10VA Amps	OK, Fail		
P Xducer 10VA Amps	OK, Fail		
Vent Int Bd 10VA Volts	OK, Fail		
Vent Int Bd 10VA Amps	OK, Fail		
Vent Valves 10VA Volts	OK, Fail		
Vent Valves 10VA Amps	OK, Fail		
Acces 1 10VA Volts	OK, Fail		

OK, Fail

Turn the ComWheel to select the second page.

Acces 1 10VA Amps

Push the ComWheel to return focus to the Power Diagnostics selection menu.

Power Diagnostics
Power Controller
Display Unit Power
Anes Control Board
Mixer
Vent Interface Bd
-> Previous Menu

### Anesthesia Control Board Power Readings, page 2

Label	Value Format
Gas Unit 10VA Volts	OK, Fail
Gas Unit 10VA Amps	OK, Fail
Mixer 10VA Volts	OK, Fail
Mixer 10VA Amps	OK, Fail
Alt 02 10VA Volts	OK, Fail
Alt O2 10VA Amps	OK, Fail
Periph 1 10VA Volts	OK, Fail
Periph 1 10VA Amps	OK, Fail
Periph 2 10VA Volts	OK, Fail
Periph 2 10VA Amps	OK, Fail

Turn the ComWheel to return to the first page.

Push the ComWheel to return focus to the Power Diagnostics selection menu.

(HPDU) 8b-10 09/07 1009-0357-000

### 8b.3.5 Electronic Mixer Power Diagnostics

Selecting *Mixer* brings up the Electronic Mixer Power Diagnostics page.

# Power Diagnostics Power Controller Display Unit Power Anes Control Board Mixer Vent Interface Bd -> Previous Menu

### **Mixer Power Readings**

Label	Value Format	Units	Normal range
From Anes Cntrl Bd Mixer 10VA Volts	OK, Fail		
Board Supplies			
12.5 V	XX.XX	Vdc	11.80 to 13.00
5.5V	X.XXX	Vdc	5.390 to 5.610
3.3V CPU	X.XXX	Vdc	3.220 to 3.380
2.5V ADC Ref	X.XXX	Vdc	2.470 to 2.530

### 8b.3.6 Ventilator Interface Board Power Diagnostics

Selecting **Vent Interface Board** brings up the Ventilator Interface Board Power Diagnostics page.

Power Diagnostics
Power Controller
Display Unit Power
Anes Control Board
Mixer
Vent Interface Bd
-> Previous Menu

### **Vent Interface Bd Power Readings**

Label	Value Format	Units	Normal range
From Anes Cntrl Bd Vent Int Bd 10VA Volts Vent Valves 10VA Volts	OK, Fail OK, Fail		
Board Supplies			
Vent Int Bd 10VA Volts	XX.XX	Vdc	11.30 to 13.13
Vent Valves 10VA Volts	XX.XX	Vdc	11.30 to 13.13
3.2Vdc (12bit Vref)	X.XXX XXXX	Vdc Counts	3.179 to 3.221
1.22Vdc (10bit Vref)	X.XXX XXXX	Vdc Counts	1.074 to 1.367
+6.0Vdc	X.XXX	Vdc	5.510 to 6.500
-6.0Vdc	-X.XXX	Vdc	-6.720 to -5.280

(HPDU) 8b-12 09/07 1009-0357-000

### **8b.4 Ventilation Diagnostics**

The service application provides several pages for ventilator diagnostics.

Selecting **Vent Diagnostics** on the Special Functions menu brings up the following menu selections in the left-hand frame and the instructions in the right-hand frame:

Special Functions
Power Diagnostics
Vent Diagnostics
Gas Diagnostics
Mixer Service
Display Diagnostics
Compatibility Table
System Download Log
CF Card Install Log
View Install Errors
-> Main Menu

Vent Diagnostics	cs Ventilator Diagnostics		
Vent Status Vent Flow and Pressure	Vent Status shows the settings of breathing circuit switches.		
Toggle Gas Inlet VIv Toggle Circuit	Vent Flow and Pressure shows vent sensor readings.		
Set Flow Valve -> Previous Menu	Other menu commands let you control ventilator valves, and send gas to either the circle or the non-circle circuit.		

Vent Diagnostics	Action when selected
Vent Status	Selecting <b>Vent Status</b> displays the Ventilation Status page.
Vent Flow and Pressure	Selecting <b>Vent Flow and Pressure</b> displays the Ventilation Flow and Pressure page.
Toggle Gas Inlet VIv	Selecting <b>Toggle Gas Inlet VIv</b> alternately sets the Gas Inlet Valve to the Open or Closed position and sets the right-hand screen to the Ventilator Flow and Pressure page.
Toggle Circuit	Selecting <b>Toggle Circuit</b> toggles between commanding Circle or Non-Circle for the SCGO valve and sets the right-hand screen to the Ventilator Flow and Pressure page (does not apply to machines with ACGO).
Set Flow Valve	Selecting <b>Set Flow Valve</b> sets the right-hand screen to the Ventilator Flow and Pressure page and allows the user to increase or decrease the DAC count to the flow valve using the ComWheel.
->Previous Menu	Selecting <b>Previous Menu</b> closes the Gas Inlet Valve and the Flow Valve and returns to the Special Functions menu.

The actions of the Gas Inlet, the Flow Valve Control, and the Circuit selection are reflected on both the Ventilation Status page and the Ventilation Flow and Pressure page.

The Gas Inlet Valve must be in the Open position in order for gas to flow.

### 8b.4.1 Vent Status

Selecting **Vent Status** brings up the Ventilation Status page.

# Vent Diagnostics Vent Status Vent Flow and Pressure Toggle Gas Inlet VIv Toggle Circuit Set Flow Valve -> Previous Menu

### **Ventilator Status**

Menu Item	Value
Gas Inlet Valve Setting	Open or Closed
Gas Inlet Valve Feedback	Open or Closed, or Fault
Circuit Setting	Circle or Non-Circle
Circuit Feedback	Circle, Non-Circle, or Fault
Vent Drive Gas	Air or O2
ABS Installed	Yes or No
Flush Valve	Not Pressed or Pressed
02 Cell Status	Connected or None
Bag/Vent Switch	Bag or Vent
ACGO/SCGO Configuration	ACGO or SCGO
Over Pressure Circuit	OK or High Pressure

(HPDU) 8b-14 09/07 1009-0357-000

# 8b.4.2 Vent Flow and Pressure

Selecting **Vent Flow and Pressure** brings up the Ventilation Flow and Pressure page.

### **Ventilator Flow and Pressure**

Vent Diagnostics
Vent Status
Vent Flow and Pressure
Toggle Gas Inlet VIv
Toggle Circuit
Set Flow Valve
-> Previous Menu

Menu Item	Value	Units	Value	Units
Gas Inlet Valve Setting	Open or Clos	ed		
Gas Inlet Valve Feedback	Open or Clos	Open or Closed, or Fault		
Circuit Setting	Circle or Non	-Circle		
Circuit Feedback	Circle, Non-C	ircle, or Fault		
Inspiratory Flow	XXX.X	I/min	XXXX	Counts
Expiratory Flow	XXX.X	I/min	XXXX	Counts
Airway Pressure	XXX	cmH20	XXXX	Counts
Manifold Pressure	XXX	cmH20	XXXX	Counts
02 Cell	XXX	%	XXXX	Counts
ADC Ref Voltage	X.XXX	Vdc		
Flow Valve Setting	XXX.X	I/min	XXXX	Counts
Flow Valve Feedback	XXXX	mV	XXXX	Counts
Flow Valve Current	XXXX	mA	XXXX	Counts
Flow Valve DAC Setting	XXXX	Counts		
$\uparrow \uparrow \uparrow$				
Turn ComWheel to adjust flow Press the ComWheel to change flow setting				
(appears when Set Flow Valve is selected on the Vent Diagnostics menu)				

### **8b.5 Gas Diagnostics**

Selecting **Gas Diagnostics** on the Special Functions menu brings up the following menu selections in the left-hand frame and the instructions in the right-hand frame.

<b>Special Functions</b>
Power Diagnostics
Vent Diagnostics
Gas Diagnostics
Mixer Service
Display Diagnostics
Compatibility Table
System Download Log
CF Card Install Log
View Install Errors
-> Main Menu

Gas Diagnostics	Gas Diagnostics	
Gas Supplies		
Mixer Output	To troubleshoot a problem with gas supplies or flows, select Gas Supplies. Then work down through other menu	
Mixer Tests and Pres	selections.	
Mixer Temperatures	To turn on gas flows, select Set O2, N2O, or Air Flow.	
Set 100% 02 Flow	Available flows are 0.2, 0.5, 5, and 10 l/min and OFF.	
Set 100% N20 Flow	Turning ON a different gas will turn OFF any other gas.	
Set Air Flow	All gas flows stop if you go back to the Previous Menu.	
Breathing Sys Leak		
-> Previous Menu		
Settings		
O2 flow: XXX		
N20 flow: XXX		
Air flow: XXX		

Selecting any of the first four items in the left-hand frame, brings up the corresponding diagnostic page.

Press the ComWheel to return focus to the right-hand frame to make another selection.

### **Start gas flow**

Selecting one of the three gas flows brings up a page on which you can set one of the following flow values for the selected gas:

• 0.2 I/min, 0.5 I/min, 5.0 I/min, 10.0 I/min, 0FF.

Press the ComWheel to start the selected flow.

Press the ComWheel again to return to the flow page to observe the Airway Pressure reading. Or select one of the four diagnostics pages to view the displayed conditions.

Selecting Breathing Sys Leak, opens the GIV valve and sets a flow through the inspiratory flow valve. While on the Breathing Sys Leak page, you can set 02 flow to the breathing system and observe the Airway Pressure reading.

### Stop gas flow

Gas will continue to flow at the set rate until you:

- set a different flow for the same gas,
- set a flow for a different gas (automatically turns the previous gas flow OFF),
- set the current gas flow to OFF (all flows OFF),
- exit Gas Diagnostics to the Special Functions menu.

(HPDU) 8b-16 09/07 1009-0357-000

### 8b.5.1 Setting Gas Flow

Gas Diagnostics

Gas Supplies

Mixer Output

Mixer Tests and Pres

Mixer Temperatures

Set 100% 02 Flow

Set 100% N20 Flow

Set Air Flow

Breathing Sys Leak

-> Previous Menu

Settings
O2 flow: XXX
N2O flow: XXX
Air flow: XXX

Selecting **Set 100% 02 Flow** brings up the Set Gas Flow page on which you can select one of the following flow values (or OFF if previously set to flow):

- 0.2 l/min
- 0.5 l/min
- 5.0 l/min
- 10.0 l/min
- OFF

Selecting **Set 100% N20 Flow** or **Set Air Flow** brings up a comparable page for setting the above flows for the selected gas. Selecting a flow for a new gas turns the previous gas flow to OFF.

### **Set Gas Flow**

	Value	Units	
Set 02 Flow	X.X	I/min	
	<b>ተ</b> ተተ		
Turn the ComWheel to select the gas flow. Press the ComWheel to change flow setting.			
Airway pressure:	XXX	cmH2O	

After setting a gas flow, push the ComWheel again to return to the Flow Setting page to observe the Airway Pressure reading, or select one of the four diagnostics pages to view the displayed conditions.

All gas flow stops when you exit Gas Diagnostics to the Special Functions menu.

# 8b.5.2 Gas Supplies

Each gas supply shows the derived pressure in kPa and psi along with the raw voltage from the pressure transducer.

If a supply module for a gas is not installed, or if the transducer is disconnected, the supply shows dashed lines.

If a supply module is installed but no supply is connected, the supply shows approximately 0.500 Vdc.

The remaining items show the state of the Electronic Mixer selector valves.

### **Gas Supplies**

Gas Diagnostics		
Gas Supplies		
Mixer Output		
Mixer Tests and Pres		
Mixer Temperatures		
Set 100% 02 Flow		
Set 100% N20 Flow		
Set Air Flow		
Breathing Sys Leak		
-> Previous Menu		

J	
02 flow: XXX	
120 flow: XXX	

Settings

۸i۶	fla		VVV
НI	IIC	w.	XXX

Lab	el	psi	kPa	Vdc
02 (	Cylinder 1	XXXX	XXXXX	XX.XXX
02 (	Cylinder 2	XXXX	XXXXX	XX.XXX
Air C	Cylinder	XXXX	XXXXX	XX.XXX
N20	Cylinder	XXXX	XXXXX	XX.XXX
02 F	Pipeline	XXXX	XXXXX	XX.XXX
Air P	ripeline	XXXX	XXXXX	XX.XXX
N20	Pipeline	XXXX	XXXXX	XX.XXX
02.5	Select Valve	Open	Flow ON	
Air S	select Valve	Closed	Flow OFF	
N20	Select Valve	Closed	Flow OFF	
Alt C	02 Button	Not pressed		

Note: dashed pressure means supply is not installed or transducer is disconnected.

(HPDU) 8b-18 09/07 1009-0357-000

### 8b.5.3 Mixer Output

Selecting *Mixer Output* brings up the Mixer Output Gas Diagnostics page.

This data comes from the Anesthesia Control board. The Flow Verify signals are rough calculations of the mixer flow based on pressure drop and temperature. The ADC reference voltage is used to convert flow signals.

### **Mixer Output**

Gas Diagnostics
Gas Supplies
Mixer Output
Mixer Tests and Pres
Mixer Temperatures
Set 100% 02 Flow
Set 100% N20 Flow
Set Air Flow
Breathing Sys Leak
-> Previous Menu

Settings	
O2 flow: XXX	
N20 flow: XXX	
Air flow: XXX	

Label	Value	Units	
Measured 02 Flow	XX.XX	I/min	
Meas. Balance Gas Flow	XX.XX	I/min	
inleds. Edidlice Gas Flow	^^.^^	1/111111	
O2 Flow Verify	XX.XX	l/min	
O2 Flow Signal	X.XXX	Vdc	
O2 Prop Valve Drive	XXXX	mA	
Balance Gas ID	None, Air, N20	)	
Balance Flow Verify	XX.XX	I/min	
Balance Flow Signal	X.XXX	Vdc	
Balance Prop Valve Drive	XXXX	mA	
O2 Select Valve	Open		
Air Select Valve	Closed		
N20 Select Valve	Closed		
ADC Ref Voltage	X.XXX	Vdc	

# 8b.5.4 Mixer Tests and Pressure

Selecting *Mixer Tests and Pres* brings up the Mixer Tests and Pres Gas Delivery Diagnostics page.

### **Mixer Tests and Pressures**

# Gas Diagnostics Gas Supplies Mixer Output Mixer Tests and Pres Mixer Temperatures Set 100% 02 Flow Set 100% N20 Flow Set Air Flow Breathing Sys Leak -> Previous Menu

Settings

02 flow: XXX

N20 flow: XXX

Air flow: XXX

Last Power-Up Tests			
O2 Proportional Valve Leak	rtional Valve Leak  Not done. No supply pressure  Not done. Selector valve incorrect sta  Pass  Fail. Selector valve leaks  Fail. Proportional valve leaks		orrect state
Alt O2 Valve Leak	Pass; Fail		
Balance Gas Prop Valve Leak	Not done. No supply pressure Not done. Selector valve incorrect state Pass Fail. Selector valve leaks Fail. Proportional valve leaks		
O2 Flow Test	Not done. No supply pressure; Not done. Selector valve incorrect state; Pass; Fail, 3 I/min test; Fail, 10 I/min test		
Balance Flow Test  Not done. No supply pressure; Not done. Selector valve incompass; Fail, 3 I/min test; Fail, 10 I/min test			
Balance Gas ID	None; Air; N2O		
Pressure Data	Value	Units	kPa
O2 Pressure (P1)	XX.XX	psi	XXX.XX
02 Pres Cal (P1)	X.XXX	Vdc	
Balance Pressure (P2)	XX.XX	psi	XXX.XX
Balance Pres Cal (P2)	X.XXX	Vdc	
Mixer Output Pres (P3)	XX.XX	psi	XXX.XX
Mixer Output Pres Cal (P3)	X.XXX	Vdc	
ADC Ref Voltage	X.XXX	Vdc	

(HPDU) 8b-20 09/07 1009-0357-000

### **8b.5.5 Mixer Temperature**

Gas Diagnostics

Gas Supplies

Mixer Output

Mixer Tests and Pres

Mixer Temperatures

Set 100% 02 Flow

Set 100% N20 Flow

Set Air Flow

Breathing Sys Leak

-> Previous Menu

Settings

O2 flow: XXX

N20 flow: XXX

Air flow: XXX

Selecting *Mixer Temperatures* brings up the Mixer Temperatures Gas Diagnostics page.

### **Mixer Temperatures**

Sensor Data	Value	Units
O2 Temp (T1)	XX.X	Deg C
O2 Temp Volts (T1)	X.XXX	Vdc
Balance Temp (T2)	XX.X	Deg C
Balance Temp Volts (T2)	X.XXX	Vdc
ADC Ref Voltage	X.XXX	Vdc

### 8b.5.6 Breathing System Leak Test

Selecting **Breathing Sys Leak** brings up the System Leak Test page.

Entering the System Leak Test page, opens the Gas Inlet Valve and sets up a constant flow through the Inspiratory Flow Valve.

On the Breathing Sys Leak page you can set an O2 flow through the Gas Mixer and view the Airway Pressure.

### **System Leak Test**

Gas Diagnostics

Gas Supplies

Mixer Output

Mixer Tests and Pres

Mixer Temperatures

Set 100% 02 Flow

Set 100% N20 Flow

Set Air Flow

Breathing Sys Leak

Settings

O2 flow: XXX

N2O flow: XXX

Air flow: XXX

-> Previous Menu

This is the Breathing System Leak Test. Refer to the Technical Reference Manual for this procedure.		
	Value	Units
Set 02 Flow	X.X	I/min
	$\uparrow\uparrow\uparrow$	
Turn the ComWheel to select the gas flow.  Press the ComWheel to change flow setting.		
NOTE: Flow Valve and Ga	as Inlet Valve are OPI	EN.
Airway pressure:	XXX	cmH20

All gas flow stops when you exit Gas Diagnostics to the Special Functions menu.

**Note** This procedure is detailed in Test 6, "Testing the bellows module and the Bag/Vent switch" on page 7-21.

(HPDU) 8b-22 09/07 1009-0357-000

### **8b.6 Mixer Service Menu**

Selecting *Mixer Service* brings up the Mixer Service Functions screen. This screen includes and automatic and a manual leak test of the Mixer's balance gas inlet check valves.

### **Automatic**

The **BalGas Auto Leak** test does not require disassembly of the system. It is a sensitive test that will Pass check valve that have leak rates within specifications; however, it may Fail some Mixers with marginal but acceptable leak rates.

- Mixers that Pass the auto leak test do not require further testing and can be left in service.
- Mixers that Fail the auto leak test should be further tested using the manual leak test.

### Manual

The manual BalGas Leak Test requires disassembly of the system.

- Mixers that Pass the manual leak test have acceptable leak rates and can be left in service.
- Mixers that Fail the manual leak test should be replaced.

Special Functions
Power Diagnostics
Vent Diagnostics
Gas Diagnostics
Mixer Service
Display Diagnostics
Compatibility Table
System Download Log
CF Card Install Log
View Install Errors

Mixer Service	Mixer Service Functions
BalGas Auto Leak	
BalGas Leak Test	To troubleshoot mixer problems.
-> Previous Menu	

### **Automatic leak test**

-> Main Menu

Selecting **BalGas Auto Leak** brings up the Service Instructions page for the automatic mixer leak test.

Service command	Service Instructions
Start Test -> Previous Menu	This will perform an Automatic Leak Test of the Balance Gas Check Valve

- 1. Connect an O<sub>2</sub> supply.
- 2. Select **Start Test** to perform the automatic leak test.
  - If Pass, balance gas check valve leak rate is acceptable.
  - If Fail, verify leak rate using the following manual leak test.

### Manual leak test

Selecting **BalGas Leak Test** brings up the Service Instructions page for the manual mixer leak test.

Service command	Service Instructions
Start Test	
-> Previous Menu	This will perform a Manual Leak Test of the Balance Gas Check Valve

Select **Start Test** to perform the manual leak test as follows.

Selecting Start Test brings up the manual leak test setup instructions.

Service command
Start Test
-> Previous Menu

### **Service Instructions**

This will perform a Manual Leak Test

- 1) Disconnect Oxygen and Bal Gas Supplies.
- 2) Remove rear cover (Avance) or dashboard (Aisys).
- 3) Disconnect mixer outlet tube (MIXER VAP IN) at the inlet to the vaporizer manifold (Avance), or Electronic Vaporizer (Aisys), and plug the Mixer outlet (1/4 Legris).
- 4) Disconnect the Mixer Alt O2 inlet tube at the On/Standby switch (SW4 ALT O2 IN).
- 5) Connect a pressurization device with gauge to the Mixer Alt-O2 inlet tubing. Minimize the tubing length, as additional volumes will affect the leak rate.

Push ComWheel to 'Start Test'.

### **Service Instructions**

- 1) Slowly pressurize the Mixer to 400 mmHg (over a 5 second period), as read on the test device.
- 2) The Pressure shown on the test gauge should not decay to ambient (zero) in less than 10 seconds.
- 3) Select 'End Test' when done.

Push ComWheel to 'End Test'.

### **Service Instructions**

- 1) Remove all test fixtures.
- 2) Remove plugs and reconnect all tubing disconnected at the start of this test.
- 3) Re-attach Oxygen supply and activate 'Confirm'.

Push ComWheel to 'Confirm'.

(HPDU) 8b-24 09/07 1009-0357-000

## **8b.7 Display Diagnostics**

Selecting **Display Diagnostics** brings up the Display Diagnostics menu.

Special Functions
Power Diagnostics
Vent Diagnostics
Gas Diagnostics
Mixer Service
Display Diagnostics
Compatibility Table
System Download Log
CF Card Install Log
View Install Errors
-> Main Menu

Display Diagnostics		
Test LEDs	Board part #:	XXX-XXX-XXX B XXX
Test Speaker	Serial #:	ABCXXXXX
Test Backlight 1	BIOS Version:	XX.XX
1000 Baokingine 1	FPC Version:	XX.XX
Test Backlight 2	MBC Version	X
Test Soft Keys	MAC Address:	XX XX XX XX XX XX
Test Keys and Battery	Internal CF disk:	
Test LCD Pixels	Card geometry:	
	CPU Temperatu	ure = XX.X deg C
-> Previous Menu	CPU Fan Spee	ed = XXXX RPM

Display Diagnostics	Action when selected	
Test LEDs	Selecting <b>Test LEDs</b> causes the red and yellow LEDS next to the <b>Silence Alarms</b> key to flash 5 times.	
Test Speaker	Selecting <b>Test Speaker</b> causes the speaker to sound for 2 seconds.	
Test Backlight 1	Selecting <b>Test Backlight 1</b> turns backlight 2 off for 4 seconds. "If screen goes blank during test then a backlight is out."	
Test Backlight 2	Selecting <b>Test Backlight 2</b> turns backlight 1 off for 4 seconds. "If screen goes blank during test then a backlight is out."	
Test Soft Keys	On the <b>Test Soft Keys</b> screen, pressing a softkey will cause a highlighted square to appear next to the soft key. A second press of the softkey removes the square.	
Test Keys and Battery	Refer to section 8b.7.1.	
Test LCD Pixels	On the <b>Test LCD Pixels</b> screen, press the ComWheel to sequence through the color screens. The first press results in a "blank" screen.	
-> Previous Menu	Selecting <b>Previous Menu</b> returns to the Special Functions screen.	

# 8b.7.1 Test Keys and Battery

**Note**: System Software version 4.X only includes the **Test Keys** function. System Software version 5.X includes the combined **Test Keys and Batery** function.

Selecting **Test Keys and Battery** brings up an expanded test keys screen which includes the battery test.

### **Display Diagnostics**

Test LEDs

Test Speaker

Test Backlight 1

Test Backlight 2

Test Soft Keys

### Test Keys and Battery

Test LCD Pixels

-> Previous Menu

### **Test Soft Keys and Battery**

Press each softkey to display a box.

Press again to clear it.

To start the battery test, shut off the AC mains. Then wait for the end of the 2-minute countdown.

(You can test the keys while you wait.)

After turning off the AC main inlet switch, the test begins automatically in a few seconds.

### **Test Soft Keys and Battery**

Press each softkey to display a box.

Press again to clear it.

Countdown = XXX sec

Battery 1 (right) = XX.XX Vdc

Battery 2 (left) = XX.XX Vdc

Battery Current = -X.XX A

Battery test PASSED/FAILED

Restore AC mains power.

Note

If the Avance system is operating under battery power, the test begins as soon as you make the **Test Keys and Battery** selection, bypassing the instructions.

(HPDU) 8b-26 09/07 1009-0357-000

### **8b.8 Compatibility Table**

Selecting **Compatibility Table** brings up a list of the current software components that were last downloaded on to the system. In essence, it is the latest listing that appears in the Revision Log (*Refer to section 8b.8.1*). The Compatibility Table allows you to view the current log directly without having to scroll to it in the Revision Log.

### **Special Functions**

Power Diagnostics

Vent Diagnostics

Gas Diagnostics

Mixer Service

**Display Diagnostics** 

### Compatibility Table

System Download Log CF Card Install Log

View Install Errors

-> Main Menu

```
# System s/n ABCDXXXXXXXXX loaded by Download XXX on (day) (date) (time)
# using HPDU Softwr Loader & Svc Ver XX.XX from card [XXXX] to [XXXX]
System Version XX.XX
Avance FPC, *
                                   (Software Level), (File Name) (# n/a ) Front Panel CN
Avance PSC, (Stock Number) (RevX),
                                   (Software Level), (File Name) (Serial #) Power Controll
Avance ACB, (Stock Number) (Rev X),
                                   (Software Level), (File Name) (Serial #) AnesControl B
Avance VNT, (Stock Number) (Rev X),
                                   (Software Level), (File Name) (Serial #) Vent Intface B
Avance MXR, (Stock Number) (RevX), (Software Level), (File Name) (Serial #) Electronic Mix
Avance DUB, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Dsply Unit BIO
Avance DUA, (Stock Number) (Rev X), (Software Level), (File Name) (# n/a) Dsply Unit App
                 n/a ) , (Software Level), (File Name) (# n/a ) ModBus Control
Avance MBC, (
Avance MHB, (
                 n/a
                                   (Software Level), (File Name) (# n/a ) Dsply Unit Fon
                       )
Avance HGG, (
                                   (Software Level), (File Name) (# n/a ) Dsply Unit Fon
                 n/a )
    < End of File>
```

# 8b.8.1 System Download Log

Whenever a Software Download is completed, the specific software download is recorded

- in the System Download Log that resides on the system (Display Unit)
- and in the CF Card Install Log that resides on the Compact Flash Card.

Selecting **System Download Log** brings up the Revision Log for the system. The log includes chronological entries for every Software Download that was completed to the system.

Note: To view currently downloaded system software, scroll to last entry in the log (<End of File>) or view the Compatibility Table, instead.

```
# System s/n ABCDXXXXXXXX loaded by Download XXX on (day) (date) (time)
# using HPDU Softwr Loader & Svc Ver XX.XX from card [XXXX] to [XXXX]
System Version XX.XX
Avance FPC, *
                                , (Software Level), (File Name) (# n/a ) Front Panel CN
Avance PSC, (Stock Number) (RevX), (Software Level), (File Name) (Serial #) Power Controll
Avance ACB, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) AnesControl B
Avance VNT, (Stock Number) (RevX), (Software Level), (File Name) (Serial #) Vent Intface B
Avance MXR, (Stock Number) (RevX), (Software Level), (File Name) (Serial #) Electronic Mix
Avance DUB, (Stock Number) (Rev X), (Software Level), (File Name) (Serial #) Dsply Unit BIO
Avance DUA, (Stock Number) (Rev X), (Software Level), (File Name) (# n/a ) Dsply Unit App
Avance MBC. (
                 n/a ) , (Software Level), (File Name) (# n/a ) ModBus Control
Avance MHB, (
                 n/a
                                  (Software Level), (File Name) (# n/a ) Dsply Unit Fon
                       ) ,
Avance HGG, (
                                   (Software Level), (File Name) (# n/a ) Dsply Unit Fon
                 n/a )
```

### **Note**

The Stock Number listed is for the board assembly and may not represent an orderable service item. Refer to the parts lists in Section 10 for service level stock numbers.

The Front Panel Control (FPC), Display Unit Application (DUA), the ModBus Control (MBC), the two font files (MHB and HGG) reside, along with the Display Unit BIOS (DUB), on the Display Unit CPU board.

# 8b.8.2 CF Card Install Log

Selecting **CF Card Install Log** brings up the CF (CompactFlash) Card Install Log for the software download card. The log includes chronological entries for every Software Download that was completed with the card.

(HPDU) 8b-28 09/07 1009-0357-000

# 8b.8.3 View Install Errors

If Software Download detects an incompatible subsystem, an error message noting the incompatible subsystem is recorded on the CF card.

### ERROR!! INCOMPATIBLE SOFTWARE.

The software version on the CF card is not compatible with the installed XXX subsystem.

Installed part #: (Stock Number) (Rev X), swver XX.XX CF card part #: (Stock Number) (Rev X), swver XX.XX

Note: The Install Errors log includes two "screen dumps" for each occurrence of an error. The last "screen dump" includes the error message at the point where the incompatibility was detected. You can scroll up to the next "screen dump" in sequence to view the completed downloads.

Notes

(HPDU) 8b-30 09/07 1009-0357-000

# **9 Repair Procedures**

In this section	9.1 Circuit board replacement precautions	9-4
	9.2 How to bleed gas pressure from the machine	9-5
	9.3 How to remove the rear panels	9-6
	9.3.1 To remove the rear upper panel	9-6
	9.3.2 To remove the lower access panels	
	9.4 How to remove the tabletop	9-7
	9.5 Servicing the Display Unit (DU)	9-8
	9.5.1 Remove the Display Unit	9-8
	9.5.2 Disassemble the Display Unit	
	9.5.3 To replace the CPU board	9-10
	9.5.4 To replace the LCD display	9-11
	9.5.5 To replace the backlights	9-13
	9.5.6 To replace the Inverters	9-14
	9.5.7 To replace the front enclosure or components	9-15
	9.6 Servicing the High Performance Display Unit (HPDU)	9-17
	9.6.1 Remove the Display Unit	9-17
	9.6.2 Disassemble the Display Unit	9-18
	9.6.3 CPU Fan	9-18
	9.6.4 To replace the CPU board	9-19
	9.6.5 To replace the LCD display	9-20
	9.6.6 To replace the backlights	9-21
	9.6.7 To replace the Inverters	9-22
	9.6.8 To replace the front enclosure or components	9-23
	9.7 Replacing the Display and MGAS cables	9-25
	9.7.1 Remove the MGAS oxygen partition	9-25
	9.8 Servicing the lower electrical enclosure components	9-26
	9.8.1 Power Controller board (original)	9-26
	9.8.2 Power Controller board and Universal Power Supply	9-27
	9.8.3 Anesthesia Control board	9-28
	9.8.4 Backup batteries	9-29
	9.8.5 Fan	9-30
	9.8.6 Display Connector board	9-30
	9.9 Servicing the pan electrical enclosure components	9-31
	9.9.1 Electronic Gas Mixer assembly	9-31
	9.9.2 Ventilator Interface board	9-32
	9.9.3 Filter board	9-33
	9.9.4 Pan Connector board	9-34
	9 9 5 Pan enclosure fan	9-34

1009-0357-000 09/07 9-1

9.10 Servicing the Vent Engine	. 9-35
9.10.1 To remove the Vent Engine	9-36
9.10.2 Replacing Vent Engine components	. 9-37
9.10.3 Replacing GIV components	. 9-38
9.11 Servicing the pipeline inlet manifold components	. 9-39
9.11.1 Replace pipeline inlet filter	. 9-39
9.11.2 Replace pipeline inlet check valve	
9.11.3 Replace the inlet manifold	
9.12 Service the cylinder supply modules	. 9-41
9.12.1 Replace primary regulator module (complete replacement)	
9.12.2 Replace cylinder inlet filter	
9.12.3 Replace cylinder check valve	
9.13 Replace gas-supply pressure transducers	
9.14 Service vaporizer manifold parts	
9.14.1 Repair manifold port valve	
9.14.2 Checkout procedure for manifold port valve	
9.14.3 Replace vaporizer manifold check valve	
9.14.4 Replace vaporizer pressure relief valve	
9.14.5 Replace vaporizer manifold	
4.15 Replace ACGO selector switch	
9.16 Clean or replace ACGO port flapper valve	
9.17 Replace the APL valve	
9.18 Replace the bag support arm	. 9-54
9.18.1 Servicing the bag support arm	. 9-55
9.18.2 Replace friction pad in lower bag arm assembly	
9.18.3 Replace bag port housing	
9.19 Replace system switch assembly	
9.20 Replace Alt O2 components	. 9-60
9.21 Replace auxiliary O <sub>2</sub> flowmeter	. 9-61
9.22 Replace the suction regulator	. 9-62
9.23 Replace task light components	. 9-63
9.23.1 To replace the task light switch	. 9-63
9.23.2 To replace the upper task light	. 9-63
9.23.3 To replace the lower task light	. 9-64
9.24 Replace ABS breathing system components	. 9-65
9.24.1 Replace Bag/Vent switch assembly	. 9-65
9.24.2 Replace bellows base latch assembly	. 9-66
9.24.3 EZchange Canister spring replacement	
9.25 Replace casters	
9.26 Reconfigure sample gas return line	. 9-69
9 27 Change drive gas	9-70

9-2 09/07 1009-0357-000

### **⚠ WARNING**

### To prevent fires:

- Use lubricants approved for anesthesia or O<sub>2</sub> equipment, such as Krytox.
- Do not use lubricants that contain oil or grease; they burn or explode in high
   O<sub>2</sub> concentrations.
- All covers used on the system must be made from antistatic (conductive) materials. Static electricity can cause fires.
- Obey infection control and safety procedures. Used equipment may contain blood and body fluids.
- A movable part or a removable component may present a pinch or a crush hazard. Use care when moving or replacing system parts and components.
- Some internal parts have sharp edges and can cause cuts or abrasions. Use care when servicing internal components.
- After repairs are completed, always perform the checkout procedure. Refer to Section 3 of this manual.

### **⚠** CAUTION

Electrostatic discharge through circuit boards may damage the components on them. Wear a static control wrist strap before touching the circuit boards. Handle all circuit boards by their non-conductive edges. Use anti-static containers when transporting them.

1009-0357-000 09/07 9-3

### 9.1 Circuit board replacement precautions

The Aisys anesthesia system has processors on several boards. On three of these boards, information such as the machine serial number and optional ventilation modes (PCV, SIMV, and PSVPro) are stored redundantly.

During power-up, the machine serial number and installed options information stored on the boards are compared. If one board differs, information from the two agreeing boards will be written to the new board. If three boards differ (in the case of two boards replaced) the system defaults to "NO OPTIONS" and default machine serial number.

To retain the installed options, install only one replacement board at a time.

If multiple boards are to be installed, install the first board, load software on the new board, and power-up the machine in normal mode. Repeat this procedure for each board installation.

The following table lists the actions required after replacing printed circuit boards:

Board Name (Short Name)	Required Action After Installation
Display Unit CPU (DU CPU)	Load Software (see <b>Note</b> ). Check / Re-Configure Machine Configurations. Affix the new Key Code and Board ID Label to Vent Casting. Preoperative Checkout.
Anesthesia Control Board (ACB)	Load Software. Check / Re-Configure Machine Configurations. User Calibrations (O2 Cell, Flow Sensor, etc.). Gas Transducer Zero. All Ventilator Calibrations. Preoperative Checkout.
Power Controller (PCB)	Load Software. Preoperative Checkout.
Ventilator Interface Board (VIB)	Load Software. User Calibrations (O2 Cell, Flow Sensor, etc.). All Ventilator Calibrations. Preoperative Checkout.
Electronic Mixer (Mixer)	Load Software. Zero Mixer Pressure Sensors. Preoperative Checkout.
All Others	Preoperative Checkout.
<b>Note</b> : Flash software starts loading immediately when the Download Application first boots.	

**Note**: Flash software starts loading immediately when the Download Application first boots. Do not interrupt the Flash download. Allow the download to complete before proceeding.

9-4 09/07 1009-0357-000

### 9.2 How to bleed gas pressure from the machine

Before disconnecting pneumatic fittings, bleed all gas pressure from the machine.

- 1. Close all cylinder valves and disconnect all pipeline supplies from the source.
- 2. Set the system switch to On.
- 3. Ensure that all cylinder and pipeline pressures read zero.
- 4. Establish a flow for the affected gas to bleed down the pressure.
- 5. Set the system switch to Standby.

1009-0357-000 09/07 9-5

### 9.3 How to remove the rear panels

You must remove the rear upper panel to repair or replace many of the machine's components. To access the rear electronics enclosure, you must remove the lower access panels.

### 9.3.1 To remove the rear upper panel

- 1. Bleed all gas pressure from the machine (Section 9.2).
- 2. Ensure that all cylinder and pipeline pressures read zero before proceeding.
- 3. Disconnect all electrical cables.
- 4. To remove the rear panel, fully loosen the three captive screws that hold the panel in place. Remove the panel.
  - If the machine includes integrated suction, disconnect the two tube fittings from the overflow safety trap manifold.
  - If the machine includes electrical outlets, lower the panel and place it so that it does not stress the power cable.

# 9.3.2 To remove the lower access panels

- 1. Disconnect the power cord from the AC mains supply.
- 2. Bleed all gas pressure from the machine (Section 9.2).
- 3. Ensure that all cylinder and pipeline pressures read zero before proceeding.
- 4. If present, remove the inboard cylinders.
- 5. Remove the small upper access panel (A) to access the display connector board.




6. Remove the large lower access panel (B) to access the electrical enclosure.





(original)



9-6 09/07 1009-0357-000

### 9.4 How to remove the tabletop

The tabletop is held in place with five captive screws along the periphery of the pan assembly (accessed from below the rim of the tabletop).

- One screw (A) is in a deep recess at the right-rear corner of the tabletop.
- Two screws (**B**) are at the front of the tabletop: one screw is at the right corner of the tabletop, one is near the O<sub>2</sub> Flush button.
- To access the remaining two screws (**C**), you must remove the ABS: one screw is at the left corner of the tabletop, one is near the APL Valve.



1009-0357-000 09/07 9-7

### 9.5 Servicing the Display Unit (DU)


Note

To service the HPDU, refer to Section 9.6.

The item numbers appearing in parenthesis in this section refer to items in the parts list in Section 10.7.

The fan filter (26) and the access door (27) to the PCMCIA interface can be replaced with the Display Unit in place.

To replace the filter, slide the filter capsule (25) downward to remove it from the Display Unit.




To service other components of the Display Unit, you must first remove the Display Unit from the machine.

# 9.5.1 Remove the Display Unit

The Display Unit is held in place with four screws: three screws fit into keyhole slots and do not have to be fully removed; the fourth screw secures the Display Unit to the front bezel and must be removed to free the Display Unit from the machine.

- 1. Open the rear access door directly behind the Display Unit.
- 2. Loosen (L) the top two and the bottom-right mounting screws.
- 3. Remove (R) the lower-left mounting screw.



- 4. Place a protective pad on the tabletop.
- 5. From the front of the machine, slide the Display Unit toward the center of the machine to free it from the keyhole slots.
- 6. Lower the Display Unit face down on the protective pad.
- 7. Remove the cables from the rear connector panel.


9-8 09/07 1009-0357-000

# 9.5.2 Disassemble the Display Unit

Place the Display Unit face down on an anti-static pad. Before removing the rear enclosure, ensure that the release tabs on the PCMCIA frame are fully depressed.

- 1. Loosen (L) the four captive screws at each corner of the rear enclosure.
- 2. Lift the rear enclosure slightly and pivot it away from the lower enclosure at the bottom side of the Display Unit.





At this point, you can replace the following items (the item numbers refer to the parts list in Section 10.7):


- the fan (24)
- the connector panel assembly (21)
- the encoder assembly (12)
- the IRDA board (\*) or IRDA board harness (\*\*)
   (First production DUs included an IRDA board that is no longer used. System software does not support the IRDA board, which is now obsolete.)
- the battery (5)
- the **speaker** (20) To access the mounting screws for the speaker, you must first remove the ten screws that hold the mounting plate to the front enclosure so that you can raise the bottom edge of the assembly slightly Refer to section 9.5.4.)
- the **rear enclosure** (1) You can transfer the captive screws to the new enclosure. However, the **gasket** (28) is held in place with adhesive. When replacing the rear enclosure, also include a new gasket.

To replace the remaining items requires further disassembly.

1009-0357-000 09/07 9-9

# 9.5.3 To replace the CPU board


- 1. Remove the **connector panel assembly** (21) two screws.
- 2. Disconnect the following cables:
  - Inverter harnesses (A)
  - Membrane switch flex-cable at ZIF (zero insertion force) connector (B)
  - Speaker cable (C)
  - Encoder assembly cable (D)
  - If present, IRDA board cable (E)
  - Membrane switch flex-cable at ZIF (zero insertion force) connector (F)
  - Fan cable (G)
  - LCD cable (H)



- 3. Remove the remaining four screws (I) that hold the CPU board to the mounting plate.
- 4. Remove the CPU board from the mounting plate.
- 5. If you are **replacing the PCMCIA frame** (3) on an existing CPU board
  (remove four screws on back of CPU
  board), you must also apply new
  gaskets (4) to the frame. Align the
  ends of the gaskets with the top
  edge of the frame.
- 6. Transfer the battery (5) to the new
- CPU board.

7. Reassemble in reverse order.


- 8. Download latest software (Section 8a.6).
- 9. Reconfigure the Machine Configurations (Install/Setups).
- 10. Affix a new Key Code and Board ID label to the Vent casting.




9-10 09/07 1009-0357-000

### 9.5.4 To replace the LCD display


- 1. Disconnect the following cables:
  - Membrane switch flex-cable at ZIF (zero insertion force) connector (A)
  - Encoder assembly cable (B)
  - If present, IRDA Interface cable (C) remove IRDA Interface board (\*)
  - Membrane switch flex-cable at ZIF (zero insertion force) connector (**D**)
  - Fan cable (E)
  - LCD cable (F)
- 2. Remove the ten screws (circled) that hold the mounting plate to the front enclosure.




- 3. Remove the mounting plate assembly from the front enclosure.
- 4. Disconnect the backlight harnesses (**G**) from the inverter boards.
- 5. Slide the grommet (**H**) out of the mounting plate slot (transfer to new LCD).
- 6. Remove the four screws (circled) that hold the LCD to the mounting plate.



- 7. Lift the left side of the LCD display slightly away from the mounting plate to pull some of the display ribbon cable (I) to the top side of the plate. Flip the LCD over to the left of the assembly.
- 8. Disconnect the display ribbon cable (J).





9. Reassemble in reverse order.

**Note**: When replacing the LCD, pull the excess ribbon cable to the bottom side of the plate as you lower the LCD on to the plate. For the backlight harness grommet (**H**), ensure that the slit in the grommet faces toward the inside of the keyhole.

9-12 09/07 1009-0357-000

#### 9.5.5 To replace the backlights

#### **Note**

Early production Avance Display Units (DU) included a DG31 LCD display. Subsequent display units (and the HPDU) include a DG41 LCD display. The backlight in the DG31 display is not compatible with the DG41 display.

It is impossible to determine which LCD is installed without disassembling the Display Unit. The serial numbers reference below are approximate, some machines near the cutoffs may have either the DG31 or the DG41 installed.

- DG31 Backlight Kit (Avance serial numbers less than ABNG00380)
- DG41 Backlight Kit (Avance serial numbers greater than ABNH00400)


To replace the backlights in a DG41 display, refer to Section 9.6.6.


To replace the backlights in a DG31 display, continue below.

Note: When replacing a backlight or a backlight inverter, you must replace both inverters and the backlight assembly found in the Backlight Kit.

The backlight replacement kit for a DG31 display includes a backlight assembly (with two backlights) and two inverters with mounting hardware. To replace the backlight assembly follow the procedure in Section 9.5.4 to gain access to the assembly. To replace the inverters, follow the procedure in the next section.

- 1. Remove the one screw (**K**) that holds the backlight assembly to the LCD.
- 2. Slide the backlight assembly to the left to free it from the retaining tabs and then lift it out of the holder.





- 3. Transfer the grommet (**H**) to the new backlight assembly.
- 4. Reassemble in reverse order.

### 9.5.6 To replace the Inverters


Note: When replacing a backlight or a backlight inverter, you must replace both inverters and the backlight assembly found in the Backlight Kit.

The Display Unit includes two inverters (one for each backlight).


The inverters "sandwich" the mounting plate and use it as a heatsink. Follow the procedure in Section 9.5.4 to gain access to the inverters. Replace one inverter at a time.

- 1. Disconnect the backlight cable (A) from the inverter.
- 2. Remove the two Nylon screws (B) that hold the inverter to the backplate.
- 3. Slide the inverter out of the sleeve and disconnect it from the CPU harness (C).
- 4. Reassemble in reverse order.

#### The "front" inverter



#### The "rear" inverter



9-14 09/07 1009-0357-000


#### 9.5.7 To replace the front enclosure or components

Disassemble the Display Unit following procedures in the previous sections to the point where you have removed the mounting plate assembly from the front enclosure.

If you are replacing the front enclosure, you can transfer the encoder (12) assembly to the new enclosure; but, you must build up the replacement enclosure with:

- a new window (11)
- new membrane switches right-side (14), lower (15), left-side spacer (19)
- new keypads right-side (16), lower (17), left-side blank (18)
- new EMC gasket (10)

If you are replacing a keypad or a membrane switch, you must replace both items.




#### To replace the window


- 1. Place the front enclosure face up on a flat surface.
- 2. Press down on one corner of the window to free it from the enclosure.
- 3. Work your way around the window until you can get a hold of it from the back.
- 4. Slowly pry the window from the enclosure.
- 5. Place the front enclosure face down on a flat surface, taking care not to damage the encoder.
- 6. Remove any remaining residue from the mounting area; clean with isopropyl alcohol.
- 7. Remove the inside protective material from the front of the window.
- 8. Peel the front outside frame of the release liner.
- 9. Lower the window straight down in the enclose, noting the notch in the window and the matching nub on the enclosure.
- 10. Before seating the window, position it in contact with the bottom and right sides of the frame (see arrows) so that the larger gap between the window and the enclosure is at the top and left edges (as viewed from behind).
- 11. Remove the protective film from the back side of the window.

### To replace a membrane switch and keypad

- 1. Remove the screw that attaches the grounding strap to the enclosure.
- 2. Pry the membrane switch and keypad from the enclosure.
- 3. Remove any remaining residue from the mounting area; clean with isopropyl alcohol.
- 4. Remove the backing from the membrane. Be sure to remove the small backing below the flex cable. Be careful not to allow the ribbon cable to adhere to the backing.



5. Insert the flex cable and ground strap through the slot in the enclosure. Ensure that all of the ground strap passes through the slot an does not remain folded over under the membrane.




- 6. Carefully lower the membrane straight down to the enclosure. Seat the membrane in place.
- 7. Remove the backing from the keypad and install it over the membrane switches.
- 8. Attach the ground strap to the enclosure.

### To install the EMI gasket

To fully seal the Display Unit enclosure, you will need approximately 2.3 meters of EMC gasket (10). Cut the gasket into five strips shown below.

Insert a continuous length of gasket in the outside grove of the enclosure (sparingly apply "Super Glue Gel" to the channels near the corners before installing the gasket).

Insert individual lengths of gasket in the inside grove around the window (sparingly apply "Super Glue Gel" to the channels near the corners before installing the gasket).




9-16 09/07 1009-0357-000

#### 9.6 Servicing the High Performance Display Unit (HPDU)

Note


To service the DU, refer to Section 9.5.

The item numbers appearing in parenthesis in this section refer to items in the parts list in Section 10.8.



The fan filter (26) and the access door (3) to the PCMCIA interface can be replaced with the Display Unit in place.

To replace the filter, slide the filter capsule (25) downward to remove it from the Display Unit.




To service other components of the Display Unit, you must first remove the Display Unit from the machine.

### 9.6.1 Remove the Display Unit

The Display Unit is held in place with four screws: three screws fit into keyhole slots and do not have to be fully removed; the fourth screw secures the Display Unit to the front bezel and must be removed to free the Display Unit from the machine.

- 1. Open the rear access door directly behind the Display Unit.
- 2. Loosen (L) the top two and the bottom-right mounting screws.
- 3. Remove (R) the lower-left mounting screw.




- 4. Place a protective pad on the tabletop.
- 5. From the front of the machine, slide the Display Unit toward the center of the machine to free it from the keyhole slots.
- 6. Lower the Display Unit face down on the protective pad.
- 7. Remove the cables from the rear connector panel.

9.6.2 Disassemble the Display Unit

Place the Display Unit face down on an anti-static pad.

- 1. Loosen (L) the four captive screws at each corner of the rear enclosure.
- 2. Lift the rear enclosure slightly and pivot it away from the lower enclosure at the bottom side of the Display Unit.





At this point, you can replace the following items (The item numbers refer to the parts list in Section 10.8):

- the internal Compact Flash card (27)
- the external Compact Flash Kit (28)
- the fan (24) for the HPDU this is a 12-volt fan
- the connector panel assembly (21)
- the encoder assembly (12)
- the battery (5)
- the **speaker** (20) (To access the mounting screws for the speaker, you must first remove the ten screws that hold the mounting plate to the front enclosure so that you can raise the bottom edge of the assembly slightly Refer to section 9.6.5.)
- the rear enclosure (1) You can transfer the captive screws to the new enclosure.
   However, the gasket (2) is held in place with adhesive. When replacing the rear enclosure, also include a new gasket.

To replace the remaining items requires further disassembly.

#### 9.6.3 CPU Fan


**CAUTION**: Do not remove the heatsink from the CPU board.

- Note the orientation of the fan harness.
   Disconnect the fan harness from the CPU board.
- 2. Remove the fan (29) from the CPU heatsink, leaving the heatsink in place.
- 3. Remove the heatsink (discard) from the replacement fan assembly.
- 4. Noting the orientation of the fan harness, secure the fan to the CPU heatsink.
- 5. Connect the fan harness to the CPU board.

9-18 09/07 1009-0357-000

### 9.6.4 To replace the CPU board


- 1. Remove the **connector panel assembly** (21) two screws.
- 2. Disconnect the following cables:
  - Inverter harnesses (A)
  - Right membrane switch flex-cable at ZIF (zero insertion force) connector (B)
  - Speaker cable (C)
  - Encoder assembly cable (D)
  - Lower membrane switch flex-cable at ZIF (zero insertion force) connector (E)
  - Fan cable (F)
  - LCD cable (G)
  - Left membrane switch flex-cable at ZIF (zero insertion force) connector (H)




- 3. Remove the remaining four screws (I) that hold the CPU board to the mounting plate.
- 4. Remove the CPU board from the mounting plate.
- 5. Transfer the battery (5) to the new CPU board.
- 6. Reassemble in reverse order.
- 7. Download the latest software (Section 8a.6).
- 8. Reconfigure the Machine Configurations (Install/Setups).
- 9. Affix a new Key Code and Board ID label to the Vent casting.

### 9.6.5 To replace the LCD display

- 1. Disconnect the following cables:
  - Right membrane switch flex-cable at ZIF (zero insertion force) connector (A)
  - Encoder assembly cable (B)
  - Lower membrane switch flex-cable at ZIF (zero insertion force) connector (C)
  - Fan cable (**D**)
  - LCD cable (E)
  - Left membrane switch flex-cable at ZIF (zero insertion force) connector (F)
- 2. Remove the ten screws (circled) that hold the mounting plate to the front enclosure.




- 3. Remove the mounting plate assembly from the front enclosure.
- 4. Disconnect the backlight harnesses (**G**) from the inverter boards.
- 5. Slide the grommet (**H**) out of the mounting plate slot (transfer to new LCD).
- 6. Remove the four screws (circled) that hold the LCD to the mounting plate.



9-20 09/07 1009-0357-000

- 7. Lift the left side of the LCD display slightly away from the mounting plate to pull some of the display ribbon cable (I) to the top side of the plate. Flip the LCD over to the left of the assembly.
- 8. Disconnect the display ribbon cable (J).



9. Reassemble in reverse order.


**Note**: When replacing the LCD, pull the excess ribbon cable to the bottom side of the plate as you lower the LCD on to the plate. For the backlight harness grommet (**H**), ensure that the slit in the grommet faces toward the inside of the keyhole.


### 9.6.6 To replace the backlights

Note: When replacing a backlight or a backlight inverter, you must replace both inverters and the backlight assembly found in the Backlight Kit.

The backlight replacement kit includes a backlight assembly (with two backlights) and two inverters with mounting hardware. To replace the backlight assembly follow the procedure in Section 9.5.4 to gain access to the assembly. To replace the inverters, follow the procedure in the next section.

- 1. Remove the one screw (**K**) that holds the backlight assembly to the LCD.
- 2. Slide the backlight assembly to the left to free it from the retaining tabs and then lift it out of the holder.





- 3. Transfer the grommet (**H**) to the new backlight assembly.
- 4. Reassemble in reverse order.

### 9.6.7 To replace the Inverters


Note: When replacing a backlight or a backlight inverter, you must replace both inverters and the backlight assembly found in the Backlight Kit.

The Display Unit includes two inverters (one for each backlight).


The inverters "sandwich" the mounting plate and use it as a heatsink. Follow the procedure in Section 9.6.5 to gain access to the inverters. Replace one inverter at a time.

- 1. Disconnect the backlight cable (A) from the inverter.
- 2. Remove the two Nylon screws (B) that hold the inverter to the backplate.
- 3. Slide the inverter out of the sleeve and disconnect it from the CPU harness (C).
- 4. Reassemble in reverse order.

#### The "front" inverter



#### The "rear" inverter



9-22 09/07 1009-0357-000


#### 9.6.8 To replace the front enclosure or components

Disassemble the Display Unit following procedures in the previous sections to the point where you have removed the mounting plate assembly from the front enclosure.

If you are replacing the front enclosure, you can transfer the encoder (12) assembly to the new enclosure; but, you must build up the replacement enclosure with:

- a new window (11)
- new membrane switches right-side (14), lower (15), left-side spacer (16)
- new keypads right-side (17), lower (18), left-side blank (19)
- new EMC gasket (10)

If you are replacing a keypad or a membrane switch, you must replace both items.




#### To replace the window


- 1. Place the front enclosure face up on a flat surface.
- 2. Press down on one corner of the window to free it from the enclosure.
- 3. Work your way around the window until you can get a hold of it from the back.
- 4. Slowly pry the window from the enclosure.
- 5. Place the front enclosure face down on a flat surface, taking care not to damage the encoder.
- 6. Remove any remaining residue from the mounting area; clean with isopropyl alcohol.
- 7. Remove the inside protective material from the front of the window.
- 8. Peel the front outside frame of the release liner.
- 9. Lower the window straight down in the enclose, noting the notch in the window and the matching nub on the enclosure.
- 10. Before seating the window, position it in contact with the bottom and right sides of the frame (see arrows) so that the larger gap between the window and the enclosure is at the top and left edges (as viewed from behind).
- 11. Remove the protective film from the back side of the window.

### To replace a membrane switch and keypad

- 1. Remove the screw that attaches the grounding strap to the enclosure.
- 2. Pry the membrane switch and keypad from the enclosure.
- 3. Remove any remaining residue from the mounting area; clean with isopropyl alcohol.
- 4. Remove the backing from the membrane. Be sure to remove the small backing below the flex cable. Be careful not to allow the ribbon cable to adhere to the backing.



5. Insert the flex cable and ground strap through the slot in the enclosure. Ensure that all of the ground strap passes through the slot an does not remain folded over under the membrane.




- 6. Carefully lower the membrane straight down to the enclosure. Seat the membrane in place.
- 7. Remove the backing from the keypad and install it over the membrane switches.
- 8. Attach the ground strap to the enclosure.

### To install the EMI gasket

To fully seal the Display Unit enclosure, you will need approximately 2.3 meters of EMC gasket (10). Cut the gasket into five strips shown below.

Insert a continuous length of gasket in the outside grove of the enclosure (sparingly apply "Super Glue Gel" to the channels near the corners before installing the gasket).

Insert individual lengths of gasket in the inside grove around the window (sparingly apply "Super Glue Gel" to the channels near the corners before installing the gasket).




9-24 09/07 1009-0357-000

#### 9.7 Replacing the Display and MGAS cables

To access the connections at the Display Connector board, remove the small access panel at the rear of the machine (Section 9.3.2).

The top side of the Display Connector board accepts the following cables:

- Airway Module (MGAS) Power Supply board (1).
- System Signal Interface to Display Unit (2).
- System Power Interface to Display Unit (3).




To replace any of these cables, you must remove the oxygen partition that surrounds the MGAS enclosure and MGAS power supply.


## 9.7.1 Remove the MGAS oxygen partition

- 1. Remove the rear panel (Section 9.3.1).
- 2. Remove the MGAS (Airway) module.
- Remove the MGAS module guide: three screws from outside of machine, two screws at MGAS power supply box inside machine (Section 10.22.17).
- 4. Remove the Display Unit (Section 9.5.1).
- 5. From the front of the machine, remove the 11 screws that hold the oxygen partition to the front bezel. Do not remove the two screw (**X**) that hold the front bezel to the vertical support of the vaporizer manifold. When replacing the partition, position it over the two locator pins (**P**) at the top of the bezel opening.
- 6. Remove the foam plug from under the pan.
- 7. Replace the cables as necessary.
- 8. When replacing the partition, ensure that cable are positioned in their respective retaining grooves and that they are long enough to reach the connectors without stressing the cable. Be careful not to pinch the task light harness.
- 9. Reassemble in reverse order.

Front View



**Rear View** 



#### 9.8 Servicing the lower electrical enclosure components



9.8.1 Power Controller board (original)


The lower electrical enclosure includes the following components (Section 10.9:

- the Power Controller board
- the Anesthesia Control board
- the Display Connector board
- the backup batteries and the lower enclosure fan.

To replace these components, remove the large access panel at the rear of the machine (Section 9.3). To replace the Display Connector board, also remove the small access panel.

The replacement Power Controller board includes the mounting plate.

- 1. Disconnect the cables coming from the following components:
  - the Display Connector board (A),
  - the Anesthesia Control board (B),
  - the batteries (C),
  - the fan (**D**),
  - the line filter (E).




- 2. Loosen the four screws (circled) that hold the Power Controller mounting plate.
- 3. Lift the Power Controller assembly slightly to release it from the keyhole slots.
- 4. To replace the Power Controller assembly, reassemble in reverse order.
- 5. Load Software.

9-26 09/07 1009-0357-000

#### 9.8.2 Power Controller board and Universal Power Supply

#### To replace the Power Controller board:

- 1. Disconnect the cables coming from the following components:
  - the Display Connector board (A),
  - the Anesthesia Control board (B),
  - the Power Supply (C).
  - the batteries (D),
  - the fan (E),



- 2. Remove the four lower screws (**circled**) that hold the Power Controller board to the mounting plate.
- 3. Loosen the two upper screws (keyhole) and remove the Power Controller board.
- 4. To replace the Power Controller board, reassemble in reverse order.
- 5. Load Software.

### To replace the power supply (U-Frame) or to provide access to the Anesthesia Control board:

- 6. Disconnect all the cables listed above, except the cable between the Power Supply (**G**) and the Power Controller board (**C**).
- 7. Loosen the four screws (at each corner) that hold the PCB/PS mounting plate.
- 8. Lift the PCB/PS assembly slightly to release it from the keyhole slots.

#### To replace the power supply:

- 1. Disconnect the output cable (**G**) from the power supply.
  - Note position of control connector when reconnecting harness.
- 2. From the back of the mounting plate, remove the four screws (**circled**) that hold the power supply to the mounting plate.
- 3. To replace the Power Supply, reassemble in reverse order.


### 9.8.3 Anesthesia Control board

To replace the Anesthesia Control board, first remove

- the Power Controller board (Section 9.8.1)
- or the PCB/PS assembly (Section 9.8.2).

Then, follow the procedure below:

- 1. Disconnect the cables coming from the following components:
  - the large ribbon cable from the Pan Connector board (A),
  - the harness from the Display Connector board (B),
  - the small ribbon cable from the Display Connector board (C),
  - the harness from the Power Controller board (**D**),
  - the harness from the Pan Connector board (E).



- Loosen the six screws (circled) that hold the Anesthesia Control board to the enclosure.
- 3. Lift the Anesthesia Control board slightly to release it from the keyhole slots.
- 4. To replace the Anesthesia Control board, reassemble in reverse order.
- 5. Load Software.
- 6. Check/Reconfigure the Machine Configurations (Install/Setups).
- 7. Perform User Calibrations (Section 4a.5.1 or Section 4b.5.2).
- 8. Zero Gas Transducers (Section 4a.5.7 or Section 4b.5.8).
- 9. Perform all Ventilator Calibrations (Section 5.4).


9-28 09/07 1009-0357-000

#### **9.8.4 Backup batteries**


#### To remove the batteries

1. Disconnect the battery cable from the Power Controller board (A).

(original)



(PCB/U-FramePS)



- 2. Loosen the two nuts (circled) that hold the battery restrainer to the enclosure.
- 3. To remove the restrainer, slide it forward (toward you) and then to the left.
- 4. Remove the battery pack from the machine.
- 5. Set the batteries upright (terminals up) and slide them out of the tray.
- 6. Remove the harness from the batteries.






#### To replace the batteries

- 1. Gently bend the tabs up, just enough for clearance to install the harness.
- 2. Place the batteries side by side and install the harness.
- 3. Slide the tray over the side of the batteries.
- 4. Guide the harness into the slot on the edge of the tray.
- 5. Position the battery tray assembly upright.
- 6. To replace the battery tray assembly, reassemble in reverse order.
- 7. Allow the batteries to charge.
- 8. Recycle old batteries in same packaging according to local requirements.






#### 9.8.5 Fan

The fan is mounted to the side of the enclosure and draws air into the enclosure through the filter on the AC Inlet assembly. For easier access, temporarily remove the backup batteries (Section 9.8.4).


- 1. Disconnect the fan harness from the Power Controller board (**Fan1**).
- 2. Remove the two screws that hold the fan to the enclosure.
- 3. To replace the fan, reassemble in reverse order.
  - Ensure that the fan is oriented with the flow direction arrow pointing toward the inside of the enclosure.



### 9.8.6 Display Connector board

To replace the Display Connector board, first remove the Power Controller board (Section 9.8.1). Then, follow the procedure below:

- 1. Disconnect the cables from the top of the Display Connector board:
  - Airway Module (MGAS) Power Supply board (A).
  - System Signal Interface to Display Unit (B).
  - System Power Interface to Display Unit (C).
- 2. Disconnect the cables from the top of the Display Connector board:
  - Harness from Power Controller board (D).
  - Harness from Anesthesia Control board (MGAS power) connector (E).
  - Ribbon cable from Anesthesia Control board (signal) connector (F).



- 3. Loosen the six screws (**circled**) that hold the Display Connector board to the enclosure.
- 4. Slide the Display Connector board slightly rearward to release it from the keyhole slots.
- 5. To replace the Anesthesia Control board, reassemble in reverse order.

9-30 09/07 1009-0357-000

#### 9.9 Servicing the pan electrical enclosure components



The pan electrical enclosure includes the following components (Section 10.10):


- the Electronic Gas Mixer assembly
- the Ventilator Interface board
- the Filter board
- the Pan Connector board
- the pan enclosure fan

To replace these components, remove the tabletop (Section 9.4) and the pan enclosure cover.

#### 9.9.1 Electronic Gas Mixer assembly

The following procedure describes how to replace the Electronic Gas Mixer assembly.


- 1. Bleed all gas pressure from the machine (Section 9.2).
- 2. Ensure that all cylinder and pipeline pressures read zero before proceeding.
- 3. Disconnect the inlet tubing elbow fittings from the manifold. If the machine does not include  $N_2O$ , transfer the plug from the  $N_2O$  inlet to the replacement assembly.
- 4. Disconnect outlet tubing from the elbow fitting (A).



- 5. Disconnect the ribbon cable from the Pan Connector board (B).
- 6. Disconnect the fan harness (C).
- 7. Remove the two screws (**D**) that hold the manifold to the enclosure.
- 8. Remove the mounting screw (**E**) at the front edge of the main circuit board.
- 9. To replace the Electronic Gas Mixer assembly, reassemble in reverse order.
- 10.Load Software.
- 11. Zero Mixer Pressure Transducers (Section 4a.5.9 or Section 4b.5.10).

### 9.9.2 Ventilator Interface board

- 1. Disconnect the white and black inline tubing fittings from the Inspiratory pressure transducer (A).
- 2. Disconnect the blue and yellow inline tubing fittings from the Expiratory pressure transducer (**B**).
- 3. Disconnect the white inline tube fitting from the Manifold pressure transducer and the black inline tube fitting from the Airway pressure transducer (**C**).



- 4. Disconnect the harness from the Filter board (D).
- 5. Disconnect the harness from the Pan Connector board (E).
- 6. Remove the four mounting screws, one at each corner, that hold the board to the enclosure.
- 7. To replace the Ventilator Interface board, reassemble in reverse order.
  - Ensure that the tubing fittings are connected to like color fittings and that the tubing will not kink when the cover is replaced.
  - Transfer the grommets to the Manifold and Airway pressure tubing (F).
- 8. Load Software.
- 9. Perform User Calibrations (Section 4a.5.1 or Section 4b.5.2).
- 10. Perform all Ventilator Calibrations (Section 5.4).


9-32 09/07 1009-0357-000

#### 9.9.3 Filter board


Before replacing the Filter board, you must apply an EMI gasket around the edges of the board that face the enclosure. Refer to section Section 10.10 for required parts.

You can replace the Filter board without removing the enclosure cover, if you do not have to remove cover for other reasons.

- 1. Cut the gasket pieces to length and apply them to the board as shown:
  - Do not block the mounting holes.
  - Keep the corner gaps to a minimum.



- 2. Disconnect all cables and harnesses from the front-side of the board:
  - Transducer cables from the pipeline and cylinder gas supplies (A)
  - Harness from system switch (B)
  - Harness from SCGO/ACGO (C)
  - Cable from flow sensors (D)
  - Cable from Vent Engine board (E)
  - Harness from ABS switches (F)
  - Harness from task lights (G)
  - Harness from Alt O<sub>2</sub> (H)




- Remove the six screws that hold the Filter board to the enclosure. (If the enclosure cover is removed, disconnect the three harnesses from the back-side of the board before removing the mounting screws.)
- 4. If the cover is in place, pull the Filter board away from the enclosure and disconnect the three harnesses from the back-side of the board.
- 5. To replace the Filter board, reassemble in reverse order.

#### 9.9.4 Pan Connector board

The Pan Connector board is a wiring interface between the lower electrical enclosure components and the pan electrical enclosure components (refer to Figure 11-11).


- 1. Disconnect the wiring from the top-side of the Pan Connector board:
  - Two harness connectors coming from the Filter board.
  - One harness connector coming from the Ventilator Interface board.
  - One ribbon cable connector coming from the Electronic Gas Mixer.
- 2. Remove the four screws that hold the Pan Connector board to the enclosure.
- 3. Lift the Pan Connector board away from the enclosure and disconnect the large ribbon cable and the harness from the under-side of the board.
- 4. To replace the Pan Connector board, reconnect the large ribbon cable and the harness to the under-side of the board.
- 5. Lower the Pan Connector board over the gasketed opening in the pan enclosure.
  - Keep the extra length of ribbon cable in the pan area.
  - Fold the ribbon cable under the enclosure toward the front of the machine.
- 6. Secure the board to the enclosure.
- 7. Reconnect the harnesses and ribbon cable to the top-side of the board.
- 8. Reassemble in reverse order.



### 9.9.5 Pan enclosure fan


The fan is mounted to the back-side of the enclosure and draws air into the enclosure.

- 1. Disconnect the fan harness from the Mixer board (A).
- 2. Remove the two screws that hold the fan to the enclosure.
- 3. To replace the fan, reassemble in reverse order.
  - Ensure that the fan is oriented with the flow direction arrow pointing toward the inside of the enclosure.



9-34 09/07 1009-0357-000

### 9.10 Servicing the Vent Engine




The Vent Engine is found in a housing located below the breathing system bellows assembly.

The Vent Engine includes the following subassemblies.

- Vent Engine Connector board (1)
- Gas Inlet Valve Assembly (2)
- Inlet Filter (3) located under the gas inlet valve
- Inlet Valve Solenoid (4)
- Drive Gas Regulator (5)
- Flow Control Valve (6)
- Reservoir (7)
- Drive Gas Check Valve (8)
- Interface Manifold (9)


To replace any of the Vent Engine components, you must first remove the Vent Engine from the housing (refer to Section 9.10.1).



### 9.10.1 To remove the Vent Engine

- 1. Disconnect pipeline supplies; close cylinder valves; bleed off pressure.
- 2. Remove the ABS breathing system.
- 3. Remove the Exhalation valve.
- 4. Remove the scavenging downtube.
- 5. Loosen the five captive screws (A) that hold the Vent Engine cover to the housing. Raise the cover to access the Vent Engine.






- 6. Disconnect the Vent Engine harness (B).
- 7. Disconnect the white tube-coupler (**c**) inline with tube to manifold pressure transducer on the Ventilator Interface Board.
- 8. If present, disconnect the black tube-coupler (**D**), inline with tube to AGSS flow indicator.
- 9. Disconnect the drive gas hose (E).
- 10. Loosen the three captive screws (F) that hold the engine manifold to the housing.
- 11. Lift the Vent Engine out of the housing.
- 12. To replace the Vent Engine, reassemble in reverse order.

9-36 09/07 1009-0357-000

#### 9.10.2 Replacing Vent Engine components

Refer to Section 6 for Vent Engine components that are to be serviced under regular maintenance. Most of the components on the Vent Engine can be replaced by removing the mounting screws and reusing them to secure the replacement part.



**Gas Inlet Valve (A)** Inspect the two o-rings that seal it to the manifold. Replace as necessary. To replace GIV shuttle valve components, refer to Section 9.10.3.

**Inlet Filter (B)** Install the filter with the smooth side facing up. Inspect the o-ring. Replace as necessary.

**Regulator (C)** Inspect the two o-rings that seal it to the manifold. Replace as necessary. Perform the Drive Gas Regulator calibration in Section 5.3.

Interface Manifold (D) Inspect the two o-rings that seal it to the manifold. Replace as necessary. Lubricate o-rings sparingly with Krytox.

**Drive Gas Check Valve** Inspect the o-ring that seals it to the manifold. Replace as necessary.

(E) Clean the seat on the manifold and the seal on Drive Gas Check Valve with isopropyl alcohol.

**Inspiratory Flow Valve** Note orientation of the flow valve. Inspect the two o-rings that seal it to the manifold.

(F) Replace as necessary.

Perform the Inspiratory Flow Valve calibration in Section 5.4.3.


**Reservoir (G)** Inspect the two o-rings: reservoir to manifold, reservoir to screw head. Replace as necessary.

Inlet Valve Solenoid (H) Inspect seal between solenoid and GIV body.

Replace as necessary (included with solenoid).

### 9.10.3 Replacing GIV components

Lubricate items marked with an asterisk (\*) sparingly with Krytox.



- 1. Remove the retaining ring (A) and the GIV cap (B).
- 2. Use pneumatic pressure to remove the shuttle. Cover the shuttle with a cloth and briefly apply pressure (connect the drive gas hose or use pipeline pressure) through the drive gas inlet.
- 3. Remove the upper o-ring (**C**) and the lower o-rings (**D**).
- 4. Install the lower o-ring (**D**\*).
- 5. Lubricate the shuttle (**E**) at the three areas (\*) shown: the circumference of the shuttle where the upper and lower u-cup seals are placed and the body part of the shuttle that slides along the lower o-ring.
- 6. Install the lower u-cup seal  $(F^*)$  and the upper u-cup seal  $(G^*)$  on the shuttle.
- 7. Press the shuttle assembly into the GIV manifold.
- 8. Install the upper o-ring (C\*).
- 9. Install the cap (B) and the retaining ring (A).
- 10. Reassemble in reverse order.

9-38 09/07 1009-0357-000

#### 9.11 Servicing the pipeline inlet manifold components

The pipeline inlet filter and the inlet check valve can be replaced without removing the pipeline manifold from the machine. To replace the pressure transducer, you have to remove the manifold.

### 9.11.1 Replace pipeline inlet filter


- 1. Remove the pipeline inlet fitting.
- 2. Pull the pipeline inlet filter out of the fitting. The o-ring should come out with the filter.




3. Install the new pipeline inlet filter in the pipeline inlet fitting. The new filter comes with an o-ring.

## 9.11.2 Replace pipeline inlet check valve

- 1. Remove the rear panel (Section 9.3).
- 2. Remove the pipeline inlet fitting.
- 3. The Air and O<sub>2</sub> pipeline manifolds include a drive gas connection at the back of the manifold. Remove the drive gas tube or plug to access the check valve.
- 4. From the back of the pipeline manifold, use a thin tool to push out the check valve. (For an N<sub>2</sub>O manifold, you will have to carefully apply pressure at the outlet of the manifold — with a syringe for example — to gently force the check valve out of the manifold).




5. Push the new check valve into the opening, using the same thin tool. The new check valve includes an o-ring — orient it toward the pipeline inlet. **Note:** Make sure to push the new check valve all the way back into the opening until it bottoms out on the shoulder.




6. Install the pipeline inlet fitting.

### 9.11.3 Replace the inlet manifold

- 1. Remove the rear panel (Section 9.3).
- 2. Disconnect the tubing from the manifold outlet(s).
- 3. Remove the two screw that hold the manifold to the side extrution.





- 4. Transfer the following item to the replacement manifold or install new as required.
  - pipeline check valve (1)
  - inlet filter (2)
  - inlet fitting (3) and o-ring (4)
  - relief valve (5)
- 5. Transfer the pressure transducer to the new supply (Section 9.13).
  - Ensure the o-ring is in place.
  - Install the transducer.
- 6. To reassemble, perform the previous steps in reverse order.
- 7. Perform the checkout procedure (Section 3).

9-40 09/07 1009-0357-000


#### 9.12 Service the cylinder supply modules

#### **⚠ WARNING**

Be careful not to expose internal components to grease or oil (except Krytox or equivalent).

#### 9.12.1 Replace primary regulator module (complete replacement)

- 1. Bleed all gas pressure from the machine (Section 9.2).
- 2. Ensure that all cylinder and pipeline pressures are at zero before proceeding.
- 3. Remove the rear panel (Section 9.3).
- 4. Disconnect the output tube fitting.
- 5. Remove the three mounting screws and lockwashers.
- 6. Remove the elbow fitting from the replacement gas supply.
- 7. Transfer the pressure transducer to the new supply (Section 9.13).
  - Remove any teflon tape remnants from the transducer mounting threads (transducer and module).
  - Apply 1-1/4 turns of new teflon tape around the treads. Verify that the first few threads are free of tape.
  - Install the transducer.
- 8. To reassemble, perform the previous steps in reverse order.
  - Pull on the cylinder output fitting to ensure it is locked in place.
- 9. Check the output of the regulator BEFORE you install the rear panel. Adjust if necessary (Section 5.1).
- 10. Perform the checkout procedure (Section 3).



### 9.12.2 Replace cylinder inlet filter


- 1. Open the cylinder yokes.
- 2. Remove the inlet adapter from the cylinder yoke, using a 4 mm hex wrench.

**Note:** A brass retaining ring keeps the filter inside the inlet adapter.

3. Thread a 6-mm screw (two turns only) into the brass retaining ring and pull it out.

#### **⚠** CAUTION

Be careful not to crush the filter. Do not thread in the screw more than two full turns.



- 4. Remove the filter.
- 5. Install the new filter and brass retaining ring.
- 6. Install the inlet adapter in the cylinder yoke.
- 7. Perform the checkout procedure (Section 3).

## 9.12.3 Replace cylinder check valve

The cylinder check valve is not a replaceable item. If the check valve is defective, you must replace the complete cylinder supply module.

9-42 09/07 1009-0357-000

#### 9.13 Replace gas-supply pressure transducers

The gas-supply pressure transducer includes an integral cable that connects to the Filter board on the pan enclosure. The transducer itself is mounted directly to the supply module. To replace a pressure transducer (pipeline or cylinder) you have to remove the module from the machine.


- 1. To access the Filter board, remove the tabletop (Section 9.4).
- 2. Disconnect the transducer cable from the Filter board.
- 3. Remove the supply module to access transducer.
  - For cylinder supplies, refer to Section 9.12.
  - For pipeline supplies, refer to Section 9.11.
- 4. Remove the transducer from the module.
- 5. Install the new transducer.
  - For pipeline transducers:
    - Be sure that an o-ring is in place.
  - For cylinder transducers:
    - Remove any teflon tape remnants from the module.
    - Apply 1-1/4 turns of teflon tape around the treads of the transducer. Verify that the first few threads are free of tape.
    - Install the transducer.
- 6. To reassemble, perform the previous steps in reverse order.
- 7. Perform the checkout procedure (Section 3).



#### 9.14 Service vaporizer manifold parts


## 9.14.1 Repair manifold port valve

- 1. Set the system switch to Standby.
- 2. Remove the vaporizers from the vaporizer manifold.
- 3. Using a 14-mm wrench, carefully remove the valve nipple (threaded).



4. Disassemble as necessary to replace parts. The following illustration shows the parts.

**Note:** The port valve replacement kit includes the valve cartridge assembly and the seal. The kit does not include o-rings.



- 5. When installing a new valve cartridge assembly into the vaporizer manifold, put a light coat of Krytox on the bottom portion of the cartridge. The bottom portion of the cartridge is defined as the brass surface that is inserted in the lower spring. **Note:** Do not apply Krytox to the valve seal.
- 6. Verify that the parts are free of dust and dirt.
- 7. To reassemble, perform the previous steps in reverse order.
- 8. Complete the port valve checkout procedure described below (Section 9.14.2).

9-44 09/07 1009-0357-000

# 9.14.2 Checkout procedure for manifold port valve

Use the Vaporizer Manifold Valve Test Tool to perform the checkout procedure for the manifold port valve. This tool and test procedure are intended for use only when the valve cartridge assembly is replaced.


#### Note

This replacement and test procedure is a service action and is not part of the maintenance program.


- 1. Set the system switch to Standby.
- 2. After replacing the valve cartridge assembly, remove the vaporizer port o-ring.
- 3. Attach the valve tester to the top of the valve by sliding the bottom of the tester onto the o-ring groove.
- 4. Tighten the tester screw down onto the valve until the screw bottoms out on the top of the valve. The tester o-ring should create a seal with the top of the valve.
- 5. With the Inspiratory Flow Sensor open to atmosphere (do no plug the right-hand port), verify that the system passes the Low P leak check in the System Checkout (Section 3.4.1).
- 6. Remove the valve tester.
- 7. Reassemble the vaporizer port o-ring.
- 8. Conduct a negative low-pressure leak test on the system.
- 9. Perform system "All checks" (Section 3.4).

## 9.14.3 Replace vaporizer manifold check valve

- 1. Set the system switch to Standby.
- 2. Remove the vaporizers from the vaporizer manifold.
- 3. Remove the upper rear panel.
- 4. Disconnect the tubing from the valve block.




- 5. Remove the valve block.
  - To access the left-hand mounting screw (**A**), remove the right (viewed from front) side panel (**B**).
  - The right-side extrusion include an access hole (**C**) for removing the left-hand mounting screw.
  - **Note:** For early production machines that do not have an access hole in the extrusion, you must remove the vaporizer manifold to remove the valve block.



9-46 09/07 1009-0357-000

**Note** The valve body, o-ring, and flapper do not come out with the block. They stay intact at the bottom of the vaporizer manifold.


6. Pull the flapper out of the valve body.



- 7. Using a hex wrench, put the wrench through one of the holes of the valve body and pull down to remove the valve body and o-ring.
- 8. Verify that parts are free of dust and dirt.
- 9. Replace the flapper by inserting the flapper stem and gently pulling the stem until the flapper secures to the valve body.
- 10. Lightly lubricate the o-ring with Krytox.
- 11. Place the lubricated o-ring on the valve body port at the bottom of the manifold.
- 12. Gently install the valve body in the manifold:
  - Check that the o-ring makes a good seal between the manifold and the valve body.
  - Check that the flapper valve makes solid contact with the valve body.
- 13. Install the valve block.
- 14. Reconnect the tubing to the valve block. Pull on the tube to ensure that it is locked in the fitting.
- 15. Install the vaporizer front panel.
- 16. Perform the checkout procedure (Section 3).

# 9.14.4 Replace vaporizer pressure relief valve

- 1. Set the system switch to Standby.
- 2. Remove the vaporizers from the vaporizer manifold.
- 3. Remove the upper rear panel (Section 9.3).
- 4. Using a 13mm open ended wrench, remove the vaporizer pressure relief valve by turning counterclockwise.




- 5. Verify that the parts are free of dust and dirt.
- 6. Install a new vaporizer pressure relief valve.
- 7. To reassemble, perform the previous steps in reverse order.
- 8. Perform the checkout procedure (Section 3).

9-48 09/07 1009-0357-000

## 9.14.5 Replace vaporizer manifold

- 1. Remove the upper rear panel (Section 9.3).
- 2. Remove the Display Unit.
- 3. Remove the right side panel (A).
- 4. From the front of the machine, remove the two screws (**B**) that hold the front bezel to the vertical support.
- 5. From the back of the machine, remove the two screws (**C**) that hold the vaporizer manifold vertical support to the horizontal bracket.
- 6. From the back of the machine, remove the two screws (**D**) that hold the vertical support to the vaporizer manifold.
- 7. Remove the vertical support from the machine.
- 8. While holding the vaporizer manifold, remove the two screws (**E**) at the right-hand extrusion to release the manifold.




- 9. Install the new vaporizer manifold in reverse order. Do not fully tighten the screws until they are all in place.
  - Attach the new manifold to the right-hand extrusion (E).
  - Attach the vertical support to the vaporizer manifold (**D**).
  - Attach the bottom of the vertical support to the horizontal bracket (C).
  - Attach the vertical support to the front bezel (B).
- 10. Tighten the mounting screws in the following order: E, D, C, B.
- 11. Reassemble the machine.
- 12. Perform the checkout procedure (Section 3).

## 4.15 Replace ACGO selector switch

#### Removal

- 1. Remove the tabletop (Section 9.4).
- 2. Clip the tie wraps (1) from the outlet barb fittings at the side of the switch.



- 3. Disconnect the fresh gas (2) and flush (3) tubes at the back of the switch.
- 4. Disconnect the wires from the ACGO mode microswitch (4) at the back of the selector switch.
- 5. Disconnect the wires from the flush pressure switch (5) on top of the selector switch.
- 6. Set the ACGO selector switch to ABS.
- 7. Back out the selector switch mounting screws (6) until the tips are flush with the face of the mounting casting.
- 8. While pushing the selector knob toward the machine and holding it steady, push the valve body toward the knob and rotate it counterclockwise to separate the valve body from the knob assembly.
- 9. Remove the knob assembly and protective shroud from the machine.
- 10. Remove the valve from the silicone output tubes.

#### Replacement

- 1. Remove the knob assembly from the valve body.
- 2. Back out the selector switch mounting screws until the tips are flush with the face of the mounting casting.
- 3. Guide the outlet fittings of the valve body into their respective silicone tubes.
- Hold the selector knob with the indicator mark facing down. Turn the chrome collar to its maximum counterclockwise position (as viewed from the front).





9-50 09/07 1009-0357-000

- Place the shroud over the knob and guide the assembly into the pan opening.
- 6. Ensure that the indicators on the shroud align with label on the pan and the alignment tab mates with the alignment hole in the pan.
- 7. While holding the knob assembly steady against the pan, place the valve assembly over the knob actuator. Using moderate force press the two assemblies together. The knob should rotate to the ACGO position.
- 8. While continuing to force the assemblies together, rotate the knob assembly to the ABS position. The assemblies should snap into place.
- 9. Verify proper alignment of the knob with the setting indicators. Tighten the mounting screws evenly to secure the switch assembly to the pan.
- 10. Secure the outlet tubing with tie wraps.
- 11. Connect the fresh gas and flush gas tubing. Pull on the tubing to ensure that it is locked in the fitting.
- 12. Reconnect the wires to the ACGO mode microswitch at the back of the valve (top two terminals).
- 13. Reconnect the wires to the flush pressure switch at the top of the valve (upper and lower terminals).
- 14. Replace the tabletop.

#### **Test procedure**


- 1. Confirm that flush flow and 5 L/min fresh gas flow are diverted to the ACGO port and the ABS in the respective knob positions.
- 2. Confirm that the ventilator display indicates ACGO mode when the valve is set in the ACGO position.
- 3. Test the function of the flush pressure switch:
  - DU "Vent Diagnostics Status" Section 8a.3.1)
  - HPDU (4.X or 5.X) "Special Functions Vent Diagnostics Vent Status" Section 8b.4.1)
  - HPDU (5.X) PC Service App "Vent Subsystem Vent Status" Section 12.9.1)
- 4. Perform the low-pressure leak test (Section 3.4.1).
- 5. Perform the checkout procedure (Section 3).







## 9.16 Clean or replace ACGO port flapper valve

- 1. Remove the tabletop (Section 9.4).
- 2. Remove the ACGO cap mounting screws.
- 3. Remove the cap.
- Examine the flapper and disk for obstructions or debris. Clean with isopropyl alcohol if necessary; retest.
- 5. If leak persists, replace the flapper.
  - Remove the flapper from the check valve disk.
  - Clean the new flapper with isopropyl alcohol.
  - Apply a drop of isopropyl alcohol to the center hub of the new flapper.
  - Before the alcohol evaporates, align the center hub of the new flapper with the center hole of the check valve disc.
  - While pressing the flapper against the disc, use your fingernail to help pull the hub through the disc from the other side.








- 6. Lubricate the o-ring sparingly with Krytox (do not get Krytox on the flapper).
- 7. Insert the flapper assembly into the ACGO outlet with the flapper up.
- 8. Replace the cap.


9-52 09/07 1009-0357-000

## 9.17 Replace the APL valve

- 1. Remove the ABS breathing system.
- 2. The APL valve is held in place with a spring and a retainer (**A**) that snaps into a recess in the lower body of the APL valve. To release the retainer, place an appropriately sized straight blade screwdriver into the housing cutout (**B**). Twist the screwdriver to release the retainer.



- 3. Place the new APL valve into position with the setting indicator facing to the front of the machine.
- 4. Place the spring into the retainer.
- 5. While holding the APL valve tight to the housing, snap the spring and retainer onto the valve body from below.
- 6. Reinstall the ABS breathing system.
- 7. Perform the checkout procedure (Section 3).



### 9.18 Replace the bag support arm

1. Remove the ABS breathing system from the machine.

Note: To help prevent the bag arm mounting from loosening, current production machines use socket head screws and flat washers to secure the bag arm to the casting instead of Posidriv screws and lockwashers. Refer to the parts section for stock numbers.

- 2. From the underside of the casting, remove the hardware (**A**) that holds the arm in place.
  - If either of the pins (see below) remain in the casting, remove them from the casting.



- 3. Install the new bag support arm assembly.
  - Position the bag arm over mounting pattern of 4 small holes in the support casting. The arm should extend towards the front of the machine. Align the two pins
     (B) extending from the base of the bag arm assembly, with two of the small holes in the casting that are in line with the APL valve.




- Lower the bag arm, pushing the two pins into the holes.
- From the underside of the casting, secure the bag arm with two M3x20 socket head screws and flat washers.
- 4. Test the force required to swing the bag arm from side to side and adjust if necessary.

**Note**: The adjustment nut is initially set so that 5-mm of exposed thread extends from the adjusting nut. With use, the force required to move the arm increases and may require readjustment.

The force is adjusted by turning the lock nut (8-mm socket) which is accessible from underneath the support casting. Turn clockwise to increase the force and counterclockwise to reduce the force.

- Swing the bag arm sideways through the 90 degree arc permitted by its internal stop.
- Adjust to just enough friction to prevent the bag arm from swinging sideways as
  the bag height is being changed. The bag arm height is changed by squeezing the
  lock release lever (C) at the free end of the bag arm and rotating it to the desired
  position.
- 5. Replace the ABS breathing system.




9-54 09/07 1009-0357-000


# 9.18.1 Servicing the bag support arm

Service parts for the bag support arm include the upper and lower assemblies.


To replace either assembly:

- 1. Remove the bag support arm from the machine (Section 9.18).
- To separate the upper assembly from the lower assembly, use a small (2.5-mm) pin punch from the bottom to drive the dowel pin up and out.
- 3. To assemble the bag arm, apply a light coat of Krytox to the area of the upper arm (**A**) that extends into the lower arm (including the dowel pin groove).





- 4. Insert the upper assembly into the lower assembly. Align the surface (**B**) of the upper assembly with the surface (**C**) of the lower assembly.
- 5. Insert the dowel pin into the hole (from the top side as shown). Drive the dowel pin into the bag arm until it is flush with the top surface.




### 9.18.2 Replace friction pad in lower bag arm assembly


- 1. Remove the ABS breathing system from the machine.
- 2. Using an 8-mm socket, remove the nut (**A**), shoulder washer, and spring from the lower assembly.
- 3. Lift the bag support arm off of the swivel post. Remove the old friction pad.
- 4. Wipe any residue and friction particles from the post.
- Insert a new friction pad into the base. Keep approximately
   1 mm of space between the end of the pad and the bottom of the base.

**Note**: Align the friction pad gap with the seam (**B**) in the base. Position the retaining screw so the pin (**C**) at the base is perpendicular to the seam.

- 6. With the bag support arm facing forward, place the base of the arm over the swivel post. Ensure that the slot in the base (**D**) engages the tab (**E**) on the swivel post.
- Replace the spring, shoulder washer and nut. Tighten the nut until 5 mm of thread extends beyond the nut.
- 8. Follow the procedure in Section 9.18 to adjust the force required to swing the bag arm from side to side.









9-56 09/07 1009-0357-000


## 9.18.3 Replace bag port housing

 Remove the bag support arm cover (A) — screw and lockwasher from below.




- 2. Remove nut (**B**) to remove the release lever (**C**).
- 3. Remove the retaining ring (**D**).
- 4. Slide the bag port housing (**E**) off the end of the bag support arm.
- 5. Before installing the new bag port housing, clean and lubricate sparingly with Krytox the exposed metal end (**F**) and the guide slot (**G**) of the bag support arm.






- 6. Slide the new bag port housing onto the bag arm. Secure it with the retaining ring.
- 7. Lubricate sparingly with Krytox the pivot boss (**H**) before replacing the release lever.
- 8. After replacing the release lever, adjust the mounting nut so that a 2-mm gap remains between the lever and housing when the release lever is fully depressed.
- 9. Replace the bag arm cover.



## 9.19 Replace system switch assembly


- 1. Bleed all gas pressure from the machine (Section 9.2).
- 2. Ensure that all cylinder and pipeline pressures read zero before proceeding.
- 3. Remove the tabletop (Section 9.4).
- 4. Remove the gauge panel mounting screws and move the panel forward to access the system switch.



- 5. Disconnect the wires from the electrical switch.
- 6. Back out the system switch mounting screws just enough to allow the knob collar to be released.
- 7. While holding the switch assembly, push in the knob and turn it counterclockwise.
- 8. Pull the knob and collar out from the front and remove the switch assembly.

9-58 09/07 1009-0357-000


- 9. Install the replacement switch assembly:
  - a. Transfer the 8-mm plugs from the old system switch to the new system switch on the pneumatic module (pull on the plug to ensure that it is locked into the module).
  - b. Turn back the system switch mounting screws until their tips recede.
  - c. Orient the switch assembly with the plugged fittings toward the right and the  $\rm O_2$  fittings toward the left.
  - d. Install the switch assembly through the gauge panel.
  - e. Push the knob collar in with the indicator up and turn it clockwise until it locks.
  - f. Tighten the mounting screws. Make sure that the top edge of the switch assembly is parallel to the top edge of the gauge panel.
  - g. Loosen the two outside screws on the electrical module.
  - h. Insert the wires in the electrical module and tighten the screws.
  - i. Pull the wires on the electrical module to ensure that there is a good connection.
  - Transfer the tubing from the old system switch to the new system switch on the pneumatic module (pull on the tubing to ensure that it is locked into the module).



- 10. Test the replacement switch assembly:
  - a. Connect an  $O_2$  supply.
  - b. Connect the power cable to an electrical outlet.
  - c. Set the system switch to On.
  - d. Make sure that the display comes On.
  - e. Select Alt O2 flow.
  - f. Increase the Alt  $O_2$  flow. Make sure that gas flows.
  - g. Make sure that you do not feel or hear any leaks.
  - h. Set the system switch to Standby.
  - i. Make sure all gas flow stops and the display turns Off.
- 11. Reinstall the gauge panel and the tabletop.
- 12. Perform the checkout procedure (Section 3).

## 9.20 Replace Alt 0<sub>2</sub> components

- 1. Bleed all gas pressure from the machine (Section 9.2).
- 2. Ensure that all cylinder and pipeline pressures read zero before proceeding.
- 3. Remove the tabletop (Section 9.4).
- 4. Remove the gauge panel mounting screws and move the panel forward to access the Alt O2 components.



#### Alt 0<sub>2</sub> Flowmeter (A)

Disconnect the tubing from the flowmeter.

Remove the four screws that hold the flowmeter mounting bracket to the front panel. Transfer the mounting bracket to the new flowmeter.

#### **Needle Valve Assembly**

Loosen the set screw that holds the knob to the needle valve; remove knob.

(B) Disconnect the tubing from the needle valve assembly.

Remove the four screws that hold the needle valve assembly to the front panel.

Transfer the mounting plate to the new needle valve assembly.

#### Alt 0<sub>2</sub> Switch (C)


Disconnect the switch harness.

When replacing the switch, face the tab on the washer toward the switch body (tab not used for positioning).

9-60 09/07 1009-0357-000

## 9.21 Replace auxiliary 02 flowmeter


- 1. Bleed all gas pressure from the machine (Section 9.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the tabletop (Section 9.4).
- 4. Remove the adjustment knob from the flowmeter; pull forward.
- 5. Remove the gauge panel mounting screws and move the panel forward to access the flowmeter.




- 6. Disconnect the inlet tube fitting (A).
- 7. Disconnect the tube (B) from the outlet fitting.
- 8. Remove the four screws (C) that hold the flowmeter mounting bracket to the front panel.
- 9. Transfer the mounting bracket to the new flowmeter.
- 10. Reassemble in reverse order.
- 11. Perform the checkout procedure (Section 3).

## 9.22 Replace the suction regulator

- 1. Lower the upper rear panel (Section 9.3).
- 2. Disconnect the white vacuum (**A**) and black suction (**B**) fittings from the rear panel. Do not remove the tubing from the regulator.
- 3. If you are replacing a Venturi Drive regulator, disconnect the tube (**C**) from the control port of the regulator assembly.





- 4. Disconnect the two screws (**D**) that hold the regulator assembly to the mounting bracket.
- 5. Remove the regulator assembly from the front panel.
- 6. Transfer the tubing to the new regulator:
  - Attach the vacuum source tube (white fitting) to the lower connector (E).
  - Attach the suction tube (black fitting) to the upper connector (F).
- 7. Guide the tubes into the front panel opening.
- 8. While holding the regulator assembly against the front panel, attach the retaining bracket to the regulator. Tighten the screws to secure the regulator assembly.
- 9. If applicable, attach the control port tube to the regulator assembly (**C**).
- 10. Attach the vacuum and suction fitting to the rear panel manifold.
- 11. Replace the rear panel.
- 12. Perform the checkout procedure (Section 3).


9-62 09/07 1009-0357-000

## 9.23 Replace task light components

The Avance machine includes two task lights (upper and lower) that are controlled by a common switch. The task light switch and the upper tack light are accessible from the front of the machine. To service the lower task light you have to remove the rear panel.

### 9.23.1 To replace the task light switch

1. Remove the four screws (A) that hold the task-light lens to the upper shelf.



- 2. Using a small needle-nose pliers, disconnect the switch harness from the task-light circuit board connector (**B**).
- 3. Remove the two screws (**C**) that hold the switch retainer plate to the upper shelf.
- 4. Transfer the switch retainer plate to the new switch, counter-sunk side to the outside.
- 5. Mount the switch to the upper shelf.
- 6. Remount the task-light assembly. Ensure that the switch harness and the task-light harness wires are positioned in their respective recesses and are not pinched under the task-light lens.


### 9.23.2 To replace the upper task light

- 1. Using a small needle-nose pliers,
  - disconnect the switch harness from the task light circuit board connector (B).
  - disconnect the task-light harness from the task light circuit board connector (D).
- 2. Slide the task-light circuit board out of the lens.
- 3. Slide the new task light into the lens, ensuring that the connectors are aligned with the lens cutouts.
- 4. Plug the task-light harness and the switch harness into their respective connectors on the task-light circuit board. Use a small screwdriver to push the connectors securely into place.
- Remount the task-light assembly. Ensure that the switch harness and the task-light harness wires are positioned in their respective recesses and are not pinched under the task-light lens.

### 9.23.3 To replace the lower task light

1. Remove the rear panel (Section 9.3.1).

The lower task light is located directly below the MGAS partition. It slides into two slots in the front bezel (**A**) and is retained with a small moveable bracket (**B**).



- 2. Rotate the retaining bracket counterclockwise to free the tack light assembly.
- 3. Disconnect the harness from the task light assembly and connect it to the new task light.
- 4. Place the task light into the slots (lights pointing down).

Note

If the machine includes an integrated suction regulator, it may be too cramped for you to place the task light. Use long-nose pliers or a similar tool to guide the task light in place, being careful not to damage the lights. Or you can temporarily remove the suction regulator to gain more room.


- 5. Rotate the bracket in place to retain the task light.
- 6. Reassemble in reverse order.

9-64 09/07 1009-0357-000


## 9.24 Replace ABS breathing system components

### 9.24.1 Replace Bag/Vent switch assembly

- 1. Remove the ABS breathing system.
- 2. From the underside, remove the bellows base manifold (**A**) and fully loosen the two captive screws (**B**) at the bag port side of the APL/BTV manifold.



3. From the topside, rotate the Bag/Vent switch cartridge counterclockwise until the Bag/Vent switch outlet port (**C**) clears the bellows housing.




- 4. Lift out the Bag/Vent switch cartridge from the housing.
- 5. Replace the Bag/Vent switch cartridge in reverse order.
- 6. Reinstall the ABS breathing system.
- 7. Perform the checkout procedure (Section 3).


# 9.24.2 Replace bellows base latch assembly

To replace the latch assembly, you must disassemble the bellows base assembly to the point where you can remove the guide (**A**) and latch assembly (**B**) as a unit.

- 1. Remove the Bag/Vent switch cartridge (Section 9.24.1).
- 2. Remove the two remaining screws (**C**) that hold the APL/BTV manifold to the bellows base assembly. Remove the APL/BTV manifol.



- To remove the guide/latch assembly, remove two mounting screws (D) from the
  underside. Remove two additional mounting screws from the topside. Remove the
  guide/latch assembly from the bellows base assembly.
- 4. Separate the latch assembly from the guide assembly.
- 5. To install the new latch assembly, put the spring (**E**) into place in the guide assembly (long leg down).
- 6. Place the latch assembly on the guide assembly so that the latch engages the short leg of the spring. Secure the latch assembly (**F**) to the guide assembly.





- 7. Mount the guide/latch assembly into the bellows base assembly.
  - Extend the latch (**G**) while placing the assembly into the base.
- 8. Reassemble the breathing system in reverse order.
- 9. Perform the checkout procedure (Section 3).

9-66 09/07 1009-0357-000

#### 9.24.3 EZchange Canister spring replacement

- 1. Detach the EZchange module from the machine.
- 2. Remove the two M3 screws (A) that hold the module cover (B); set the cover aside.





- 3. Remove the two M3 shoulder screws (**C**) that fasten the canister latch lever (**D**).
- 4. Remove the latch lever, the switch actuator lever (**E**) and the spring; discard the spring.
- 5. Place the new spring on the module (as shown below). Position the switch actuator lever over the spring. Ensure the spring hooks are fully engaged into the posts on the manifold and the actuating lever.





- 6. Clean any residual Loctite debris from the M3 shoulder screws removed in Step 3.
- 7. Place the canister latch lever in position. Apply Loctite 242 to the threads of the two M3 shoulder screw threads and secure the canister latch level.
- 8. Check the switch actuator lever to ensure free movement. If sticking is observed, loosen the M3 shoulder screw approximately 1/8 of a turn until free movement of the switch actuator lever is observed.
- 9. Install the module cover.
- 10. Install the EZchange module.
- 11. Verify that the following message appears on the screen when the absorber canister is released.
  - 'No CO2 absorption' for Aespire machines
  - 'CO2 Absorber Out of Circuit' for Avance and Aisys machines
- 12. Perform the Preoperative Checkout Procedure (refer to the User's Reference manual).

## 9.25 Replace casters

#### **⚠ WARNING**

Replacing a caster requires at least two people to maneuver and tip the machine. Personal injury and/or machine damage is possible if one person attempts this procedure alone.

1. Disconnect all pipeline hoses from the wall and the machine, close all gas cylinders, unplug the power cord, and set the system switch to standby.

#### **⚠** CAUTION

Remove the vaporizers before tipping the machine. If a vaporizer is inverted, it must be set to 5% and purged for 30 minutes with a 5 L/min flow. The interlock system prevents purging more than one vaporizer at a time.

2. Remove the absorber, the vaporizers, gas cylinders, drawers and all auxiliary equipment.

#### **⚠** CAUTION

To prevent damage, do not tip the Avance machine more than 10 degrees from vertical.

3. Block the opposite wheels; then, block up the machine until there is enough room to remove the defective caster.

To block up the machine, tip and slide blocks under the caster base. Raise both sides evenly until the unit is high enough to remove the caster.

- The casters are threaded into the base and held with a Loctite compound.
   Remove the caster with an appropriately sized open-end wrench.
- 5. If required, clean the threads of the new caster with denatured alcohol.
- Apply Loctite 242 to the threads of the new caster. Install the caster securely into place.
- 7. Make sure the caster turns freely.
- 8. Carefully lower the machine to the floor.
- 9. Perform the checkout procedure (Section 3).



9-68 09/07 1009-0357-000

## 9.26 Reconfigure sample gas return line


**Note** In the U.S., it is not permitted to return sample gas to the breathing circuit.

Outside the U.S., consult with hospital guidelines on filtering of re-breathed gas before allowing the sample gas return line to be reconfigured.

#### **Procedure**

Sample gas return is directed to the scavenging system as a factory default. Perform the following to reroute the sample gas back to the breathing system.

- 1. Remove the tabletop.
- 2. Port 4 (**A**) of the ABS breathing system is connected to the expiratory circuit, downstream of the expiratory check valve. As a factory default, Port 4 is plumbed with a length of tubing that is plugged (**B**) at the far end.
- 3. Remove the plug from the tube.
- 4. Find the sample return line at the left-rear corner of the pan assembly. The sample return line includes an inline connector (**C**) at the point where the sample line goes down into the vent engine housing.
- 5. Separate the scavenging tube, removing the inline connector from the portion of the tube that extends into the vent engine housing. Plug the open end of the scavenging tube (**D**) with the plug removed above.
- 6. Insert the inline connector from the sample return port into the open end of the return tube (**E**) to Port 4. Pull on the connector to ensure that it is securely connected.
- 7. Replace the tabletop.
- 8. Perform the checkout procedure (Section 3).



### 9.27 Change drive gas

#### **⚠** CAUTION

If you change the drive gas, you must also change the drive gas selection on the ventilator service setup screen. Refer to Section 4 of the ventilator Technical Reference manual.

 If the drive gas selection and the actual drive gas do not agree, volumes will not be correct.

The ventilator will alarm with the message "Low Drive Gas Press" if the selected drive gas pressure, either  $O_2$  or Air, is lost.

1. Remove the rear panel (Section 9.3).

Note:

The  ${\rm O}_2$  and Air pipeline manifolds have a drive gas connection at the back. The connection not in use is plugged.

- 2. Remove the plug from the new connection.
- 3. Disconnect the drive gas hose from the present connection.
- 4. Install the plug in this connection (pull on the plug to ensure that it is locked into the fitting).
- 5. Reroute the drive gas hose so that it does not cause kinks in other tubing.
- 6. Connect the drive gas hose to the new connection (pull on the hose connector to ensure that it is locked into the fitting).
- 7. Do a high-pressure leak test (Section 3.8).
- 8. Enter the service mode and select the correct drive gas.
- 9. Test the primary regulator. Verify that it functions within specifications now that it will be supplying drive gas to the ventilator (Section 5.1).
- 10. Perform the checkout procedure (Section 3).

9-70 09/07 1009-0357-000

## **10 Illustrated Parts**

In this section	10.1 Service tools	10-3
	10.1.1 Software tools	10-3
	10.1.2 Manifold pressure test adapter	10-3
	10.1.3 Test Devices	
	10.1.4 Lubricants and Adhesives	
	10.1.5 Test Tools	
	10.2 External components - front view	
	10.3 External components - front view references	10-7
	10.4 External Components - rear view	10-8
	10.5 AC Power cords and AC Inlet filter	10-9
	10.6 AC Inlet/Outlet Components	10-10
	10.7 Display Unit (DU)	10-12
	10.8 High Performance Display Unit (HPDU)	10-14
	10.9 Lower electronic enclosure components	10-16
	10.9.1 Anesthesia Control and Display Connector board	10-17
	10.10 Pan electronic enclosure components	10-18
	10.11 Electronic Gas Mixer	10-19
	10.12 Pipeline inlet fittings	10-20
	10.13 Cylinder Gas Supplies	10-21
	10.13.1 Cylinder inlet fittings	10-22
	10.14 Vaporizer manifold	10-23
	10.15 Vent Engine Housing	10-24
	10.16 Vent Engine	10-25
	10.16.1 Vent Engine - under side	10-26
	10.17 ABS to machine Interface Components (SCGO)	10-27
	10.18 ABS to machine Interface Components (ACGO)	10-28
	10.19 Flush Regulator and Flush Valve	10-29
	10.20 Front panel, Alt O2, and system switch	10-30
	10.21 Breathing system interface	10-31

10.22.1 APL Valve       10-32         10.22.2 Bag/Vent Switch       10-33         10.22.3 Absorber canister       10-34         10.22.4 Flow Sensor Module       10-35         10.22.5 Breathing Circuit Module       10-36         10.22.6 Exhalation valve       10-37         10.22.7 Bellows       10-38         10.22.8 Bellow base       10-39         10.22.9 Bag Arms       10-40         10.22.10 EZchange Canister system (CO2 Bypass)       10-41
10.22.3 Absorber canister       10-34         10.22.4 Flow Sensor Module       10-35         10.22.5 Breathing Circuit Module       10-36         10.22.6 Exhalation valve       10-37         10.22.7 Bellows       10-38         10.22.8 Bellow base       10-39         10.22.9 Bag Arms       10-40
10.22.4 Flow Sensor Module       10-35         10.22.5 Breathing Circuit Module       10-36         10.22.6 Exhalation valve       10-37         10.22.7 Bellows       10-38         10.22.8 Bellow base       10-39         10.22.9 Bag Arms       10-40
10.22.5 Breathing Circuit Module       10-36         10.22.6 Exhalation valve       10-37         10.22.7 Bellows       10-38         10.22.8 Bellow base       10-39         10.22.9 Bag Arms       10-40
10.22.6 Exhalation valve       10-37         10.22.7 Bellows       10-38         10.22.8 Bellow base       10-39         10.22.9 Bag Arms       10-40
10.22.7 Bellows       10-38         10.22.8 Bellow base.       10-39         10.22.9 Bag Arms       10-40
10.22.8 Bellow base.       10-39         10.22.9 Bag Arms       10-40
10.22.9 Bag Arms
10.22.10 EZchange Canister system (CO <sub>2</sub> Bypass)
<u>-</u>
10.22.11 Condenser
10.22.12 Anesthetic Gas Scavenging System — AGSS
10.22.13 Passive AGSS
10.22.14 Adjustable AGSS
10.22.15 Active AGSS
10.22.16 AGSS gauge, and sample return
10.22.17 Airway module (MGAS) components
10.23 Integrated Suction Regulator
10.23.1 Major Components (Continuous and Venturi suction)
10.23.2 Suction Control Module
10.23.3 Venturi assembly
10.24 Auxiliary O <sub>2</sub> Flowmeter
10.25 Rear panel components         10-56
10.26 Tabletop components
10.27 Right-side Components
10.28 External components - lower assembly
10.29 Drawer
10.30 Legris quick-release fittings
10.31 Vent Drive and low-pressure tubing
10.32 Tubing for use with Legris fittings
10.33 Cables and harnesses (power supply)
10.34 Cables and harnesses in lower electronic enclosure
10.34.1 Machines with original Power Controller board
10.34.2 Machines with Universal Power Supply (U-Frame) 10-70
10.35 Cables and harnesses (Filter Board interface)
10.36 Cables and harnesses in Pan enclosure
10.37 Optional Monitor Display mounts
10.38 Display arm mounting kits for optional equipment

10-2 09/07 1009-0357-000

### 10.1 Service tools


## 10.1.1 Software tools

Item	Description	Version	Stock Number
1	Avance System Software (for HPDU)	5.X 4.X	M1087631-S M1062100-S
1	Avance System Software (for DU)	3.20 3.00 2.02	M1055796 1009-6163-000 1009-6135-000
2	Compact Flash adapter, PCMCIA carrier (required only for DU)		1009-5874-000
3	Service Application, PC based (for System Software 5.X)	1.X	M1088110
3	Cable, DU serial port to PC serial port		1011-3984-000

# 10.1.2 Manifold pressure test adapter

The manifold pressure test adapter is used to tee into the manifold pressure line for the Manifold P Span calibration (Section 5.4.2).


Assemble the adapter using the parts shown.



Item	Description	Stock Number
1	Coupler, male - white	1503-3236-000
2	Tee (male barb)	1009-3011-000
3	Coupler, female - white	1503-3119-000
4	Tubing (low-pressure) 1/4 inch	1605-1001-000

## **10.1.3 Test Devices**

Item	Tool	Stock Number
1	Test flowmeter, 6–50 L/min (Suction Flow Test)	1006-8431-000



#### **Not Shown**

Low-pressure Leak Test Device	(negative pressure)	0309-1319-800
Low-pressure Leak Test Device	(positive pressure - ISO)	1001-8976-000
Low-pressure Leak Test Device	(positive pressure - BSI)	1001-8975-000
Flow test device capable of measuring 0 with an accuracy of $\pm 2\%$ of reading	Refer to section 6.7	
Vacuum test gauge capable of measuring 0 to 550 mm Hg with an accuracy of $\pm 1\%$ of reading		Refer to section 6.9
Test device capable of measuring 0-30	Refer to section 6.9	
Leakage current test device	Refer to section 3.16	
Test device capable of measuring 689 kl	Refer to section 5.1.1	


### **10.1.4 Lubricants and Adhesives**

Item	Description	Stock Number
1	Lubricant, Krytox GPL 205, 2 oz	1001-3854-000
2	Lubricant, Dow 111, 5.3 oz	6700-0074-200
3	Thread Lock, Loctite No 24221, 10 ml	0220-5017-300
4	"Super Glue Gel", Loctite 454	6812-2160-010


10-4 09/07 1009-0357-000

## **10.1.5 Test Tools**

Item	Tool	Stock Number
1	Leak Test Tool Kit, ABS breathing system	1407-7013-000
1a	Test Tool, bulkhead	1407-8500-000
1b	Plug, tapered 27x12 mm	1407-8505-000
1c	Plug, tapered 24x18 mm	1407-8506-000
1d	Test Tool, circle module (2 each)	1407-8502-000
1e	Plug, service B/S 11 mm (2 each)	1407-8504-000
<b>1</b> f	Plug, service BTV 18 mm (2 each)	1407-8503-000
2	Adapter, positive low-pressure leak test	1009-3119-000
3	PEEP/INSP Calibration Flow Orifice	1504-3016-000
4	Airway Pressure Sensing Tee	1504-3011-000
5	Vaporizer Manifold Valve Test Tool	1006-3967-000
6	Plug, stopper	2900-0001-000
Not Sho	wn	
	Tool to help disconnect tubing from Legris fittings	2900-0000-000
	Test Lung	0219-7210-300
	Leak detection fluid, Snoop	obtain locally




## 10.2 External components - front view



Item	Description	Stock Number
1	Caster, 125-mm with brake (front)	1006-3070-000
2	Caster, 125-mm no brake (rear)	1006-3071-000
3	Cover, cable channel	1009-3020-000
4	Upper shelf	1009-3022-000
5	Bolt, M6x40	0144-2131-911
6	Lockwasher, M6 internal	0144-1118-130
7	Label, S/5 Avance	1009-3206-000
8	Handle, side	1009-3033-000
9	Handle, Medirail	1009-3101-000
10	Screw, M6x12 Sems	0144-2436-106
11	Spacer	1009-3102-000
12	Screw. M6x70	0144-2131-923
13	Lockwasher M6 external	9213-0560-003
14	Screw, M6x20	0144-2131-921
15	Shim	1009-3131-000
16	Task Light PCB, Upper	1009-5855-000
17	Lens, Task Light	1011-3308-000
	Screw	0142-4254-106
18	Switch Assembly, task light	1009-5587-000
19	Plate, switch mounting retainer	1009-3143-000
	Screw	0140-6226-107
20	Task Light PCB, Lower	1009-5857-000


10-6 09/07 1009-0357-000

## **10.3 External components - front view references**



Item	Description	Section number
1	Airway module (MGAS) components	Refer to section 10.22.17
2	"AGSS gauge, and sample return"	Refer to section 10.22.16
3	Vent Engine Housing	Refer to section 10.15
4	Anesthetic Gas Scavenging System — AGSS	Refer to section 10.22.12
5	Breathing System	Refer to section 10.22
6	Drawer	Refer to section 10.29
7	Tabletop components	Refer to section 10.26
8	Display Unit (DU)	Refer to section 10.7
9	External components - lower assembly	Refer to section 10.28
10	Auxiliary O <sub>2</sub> Flowmeter	Refer to section 10.24
11	Integrated Suction Regulator	Refer to section 10.23
12	Front panel, Alt 02, and system switch	Refer to section 10.20
13	ABS to machine Interface Components (SCGO)	Refer to section 10.17
14	Vaporizer manifold	Refer to section 10.14
15	Right-side Components	Refer to section 10.27
16	Optional Monitor Display mounts	Refer to section 10.37


## **10.4 External Components - rear view**



Item	Description	Stock Number
1	AC Inlet	Refer to section 10.5
2	Pipeline Inlets Label, pipeline inlet blank	Refer to section 10.12 1009-3197-000
3	Cylinder Gas Supplies	Refer to section 10.13
4	Electrical Power Outlet	Refer to section 10.6
5	Suction items	Refer to section 10.23
6	Rear panel items	Refer to section 10.25

10-8 09/07 1009-0357-000

## 10.5 AC Power cords and AC Inlet filter




Item	Description	Stock Number
1	Power Cord	
-	100-120V~ 50-60Hz, NEMA, Japan and US	1006-3907-000
	220-240V~ 50-60Hz, AS 3112, Australia	1006-3888-000
	220-240V~ 50-60Hz, GB2099, China	M1053942
	220-240V~ 50-60Hz, BS1363, UK	1006-3884-000
	220-240V~ 50-60Hz, BS546, India and South Africa	1006-3885-000
	220-240V~ 50-60Hz, CEE 7/7, EURO and France	1001-3380-000
	220-240V~ 50-60Hz, Danish	1011-3696-000
	220-240V~ 50-60Hz, SEV 1011, Swiss	1006-3889-000
	220-240V~ 50-60Hz, NEMA, Peruvian	1006-3882-000
2	Clamp, power cord retainer (all except China) (China)	1009-3103-000 M1054192
3	Screw, M4x8 Pozidriv DIN84	1006-3178-000
4	Filter, foam	1009-3064-000
5	Retainer, filter	1009-3058-000

## 10.6 AC Inlet/Outlet Components

Item	Description	Stock Number
1	Inlet, 100/120 AC, with line filter and 15 A circuit breaker Inlet, 220/240 AC, with line filter and 8 A circuit breaker	1009-5698-000 1009-5757-000
2	Fuse, 5A - 5x20mm Fuse holder	1202-3345-000 1009-5674-000
3	Circuit board, Inrush, 100-120V Circuit board, Inrush, 220-240V	1006-3245-000 1006-3246-000
4	Filter, AC Line, 6VW1, 100-240V	1009-5690-000
5	Stud, Equal Potential, 6 mm	0208-0070-300
6	Harness, 100/120 V to Toroid Harness, 220/240 V to Toroid	1009-5752-000 1009-5753-000
7	Harness, to 100/120 V outlets Harness, to 220/240 V outlets	1009-5716-000 1009-5717-000
8	Outlet Receptacle, Australia and China, AS 3112	1001-3305-000
	Outlet Receptacle, Danish, AFSNIT 107-2-D1	1011-3910-000
	Outlet Receptacle, EURO, CEE 7/7	1202-3551-000
	Outlet Receptacle, France, CEE 7/4 Support Frame, snap in	1006-4421-000 1006-4422-000
	Outlet Receptacle, India and South Africa, BS 546	1006-3805-000
	Outlet Receptacle, Japanese	1006-3578-000
	Outlet Receptacle, NA, Nema 5-15	1006-3555-000
	Outlet Receptacle, Swiss, SEV 1011	1006-3807-000
	Outlet Receptacle, UK, BS1363	1001-3309-000
9	Circuit Breaker, 1A, Rocker	1009-5722-000
	Circuit Breaker, 2A Rocker	1009-5721-000
	Circuit Breaker, 3A Rocker	1009-5720-000
	Circuit Breaker, 4A Rocker	1009-5719-000
10 10a	Toroid, 100-240V Toroid, 100-240V - used with no outlets	1009-5692-000 1009-5758-000
11 11a	Screw, M6x70 Screw, M6x60	0144-2131-923 0144-2131-914
12	Lockwasher, M6	9213-0560-003
13	Washer	0402-1107-500
14	Cover, transformer	1009-3063-000
15	Screw, M4x8 DIN84 (for transformer cover)	1006-3178-000

10-10 09/07 1009-0357-000





## 10.7 Display Unit (DU)

Item	Description	Stock Number		
1	Enclosure, rear	1009-5673-000		
2	CPU Board, display unit (with PCMCIA frame)	1009-8289-000		
3	Frame, PCMCIA	1009-5761-000		
4	Gasket, knife edge (2 each)	1009-5804-000		
5	Battery, Lithium 3V (positive side up)	1009-5800-000	Stock Number	
6 ( <b>see</b>	Display, LCD 12-inch color (original DG31 is obsolete)		Backlight Kit	
note)	Display (DG41), LCD 12-inch color (includes backlights)	1009-5938-000	for DG 31 display	
7	Backlight Kit for DG41 display	1009-8243-000	1009-8243-000	
	(backlight assembly, 2 inverters, and hardware)			
7a	Harness, inverter	1009-5527-000		
7b	Spacer, 8mm Nylon	1009-5695-000		
7c	Insulator, tube	1009-3149-000		
8	Grommet, diagonal cut (backlight cable)	1009-3151-000		
9	Enclosure, front	1009-5672-000		
10	Gasket, EMC 1.8mm OD hollow RND (2.3 m per enclosure)	1009-5802-000		
11	Window	1009-5676-000		
12	Encoder assembly	1503-3012-000		
13	Knob, ComWheel	898794		
14*	Membrane switches, right	1009-5505-000		
15*	Membrane switches, lower	1009-5507-000		
16	Keypad, right-side (part of keypad set)	Refer to Table 1		
17	Keypad, lower (part of keypad set)	Refer to Table 1		
18	Keypad, blank (part of keypad set)	Refer to Table 1		
19*	Spacer, blank keypad	1009-5870-000		
20	Speaker assembly, 8-ohm	1605-3263-000		
21	Rear Connector Panel Assembly (with interface boards)	1009-8244-000		
22	Cable, ribbon CPU to Display	1009-5520-000		
23	Grommet	1009-3152-000		
24	Fan, 5Vdc	1504-3516-000		
25	Capsule, fan filter	896089		
26	Filter, fan	897010		
27	Door, PCMCIA	1009-5679-000		
28	Gasket, cover plate	1009-5678-000		
29**	Screw, M3x6 Sems	0140-6219-128		
30**	Screw, M4x8 Sems	0140-6226-113		
31**	Screw, M4x12 relieved body	1504-3001-000		
32**	Lockwasher, M4 external	9213-0540-003		
33**	Screw, M3x16	1504-3003-000		
34**	Lockwasher, M3 external	9213-0530-003		
35**	Screw, M3x6 Nylon	9211-1730-065		
36**	Screw, M2x16	0140-6216-100		
37 Screw, M4x12 Pan Washer HD 1009-3341-000  * When replacing a backlight or a backlight inverter, you must replace both inverters and the backlight assembly found in				

<sup>\*</sup> When replacing a backlight or a backlight inverter, you must replace both inverters and the backlight assembly found in the Backlight Kit (Item 7).

Note: Refer to Section 9.5.5 regarding backlights used in the Display Unit.

10-12 09/07 1009-0357-000

<sup>\*\*</sup> Keypads will likely be damaged during membrane switch replacement; order parts accordingly.

<sup>\*\*\*</sup> Refer to Table 2 for where used.




Table 1:		
Language	Keypad Set	
Chinese	1009-5932-000	
Czech	1009-5923-000	
Danish	1009-5931-000	
Dutch	1009-5918-000	
English	1009-5915-000	
Finnish	1009-5922-000	
French	1009-5916-000	
German	1009-5917-000	
Greek	1009-5927-000	
Hungarian	1009-5928-000	
Italian	1009-5919-000	
Japanese	1009-5929-000	
Norwegian	1009-5925-000	
Polish	1009-5924-000	
Portuguese	1009-5921-000	
Russian	1009-5930-000	
Spanish	1009-5920-000	
Swedish	1009-6140-000	
Turkish	1009-5926-000	

Tabl	e 2:
Hardware Item where used (Qty)	
Spe	aker: 29(2)
Doo	r: 29(2)
CPU	to plate: 30(4)
Rea	r connector panel assembly: 30(2)
Grou	und straps for keypads: 30(2)
Mou	inting plate to Front enclosure: 30(10)
Rea	r enclosure: 31(4), 32(4)
Fan:	33(4), 34(4)
Inve	rters: 35(4)
PCN	ICIA frame: 36(4)



### **10.8 High Performance Display Unit (HPDU)**

	Item	Description	Stock Number
		HPDU Assembly, complete - without keypads	1011-8084-000-S
	1	Enclosure, rear (HPDU)	1009-6089-000
	2	Gasket, cover plate	1009-5678-000
	3	Door	1009-5679-000
	4	CPU Board, HPDU	1009-5933-000
	5	Battery, Lithium 3V (positive side up)	1009-5800-000
	6	Display (DG41), LCD 12-inch color (includes backlights)	1009-5938-000
	7*	Backlight Kit for DG41 display	1009-8422-000
		(backlight assembly, 2 inverters, 2 insulators, and hardware)	
	7a	Harness, inverter	1009-5527-000
	7b	Spacer, 8mm Nylon	1009-5695-000
	7c	Insulator, tube	1009-3149-000
	8	Grommet, diagonal cut (backlight cable)	1009-3151-000
	9	Enclosure, front	1009-5672-000
	10	Gasket, EMC 1.8mm OD hollow RND (2.3 m per enclosure)	1009-5802-000
	11	Window	1009-5676-000
	12	Encoder assembly	1503-3012-000
	13	Knob, ComWheel	898794
	14	Membrane switches, right	1009-5505-000
	15	Membrane switches, lower	1009-5507-000
	16	Membrane switches, left	1009-5506-000
	17**	Keypad, right-side (part of keypad set)	Refer to Table 1
	18**	Keypad, lower (part of keypad set)	Refer to Table 1
	19**	Spacer, blank keypad	1009-5870-000
	20	Speaker assembly, 8-ohm	1605-3263-000
	21	Connector Panel Assembly (with interface boards) - HPDU	M1056239
	22	Cable, ribbon CPU to Display	1009-5520-000
	23	Grommet	1009-3152-000
	24	Fan, 12Vdc (HPDU)	1009-6096-000
	25	Capsule, fan filter	896089
	26	Filter, fan	897010
	27	Compact Flash Card, formatted	M1055008
	28	External Compact Flash Interface Replacement Kit	M1055008 M1056240
	20 29***	Fan, CPU - with heatsink 12V	1009-6095-000
Where Used (Not Shown)	23	i an, or o - with heatshir 12 v	1009-0090-000
Speaker (2) Door (2) –	34	Screw, M3x6 Sems	0140-6219-128
Rear enclosure (4) –	35	Screw, M4x12 relieved body	1504-3001-000
Rear enclosure (4) –	36	Lockwasher, M4 external	9213-0540-003
Fan (4) –	37	Screw, M3x16	1504-3003-000
Fan (4) –	38	Lockwasher, M3 external	9213-0530-003
Inverters (4) –	39	Screw, M3x6 Nylon	9211-1730-065
CPU to plate (4) –	39 40	Screw, M4x8 Sems	0140-6226-113
	40	OUIGW, IVITAO OCIIIO	0140-0220-113
Connector panel (2) –			
Ground straps, keypads (2) –			
Mounting plate (10) –	4.4	0 144 40 0 144 1410	4000 00 11 000
	41	Screw, M4x12 Pan Washer HD	1009-3341-000
		eplacing a backlight or a backlight inverter, you must replace bo	in inverters and the
	hacklight	assembly found in the Backlight Kit (Itom 7)	

<sup>\*</sup> When replacing a backlight or a backlight inverter, you must replace both inverters and the backlight assembly found in the Backlight Kit (Item 7).

10-14 09/07 1009-0357-000

<sup>\*\*</sup> Keypads will likely be damaged during membrane switch replacement; order parts accordingly.

<sup>\*\*\*</sup> Remove (discard) the heatsink from the replacement fan before installing it onto the CPU heatsink on the HPDU CPU board (Refer to section 9.6.3).







Table 1:		
Language	Keypad Set	
Chinese	M1080398	
Czech	M1080400	
Danish	M1080401	
Dutch	M1080404	
English	M1080395	
Finnish	M1080406	
French	M1080423	
German	M1080429	
Greek	M1080431	
Hungarian	M1080432	
Italian	M1080434	
Japanese	M1080436	
Norwegian	M1080438	
Polish	M1080458	
Portuguese	M1080440	
Russian	M1080441	
Spanish	M1080443	
Swedish	M1080445	
Turkish	M1080448	





### **10.9 Lower electronic enclosure components**






Item	Description	Stock Number
1	Power Controller board (with mounting plate) Replacement Kit	1009-8290-000 M1052795
2	Power Controller board (used with U-Frame)	1011-3572-000-S
3	Battery, sealed lead acid, 12V 4AH	1009-5682-000
4a 4b	Harness, battery Cable, Flex, Battery 4Ah Fused	1009-5557-000 M1049280
5	Tray, battery	1009-3133-000
6	Bracket, battery restraint	1009-3060-000
7a 7b	Label, battery service instructions (harness) Label, battery service Instructions (flex-cable)	1009-5530-000 M1052828
8	Power Supply, universal 150W	M1052831
9	Harness, power supply to PCB	M1049276
10	Harness, J3-PCB to J5-DCB	1009-5552-000
11	Fan (flow into enclosure)	1009-5697-000


10-16 09/07 1009-0357-000

### 10.9.1 Anesthesia Control and Display Connector board



Item	Description	Stock Number
1	Anesthesia Control board (tested)	1009-8291-000
2	Display Connector Board	1009-3005-000
3	Cable, to Display Unit system power interface	1009-5571-000
4	Cable, to Display Unit system signal interface	1009-5572-000
5	Cable, to Airway Module power supply	1009-5555-000
6	Cable, ribbon J2-ACB to J9-DCB	1009-5561-000
7	Harness, J7-ACB to J6-DCB	1009-5556-000
8	Cable, ribbon J1-ACB to underside of Pan Connector Board	1009-5549-000
9	Harness, J3-ACB to underside of Pan Connector Board	1009-5560-000
10	Harness, J4-ACB to J4-PCB	1009-5551-000

### 10.10 Pan electronic enclosure components



Item	Description	Stock Number
1	Fan (flow into enclosure)	1009-5680-000
2	Pan Connector Board	1009-3003-000
3	Gasket, Pan Connector Board	1009-5536-000
4	Harness, to J3-ACB	1009-5560-000
5	Cable, ribbon to J1-ACB (fold excess cable into pan area)	1009-5549-000
6**	Filter Board, ABS	1009-3007-000
7	Gasket, EMI 272 mm (2 each required per board)	1009-5811-000
8	Ventilator Interface Board, calibrated	1009-8236-000
9	Gas Mixer Assembly, complete	Refer to section 10.11
10	Screw, M4x6	1009-3283-000
11	Screw, M4x40	0140-6226-128
12	Lockwasher, M4 external	9213-0540-003
13	Plug, 4-mm (if no N <sub>2</sub> O)	1006-3530-000
14*	Cover, pan electronic enclosure	1009-3047-000
15*	Gasket, EMI electronic enclosure (4 required)	1009-5811-000


st The top side of the cover has rounder edges. Apply the EMI gasket to the underside (sharp edges) of the cover.

10-18 09/07 1009-0357-000

<sup>\*\*</sup> Also requires Item 7, EMI gasket, two each.

### **10.11 Electronic Gas Mixer**


#### **⚠ CAUTION** Ensure a clean environment when servicing the gas mixer.




Item	Description	Stock Number
1	Mixer Assembly - complete	1011-8000-000
2	Valve, 2-way NO (includes screws and gasket)	1009-3014-000
3	Valve, 2-way NC (includes screws and gasket)	1009-3013-000
4	Valve, proportional	1011-3560-000
5	O-ring (2 used with each proportional valve)	6027-0000-165
6	Screw, M3x16 (2 used for mounting each valve)	1504-3003-000
7	Lockwasher, M3 external	9213-0530-003
8	Valve, 3-way NC (includes screws and gasket)	1009-3346-000
9	Flex-cable, valve interface	1009-3359-000
10	Outlet check valve, replacement kit	1009-8246-000
	(includes o-ring and flapper valve)	
11	Retainer, flapper valve	1011-3516-000
12	O-ring, retainer	1011-3518-000
13	Elbow, 1/4 inch tube to 1/8 inch NPT	1011-3071-000
14	Cable, TSI interface	1011-3082-000
* Lubria	note anaringly with Krytov	

<sup>\*</sup> Lubricate sparingly with Krytox.

### 10.12 Pipeline inlet fittings






Item	Description	Stock Number
1	Pipeline inlet - 0 <sub>2</sub> fittings	
	Body, O <sub>2</sub> DISS	1006-5149-000
	Body, O <sub>2</sub> NIST	1006-5158-000
	Body, O <sub>2</sub> DIN	1006-5161-000
	Body, O <sub>2</sub> G 3/8 BSPP	1006-5170-000
	Pipeline inlet assembly O <sub>2</sub> France	1006-8363-000
	Pipeline inlet assembly 0 <sub>2</sub> Canada	1006-8360-000
	Pipeline inlet assembly 0 <sub>2</sub> Australia	1006-8396-000
1	Pipeline inlet - N20 fittings	
	Body, N <sub>2</sub> O DISS	1006-5150-000
	Body, N <sub>2</sub> O NIST	1006-5159-000
	Body, N <sub>2</sub> O DIN	1006-5162-000
	Body, N <sub>2</sub> O G 3/8 BSPP	1006-5171-000
	Pipeline inlet assembly N <sub>2</sub> O France	1006-8362-000
	Pipeline inlet assembly N <sub>2</sub> O Canada	1006-8359-000
	Pipeline inlet assembly N <sub>2</sub> O Australia	1006-8397-000
1	Pipeline inlet Air fitting	
	Body, Air DISS	1006-5151-000
	Body, Air NIST	1006-5160-000
	Body, Air DIN	1006-5163-000
	Body, Air G 3/8 BSPP	1006-5172-000
	Pipeline inlet assembly Air France (service kit)	1006-8361-000
	Pipeline inlet assembly Air Canada (service kit)	1006-8358-000
_	Pipeline inlet assembly Air Australia (service kit)	1006-8398-000
2	O-ring, bore seal	
	$\mathrm{O}_2$ and $\mathrm{N}_2\mathrm{O}$	0210-0479-300
_	Air	0210-0539-300
3	Sintered metal filter with o-ring	1006-8351-000
4	Pipeline check valve with o-ring	1505-3273-000
5	Gas Inlet Manifold (replacement) 0 <sub>2</sub>	1009-8066-000
	$N_2O$	1009-8067-000
	Air	1009-8068-000
6	Relief valve, 689/758 kPa (100/110 psi)	1011-3049-000
7	Screw, M4x20	0144-2124-218
	Lockwasher, M4	9213-0540-003
8	Transducer, pipeline pressure (includes cable)	1011-3000-000

10-20 09/07 1009-0357-000

## **10.13 Cylinder Gas Supplies**




Item	Description	Pin Index	DIN	DIN (Large Cylinder)
1	Gas supply ${\rm O_2}$	1006-3201-000	1006-3207-000	1006-3880-000
1	Gas supply N <sub>2</sub> O	1006-3225-000	1006-3208-000	1006-3881-000
1	Gas supply Air	1006-3203-000	1006-3209-000	


Item	Description	Stock Number
2	Cylinder inlets (Pin Index or DIN for external cylinder)	Refer to section 10.13.1
3	Tee handle beige	0219-3372-600
4	Spacer, gas block (2) Screw, M8 x 25 long socket head cap (2)	1001-4077-000 9211-0680-253
5	Clamp, yoke	1001-4076-000
6	Label Set, cylinder supply, O <sub>2</sub> Label Set, cylinder supply, N <sub>2</sub> O Label Set, cylinder supply, Air	1006-3854-000 1006-3855-000 1006-3856-000
7	Screw M6x25 socket head cap/ 3 per supply Lockwasher (for above screw) Internal M6	9211-0660-254 0144-1118-130
8	Transducer, cylinder pressure (includes cable) Refer to Section 9.13 for replacement procedure.	1011-3001-000

### 10.13.1 Cylinder inlet fittings

#### 1 Pin Index



#### 2 DIN (external cylinder)




Item	Description	Stock Number
1	Cylinder inlets (Pin Index)	
1a	Gasket	0210-5022-300
1b*	O-ring	9221-3013-116
1c	Adapter, inlet	1001-4075-000
1d	Filter, sintered bronze	9914-6380-000
1e	Retaining ring, filter	1001-5954-000
2	Cylinder inlets (DIN)	
2a	Screw, M8x16	0144-2140-242
2b	Sealing ring (DIN)	1009-3356-000
2c	DIN Adapter (O <sub>2</sub> )	1006-4000-000
	DIN Adapter (N <sub>2</sub> 0)	1006-4001-000
	DIN Adapter (Air)	1006-4002-000
2d	O-ring, 0.687 ID, 0.812 OD	0210-0544-300
2e	Filter, sintered bronze	9914-6380-000
2f	Sealing ring, N <sub>2</sub> O DIN Conn 11	1202-3641-000
2g	Adapter, large cylinder N <sub>2</sub> O	1006-4028-000
2h	Adapter, large cylinder O <sub>2</sub>	1006-4027-000

10-22 09/07 1009-0357-000


Lubricate sparingly with Krytox

## 10.14 Vaporizer manifold



Item	Description	Stock Number
	Manifold assembly, complete, two position	1006-8355-000
	Manifold assembly, complete, one position	1009-8065-000
1	O-ring, 0.687 inch ID 0.812 inch OD	0210-0544-300
2	Spring, compression	1006-3736-000
3	Valve kit, includes seal	1006-8373-000
4	Seal	1006-3690-000
5	O-ring, 14.3 mm ID	1102-3043-000
	(Package of 6 o-rings)	1102-3016-000
6	Nipple, vaporizer port (New Style)	1006-4215-000
7	Screw, M2.5 - 0.45x6 PAN, Pozidriv, SST	1006-3037-000
8	Spring, Dzus	1102-3056-000
9	Seat, check valve	1006-1352-000
10	O-ring 27.1 OD 21.89 mm ID	1006-3866-000
11	Flapper	0211-1451-100
12	Housing	1009-8477-000
13	Screw, M4 x 30, cap head	9211-0640-304
14	Valve, relief, 5.5 psi, 7/16-20 THD	1006-4128-000
15	Flexible tubing, $1/4$ inch, mixed gas	1001-3064-000


## **10.15 Vent Engine Housing**




Item	Description	Stock Number	Qty
1	Vent Engine Cover Plate Assy	1407-7009-000	
2	CASTING VENT ENG HOUSING	1407-3301-000	
3	TAB GUIDE BELLOWS BASE	1407-3313-000	(2)
4	SCR M3X16 POSI DR PAN HD A4 SST	1504-3003-000	(2)
5	Cap, Plug	1406-3524-000	
6	FITTING PNL MOUNT 3.18 HOSE BARB UNION	1504-3014-000	(2)
7	PLUG HOLE 15.9 DIA NYLON MICRO PLASTICS	1006-1473-000	
8	PLATE CONN VENT	1407-3321-000	
9	SCR M4X8 POZI-DR DIN84 PAN SERRATED	1006-3178-000	(3)
10	Harness, Vent Engine Board to Connector Plate	1009-5545-000	
11	BLOCK LATCHING DSUB CONN	1504-3617-000	(2)
12	SCR 4-40 X 3/8 SKT BCG HD CAP	0144-2117-206	(2)
13	CLIP-SUCTION BAG HOSE	1407-3327-000	
14	SCR M5 X 16 PAN PH HD SST	9211-8350-163	(2)
15	Lockwasher	0144-1118-220	(2)
16	Vent Engine	Refer to section 10.16	

10-24 09/07 1009-0357-000

### 10.16 Vent Engine





2 4 5

-	
Vent Engine Assembly, Service (Avance/Aisys)	1009-8216-000
Gas Inlet Valve (GIV) components	Refer to section 9.10.3
Solenoid, 3-way NO (12 Volt)	1503-3853-000
Screw, m1.6x14	1006-4730-000
Retaining ring, 34.9 mm	1500-3158-000
Cap, inlet valve	1503-5006-000
O-ring, upper Viton	9221-3032-116
Shuttle, inlet valve	1503-5018-000
U-cup, upper EDPM (fits on shuttle valve)	1503-3090-000
U-cup, lower Viton (fits on shuttle valve)	1503-3089-000
O-ring, lower Viton	1503-3108-000
Filter (under GIV), 2-micron (install course side DOWN)	1504-3708-000
Fitting, manifold pressure	1500-3116-000
Reservoir, pneumatic engine	1504-3704-000
O-ring, base, 56.87 ID x 60.43 OD	1504-3614-000
O-ring, screw head, 0.219 ID x 0.344 OD	0210-0686-300
Screw, M6x90	1504-3004-000
Flow control valve (12 Volt)	1503-3854-000
O-ring under flow control valve (2 each)	1503-3056-000
Drive gas check valve	1503-3006-000
O-ring under drive gas check	1503-3213-000
Interface Manifold	Refer to section 10.16.1
	Gas Inlet Valve (GIV) components Solenoid, 3-way NO (12 Volt) Screw, m1.6x14 Retaining ring, 34.9 mm Cap, inlet valve O-ring, upper Viton Shuttle, inlet valve U-cup, upper EDPM (fits on shuttle valve) U-cup, lower Viton (fits on shuttle valve) O-ring, lower Viton Filter (under GIV), 2-micron (install course side DOWN) Fitting, manifold pressure Reservoir, pneumatic engine O-ring, base, 56.87 ID x 60.43 OD O-ring, screw head, 0.219 ID x 0.344 OD Screw, M6x90 Flow control valve (12 Volt) O-ring under flow control valve (2 each) Drive gas check valve O-ring under drive gas check

1504-3623-000

1011-3165-000


4-pin connector 7 8 9 10

1009-0357-000 09/07 10-25

Vent Engine Connector board (not part of assembly)

Regulator, 172 kPa BCG

#### 10.16.1 Vent Engine - under side



Item	Description	Stock Number
1	Interface Manifold, pneumatic engine (with free breathing valve and mechanical overpressure valve)	1504-8505-000
	Ö-ring, 12.42 ID x 15.98 OD (2)	1006-3615-000
1a	Seat, free breathing valve	1503-3204-000
1b*	Valve, flapper	0211-1454-100
1c	O-ring	1503-3208-000
2	Fitting, 6.35-mm (1/4-inch)	1504-3621-000
3	Plug, 6.35-mm (1/4-inch)	1503-3245-000
4	Fitting, barbed	1504-3014-000
5	Manifold	1503-3843-000
6**	Gasket, manifold	1503-3845-000
7	Plate, manifold	1503-3844-000
8***	Screw, M4x8 Pozidriv PAN	1006-3178-000

If necessary, clean with alcohol before installing new; trim off flush with outside surface of seat


(refer to removed flapper).
Install gasket into manifold. Check to see that it is properly positioned.

\* Carefully install plate onto manifold making sure not to disturb the gasket.

First, start all screws. Then, torque to 1.7 N-m (15 lb-in) using sequence shown.


10-26 09/07 1009-0357-000

### 10.17 ABS to machine Interface Components (SCGO)



Item	Description	<b>Stock Number</b>
1	SCGO Selector Module, complete	1009-3098-000
1a	Flush pressure switch (includes o-ring)	1006-3972-000
1t	O-ring	1006-3213-000
10	Screws, M3x20	0144-2124-201
10	Switch, mode (CGO/SCGO), kit	1009-3282-000
16	e Valve, relief 150 cmH20	1009-3052-000
1	f Solenoid kit CGO	1009-3279-000
2	Tubing, silicone (110 mm, 100 mm)	1009-3164-000
3	Tie wrap	0203-5915-300


### 10.18 ABS to machine Interface Components (ACGO)




Item	Description	Stock Number
	Service Kit to covert SCGO Avance machine to ACGO	1009-8249-000
1	Port, ACGO body	1009-3096-000
2	Screw, M4x30	9211-0640-304
3	Lockwasher, M4	9213-0540-003
4	Cap, ACGO check valve	1009-3095-000
5	Screw, M4x8	9211-1040-069
6	Disk, ACGO check valve	1009-3062-000
7	Flapper, ACGO check valve	1009-3097-000
8*	O-ring	0210-0543-300
9	Fitting, elbow barbed	1009-3160-000
10*	O-ring	0210-0691-300
11	Screw, M3x6	9211-1030-055
12	ACGO Selector Switch, complete	1009-3099-000
	(without guard - item 13)	
12a	Flush pressure switch	1006-3972-000
12b	O-ring	1006-3213-000
12c	Screws, M3x20	0144-2124-201
13	Guard	1009-3140-000
14	Tubing, silicone	1009-3164-000
15	Tie wrap	0203-5915-300
* Lubric	ate sparingly with Krytox.	

10-28 09/07 1009-0357-000

## **10.19 Flush Regulator and Flush Valve**






Item	Description	Stock Number
1	Flush valve, without button	1006-8357-000
2	Flush Button with rod	1011-3354-000
3	Spring	1006-3186-000
4	E-ring	0203-5225-300
5	Bracket	1011-3355-000
6	Screw, M4x8	1006-3178-000
7	Screw, M4x12	0140-6226-111
8	Lockwasher, M4	9213-0540-003
9	Regulator, O <sub>2</sub> Flush	1011-3168-000

## 10.20 Front panel, Alt 02, and system switch


Table 1:	
Language	Alt 0 <sub>2</sub> Label
	-
Chinese	1009-3332-000
Czech	1009-3322-000
Danish	1009-3331-000
Dutch	1009-3317-000
English	1009-3315-000
Finnish	1009-3321-000
French	1009-3316-000
German	1009-3314-000
Greek	1009-3326-000
Hungarian	1009-3327-000
Italian	1009-3318-000
Japanese	1009-3328-000
Norwegian	1009-3324-000
Polish	1009-3323-000
Portuguese	1009-3320-000
Russian	1009-3330-000
Spanish	1009-3319-000
Turkish	1009-3325-000



Item	Description	Stock Number
1	Switch, system On/Standby	1006-8452-000
2	Harness, On/Standby system switch	1009-5542-000
3	Panel	1009-3019-000
4	Label, Alt O <sub>2</sub>	See Table 1
5	Needle Valve assembly, flow control	1011-3429-000
6	Knob (set screw not included)	1011-3472-000
	Set screw	9211-0830-053
7	Bracket	1009-3127-000
	Screw, M4x8 (bracket to needle valve - 2 each)	0144-2436-108
	Screw, M4x10 (assembly to front panel - 4 each)	1009-5534-000
8	Switch, Alt O <sub>2</sub> (includes harness)	1009-5517-000
9	Flowmeter, Alt O <sub>2</sub>	1011-3428-000
10	Bracket, flowmeter	1009-3126-000
	Screw, 10-32x3/8 (bracket to flowmeter - 2 each)	0140-6631-107
	Screw, M4x10 (assembly to front panel - 4 each)	1009-5534-000
11	LED assembly, mains green	1009-5514-000

10-30 09/07 1009-0357-000

### **10.21** Breathing system interface




	Item	Description	Stock Number	Qty
	1	Assembly, main support casting	1407-7010-000	
	2	Bolt, M6x16 flange	1009-3125-000	(5)
	3	Handle, grip	1407-3317-000	
	4	Screw, M6X16 Sems	0144-2436-109	(2)
	5	Latch, push to close	1407-3309-000	
	6	Latch, push to close w/microswitch	1407-3310-000	
	7	Screw, SKT HD CAP M3x8 SST	1006-3865-000	(2)
	8	Washer, lock external M3	9213-0530-003	(2)
	9	Port, plug circuit (ACGO)	1407-3333-000	
Note: SCGO only applies	10	Port, fresh gas (SCGO)	1407-3314-000	
to Avance and Aisys	11*	Seal, U-Cup 12.7 ID BCG 19.05 OD EPR	1407-3320-000	(4)
machines.	12	Ring, retaining 15.88 SHAFT DIA TYPE E SST	1406-3446-000	(4)
	13	Port, sample gas	1407-3318-000	(2)
	14	Connector, bulkhead O2 Cell, with harness	1009-5586-000	
	15	Ring, retaining 9.53 SHAFT DIA TYPE E SST	1406-3277-000	
	16	Switch, subminiature w/QDISC terminals	1406-3296-000	(2)
	17	Screw, M2.5 x10	1009-3153-000	(4)
	18	Bracket, BTV switch	1407-3319-000	
	19	Lever, BTV switch	1407-3325-000	
	20	Bracket, cap BTV	1407-3324-000	
	21	Screw, SKT HD CAP M3x8 SST	1006-3865-000	(2)
	22	Washer, lock external M3	9213-0530-003	(2)
	23	Bracket, bypass switch	1407-3139-000	
	24	Paddle, switch actuator	1407-3141-000	
	25	Bracket, paddle hinge	1407-3140-000	
	26**	Screw, M6x6 set cup	1007-3329-000	(2)
	* Lubri	cate sparingly with Krytox.		

<sup>\*\*</sup> Apply Loctite 242.

# **10.22 Breathing System**


#### 10.22.1 APL Valve



Item	Description	Stock Number	QTY
1	APL Valve Assy (includes items 2 through 6)	1009-8200-000	
2	SPRING CPRSN 53.14 OD 36.8 L 1.48 N/MM	1406-3328-000	
3	RETAINER SPRING APL	1407-3404-000	
4	DIAPHRAGM APL	1406-3331-000	
5	CAGE APL	1406-3333-000	
6	POPPET APL VALVE	1406-3332-000	
7	RAMP APL	1407-3400-000	
8	COVER APL	1407-3405-000	
9	MANIFOLD APL/BTV	1407-3401-000	
10	Cover, Manifold APL/BTV (with 22-mm male bag port)	1407-3402-000	
	Cover, Manifold APL/BTV (with Australian bag port - 22 mm female)	1407-3412-000	
11	WEIGHT DEAD 14CM H2O BCG ABS NEG RELIEF	1407-3406-000	
12	SEAL ABS NEG RELIEF VLV	1407-3407-000	
13	O-RING 22 ID 30 OD 4 W SI 40 DURO	1407-3104-000	
14	O-RING 88.49 ID 95.55 OD 3.53 W SILICONE 50 DURO	1407-3403-000	
15	O-RING 1.049ID 1.2550D .103W EPDM NO 121	1407-3408-000	
16	SCR M4X16 BT SKT HD SST TYPE 316	0140-6226-115	(2)
17	Lockwasher, M4 external	9213-0540-003	(2)
18	SCR THUMB M4 SHLDR 7.5 X 7	1407-3410-000	(3)
19	RING RETAINING 3.96 SFT DIA CRESCENT SST	1407-3411-000	(3)
20	SCR M4 X 40 FL HD SST PH	0140-6226-122	(2)
21	O-RING 2.9 ID 6.46 OD 1.78 W EP 70 DURO	1407-3409-000	(2)
22	SCR SEMS M4X8 BT SKT HD W/EXT L/W SST 316	0144-2436-108	(3)


10-32 09/07 1009-0357-000

#### 10.22.2 Bag/Vent Switch



Item	Description	Stock Number	QTY
	BTV Switch Cartridge	1407-7003-000	
1	COVER BTV	1407-3500-000	
2	SCR SEMS M4X8 BT SKT HD W/EXT L/W SST 316	0144-2436-108	(2)
3	O-RING 44.02 ID 51.1 OD 3.53 W SI 70 DURO	1407-3507-000	
4	SEAL, BTV	1407-3506-000	


#### 10.22.3 Absorber canister



Item	Description	Stock Number	Qty
1	Multi-Absorber canister, reusable (does not include absorbent)	1407-7004-000	
2	Cover assembly, CO <sub>2</sub> canister	1009-8240-000	
3	Foam, CO <sub>2</sub> canister (pack of 40)	1407-3201-000	
4	O-ring	1407-3204-000	
5	Canister, CO <sub>2</sub>	1407-3200-000	
	Multi-Absorber canister, disposable (white to violet; pack of 6)	8003138	
	Multi-Absorber canister, disposable (pink to white; pack of 6)	8003963	

10-34 09/07 1009-0357-000


#### 10.22.4 Flow Sensor Module



Description	Stock Number	Qty
Flow Sensor Module (*)	1407-7001-000	
Flow Sensor (plastic - moisture resistant)	1503-3858-000	
Flow Sensor (metal - autoclavable)	1503-3244-000	
Flow Port Adapter	1503-3849-000	
Cover, Flow Sensor (Avance and Aespire machines)	1407-3000-000	
Cover, Flow Sensor (Aisys machines)	1011-3283-000	
HOLDER FLOW SNSR UPPER	1407-3002-000	
HOLDER FLOW SNSR LOWER	1407-3003-000	
SCR THUMB M6X43 SST	1406-3304-000	
SCR M4 .07 X 10 SKT CAP BUTTON HEAD SST	0144-2117-718	(2)
CUFF FLOW SNSR	1407-3004-000	(2)
LATCH FLOW SNSR	1407-3001-000	
SPR TORSION FLOW SNSR LATCH	1407-3005-000	
RING TRUARC .188 SHAFT E-RING SST	0203-5225-300	
	Flow Sensor Module (*) Flow Sensor (plastic - moisture resistant) Flow Sensor (metal - autoclavable) Flow Port Adapter Cover, Flow Sensor (Avance and Aespire machines) Cover, Flow Sensor (Aisys machines) HOLDER FLOW SNSR UPPER HOLDER FLOW SNSR LOWER SCR THUMB M6X43 SST SCR M4 .07 X 10 SKT CAP BUTTON HEAD SST CUFF FLOW SNSR LATCH FLOW SNSR SPR TORSION FLOW SNSR LATCH	Flow Sensor Module (*)  Flow Sensor (plastic - moisture resistant)  Flow Sensor (metal - autoclavable)  Flow Port Adapter  Cover, Flow Sensor (Avance and Aespire machines)  Cover, Flow Sensor (Aisys machines)  HOLDER FLOW SNSR UPPER  HOLDER FLOW SNSR LOWER  SCR THUMB M6X43 SST  SCR M4 .07 X 10 SKT CAP BUTTON HEAD SST  LATCH FLOW SNSR  LATCH  1407-3005-000

<sup>\*</sup> The flow sensors are not included in the flow sensor module.

#### 10.22.5 Breathing Circuit Module




Item	Description	Stock Number	Qty
	Breathing Circuit Module *	1407-7002-000	
1	LENS CIRCUIT CHK VALVES	1407-3101-000	
2	0-RING 44.02 ID 51.1 OD 3.53 W SI 70 DURO	1407-3507-000	(2)
3	Check Valve Assembly	1406-8219-000	(2)
3a	RETAINER DISK 26.97D 12.7H 0.76T	1400-3017-000	(2)
3b	DISC CHK V RVSBL 1.025D	0210-5297-100	(2)
3c	SEAT UNIDIRECTIONAL V B/S	1406-3396-000	(2)
3d	0-RING 20.35 ID 23.90 OD 1.78W	1406-3397-000	(2)
4	PLATE CIRCUIT FLANGE	1407-3110-000	
5	SCR SEMS M4X8 BT SKT HD W/EXT L/W SST 316	0144-2436-108	(6)
6	HOOK LATCH	1407-3604-000	
7**	0-RING 22 ID 30 OD 4 W SI 40 DURO	1407-3104-000	(2)
8**	O-RING 12.37 ID 17.6 OD	1006-3968-000	
9	MANIFOLD CIRCUIT	1407-3100-000	
10	PIN CANISTER PIVOT	1407-3109-000	
11	RING TRUARC .188 SHAFT NO 5133-18H E-RING SST	0203-5225-300	(5)
12	LEVER CANISTER LATCH	1407-3102-000	
13	PIN CANISTER LEVER	1407-3108-000	
14*	O <sub>2</sub> Cell (includes o-ring)	6050-0004-110	
	O-ring, cell	1406-3466-000	
15*	Plug with o-ring (for units without circuit $O_2$ sensing)	1503-3857-000	
	O-ring, plug	1406-3466-000	
16*	Cable, O <sub>2</sub> Cell	1009-5570-000	
TI (	5 117 1 5 101 11 11 12 11 12 11 1		

<sup>\*</sup> The  ${\rm O}_2$  cell (or plug) and the cell cable are not included in the breathing circuit module.

10-36 09/07 1009-0357-000

<sup>\*\*</sup> Lubricate sparingly with Krytox.


#### 10.22.6 Exhalation valve



Item	Description	Stock Number	Qty
	Exhalation Valve assembly	1407-7005-000	
1	Base, exhalation valve	1407-3701-000	
2	Diaphragm assembly	1503-8121-000	
3	Seat, exhalation valve	1407-3704-000	
4	Cover, exhalation valve	1407-3700-000	
5	Screw, M4x16 PH PAN HD	9211-0440-163	(3)
6	0-ring, 2.9 ID 6.46 OD 1.78 W EP 70 DURO	1407-3409-000	(3)
7	Thumbscrew, M6x43 10mm head	1406-3306-000	(2)
8	0-ring, 4.47 ID x 8.03 OD 1.78 W EPR 70 DURO	1407-3703-000	(2)
9	Retainer, disk 26.97D 12.7H 0.76T SST flutter	1400-3017-000	
10	Weight, dead 10 cm H2O	1406-3572-000	
11	Seat	1406-3571-000	
12*	O-ring, OD19.16 ID15.6 EPDM DURO 70 -016	1006-3616-000	
13	Ring, retaining 19.05 SHAFT DIA	1406-3577-000	

<sup>\*</sup> Lubricate sparingly with Krytox.


#### 10.22.7 Bellows



Item	Description	Stock Number
1	Bellows housing	1500-3117-000
2	Bellows	1500-3378-000
3	Rim	1500-3351-000
4	Pressure relief valve assy	1500-3377-000
5	Latch, base	1500-3352-000
6	Seal, base	1500-3359-000
7	Base, bellows	Refer to section 10.22.8
8	Manifold, bellows base	1407-3702-000


10-38 09/07 1009-0357-000

#### **10.22.8 Bellow base**



Item	Description	Stock Number
1	Bellows Base Assy	1407-7006-000
18	a Latch Assy	1407-7007-000
2	HOOK LATCH	1407-3604-000
3	E-Ring	0203-5225-300


#### **10.22.9 Bag Arms**



Ite	m	Description	Stock Number	Qty	
1		Bag Arm Assembly (complete)	1009-8159-000	4.7	
2		Bag Arm Upper Assembly	1407-7011-000		
	2a	Cover, bag port housing	1407-3807-000		
		Screw, M3x20	0140-6719-103		
		Lockwasher, M3 internal	9213-0430-003		
	2b	Housing, bag port	1407-3806-000		
	2c	Lever, lock release	1407-3808-000		
	2d	Ring, retaining	1406-3577-000		
	2e	Nut, M3 Nyloc	0144-3536-112		
3		Bag Arm Lower Assembly	1407-7012-000		
	3a	Pad, Friction Material			
	3b	Bag Arm post subassembly			
	3с	Spring			
	3d	Washer, shoulder			
	3e	Nut, M5 Nyloc			
4		Pin, dowel 3.18 DIA 31.8 L SST	1407-3804-000		
5*		Screw, M3x20 SKT HD CAP	0144-2124-201	(2)	
6		Washer, M3 flat	0144-1003-132	(2)	
Ite	ms if ı	no Bag Arm			
7		Clip, patient tubing	1407-3810-000		
8		Washer, shoulder	1407-3814-000		
9		Lockwasher, M4 external	9213-0540-003		
10		Screw, M4x16	9211-0440-163		
Ac	Accessory items				
11		Bag arm connector, reusable (All except China)	8004459		
		Bag arm connector, reusable (China)	M1082613		
* A	* Apply Loctite 242.				

10-40 09/07 1009-0357-000

### 10.22.10 EZchange Canister system (CO<sub>2</sub> Bypass)



Item	Description	Stock Number	Qty
	EZchange Canister module	1407-7021-000	
1	Cover, Bypass Manifold	1407-3123-000	
2	Manifold, Bypass	1407-3113-000	
3	Screw, M3x8 PT PAN PH SST	0142-4254-106	(2)
4**	O-ring, 59.92 ID 66.98 OD 3.53 W SIL 50 DURO	1407-3142-000	
5	Cap, Manifold	1407-3130-000	
6	Lever, Switch Actuator	1407-3116-000	
7	Spring, Torsion Switch Actuator Lever	1407-3117-000	
8*	Screw, M3x0.5 Shoulder 4 DIA X 4 L SST	1407-3915-000	(2)
9	Lever, Canister Latch	1407-3115-000	
10	Seal, Drain	1407-3121-000	
11**	O-ring, 37.69 ID 44.75 OD 3.53 W SIL 50 DURO	1407-3129-000	
12**	O-ring, 50.39 ID 57.45 OD 3.53 W SIL 50 DURO	1407-3143-000	
13**	O-ring, OD30 ID 22 4W SIL 40 DURO	1407-3104-000	(2)
14***	Valve, Housing Assembly Bypass	1407-3126-000	
15	Screw, Thumb M4 Shoulder 7.5 X 7	1407-3410-000	
16	Ring, Retaining 3.96 Shaft DIA SST	1407-3411-000	
17	Cradle Canister	1407-3118-000	
18	Screw, M4x10 CSK SKT HD SST TYPE 316	0140-6226-119	(2)
19**	Spacer, Shoulder 6.8 DIA x4.1 L	1407-3120-000	(2)
20	Support, Cradle Pivot	1407-3119-000	
21	Screw, M4x8 Sems BT SKT HD SST 316	0144-2436-108	(3)
* Apply	Loctite 242		

Apply Loctite 242.

<sup>\*\*</sup> Lubricate sparingly with Krytox.

\*\*\* Rubber valve seats can not be removed from assembly (Item 14).


#### **10.22.11 Condenser**

Item	Description	Stock Number	Qty
	Condenser assembly (includes all Items)	1407-7026-000	
	Condenser module (Items 1 through 16)	1407-7025-000	
	Condenser (Items 17 through 26)	1407-7024-000	
1	Cover, Bypass Manifold	1407-3123-000	
2	Manifold, Condenser	1407-3114-000	
3	Screw, PT PAN PH M3X8 SST	0142-4254-106	(2)
4**	0-ring, 63.09 ID 70.15 OD	1407-3142-000	
5	Lever, Canister Latch	1407-3115-000	
6*	Screw, M3x0.5 Shoulder 4 DIA X 4 L SST	1407-3915-000	(2)
7	0-ring, 12.37 ID 17.6 OD	1006-3968-000	
8**	0-ring, 37.69 ID 44.75 OD	1407-3129-000	
9**	0-ring, 50.39 ID 57.45 OD	1407-3143-000	
10**	0-ring, 22 ID 30 0D	1407-3104-000	(2)
11	Cap, Valve Housing	1407-3125-000	
12	Screw, Thumb M4 Shoulder 7.5 X 7	1407-3410-000	
13	Ring, Retaining 3.96 Shaft DIA SST	1407-3411-000	
14*	Pin, Condenser Manifold	1407-3131-000	(2)
15	Support, Cradle Pivot	1407-3119-000	
16	Screw, M4x8 Sems BT SKT HD	0144-2436-108	(3)
17	Tube Assembly	1407-3133-000-S	
18	Washer (part of Item 17)		(2)
19	Fork (part of Item 17)		
20	Screw, M5.28x20 (part of Item 17)		
21	Reservoir, Condenser	1407-3137-000	
22	Seal, Condenser Reservoir	1407-3136-000	
23	Spring, Compression Drain Button	1407-3135-000	
24	Button, Drain	1407-3134-000	
25	Cover, Condenser	1407-3138-000	
26	Guard	1407-3145-000	

<sup>\*</sup> Apply Loctite 242.

10-42 09/07 1009-0357-000

<sup>\*\*</sup> Lubricate sparingly with Krytox.




#### 10.22.12 Anesthetic Gas Scavenging System — AGSS

**10.22.13 Passive AGSS** Items 1 through 12 are included in all AGSS kits.

Item	Description, Common Parts	Stock Number	Qty
1	Seal, Receiver Body	1407-3901-000	
2	Reservoir	1407-3903-000	
3	Seal and scavenging down-tube	1407-3904-000	
4	Thumbscrew, M6x28.5	1406-3305-000	
5	O-ring, 4.42 ID, 9.65 OD	1407-3923-000	(2)
6	Thumbscrew, M6x43	1406-3304-000	
7	Valve, unidirectional (negative pressure relief)	1406-8219-000	
7a	Seat, Valve, Negative Pressure	1406-3396-000	
7b	Retainer, disc	1400-3017-000	
*7c	0-ring, 20.35 ID, 23.90 OD	1406-3397-000	
7d	Disc, check-valve	0210-5297-100	
8*	O-ring, 22 ID, 30 OD silicone	1407-3104-000	(2)
9*	O-ring, 21.95 ID, 25.51 OD	1406-3558-000	
10	Screw, M4x8	9211-0640-083	(2)
11	Cap, 3.18 Barb, Silicone	1406-3524-000	
12	Adapter, auxiliary inlet, 30-mm male to 30-mm male	M1003134	
13	Adapter, auxiliary inlet, 30-mm male to 19-mm male	M1003947	
Passive	AGSS Specific Parts		
14	Receiver, Passive/Adjustable	1407-3908-000	
15	Plug Assembly, tethered	1407-3909-000	(2)
16	Screw, shoulder M3	1407-3915-000	
17	Connector, 30-mm ISO, Male	1406-3555-000	
18	Adapter, scavenging, 30-mm female to 19-mm male	1500-3376-000	(5 pack)

<sup>\*</sup> Lubricate sparingly with Krytox

10-44 09/07 1009-0357-000




# 10.22.14 Adjustable AGSS

Items 1 through 12 are included in all AGSS kits.

Item	Description, Common Parts	Stock Number	Qty
1	Seal, Receiver Body	1407-3901-000	
2	Reservoir	1407-3903-000	
3	Seal and scavenging down-tube	1407-3904-000	
4	Thumbscrew, M6x28.5	1406-3305-000	
5	O-ring, 4.42 ID, 9.65 OD	1407-3923-000	(2)
6	Thumbscrew, M6x43	1406-3304-000	
7	Valve, unidirectional (negative pressure relief)	1406-8219-000	
7a	Seat, Valve, Negative Pressure	1406-3396-000	
7b	Retainer, disc	1400-3017-000	
*7c	0-ring, 20.35 ID, 23.90 OD	1406-3397-000	
7d	Disc, check-valve	0210-5297-100	
8*	O-ring, 22 ID, 30 OD silicone	1407-3104-000	(2)
9*	O-ring, 21.95 ID, 25.51 OD	1406-3558-000	
10	Screw, M4x8	9211-0640-083	(2)
11	Cap, 3.18 Barb, Silicone	1406-3524-000	
12	Adapter, auxiliary inlet, 30-mm male to 30-mm male	M1003134	
13	Adapter, auxiliary inlet, 30-mm male to 19-mm male	M1003947	
Adjusta	ble AGSS Specific Parts		
14	Receiver, Passive/Adjustable	1407-3908-000	
15	Plug Assembly, tethered	1407-3909-000	
16	Screw, shoulder M3	1407-3915-000	
17	Needle Valve Assembly (with DISS EVAC connector)	1407-3918-000	
18	Bag with 30 mm male connector	8004460	

<sup>\*</sup> Lubricate sparingly with Krytox

10-46 09/07 1009-0357-000




#### **10.22.15 Active AGSS** Items 1 through 12 are included in all AGSS kits.


Item	Description, Common Parts	Stock Number	Qty
1	Seal, Receiver Body	1407-3901-000	
2	Reservoir	1407-3903-000	
3	Seal and scavenging down-tube	1407-3904-000	
4	Thumbscrew, M6x28.5	1406-3305-000	
5	0-ring, 4.42 ID, 9.65 OD	1407-3923-000	(2)
6	Thumbscrew, M6x43	1406-3304-000	
7	Valve, unidirectional (negative pressure relief)	1406-8219-000	
7a	Seat, Valve, Negative Pressure	1406-3396-000	
7b	Retainer, disc	1400-3017-000	
*7c	0-ring, 20.35 ID, 23.90 OD	1406-3397-000	
7d	Disc, check-valve	0210-5297-100	
8*	O-ring, 22 ID, 30 OD silicone	1407-3104-000	(2)
9*	O-ring, 21.95 ID, 25.51 OD	1406-3558-000	
10	Screw, M4x8	9211-0640-083	(2)
11	Cap, 3.18 Barb, Silicone	1406-3524-000	
12	Adapter, auxiliary inlet, 30-mm male to 30-mm male	M1003134	
13	Adapter, auxiliary inlet, 30-mm male to 19-mm male	M1003947	
Active A	GSS Specific Parts		
14	Receiver, with air brake	1407-3900-000	
15	Seal, for filter and orifice	1407-3902-000	(2)
16	Filter	1406-3521-000	
A .45 II	ligh Flour Conneille Doub		
	ligh Flow Specific Parts	1406 2557 000	
17a	Connector, high flow M30 thread	1406-3557-000	
18	Orifice, high flow	1407-3920-000	
Active L	ow Flow with EVAC connector Specific Parts		
17b	Connector, low flow EVAC	1406-3597-000	
18	Orifice, low flow	1407-3919-000	
Active L	ow Flow with 25 mm connector Specific Parts		
17c	Connector, low flow 25 mm	1406-3573-000	
18	Orifice, low flow	1407-3919-000	
Active L	ow Flow with 12.7 mm hose barb connector Specific Pa	ırts	
17d	Connector, low flow 12.7 mm (1/2 inch)	1406-3574-000	
18	-none-		
	TI 111 00 100 1 1 1 1 7 7 7 7 1		
	ow Flow with 30 mm ISO male connector Specific Parts		
17e	Connector, 30 mm ISO, Male	1406-3555-000	
18	Orifice, low flow	1407-3919-000	

<sup>\*</sup> Lubricate sparingly with Krytox

10-48 09/07 1009-0357-000



### 10.22.16 AGSS gauge, and sample return








Item	Description	Stock Number
1	Flowtube, AGSS	1406-3560-000
2	Clip, AGSS flowtube	1009-3181-000
3	Label, flow indicator AGSS Label, flow indicator AGFS (for German variant) Label, blank (for machines without flow indicator)	1406-3527-000 1009-3301-000 1009-3241-000
4	Label, AGFS (for German variant)	1009-3300-000
5	Coupling, Colder (Kit includes mounting nut)	1009-8321-000


10-50 09/07 1009-0357-000


### 10.22.17 Airway module (MGAS) components















Item	Description	Stock Number
1	Guide, MGAS module	1009-3072-000
2	Screw, M4x8	1006-3178-000
3	Cable, MGAS to Connector board	1009-5555-000
4	Screw, M4x12 THD FORMING	1009-3109-000
5	MGAS monitoring board	1009-5573-000
6	Box, MGAS power supply	1009-3092-000
7*	Screw, 4-40	0140-6617-103
8	Lockwasher, #4 split	0144-1104-331
9*	Standoff	1504-3007-000
10	Cover, MGAS power supply	1009-3093-000
11	Screw, M4x8	0140-6226-113

<sup>\*</sup> Apply Loctite 242.

## **10.23 Integrated Suction Regulator**


#### 10.23.1 Major Components (Continuous and Venturi suction)



Item	Description	Stock Number		
1	Suction Control Module	Refer to section 10.23.2		
2	Venturi Assembly	Refer to section 10.23.3		
3	Cover, blank (if no Suction)	1009-3271-000		
4	Bracket, blank cover mounting	1009-3270-000		
5	Screw, M4x10 self-tapping	1009-5534-000		
6	Manifold	1009-3123-000		
7	Screw, M5x20 BHSCS PT THD FORMING	1009-3384-000		
8a*	Connector, NIST	1011-3524-000		
8b*	Connector, Barb	0221-0702-300		
8c*	Connector, Air Liquide	1009-8292-000		
8d	Muffler, for Venturi Drive	1011-3511-000		
9	Coupling, Colder insert metal	1009-3135-000		
10	Coupling, Colder body black	1009-3373-000		
11	Coupling, Colder insert black	1009-3374-000		
12	Coupling, Colder body white	1009-3371-000		
13	Coupling, Colder insert white	1009-3372-000		
14	Tubing, Tygon	Refer to section 10.31		
15	Fitting, barb to 8-mm Legris	1009-3137-000		
16	Cap, white	1009-3385-000		
17	Overflow Safety Trap	6700-0647-800		
* Apply Teflon tape to threads (not 8d)				

10-52 09/07 1009-0357-000

#### **10.23.2 Suction Control Module**




Item	Description	Stock Number
1	Gauge, 760 mmHg	1009-3227-000
	Gauge, 1 Bar	1009-3228-000
1a	O-ring, Gauge (included with gauge assy, 2ea. required)	6700-0133-500
2	Control panel assembly, with suction regulator knob and mode control knob	1009-3213-000
3	Regulator Module (plugs into manifold assembly)	6700-1225-800
За	O-ring, Regulator Module, Large (included with regulator module)	6700-0136-500
3b	O-ring, Regulator Module, Stem (included with regulator module)	0210-0527-300
4	Manifold Assembly, without Gauge and Regulator Module	1009-3277-000
5	Screw, #6 - 2 inch	1009-3340-000
6	Mounting bracket	1009-3255-000
7	Screw, #6 - 1 inch	1009-3339-000
8	Filter	0206-5159-300
9	Pilot valve adapter assembly (includes plunger, jam nut, and valve assembly)	1009-3278-000
10	Cap, white	1009-3192-000

<sup>\*</sup> Lubricate the regulator module o-rings and the mating bore of the manifold sparingly with Dow 111 lubricant.

<sup>\*\*</sup> Drop the plunger (**9a**), round end first, into the manifold. Thread the pilot valve into the manifold body. Set the mode switch to raise the plunger. Adjust the pilot valve (**9b**) so that the plunger actuates the pilot valve approximately half of its travel. Tighten the jam nut (**9c**).


### 10.23.3 Venturi assembly





Item	Description	Stock Number
1	C-clip retainer, Truarc	1500-3158-000
2	Elbow fitting, 4-mm Legris	1006-3663-000
3	Cap	1011-5002-000
4	Spoppet	1011-5001-000
5	Seal, u-cup large	1503-3090-000
6	Orifice	1011-3508-000
7	Screen, 150 mesh monel	1001-3808-000
8	Seal, u-cup small	1503-3089-000
9	Body	1011-5000-000
10	Venturi	1011-3509-000
11	Elbow fitting, 8-mm Legris	1011-3510-000
12	O-ring, large	9221-3032-116
13	O-ring, small	1503-3108-000
14	Check valve	1011-8002-000
15	Bracket, Venturi mounting	1009-3182-000
16	Screw, M5x20 BHSCS PT THD FORMING	1009-3384-000
17	Cable tie	0203-5915-300
18	Tubing, Tygon	Refer to section 10.31

10-54 09/07 1009-0357-000

# 10.24 Auxiliary 0<sub>2</sub> Flowmeter








Item	Description	Stock Number
1	Flowmeter, 1-10 L/min, Complete with fittings installed Flowmeter, 1-10 L/min, without fittings	1006-8424-000 1006-3841-000
2	Knob, gray	1011-3471-000
3	Set Screw	9211-0830-053
4	Nipple, Panel-Mount, Auxiliary O <sub>2</sub> Outlet	1006-5177-000
5	Label, blank (if no Auxiliary O <sub>2</sub> )	1009-3243-000
6*	Nut, M12x1.75, SST	0144-3132-140
7**	Flowmeter Fitting, 1/8 NPTM straight adapter Flowmeter Fitting, 1/8 NPTM elbow adapter	0204-8877-300 0204-8788-300
8**	Flowmeter Fitting Assembly, 6-mm Tubing Adapter	1006-8423-000
9	Plate, Flowmeter Mounting	1009-3126-000
10	Screw, 10-32 x 3/8	0140-6631-107
11	Screw, M4x10 self-tapping	1009-5534-000
12	Tubing (low-pressure) 250 mm - 1/4 inch	1605-1001-000
13	Plate, blank (if no Auxiliary O <sub>2</sub> )	1009-3128-000

<sup>\*</sup> Apply Loctite 242.
\*\* Apply Teflon tape.

## **10.25** Rear panel components




Item	Description	Stock Number
1	Cover, rear upper	1009-3073-000
2	Cap, hose reel	1009-3075-000
3	Screw, M5.5x20	1009-3384-000
4	Strap, hook/loop	1009-3233-000
5	Screw, M6x1.0 captive	1009-3114-000
6	Door, access (not functional for Aespire)	1009-3074-000
7	Screw, M4x12	1009-3109-000
8	Spring, cantilever	1009-3124-000
9	Screw, M3x8	0142-4254-106
10	Cover, trap bottle (if no internal suction)	1009-3173-000
11	Cover, regulator yoke (if no regulator)	1009-3121-000
12	Plate, clip cover	1009-3185-000
13	Wrench, DIN cylinder (without cable)	1202-3651-000
14	Wrench, pin index cylinder (with cable)	0219-3415-800
15	Cable	1010-3049-000
16	Ferrule, cylinder wrench cable retainer	1001-3708-000
17*	Handle, P-grip	1009-3343-000

<sup>\*</sup> Clean mounting surface with isopropyl alcohol.

10-56 09/07 1009-0357-000


## **10.26 Tabletop components**



Item	Description	Stock Number
1	Tabletop, work surface	1009-3029-000
2	Screw, relieved M4x12	1504-3001-000
3	Washer, retainer	1009-3178-000
4a	Screw, relieved M4x16	1011-3980-000
4b	Washer, M5 flat	1006-1459-000
5	Clip (used with bag arm)	1009-3142-000
6	Clip (used with bag on tube)	1009-3139-000
7	Window, check-valve	1009-3088-000
8	Palnut	1009-3090-000
9	Hook, breathing circuit	1009-3086-000
10	Bolt, shoulder	1009-3172-000
11	Washer, wave	1009-3035-000
12	Washer, Nylon	1009-3150-000

# **10.27 Right-side Components**






ı	Item	Description	Stock Number
	1	Extrusion cover	1009-3021-000
:	2	Screw, M6x20	0144-2131-921
;	3	Lockwasher, M6 internal	0144-1118-130
	4	Dovetail, RH upright	1009-3129-000
ļ	5	Screw, self tapping	1009-5534-000
(	6	Cover, pipeline inlet	1009-3091-000
	7	Screw, M4x8	1006-3178-000


10-58 09/07 1009-0357-000

## **10.28 External components - lower assembly**



Item	Description	Stock Number
1	Panel, access	1009-3059-000
2	Panel, service	1009-3141-000
3	Screw, M4x8	1006-3178-000
4	Thumbscrew	1406-3304-000
5	Ring, retaining	1406-3319-000
6	Cover, scavenger reservoir	1009-3027-000
7	Bracket, suction reservoir	1009-3107-000
8	Screw, M4x16	9211-0440-163
9	Lockwasher, M4 external	9213-0540-003
10	Clip, suction bag hose	1407-3327-000
11	Screw, M5x16 PAN HD	9211-8350-163
12	Lockwasher, M5 external;	0144-1118-220
13	Bumper, absorber	1009-3105-000
14	Label, CO <sub>2</sub> canister	1011-3946-000

### **10.29 Drawer**



Item	Description	Stock Number
1	Slide, drawer	1009-3084-000
2	Screw, M4x8 Nyloc	1009-3183-000
3	Drawer, body	1009-3078-000
4	Drawer Front, lower (with down arrow on back)	1009-3032-000
5	Drawer Front, upper (with up arrow on back)	1009-3031-000
6	Screw, M4x12	1009-3109-000

10-60 09/07 1009-0357-000

## 10.30 Legris quick-release fittings

Item	Description	Stock Number
1	Tees – (tube/tube/tube) $4 \text{ mm (N}_2\text{O})$ $6 \text{ mm (O}_2\text{)}$ $8 \text{ mm (Air)}$ $8 \text{ mm/6 mm/8 mm (SCGO pilot)}$ $3/16 \text{ inch (CO}_2 \text{ and Heliox)}$	1202-3653-000 1006-3544-000 1006-3545-000 1009-3297-000 0213-4727-300
2	<b>Tees — (tube/tube/standpipe)</b> 6 mm (O <sub>2</sub> ) 8 mm (Air - Drive gas)	1006-3862-000 1009-3370-000
3	Elbow — (tube/standpipe) 4 mm (N <sub>2</sub> 0) 6 mm (O <sub>2</sub> ) 8 mm (Air) 1/4 inch (mixed gas) 1/4 inch (45° - mixed gas)	1006-3533-000 1006-3534-000 1006-3535-000 1006-3737-000 1009-3368-000
4	Elbow – (tube/tube) 1/4 inch (mixed gas) $4 \text{ mm } (N_20)$ $6 \text{ mm } (O_2)$	1202-3804-000 1009-3040-000 1009-3041=000
5	Y 6 mm (O <sub>2</sub> ) 8 mm (Air) 8 mm Y with tailpiece 1/4 inch (mixed gas)	1009-3043-000 1009-3044-000 1009-3360-000 1006-3065-000
6	Plug $4 \text{ mm } (N_2 0)$ $6 \text{ mm } (O_2)$ $8 \text{ mm } (Air)$ $3/16 \text{ inch } (CO_2 \text{ and Heliox})$	1006-3530-000 1006-3531-000 1006-3532-000 1006-3835-000

**Note:** Not every fitting is used in all machines.


### 10.31 Vent Drive and low-pressure tubing

Item	Description		Length — Size	Stock Number
1	Coupler, female - blac	k		1503-3128-000
2	Coupler, male - black			1503-3237-000
3	Coupler, female - white			1503-3119-000
4	Coupler, male - white			1503-3236-000
5	Coupler, female - yello	W		1503-3132-000
6	Coupler, male - yellow			1407-3330-000
7	Tee (male barb)			1009-3011-000
8	Cap, plug			1406-3524-000
9	Fitting, coupler barb e	nds		1009-3077-000
10	Plug, 4-mm			1006-3530-000
11	Tee (8mm/6mm/8mm	1)		1009-3297-000
12	Tubing (silicone)		72 mm - 3/8 inch	1009-3164-000
13	Tubing (silicone)		42 mm - 3/8 inch	1009-3164-000
14	Tubing (silicone)		100 mm - 3/8 inch	1009-3164-000
15	Tubing (silicone)		110 mm - 3/8 inch	1009-3164-000
	Tube Markings (facto	ry build only)	Length — Size	
120	Aux O2 OUT	(low-pressure)	250 mm - 1/4 inch	1605-1001-000
127	RGM return	(low-pressure)	750 mm - 1/4 inch	1605-1001-000
128	unmarked	(low-pressure)	300 mm - 1/4 inch	1605-1001-000
129	unmarked	(low-pressure)	, 151 mm - 1/4 inch	1605-1001-000
130	AGSS flowtube	(low-pressure)	750 mm - 1/4 inch	1605-1001-000
134	unmarked	(low-pressure)	25 mm - 1/4 inch	1605-1001-000
135	unmarked	(low-pressure)	50 mm - 1/4 inch	1605-1001-000
136	RGM to Scavenge	(low-pressure)	200 mm - 1/4 inch	1605-1001-000
137	RGM to Circuit	(low-pressure)	300 mm - 1/4 inch	1605-1001-000
144	Venturi Pilot	(black)	330 mm - 4 mm	1009-3363-000
145	Venturi Drive	(black)	300 mm - 8 mm	1009-3296-000
146	unmarked	Tygon	260 mm - 1/2 inch	6700-0005-300
148	unmarked	Tygon	465 mm - 1/2 inch	6700-0005-300
149	unmarked		40 mm - 8 mm	1001-3063-000
150	unmarked	Tygon	180 mm - 1/2 inch	6700-0005-300
225	VENT DRIVE	(black)	440 mm - 8 mm	1009-3296-000
227	SCGO PILOT	(black)	320 mm - 6 mm	1009-3295-000
230	MANIFOLD PRESS	(low-pressure)	300 mm - 1/4 inch	1605-1001-000
240	PAW	(low-pressure)	500 mm - 1/4 inch	1605-1001-000

<sup>\*</sup> Refer to section 10.32

10-62 09/07 1009-0357-000

<sup>\*\*</sup> Sample gas return is directed to the scavenging system as a factory default. A qualified service representative can reroute the sample gas back to the breathing system (refer to Section 9.26).



### 10.32 Tubing for use with Legris fittings


Except for the Tygon tubing (Items 147 and 150), this tubing is a flexible, Nylon-type tubing for use with quick-release fittings.

Item	Description	Length — Size	Stock Number
	Tube Markings (factory build only)		
206	N20 PLINE	230 mm - 4 mm	1001-3060-000
208*	N20 CYL	330 mm - 4 mm	1001-3060-000
209	N20 PLINE - MIXER	430 mm - 4 mm	1001-3060-000
114	REGULATED 02	400 mm - 6 mm	1001-3062-000
116	AUX 02	250 mm - 6 mm	1001-3062-000
117	02 PLINE - REG IN	330 mm - 6 mm	1001-3062-000
123**	unmarked	175 mm - 6 mm	1001-3062-000
213	REGULATED 02	560 mm - 6 mm	1001-3062-000
214	REGULATED 02	270 mm - 6 mm	1001-3062-000
215	REGULATED 02	300 mm - 6 mm	1001-3062-000
217	02 PLINE - MIXER	580 mm - 6 mm	1001-3062-000
219	02 CYL - 02 PLINE	215 mm - 6 mm	1001-3062-000
112	AIR CYL- AIR PLINE	270 mm - 8 mm	1001-3063-000
211	AIR PLINE - MIXER	460 mm - 8 mm	1001-3063-000
104	VAP OUT- ACGO	840 mm - 1/4 inch	1001-3064-000
147	unmarked (Tygon)	290 mm - 1/2 inch	6700-0005-300
150	unmarked (Tygon)	180 mm - 1/2 inch	6700-0005-300
203	FLUSH VLV-SCGO	110 mm - 1/4 inch	1001-3064-000
204	FLUSH VLV-ACGO	125 mm - 1/4 inch	1001-3064-000
205	MIXER - VAP IN	600 mm - 1/4 inch	1001-3064-000

 $<sup>^{\</sup>ast}$  With an N $_2$ O cylinder supply, Item 208 connects the N $_2$ O cylinder supply to the Tee connector between Items 206 and 209.

10-64 09/07 1009-0357-000


<sup>\*\*</sup> With two inboard  $O_2$  cylinder supplies, Item 123 connects the second  $O_2$  cylinder supply to the first  $O_2$  cylinder supply.



## **10.33** Cables and harnesses (power supply)

Item	Description	Stock Number
1	Power Cord	Refer to section 10.5
2	Harness, 100/120 V to Toroid Harness, 220/240 V to Toroid	1009-5752-000 1009-5753-000
3	Harness, to 100/120 V outlets Harness, to 220/240 V outlets	1009-5716-000 1009-5717-000
4	Harness, fuse block to AC line filter	1009-5754-000
5	Harness, AC line filter to Power Controller Board	1009-5751-000
6	Harness, AC line filter to Power Supply (U-Frame)	M1052276

10-66 09/07 1009-0357-000




### 10.34 Cables and harnesses in lower electronic enclosure

#### 10.34.1 Machines with original Power Controller board

Item	Description	Stock Number
1	Cable, to MGAS (Airway module) power supply	1009-5555-000
2	Cable, to Display Unit system power interface	1009-5571-000
3	Cable, to Display Unit system signal interface	1009-5572-000
4	Harness, J3-PCB to J5-DCB	1009-5552-000
5	Cable, ribbon J2-ACB to J9-DCB	1009-5561-000
6	Harness, J7-ACB to J6-DCB	1009-5556-000
7	Cable, ribbon J1-ACB to underside of Pan Connector Board	1009-5549-000
8	Harness, J3-ACB to underside of Pan Connector Board	1009-5560-000
9	Harness, J4-ACB to J4-PCB	1009-5551-000
10	Harness, battery	1009-5557-000


10-68 09/07 1009-0357-000



#### 10.34.2 Machines with Universal Power Supply (U-Frame)

Item	Description	Stock Number
1	Cable, to MGAS (Airway module) power supply	1009-5555-000
2	Cable, to Display Unit system power interface	1009-5571-000
3	Cable, to Display Unit system signal interface	1009-5572-000
4	Harness, J3-PCB to J5-DCB	1009-5552-000
5	Cable, ribbon J2-ACB to J9-DCB	1009-5561-000
6	Harness, J7-ACB to J6-DCB	1009-5556-000
7	Cable, ribbon J1-ACB to underside of Pan Connector Board	1009-5549-000
8	Harness, J3-ACB to underside of Pan Connector Board	1009-5560-000
9	Harness, J4-ACB to J4-PCB	1009-5551-000
10	Harness, (flex-cable) battery	M1049280
11	Harness, Power Supply to PCB	M1049276


10-70 09/07 1009-0357-000



# 10.35 Cables and harnesses (Filter Board interface)

Item	Description	Stock Number
1	Harness, On/Standby system switch	1009-5542-000
2	LED assembly, mains green	1009-5514-000
3	Harness, Filter Board to On/Standby switch and LED	1009-5538-000
4	Harness, Filter Board to SCGO/ACGO	1009-5528-000
5	Harness, Filter Board to ABS flow sensors (includes tubing)	1009-8223-000
6	Cable, Filter Board to Vent Engine harness connector	1009-5521-000
7	Harness, Filter Board to $\mathrm{O}_2$ Cell and ABS switches	1009-5531-000
8	Harness, Filter Board to Task Lights	1009-5533-000
9	Harness, Filter Board to Alternate O <sub>2</sub> switch	1009-5532-000
10	Switch, Alt O <sub>2</sub> (includes harness)	1009-5517-000
11	Harness, Task Light	1009-5853-000
12	Harness, Lower Task Light	1009-5854-000
13	Harness, Vent Engine Board	1009-5545-000
14	Harness, Bag/Vent switch to Filter Board harness	1009-5585-000
15	Harness, Canister Release switch (CO <sub>2</sub> Bypass)	1407-3144-000
16	Harness, O <sub>2</sub> Cell to Filter Board harness	1009-5586-000
17	Transducer, pipeline pressure (includes cable)	1011-3000-000
18	Transducer, cylinder pressure (includes cable)	1011-3001-000
19	Harness, ACGO switch to Filter Board harness	1009-5872-000


10-72 09/07 1009-0357-000




### 10.36 Cables and harnesses in Pan enclosure

Item	Description	Stock Number
1	Cable, ribbon J1-ACB to underside of Pan Connector board	1009-5549-000
2	Harness, J3-ACB to underside of Pan Connector board	1009-5560-000
3	Harness, Pan Connector board to Filter board, 24 POSN	1009-5544-000
4	Harness, Pan Connector board to Filter board, 12 POSN	1009-5543-000
5	Harness, Filter board to VIB	1009-5546-000
6	Harness, Pan Connector board to VIB	1009-5547-000
7	Cable, ribbon, Pan Connector board to Mixer	1009-5550-000

10-74 09/07 1009-0357-000




# **10.37 Optional Monitor Display mounts**





Item	Description	Stock Number
1	Arm (long), pivoting display mount	1009-3262-000
2	Cable Guide Kit, for long arm	1009-8473-000
3	Extrusion, upper dovetail	1009-3113-000
4	Screw, M6x20	0144-2131-925

10-76 09/07 1009-0357-000

## 10.38 Display arm mounting kits for optional equipment







Item	Description	Stock Number
1	Cardiocap 5 mount	1009-3265-000
2	S/5 Flat Panel mount	1009-3266-000
3	Spacelabs Flat Panel mount	1009-3267-000
4	Spacelabs PC Scout mount	1009-3268-000

Notes

10-78 09/07 1009-0357-000

# **11 Schematics and Diagrams**

	Circuit boards are available only as complete assemblies.
Figure 11-1	System circuit diagram
Figure 11-2	Gas scavenging circuits
Figure 11-3	Pneumatic circuit diagram
Figure 11-4	Cabling block diagram
Figure 11-5	Cabling block diagram (U-Frame and HPDU)
Figure 11-6	System block diagram - original DU and Power Controller (sheet 1a of 2)
Figure 11-7	System block diagram - with HPDU and U-Frame PCB (sheet 1b of 2)
Figure 11-8	System block diagram (sheet 2 of 2)
Figure 11-9	Wiring harnesses (1 of 2)
Figure 11-10	Wiring harnesses (2 of 2)
Figure 11-11	Electrical cabling block diagram (Pan Enclosure)
Figure 11-12	Electrical cabling (Lower Electrical Enclosure)
Figure 11-13	Tubing
Figure 11-14	Schematic, AC Inlet module; 100-120 V (with isolated outlets)
Figure 11-15	Schematic. AC Inlet module: 100–120 V (no outlets)

In this section Schematics are subject to change without notice.

1009-0357-000 09/07





Figure 11-1 • System circuit diagram

09/07 1009-0357-000



#### **Key to Symbols**

 $\times$  = Plugged port (1/8 inch) for sample gas return.

Plugged port (30 mm) for auxiliary breathing system scavenging.

O = Open port (30 mm) for auxiliary breathing system scavenging.

\* = Zero to 10 l/min drive gas; zero to 10 l/min patient and fresh gas; zero to 20 l/min total typical flow.

**Note**: Active AGSS systems with a 12.7 mm connector do not include the Flow Orifice and the Flow Indicator.

Figure 11-2 • Gas scavenging circuits

1009-0357-000 09/07

#### Key to Numbered Components

- 1. Pipeline inlet
- 2. Pipeline pressure transducer
- 3. High-pressure relief valve (758 kPa / 110 psi)\*
- 4. Supply connections for the ventilator and pilot pressure for SCGO
  - a.  $O_2$  drive gas
  - b. Air drive gas
- 5. Venturi suction supply connection
  - a.  $O_2$  drive gas
  - b. Air drive gas
- 6. Cylinder inlet
- 7. Cylinder pressure transducer
- 8. Primary regulator (cylinder pressure)
- 9. Test port (primary regulator)
- 10. System switch
- 11. Selector valve
  - $a = O_2$ ; b = Air;  $c = N_2O$
- 12. Flow controller
  - $a = O_2$ ; b = balance gas
- 13. Alternate 0<sub>2</sub> disable valve
- 14. Hot-wire anemometer
  - $a = O_2$  flow sensor channel; b = balance gas flow sensor channel
- 15. Vaporizer port valve
- 16. Vaporizer
- 17. Low-pressure relief valve (38 kPa / 5.5 psi)\*
- 18. O<sub>2</sub> flush and auxiliary flowmeter regulator (241 kPa / 35 psi)\*
- 19. 0<sub>2</sub> Flush valve
- 20. Pressure switch (used with the ventilator)
- 21. Breathing system pressure relief valve (SCGO only  $-150 \text{ cmH}_2\text{O}$ )\*
- 22. To Port 3 of ABS interface (circle)
- 23. For SCGO, to Port 2 of ABS interface (non-circle Inspiratory port) For ACGO, to external 22-mm ACGO connector
- 24. Auxiliary O<sub>2</sub> flowmeter (optional)
- \* Approximate values

#### **Key to Symbols**

- ├─ ├── Pneumatic Connection
- ← Filter
- Direction of Flow
- ♦ Check Valve

Avance




Figure 11-3 • Pneumatic circuit diagram

11-4 09/07 1009-0357-000




Figure 11-4 • Cabling block diagram

1009-0357-000 09/07




Figure 11-5 • Cabling block diagram (U-Frame and HPDU)

11-6 09/07 1009-0357-000




Figure 11-6 • System block diagram - original DU and Power Controller (sheet 1a of 2)





Figure 11-7 • System block diagram - with HPDU and U-Frame PCB (sheet 1b of 2)

11-8 09/07 1009-0357-000







Figure 11-8 • System block diagram (sheet 2 of 2)



ABS/FB = ABS/Filter Board

PCB = Power Controller Board

SCGO = Switched Common Gas Outlet


VECB = Vent Engine Connector Board

PpIn = Pressure Transducer Pipeline

Pcyl = Pressure Transducer Cylinder

Figure 11-9 • Wiring harnesses (1 of 2)

09/07 1009-0357-000



ABS/FB = ABS/Filter Board
PCB = Power Controller Board
SCGO = Switched Common Gas Outlet
VECB = Vent Engine Connector Board

PpIn = Pressure Transducer Pipeline

Pcyl = Pressure Transducer Cylinder

----AC Inlet Fuses 0

Figure 11-10 • Wiring harnesses (2 of 2)

1009-0357-000 09/07

ABS/FB = ABS/Filter Board

ACB = Anesthesia Control board

DCB = Display Connector board

DU = Display Unit

MGAS = Compact Airway Module

MGAS/PS = MGAS Power Supply

PanCB = Pan Connector board

PCB = Power Controller board

VIB = Ventilator Interface Board





Figure 11-11 • Electrical cabling block diagram (Pan Enclosure)

11-12 09/07 1009-0357-000



ACB = Anesthesia Control board

DCB = Display Connector board

DU = Display Unit

HPDU = High Performance Display Unit

MGAS = Compact Airway Module


MGAS/PS = MGAS Power Supply

PanCB = Pan Connector board

PCB = Power Controller board

U-Frame = Universal Power Supply

Figure 11-12 • Electrical cabling (Lower Electrical Enclosure)



1009-0357-000 09/07




Figure 11-13 • Tubing

11-14 09/07 1009-0357-000




Figure 11-14 • Schematic, AC Inlet module; 100–120 V (with isolated outlets)

1009-0357-000 09/07




Figure 11-15 • Schematic, AC Inlet module; 100–120 V (no outlets)

11-16 09/07 1009-0357-000




Figure 11-16 • Schematic, AC Inlet module; 220–240 V (with isolated outlets)

1009-0357-000 09/07




Figure 11-17 • Schematic, AC Inlet module; 220–240 V (no outlets)

11-18 09/07 1009-0357-000

# **12 Service Application**

In this section	12.1 Avance Service Application (PC based)	12-2
	12.1.1 PC Requirements	12-2
	12.2 Startup screen – System Status	
	12.3 System Schematics	12-4
	12.3.1 Power Schematic	12-4
	12.3.2 Gas Delivery Schematic	
	12.3.3 Vent Schematic	
	12.4 Menu Items	12-7
	12.5 File menu	12-8
	12.5.1 File — Preferences	12-8
	12.6 Tools menu	12-9
	12.6.1 Tools — Communication Status	12-9
	12.6.2 Tools — System Calibrations	12-10
	12.6.3 Tools — Transfer Logs	12-11
	12.7 Power Diagnostics menu	12-12
	12.7.1 Power Diagnostics — Power Board	12-12
	12.7.2 Power Diagnostics — Anesthesia Control Board Power	12-13
	12.7.3 Power Diagnostics — Mixer Board Power	12-14
	12.7.4 Power Diagnostics — Vent Interface Board Power	12-15
	12.7.5 Power Diagnostics — Display Unit Power	12-16
	12.8 Gas Delivery Subsystem menu	12-17
	12.8.1 Gas Delivery Subsystem — Gas Supply Status	12-17
	12.8.2 Gas Delivery Subsystem — Mixer Output	12-18
	12.8.3 Gas Delivery Subsystem — Mixer Pressure and Temperature $\ldots$	12-19
	12.8.4 Gas Delivery Subsystem — Gas Delivery Status	12-20
	$12.8.5~\text{Gas Delivery Subsystem} - \text{Mixer Post/Checkout Test Results} \dots.$	12-21
	12.8.6 Gas Delivery Subsystem — Perform Mixer Checkout Tests	
	12.9 Vent Subsystem menu	12-24
	12.9.1 Vent Subsystem — Vent Status	12-24
	12.9.2 Vent Subsystem — Vent Flow and Pressure	12-25
	12.10 Window menu	12-26
	12.11 Help menu	12-26

# 12.1 Avance Service Application (PC based)

This section documents the Avance Service Application that runs on a Windows based computer and communicates with the high performance display unit (HPDU). It is compatible with Avance system software 5.0 or greater.

To enable communication with the Service Application, the Avance system must be in the Install/Service mode (or in a failed state provided the display unit is able to communicate).

The application can be used to diagnose electronically detectable failures in an Avance system.

#### **Note**

This program is for machine diagnosis. It cannot be used for machine checkout or acceptance tests.

# 12.1.1 PC Requirements

Minimum requirements to run the application include:

- Personal computer using a Pentium 600 or higher microprocessor
- Windows 2000/XP
- 1024 by 768 resolution (or higher) video adapter
- Minimum of 128 MB of RAM, 256 MB recommended
- About 150 MB free hard disk space
- Microsoft-compatible mouse or equivalent device
- Serial Port or USB port with an RS-232 adapter

The PC used should meet the GE laptop standards. Screen resolution should be set to 96 dpi.

#### **Port Setup**

The Service Application communicates with the system through the serial port on the display unit (refer to section 2.4).

By default, the communication is channeled through the COM1 port on the PC.

**Note**: Ensure that no other application (such as PDA hot-sync) is using this port while the Service Application is running.

If other ports are available, you can select an available port in "Preferences" on the File menu (Section 12.5.1).

#### **Startup Screen**


Launching the Service Application opens the System Status screen. The startup screen establishes proper communications with the system.

12-2 09/07 1009-0357-000

### 12.2 Startup screen — System Status

If proper communication is established, the System Status screen displays the software and hardware revisions of various subsystems in the tested machine.

If there is a communication problem with the system, an error message is displayed and the connection is not completed.



A series of indicator lights at the bottom of the screen give an overall assessment of the system.

- a green bullet indicates proper operation of the subsystem.
- a red X indicates a failed condition in the subsystem.


#### **Schematics**

The System Status screen includes three Softkeys that provide direct access to the schematical representation of a subsystem as shown in the following sections.

- The Power Schematic indicates the condition of power supplies throughout the system (Section 12.3.1).
- The Gas Delivery Schematic (Section 12.3.2)
   and the Vent Schematic (Section 12.3.3)
   include control devices that allow direct manipulation of components in the
   subsystem and displays the resulting output values for select downstream
   components.


# **12.3 System Schematics**

### 12.3.1 Power Schematic



12-4 09/07 1009-0357-000


### 12.3.2 Gas Delivery Schematic



The Softkeys at the bottom of the Gas Delivery schematic bring up related diagnostic screens that are also accessible from the Gas Delivery Subsystem menu.

- Gas Supplies (Section 12.8.1)
- Mixer Output (Section 12.8.2)
- Mixer Press/Temp (Section 12.8.5)

#### 12.3.3 Vent Schematic



The Softkeys at the top and bottom of the Vent schematic bring up related diagnostic screens that are also accessible from the Vent Subsystem menu.

- Vent Flow/Press (Section 12.9.2)
- **Vent Status** (Section 12.9.1)

12-6 09/07 1009-0357-000

#### 12.4 Menu Items

In addition to the schematic representations, the Service Application provides access to diagnostic screens through the following menu item structure.

File

• Preferences (Section 12.5.1)

(Section 12.5)

Exit

**Tools** 

Communication Status (Section 12.6.1)

(Section 12.6)

• System Calibrations (Section 12.6.2)

• Transfer Logs (Section 12.6.3)

**Power Diagnostics** 

• Power Schematic (Section 12.3)

(Section 12.7)

• Power Board (Section 12.7.1)

Anesthesia Control Board Power (Section 12.7.2)

• Mixer Board Power (Section 12.7.3)

Vent Interface Board Power (Section 12.7.4)

• Display Unit Power (Section 12.7.5)

**Gas Delivery Subsystem** 

• Gas Delivery Schematic (Section 12.3.2)

(Section 12.8)

• Gas Supply Status (Section 12.8.1)

Mixer Output (Section 12.8.2)

Mixer Pressure and Temperature (Section 12.8.3)

• Gas Delivery Status (Section 12.8.4)

Mixer Post/Checkout Test Results (Section 12.8.5)

Perform Mixer Tests (Section 12.8.6)

**Vent Subsystem** 

Vent Schematic (Section 12.3.3)

(Section 12.9)

Vent Status (Section 12.9.1)

• Vent Flow and Pressure (Section 12.9.2)

Window

• standard "Window" manipulation items

(Section 12.10)

• (list of all open windows)

Help

About

(Section 12.11)

### 12.5 File menu

The File menu includes the following menu items:

- Preferences (Section 12.5.1)
- Exit (quits the application)

#### 12.5.1 File — Preferences

Selections on this screen affect the format of applicable values displayed on several of the diagnostic screens.

Unit Selection	
Label	Value
Gas Supply Pressure	kPa psi bar
Airway Pressure Units	cmH2O kPa hPa mmHg mbar
Temperature Units	Celsius Fahrenheit
Gas Color Code	ANSI ISO

#### **Serial Port Selection**

If the COM1 port is in use by another application, use the **COM Port** drop-down list on the **Serial Port Selection** screen to select an alternate COM port to use with the Service Application.

Port Selection	
COM Port:	COM1
	COM6

12-8 09/07 1009-0357-000

### 12.6 Tools menu

The Tools menu includes the following menu items:

- Communication Status (Section 12.6.1)
- System Calibration (Section 12.6.2)
- Transfer Logs (Section 12.6.3)

### 12.6.1 Tools — Communication Status

Communications Status	
Label	Value
Anesthesia Control Board	OK, Fail,
Mixer Board	OK, Fail,
Ventilator Interface Board	OK, Fail,
Power Controller Board	OK, Fail,
Controls (Front Panel) Board	OK, Fail,
Gas Module (MGas)	OK, Fail, not detected,
Link Status	No Link


### 12.6.2 Tools — System Calibrations

The System Calibrations screen displays the most recent date that the system passed a User or Service calibration or test.



12-10 09/07 1009-0357-000

### 12.6.3 Tools — Transfer Logs



# 12.7 Power Diagnostics menu

The Power Diagnostics menu includes the following menu items:

- Power Schematic (Section 12.3)
- Power Board (Section 12.7.1)
- Anesthesia Control Board Power (Section 12.7.2)
- Mixer Board Power (Section 12.7.3)
- Vent Interface Board Power (Section 12.7.4)
- Display Unit Power (Section 12.7.5)

## 12.7.1 Power Diagnostics — Power Board

Label	Value Format	Units	Range
12Vdc Supply	Supply XX.XX		11.70 to 12.30
3.3Vdc Supply	X.XXX	Vdc	3.201 to 3.399
1.5 Vdc Supply	X.XXX	Vdc	1.450 to 1.550
Battery Connected	Yes, No		
Calc Battery Time	XX	Min	00 to 30
Battery 1 Volts Battery 2 Volts	XX.X < 6.0 FAIL (red) <10 T Chg (yellow) 10-15.5 (Green)	Vdc	10.00 to 16.50
Battery Current	X.XXX	Α	- 11.000 to 4.000
Battery Status	Fail Bulk Chg Over Chg Float Chg Trickle Chg On Discharge		
Date battery Tested	/		
Last Full Discharge	XX	Min	
Board Temperature	<65C OK (green) >65C <75C Warn (Yellow) >75C Fail (red)	C F	max 64 max 147
Fan Speed	Slow, Fast		
Fan 1 Voltage	XX.XX	Vdc	10.08 - 10.92
Fan 1 Status	Low, High, Fail, OK		
Fan 2 Voltage	XX.XX	Vdc	10.08 - 10.92

<sup>&</sup>lt;\*> Date battery Tested = date of Last Full Discharge (refer to section 6.10).

12-12 09/07 1009-0357-000

# 12.7.2 Power Diagnostics — Anesthesia Control Board Power

Anesthesia Cont	rol Board	Power			
Label	Value	Units	Normal range		
12 Vdc Supply	XX.XX	Vdc	11.90 to 12.90		
ADC Ref	X.XXX	Vdc	4.018 to 4.176		
Label			Value	Label	Value
Gas Select 10VA	Volts		OK, Fail	Periph1 10VA Volts	OK, Fai
Gas Select 10VA	Amps		OK, Fail	Periph1 10VA Amps	OK, Fai
Press Transducer	10VA Amp	S	OK, Fail		
				Periph2 10VA Volts	OK, Fai
Vent Interface Bo	ard 10VA V	olts/	OK, Fail	Periph2 10VA Amps	OK, Fai
Vent Interface Bo	ard 10VA A	Amps	OK, Fail		
Vent Valves 10VA	Volts		OK, Fail		
Vent Valves 10VA	Amps		OK, Fail		
Acces 1 10VA Vol	ts		OK, Fail		
Acces 1 10VA Am	ips		OK, Fail		
Gas Unit 10VA Vo	lts		OK, Fail		
Gas Unit 10VA An	nps		OK, Fail		
Mixer 10VA Volts			OK, Fail		
Mixer 10VA Amps	i		OK, Fail		
Alt 02 10VA Volts	i		OK, Fail		
Alt 02 10VA Amp	S		OK, Fail		

# 12.7.3 Power Diagnostics — Mixer Board Power

Mixer Board Power						
Label	Value Format	Units	Range			
Mixer 10VA Volts	OK, Fail					
12.5 V	XX.X	Vdc	11.8 to 13.0			
5.5V	X.XX	Vdc	5.39 to 5.61			
3.3V CPU	X.XX	Vdc	3.22 to 3.38			
2.5V ADC Ref	X.XX	Vdc	2.47 to 2.53			

12-14 09/07 1009-0357-000

# 12.7.4 Power Diagnostics — Vent Interface Board Power

Vent Interface Board Power					
Label	Value Format	Units	Range		
Vent Board 10VA Volts	OK, Fail				
Vent Valves 10VA Volts	OK, Fail				
Vent Board 12.5V	XX.XX	Vdc	11.30 to 13.13		
Vent Valves 12.5V	XX.XX	Vdc	11.30 to 13.13		
3.2 Vdc (12bit Vref)	X.XXX XXXX	Vdc Counts	3.179 to 3.221		
1.22 Vdc (10bit Vref)	X.XXX XXXX	Vdc Counts	1.074 to 1.367		
+6.0Vdc	X.XX	Vdc	5.51 to 6.50		
-6.0Vdc	-X <b>.</b> XX	Vdc	-6.72 to -5.28		

# 12.7.5 Power Diagnostics — Display Unit Power

### High Performance Display Unit (HPDU)

Display Unit Power						
Label	Value Format	Units	Range			
LCD backlight 12.5V	XX.XX	Vdc	10.50 to 14.00			
Audio 5V	X.XX	Vdc	4.50 to 5.50			
Digital 3.3V	X.XX	Vdc	2.97 to 3.63			
LCD 3.3V	X.XX	Vdc	2.97 to 3.63			
LCD Inverter B output	X.XX	Vdc	0.00 to 1.30			
LCD Inverter A output	XX.XX	Vdc	0.00 to 1.30			

12-16 09/07 1009-0357-000

# 12.8 Gas Delivery Subsystem menu

The Gas Delivery Subsystem menu includes the following menu items:

- Gas Delivery Schematic (Section 12.3.2)
- Gas Supply Status (Section 12.8.1)
- Mixer Output (Section 12.8.2)
- Mixer Pressure and Temperature (Section 12.8.3)
- Gas Delivery Status (Section 12.8.4)
- Mixer Post/Checkout Test Results (Section 12.8.5)
- Perform Mixer Tests (Section 12.8.6)

### 12.8.1 Gas Delivery Subsystem — Gas Supply Status

Gas Supply Status					
Label	Value Format	Range (kPa)	Range (psi)	Range (bar)	
O2 Cylinder 1	XXXXX.XX	0.00 to 27580.00	0.00 to 40000.14	0.00 to 275.80	
O2 Cylinder 2	XXXXX.XX	0.00 to 27580.00	0.00 to 40000.14	0.00 to 275.80	
Air Cylinder	XXXXX.XX	0.00 to 27580.00	0.00 to 40000.14	0.00 to 275.80	
N20 Cylinder	XXXX.XX	0.00 to 9805.00	0.00 to 1422.10	0.00 to 98.05	
O2 Pipeline	XXX.XX	0.00 to 697.00	0.00 to 101.09	0.00 to 6.97	
Air Pipeline	XXX.XX	0.00 to 697.00	0.00 to 101.09	0.00 to 6.97	
N20 Pipeline	XXX.XX	0.00 to 697.00	0.00 to 101.09	0.00 to 6.97	
02 Select Valve	Open, Closed (Open = connected	d to Mixer)			
Air Select Valve	Open, Closed (Open = connected	d to Mixer)			
N20 Select Valve	Open, Closed (Open = connected	d to Mixer)			
Alt O2 Valve	Open, Closed (Open = O2 bypass)				
Alt O2 Button	Not Pressed, Press	sed			
02 Flush	Not Pressed, Press	sed			
Gas Outlet Config.	Circle, SCGO, ACG	0			

# 12.8.2 Gas Delivery Subsystem — Mixer Output

Mixer Output					
Label	Value	Units	Range		
O2 Flow	XX.XX	l/min	0.00, 0.10 to 15.75		
O2 Flow Verify	XX.XX	l/min	2.00 - 18.00		
O2 Flow Signal	X.XXXX	Vdc	0.0986 to 4.0100		
O2 Prop Valve Drive	XXX	mA	0, 29 to 138		
Balance Gas ID	None, Air, N20				
Balance Flow	XX.X	l/min	0.1 - 15 for Air 0.1 - 12 for N20		
Balance Flow Verify	XX.X	l/min	2 - 18 for Air 2 - 14.4 for N2O		
Balance Flow Signal	X.XXX	Vdc	0.050 to 4.045 Vdc		
Balance Prop Valve Drive	XXX	mA	0, 29 to 138		
02 Select Valve	Open, Closed				
Air Select Valve	Open, Closed				
N20 Select Valve	Open, Closed				
ADC Ref Voltage	X.XX	Vdc	2.47 to 2.53		

12-18 09/07 1009-0357-000

# 12.8.3 Gas Delivery Subsystem — Mixer Pressure and Temperature

Mixer Pressure and Temperature						
Label	Value	Units	Range (kpa/C)	Range (psi/F)	Range (bar/C)	
02 Pressure	XX.XX		62.05 to 220.63	9.00 to 32.00	0.62 to 2.21	
02 Pressure Cal	X.XXX	Vdc	0.55 to 2.90	0.55 to 2.90	0.55 to 2.90	
Balance Pressure	XX.XX		62.05 to 220.63	9.00 to 32.00	0.62 to 2.21	
Balance Pres Cal	X.XXX	Vdc	0.55 to 2.90	0.55 to 2.90	0.55 to 2.90	
Mixer Output Pres	XX.XX		62.05 to 199.95	9.00 to 29.00	0.62 to 2.00	
Mixer Output Pres Cal	X.XXX	Vdc	0.55 to 2.90	0.55 to 2.90	0.55 to 2.90	
O2 Temp	XX.X		5.00 to 50.00	41.00 to 122.00	5.00 to 50.00	
02 Temp Volts	X.XXX	Vdc	0.25 to 3.00	0.25 to 3.00	0.25 to 3.00	
Balance Temp	XX.X		5.0 to 50.0	41.0 to 122.0	5.0 to 50.0	
Balance Temp Volts	X.XXX	Vdc	0.25 to 3.00	0.25 to 3.00	0.25 to 3.00	

# 12.8.4 Gas Delivery Subsystem — Gas Delivery Status

Gas Delivery Status				
Label	Value	Units	Range	
O2 Flow	XX.XX	l/min	0.00, 0.10 to 15.75	
Air Flow	XX.XX	l/min	0.00, 0.10 to 15.75	
N20 Flow	XX.XX	I/min	0.00, 0.10 to 12.60	

12-20 09/07 1009-0357-000

## 12.8.5 Gas Delivery Subsystem — Mixer Post/Checkout Test Results

Mixer Post/Checkout Test Results			
Label	Test Results		
Alt 02 Valve Leak	Pass not performed (if O2 Valve Leak fails) Not done. No supply pressure Not done. Selector valve incorrect state Fail. Selector valve leaks Fail. Proportional valve leaks		
02 Valve Leak	Pass Fail. O2 Bypass valve leaks		
Balance Gas Valve Leak	Pass Not done. No supply pressure Not done. Selector valve incorrect state Fail. Selector valve leaks Fail. Proportional valve leaks		
O2 Flow Test	Pass Not done. No supply pressure Not done. Selector valve incorrect state Fail, 3 LPM test failed Fail, 10 LPM test failed		
Balance Flow Test	Pass Not done. No supply pressure Not done. Selector valve incorrect state Fail, 3 LPM test failed Fail, 10 LPM test failed		

#### 12.8.6 Gas Delivery Subsystem — Perform Mixer Checkout Tests

Selecting **Perform Mixer Checkout Tests** brings up the following screen. This screen includes and automatic and a manual leak test of the Mixer's balance gas inlet check valves.



#### **Automatic**

The **BalGas Check Valve Leak Test Auto** does not require disassembly of the system. It is a sensitive test that will Pass check valve that have leak rates within specifications; however, it may Fail some Mixers with marginal but acceptable leak rates.

- Mixers that Pass the auto leak test do not require further testing and can be left in service.
- Mixers that Fail the auto leak test should be further tested using the manual leak test.

#### Manual

The BalGas Check Valve Leak Test Manual requires disassembly of the system.

- Mixers that Pass the manual leak test have acceptable leak rates and can be left in service.
- Mixers that Fail the manual leak test should be replaced.


#### **Automatic leak test**

- 1. Connect an O<sub>2</sub> supply.
- 2. Select Start Test to perform the BalGas Check Valve Leak Test Auto.
  - If Pass, balance gas check valve leak rate is acceptable.
  - If Fail, verify leak rate using the manual leak test.

12-22 09/07 1009-0357-000

#### Manual leak test

Selecting **BalGas Check Valve Leak Test Manual** brings up the manual leak test setup instructions.



Select **Start Test** to perform the manual leak test as follows.

- 1. Slowly pressurize the Mixer to 400 mmHg (over a 5-second period), as read on the test device.
- 2. The pressure shown on the test gauge should not decrese to zero in less than 10 seconds.
- 3. Select 'End Test' when done.

#### After completing the test:

- 1. Remove all test fixtures.
- 2. Re-assemble Mixer pneumtics and remove plugs from Mixer assembly.
- 3. Re-attach Oxygen supply and activate 'Confirm'.

# 12.9 Vent Subsystem menu

The Vent Subsystem menu includes the following menu items:

- Vent Schematic (Section 12.3.3)
- Vent Status (Section 12.9.1)
- Vent Flow and Pressure (Section 12.9.2)

### 12.9.1 Vent Subsystem — Vent Status

	Vent Status		
	Label	Value	
	Vent Drive Gas	Air or O2	
	ABS Installed	Installed or Not Installed	
	Flush Valve	Not Pressed or Pressed	
<*>	CO2 Bypass	Closed or Open	
	O2 Cell Status	Connected or None	
	Bag/Vent Switch	Bag or Vent	
	Circuit Feedback	Circle, Non-circle, or fault	
	Over Pressure Circuit	OK or High-Pressure	
	ACGO/SCGO Config	ACGO or SCGO	
	Gas Inlet Valve Feedback	Open or Closed	

<sup>&</sup>lt;\*> This refers to the Canister Release switch. The value defaults to Closed if the switch kit is not installed (refer to section 10.21).

12-24 09/07 1009-0357-000

# 12.9.2 Vent Subsystem — Vent Flow and Pressure

	Vent Flow and Pressure					
	Menu Item	Value	Units	Range	Counts (0-4095)	
	Inspiratory Flow	XXX.X	I/min	-120.0 to 120.0	XXXX	
	Expiratory Flow	XXX.X	I/min	-120.0 to 120.0	XXXX	
<*>	Airway Pressure	XX X	cmH <sub>2</sub> O	-20.0 to 120.0	XXXX	
<*>	Manifold Pressure	XXX	cmH <sub>2</sub> O	-20.0 to 120.0	XXXX	
	O2 Cell	XXX	%	5 to 110 %	XXXX	
	ADC Ref Voltage	X.XXX	Vdc	3.179 to 3.221		
	Flow Valve Setting	XXX.X	I/min	0.00 to 140.00	XXXX	
	Flow Valve Feedback	XXXX	mV	0 to 4095	XXXX	
	Flow Valve Current	XXXX	mA	0.0 to 102.4	XXXX	

<*>	kPA	-2.0 to 11.8
	hPA	-19.6 to 117.7
	mmHg	-14.7 to 88.3
	mBar	-19.6 to 117.7

### 12.10 Window menu

The Window menu includes the following menu items:

- Cascade
- Tile Horizontal
- Tile Vertical
- Close
- Close All
- (list of all open windows)

# 12.11 Help menu

The Help menu includes the following menu items:

About...

#### **About**



12-26 09/07 1009-0357-000